

Lipo Laser Slimming Machine

Model#WL-UU650

-----Manual-----



Contents

Chapter One	Principle and Application-----	2
	I Principle-----	2
	II Application-----	2
Chapter Two	Safety Warnings-----	3
	I Electricity & Machinery Safety-----	3
	II Contraindication Safety -----	3
	III Treatment Safety -----	3
Chapter Three	Installment and Use-----	4
	I Working Condition-----	4
	II. Technical Parameters -----	4
	III. How to install the accessories-----	4
	V. How to put the pads on Body -----	5
	VI. How to start the machine-----	5
	VII. Program setting-----	6
	VIII. Accessory List -----	6
Chapter Four	Maintenance -----	7
	I General Maintenance-----	7
	II Troubleshooting-----	8

Chapter One Principle and Application

I. Principle

Lipo Laser Machine involves the application of a highly sensitive laser to dissolve the fat deposits situated in the upper layers of the human skin. The energy emanating from the laser breaks up the fat cells to produce an oily material. This substance is either naturally eliminated through normally bodily processes or can be removed via a small incision made by the practising surgeon. This leads to a permanent loss of fat as long as a healthy diet is maintained after treatment.

It penetrates into skin surface and stimulates the fat cell membranes, changing their permeability, the fat cells reduce their overall size and intracellular fat is released then the fatty triglycerides flow out of the disrupted cell membranes and into the interstitial space, where they gradually pass through the body's natural metabolic functions with no harmful physiological effects, this process is resulting in inch loss of patients.

The subtlety of the treatment carries various advantages over the problems often associated with liposuction. Damaged blood vessels and allied physical risks are significantly reduced, and recovery time is condensed. It is also acknowledged that the laser itself is important in encouraging collagen production, thus strengthening the skin.

II. Application

Lipo laser is a non-invasive or minimally invasive procedure, depending on the amount of fat. And it is a method to get rid of the fat deposits that affect the aesthetic of the body contour and cause the cellulite "dimples" in the skin. And The laser method is for mild to moderate fat conditions and skin tightening But it cannot remove larger amounts of dis-proportioned fat.

Chapter Two Safety Warnings

I. Electricity & Machinery Safety

- 1.The machine is either 110V or 220V.So pls check the machine voltage before starting the machine.
- 2.Please Keep the voltage at stable condition.

II. Contraindication Safety

- 1.People with heart disease, High blood pressure, or configured cardiac Pacemaker.
- 2.People with acute inflammation, asthma, deep vein thrombosis, thyroncus, cancer.
- 3.People with hemorrhagic diseases, trauma or who is bleeding.
- 4.Pregnant women
- 5.Medical plastic, artificial or metal parts inside the body such as metal teeth or Silicone breast implants
- 6.People with abnormal immune system
- 7.People with Numb or insensitive to heat.
- 8.If you are feeling unwell, Prohibited in the ears, nose, eye throat.
and do not use this device in advance and tell your doctor.
- 9.Children under 12 years old are not allowed.

III. Treatment Safety

- 1.Treatment Time is less than 30mins/Body Part/Each Treatment.
- 2.The energy setting is from low to high gradually which should be subject to skin comfortableness for different people.
- 3.The patient had better wear the protective eye glasses.
- 4.If you cannot tolerant heat or unfortableness, pls stop the treatment in advance and turn to your doctor for help.

Chapter Three Installment and Use

I. Working Condition

at normal temperatures and pressures. Please do not stay in too cold, too hot, too dry and too wet.

II. Technical Parameters

Power Voltage:100-110V/60HZ,200-240V/50HZ.

Laser wave length: Laser diode AlGalnp 635-650nm

Maximum output power:60W

Power of Each Diode Laser:100mw

Total power of diode laser pads:3400mw

III. How to install the accessories with the machine.



S1,S2 is for two small pads,4pin jet.

B1,B2,B3,B4 is for four large pads,3 pin jet.

Push in the pads and install them. Click the push button and pull them out.

And plug in the power supply and power key.

V. How to put the pads on the body?

Big pads for big body part such as belly, waist and etc.

Small pads for small body part such as face.

And use the belts and adjust it a proper place.

Insert the pads in and close to the skin.



VI. How to Start the machine?

Push the red power button, then turn right the power key or push in the emergency switch. Then you can start the machine.

On the contrary, turn left the power key or push out the emergency switch, pull the red power button out, then you turn off the machine.



VII. Programs Settings.

The screen will show “MYCHWAY”

Then Click the 1th button ”Menu” and enter into main interface and then click “Menu” again, when The screen is flashing and select programs for large pads or small pads or time.

The two programs can work at the same time or seperately. Click the 6th button “ON/OFF” and t normally we suggest that the intensity should be adjusted from low to high gradually.

And recommended intensity level for big pads is 40-60% and 20-30% for small pads.

The system automatic time setting is 20 mins.



VIII. Accessory List

1xPower line(We offer 100-240V,AU/EU/UK/AU plug)

4xBig paddle with 8 diode laser//Each (Total 32 Diode laser)

2xSmall paddle with 1 led diode laser//Each (Total 2 Diode laser)

1xMetal holder for diode paddles

Bangadge x4

Chapter Four Maintenance

I. General Maintenance

1. After the treatment, pls turn off the machine, Plug out the key, Laser pads and the power supply. Use the dry towel and clean the pads and the machine.
2. Put back all the parts in original place, at normal temperatures and pressures. Please do not stay in too cold, too hot, too dry and too wet.
3. Machine needs 30mins break after 2-3hr operating.
4. Everyday we do not think that 9hr Operating.
5. If it is not non-normal use and irresistible factors We offer Two years' warranty for the units include the power, capacitor, touch-screen, and main board and inside the machine. And offer free repairment , one to one exchange or return.
6. The Pads belong to consumative electronic products, we only have three hours warranty and offer one to one exchange or return during the warranty period.
To know more, pls check our warranty card that come with the machine.

II .Troubleshooting

1.How many treatments will i need to see the results?

-You will see the results of the treatment periodically. Initially you may notice some minor changes, But normally We recommend 2 treatments per week for 8 weeks.

2.Can i have more than one area treated at one time?

-Yes if you are prepared to do more cardio vascular exercise.

3.Will the inches stay off after a treatment?

-Yes, if you follow our instruction, maintain a healthy diet and exercise routine.

4.Does the treatment hurt?

-No, it is a non-invasive treatment and here is no pain, bleed at all. But you may feel some heat because when it worked for a while.

5.Will i have to keep coming back to keep my results?

- Generally, one treatment per a single area is usually adequate to achieve the required results. However, We recommend you more than 2 months of treatment for a better results. If the results are unsatisfactory there is the opportunity for repeat treatments. Because it depends on different body.

6.What is the difference between liposuction and lipo laser?

-The procedure is similar to liposuction with some differences that may make Lipo laser a good alternative to traditional liposuction. Lipo laser essentially permanently removes your excess fat however without the suction removal of the fat. The body fat is dissolved or liquefied (turned into liquid) using a medical laser instead of suction.