

Cautions

1. Do not operate for over 30 minutes each time
2. For fastest results begin treatment with the most desired area. Do not treat the whole body in one session to avoid fatigue
3. For better results, avoid treatments during the menstrual period
4. Not suitable for patient's with heart/liver/kidney disease, cancer, varicose veins and pregnant women
5. Hyper sensitive skin .Damaged incised skin. Fresh scars
6. Overly loose skin .Severe ulcer sore
7. Severe blood vessel dilatation
8. Migraine (negative pressure treatment recommended for relaxation)

POWER SUPPLY	110V	DIMENSIONS	35.5*30*36CM
POWER CONSUMPTION	120W	WELGHT	4 KG
FREQUENCY	60Hz	ELECTRODE PLATE	1PCS
FACE PROBE HANDLE	1PCS	ELECTRODES	5PCS

Radio Frequency OPERATION MANUAL

Thank you for using this product

Does please read this manual before use and comply with the relevant matters needing attention

Principle of Habalan RF

RF Deep Pyrogenic Device emits two types of High Frequency Waves (RF: 0.3MHz or 0.5MHz) deep into subcutaneous tissues. Tissue gets vibrating, rubbing against each other, producing biological thermal energy.

This thermal energy smoothly elevates the temperature (42 degree Celsius to 43 degree Celsius) for vasodilatation. It recovers modified cellular functions, promotes blood circulation and normalizes lipolysis and hormone balance. It is also effective in accelerating collagen contraction and regeneration through stimulation of collagen layers.

RF does not stimulate the sensory nerve, motor nerve or causes any muscular shrinkage inside the body but heat up the specific inner part. The excellent human safety and its effectiveness is internationally accepted and has been proven through superior clinical results.

Need to know the RF treatment

For skin and restoration project is divided into three periods

The first stage is the collagen systolic, nursing after 1 ~ 3 days will feel tight due to contraction.


The second stage is collagen update stage, nursing ~ 1 month, 3 days after delicate degree increase, tactility is tight.

The third stage is the collagen elastic stage, nursing, 1 ~ 3 months after the skin has the effect of a tight, and become firm elasticity.

So the effect can last quite some time, if keep gold treatment time, law of intensive use, can let the skin for a long time in the third phase of an organic whole, together to present the best effect.



Panel shows



The following treatment timing is recommended:

Treatment for eye-bag and wrinkles around the eye area	1-2 minutes on each side
Treatment for face area	9-8 minutes on each side
Total time	20 minutes for the whole face care

Results

Wrinkles	Due to accelerate collagen contraction and regeneration, broken part of the skin will be lifted up, thus reducing wrinkles
Drooping Cheeks	The cheeks will be raised due to stimulation of fibroblasts, elevation of tensile strength of deep collagen tissues.
Pimples	The hyperthermia energy will normalize the sebaceous glands which control the sebum secretion and the high frequency will "killed" the bacteria that causes the infections.
Firmer skin	The hyperthermia energy helps to promote blood and lymph circulation as well as dismantle the fat for discharge. One may feels that the face has reduced in size.