



Advantages:

Freezing painless,complete hair removal,colorful touch screen,integrated photon and IPL in one,free to choose functions,intelligent interface,easy to understand,SHR(OTP) professional hair removal beauty instrument technology:optical energy+epidermal cooling.Optical energy plays a guide role in IPL technology,optical energy concentrated effect in target tissue ,it helps strengthen the target tissue absorb optical energy under the condition of low intensity.At the same time,the contact surface cooling technology on the treatment head can eliminate the thermal effect caused by optical energy over-strong,and increase the skin impedance,reduce skin to absorb optical energy,improve the treatment efficacy and safety,and it's also effective for light color target tissue.

TREATMENT ADVANTAGES

A.Import gem crystal guiding light system,the cooling performance is better,more comfortable,and more durable.

B.Israeli semiconductor cooling system gives epidermis skin minus 8 degrees cooling

C. Use American eight times atmospheric pressure light source.

D. Intelligent operation panel and user interface, intelligent software, easy to operate;

E. With OPT technology, the pulsed optical equipment not only can adjust the pulse width and pulse number, but also the pulse interval is adjustable, the most important thing is that each pulse peak value energy is consistent, clinical effect is more outstanding, safety has improved significantly and obviously relieve postoperative pain.

TREATMENT RANGE

Hair Removal: armpit hair, hairline, beard, whiskers, lip hair, body hair bikini hair or any other unwanted hair.

Treatment Advantages

1. Short treatment time: One treatment about 20-30min
2. No Injure: Noe any injure for skin
3. No pain: No anesthesia, just a little burning sensation
4. No follow-up care
5. No any bad effects after finish treatment.
6. Better effect: More than ten times the effect of other methods
7. Safety: OPT is high-tech optical pulse system with perfect combination of human biological engineering technology, changed the way past tender skin, deeply remove color, high safety coefficient.

TECHNICAL SPECIFICATIONS

Energy range: 5-50J/Square centimeter

Treatment area: 10*40mm

Wavelength: 710nm

Pulse duration: 0.1-7.0ms

Frequency: 1-10HZ

Power supply: AC220V 50HZ; AC110V 50HZ

Output power: 2000W

Dimension: 640*400*1000(mm)

Net Weight: 50kg

Atmospheric Pressure: 860hPa~1060hPa

Humidity: $\leq 80\%$

Operation

Open package ,can view the holes which connect the machine with handle like as below picture:



Connect the head with machine:



Note:Pls make sure each holes connect well,after connected can hear a sound "DA",

Then can put the pure water into machine (Do not add alkaline mineral water and no water opening machine)



Take off over flow hole, it will make sure water can be flowed when machine working;


Open Inject hole, fill water into machine ,until water come out from over flow hole.

Next connect the powr, turn on the red button(emergency button),then Red button Switch bounce.Note: after add water, do not use machine immediately,Pls make sure the machine work well,water can be flowed

interface




“Treatment account”: show total treatment times.



: Show the level of energy



: control the level of enegy,range 1-50J/cm



: control system for cooling, total 10 level



: Show whether system is working well,Respectively from left to right: pre-combustion, recharge, the water cycle, water temperature, standby mode, charging precombustion and displayed on the exclamation mark, readiness exclamation mark will disappear quickly, enter the work state, charge when light will show a exclamation mark temporarily, luminous exclamation mark end disappear. If the water cycle and water temperature are abnormal, please check the reason to turn it off immediately, please contact the manufacturer if necessary.



: Standby and prepare for two kinds of state, the default for the standby mode, switch to the ready state when starting treatment, xenon lamp

light at this time, press the handle switch light (at this time to press the lever switch, should pay attention to safety protection).



: RETURN.

SHR Fast hair removal mode interface is only one pulse number ,



Pulse width of light duration, 0.1-7.0 - ms is adjustable, Step 0.1 ms, recommend the use of 2.0 to 3.5



ms, Adjustable frequency of light speed, 1-10, step 1, it is recommended to use 2 or 3; Can be directly click pulse width and frequency, and make it into the editor, Then



through Key directly to adjust, it is recommended that the light energy to 10, and gradually increase according to the circumstance, other parameters meaning and OPT mode.

operative skills:

1, do hair removal, armpit hair and a slight tingling hairline affordable (guest) slightly burnt flavour, onset of hair removal has a slight tingling, 3-4 times can play a unique effect. Because the hair growth have 3 cycle, namely: the growing season - maturity - dormant and only in the growth period of hair can remove effectively, so the hair removal need 3-4 times.

2, sensitive skin, dark skin, skin easily appear delay response in 24 to 48 hours, not because the reaction is not obvious and optional to increase energy.

3, treatment of flare in the operating parts don't stay for long, sliding operation from down to up the same direction. Every time interval of 1 to 3 minutes, the effect is better.

5, treatment should constantly ask the guest when feeling, adjust energy. The treatment is finished for

6, cold gel apply for 20 to 30 minutes, if no discomfort can be washed off. Postoperative taboo:

The fasting, spicy food, seafood hair content, non-ferrous vegetarian foods, smoking, wine;

Do not use irritating skin care products within three days after treatment (such as soap, facial cleanser, avoid is used a function type products.

Pay attention to prevent bask in, use sunscreen (SPF20) above;

Recommend the use of pure natural class after skin care products;

If there is scabby, no peeling off by hand, their own metabolism;

6.10 days with cold water washes a face, it is prohibited to sauna, steam, after 10 days to do facial treatment.