



SR-MJ1705

Instructions

Warning

Dear users:

Thanks for choosing our company's product, in order to fully make use of the products, we sincerely suggest follow:

- ◆◆Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆Please follow the instruction to fit and operate the machine correctly.
- ◆◆Please do not remove and change any accessories of the machine.
- ◆◆It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆It forbids others to change and maintain the machine except the professional people.
- ◆◆Please contact us if there is any trouble during operation.

After Sale service

◆We offer ONE YEAR warranty. Within warranty period, all exchanges or replacement due to defects are free. (Buyer is responsible for shipping charges).

If the machine is over warranty, please also do not hesitate to contact us. we will assure you of our best service all the time.

◆How to contact us?

Phone:(+86) 189 -487-55929

Tel: 0755-86368500

Email: sales@mychway.com

This product does not belong to medical, please don't use for disease treatment!

I. Function Introduction:

1. Ultrasonic

Ultrasonic utilizes sonic vibrations; it works one million times per second. The ultrasonic wave could stimulate skin tissue and improve permeability of the membrane to promote collagen tissue, reduce wrinkles and invisible fine lines, also lift the skin for a younger and healthier look. Ultrasonic generates frictional heat which accelerates blood circulation, muscle relaxation, and decomposes cellulite to break down excess fat.

2. Skin Scrubber

With Ultrasonic Scrubber begins a new generation of aesthetic soft peeling. Using the high performance of ultrasounds, this exciting cosmetic breakthrough provides a safe, painless and perfect physiological skin regeneration. By removing the superficial scaly layers of the skin, the Scrubber stimulates the revival of the skin tissue, rendering it youthful and luminous from the first treatment. The tip of the metallic probe produces extremely rapid moments at a very low amplitude when it is exposed to ultrasound vibration. This ultrasonic emission has a spectacular effect on water molecules, which are pulverized. This phenomenon is known as micronization.

3. High Frequency

Use high-frequency current and infrared ray to diminish inflammation and ease pain. The generated ozone kills bacteria, heals and astringes the acne vulnus. It has an effect on regulating sebum secretion, improving the texture of skin, accelerating blood circulation and speeding up metabolism, anti-wrinkle and anti-aging.

4. Mole Removal Pen

By means of coagulation of high-frequency electric current, it urges the cell necrosis of the epidermis, which can thoroughly eliminate all kinds of spots on the face and the other skin surface.

5. Bio

Skin tightening after cellulite and fat breaking down; Lifting up, for example, faces, breasts, hips and etc.

II. Package list

1. Host 1set
2. BIO 1pc
3. High frequency 1pc
4. high frequency heads 3 pc
5. Ultrasonic probes 2pcs (flange head and point head)
6. Skin Scrubber 1pc
7. Spots removal pen 1pc+10 needles
8. Holder 1pc
9. Power cord 1pc

III. How to install the Accessories




IV. Treatment Introduction

Please press the power switch behind machine, then the screen will show the following menu image:



1. Ultrasonic



1.1 Operation:

1.1.1 please chose the appropriate ultrasonic head, and connect with output socket behind machine

1.1.2 please set up the treatment time, energy and mode. (please adjust the energy from low to high intensity under customers' tolerance.

1.1.3. please press the start button, then the head will work

1.2 Mode button option:

Mode A: continued vibration

Mode B: short pause vibration.

Mode C. A bit long short pause vibration.

1.3 Ultrasonic head treatment area:

1.3.1 flate ultrasonic head: used for face, body, arm, etc

1.3.2 point ultrasonic head: used for eye and mouth around, etc



1.4 Ultrasonic test way:

Please drop a bit water on the ultrasonic head, setting up high energy intensity. Then the water will vibrate.

2. Spot Removal



2.1 Operation:

Please setting up the treatment time and energy, then press start button and put the spot tips on skin.
Please adjust the energy from low to high under customer's tolerance.

3. High Frequency



3.1 Operation

3.1.1 please chose appropriate high frequency tips replacement, and connect with handle piece.

3.1.2 Please setting up the treatment time and energy

3.1.3 Please press start button and put the tips on skin.

3.1.4 Please adjust the energy from low to high under customer's tolerance.

3.2 Usage of high frequency tips:

A. Mushroom Tube for large facial areas

B. Tougue Tube for small facial areas

C. Bend Tube for spot areas



Mushroom Tube




Tougue Tube



Bend Tube

4. BIO



4.1 Operation

4.1.1 After Cleaning Skin, apply essential gel on the skin and massage skin


4.1.2 Please setting up the treatment time, energy and mode

4.1.3 Please press start button and put the tips on skin.

4.1.4 Please adjust the energy from low to high under customer's tolerance.

Note: Press button "MODE", different frequency, and could try one by one, then chose the one option customers like.

4.2 Diagram of suggested movements and areas for BIO:



5. Skin Scrubber



5.1 Operation

5.1.1 After Cleaning Skin, Apply essential gel on the skin and massage parallel with the back of the scrub head. Apply cleansing lotion on the skin for 20 seconds; gently rub the skin with the front of the scrub head at 35 to 40 degree angle

5.1.2 Please setting up the treatment time, energy and mode

5.1.3 Please press start button and put the tips on skin.

5.1.4 Please adjust the energy from low to high under customer's tolerance.

5.2 Mode button option:

Mode A: continued vibration

Mode B: short pause vibration.

Mode C: A bit long short pause vibration.

V. Attention:

Never operate this equipment on a person has or suffers from the following medical conditions described below. (The salon or operator may consider having a client sign a waiver to the effect that he/she is free of any poor health conditions.)

1. Heart Disease
2. Heart Murmur or Irregular Heart Rate
3. Has a Pacemaker or other type of heart rate balancer.
4. Pregnancy
5. Epilepsy
6. Alcoholism
7. Ulcers
8. Disease of Stomach, Intestine, Liver, or Pancreas.
9. Breast Cancer.
10. Undergoing treatment for any form of Cancer.
11. Had major surgery within the past year.
12. Asthma
13. Cold, Flu, or Upset Stomach
14. Any serious illness or disease not listed above.

VI. Do not use on these areas:

1. Sensitive Skin
2. Skin Rash
3. Damaged, Inflamed, or Infected Skin.
4. Surgical scar less than 12 months old.
5. Skin Swollen.
6. The probes can not point to Eye ball.
7. Genitals.