



MYCHUWAY

SR-MZ1805U

Instructions

Warning

Dear users:

Thanks for choosing our company's product, in order to fully make use of the products, we sincerely suggest follow:

- ◆◆ Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆ Please follow the instruction to fit and operate the machine correctly.
- ◆◆ Please do not remove and change any accessories of the machine.
- ◆◆ It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆ It forbids others to change and maintain the machine except the professional people.
- ◆◆ Please contact us if there is any trouble during operation.

After Sale service

◆ We offer ONE YEAR warranty. Within warranty period, all exchanges or replacement due to defects are free. (Buyer is responsible for shipping charges).

If the machine is over warranty, please also do not hesitate to contact us. we will assure you of our best service all the time.

◆ How to contact us?

Phone:(+86) 189 -487-55929

Tel: 0755-86368500

Email: sales@mychway.com

This product does not belong to medical, please don't use for disease treatment!


Table of Contents:

- 1. Introduction to Manual
 - 1.1 Purpose
 - 1.2 Conventions
- 2. Medical Safety.....
 - 2.1. Indications for Use.....
 - 2.2. Contraindications
 - 2.3. Cautions
- 3. Product introduction.....
- 4. Control Unit & Handpiece.....
- 5. Operation Procedures.....
 - 5.1. Ultrasound body/face probe.....
 - 5.2. Skin scrubber.....
 - 5.3. ION+/-.....
 - 5.4. Cold/hot handle.....
 - 5.5. Spray gun
- 6. Illustration of Operation.....

1. Introduction to Manual


1.1 Purpose


This instructions for use manual provides a description of the system components, its controls and displays, instructions for its operation, and other equipment information important to the user


 **Warning:** Do Not operate this machine before reading this manual thoroughly. In addition to this manual,additional clinical training may be available by the company or your local distributor.

For more information on training available please contact your local representative.

1.2 Conventions

 **Note:** Notes designate information of special interest.

 **Caution:** Cautions alert the user to precautionary steps necessary to properly operate the system.Failure to observe these cautions may void the warranty

 **Warning:** Warnings alert the user to information that is of the highest importance and vital to the safety of the patient and user

All procedures are broken down by numbered steps.Steps must be completed in the sequence they are presented

Bulleted lists indicate general information about a particular function or procedure.They do not imply a sequential procedure.

Control names are spelled as they are on the system,and they appear in **Bold text**

1.3 Cautions

Please check all handles connected well.

Please remove all the jewelry before treatment

After treatment please remember to turn off the machine

Please clean the ion handles after treatment

Sensitive skin and Contact dermatitis client dont suggest to use.

Please clean the spray bottle after treatment

2. Medical Safety

2.1. Indications for Use

The 6 in 1 ultrasound skin care machine is indicated for use as a non-invasive dermatological aesthetic treatment to:

- Body slimming,body sculpture
- Face care,lifting,whitening, moisturizing
- Skin deeply clean

2.2. Contraindications

The machine is contraindicated for use in patients with:

- Open wounds or lesions on the face and/or neck
- Severe or cystic acne on the face and/or neck

3. Product introduction

ULTRASOUND(SONO FACE/SONO BODY)

Ultrasonic utilizes sonic vibrations; it works one million times per second. The ultrasonic wave could stimulate skin tissue and improve permeability of the membrane to promoted collagen tissue, reduced wrinkle and invisible appeared fine line, also lift the skin result for younger and healthy look. Ultrasonic generate frictional heat which accelerate blood circulation, muscle relaxation, and decompose the cellulite to bread down exceed fat. This machine utilizes sonic vibrations to reduce lines, leaving the skin firm and smooth.

SKIN SCRUBBER

With Ultrasonic Scrubber begins a new generation of aesthetic soft peeling. Using the high performance of ultrasounds, this exciting cosmetic breakthrough provides a safe painless and perfect physiological skin regeneration. By removing the superficial scaly layers of the skin, the Scrubber stimulates the revival of the skin tissue rendering it youthful and luminous from the first treatment. The top of the metallic probe produces extremely rapid moments at a very low amplitude when it is exposed to ultrasound vibration. This ultrasonic emission has a spectacular effect on water molecules, which are pulverized.

IONTOPHORESIS

- Ion make skin active, help skin absorb nutrition and skin rejuvenation
- +Ion calm down skin, clean skin , shrink pores

COLD/HOT LED PHOTON

Hot temperature can help skin open pores , promote blood circulation, and help skin absorb nutrition

Cold temperature shrink pores, firmly lock moisture

HYDRO OXYGEN JET SPRAY

Nutrition or cosmetic product can be filled into the stainless steel bottle. With high pressure, the nutrition and oxygen are sprayed on the skin area which greatly enhance nutrition to be absorbed maximally, skin cleaning and skin rejuvenation.

4. Control Unit & Handpiece

Please check if all contents are included before you set- up



Accessory



Face ultrasound



Body ultrasound



+ ion -ion



Spray gun



Cold/ hot handle



Skin scrubber

5. OPERATION PROCEDURES:



Please connect all the accessories well and turn on the machine:





5.1 Sono face/body



Press choose  to adjust the working time.

Press  to decrease the energy, and press  to increase the energy.

Press  to select the model. three models A,B,C can be selected.

Press  to start or pause the operation.





Operation guide:


- Apply some ultrasonic gel on the treatment area
- Choose ultrasound working mode, Mode A~C makes the different vibration change at a certain time
- Then choose the energy level usually 3~4 is ok
- Set the working time 12~15min is ok


5.2 Skin scrubber



Press choose  to adjust the working time.

Press  to decrease the energy, and press  to increase the energy.

Press  to select the model. three models A,B,C can be selected.

Press  to start or pause the operation.

Operation guide:



Choose the proper mode: mode A is for scaling
mode B is for permeation
Mode C is for lifting



Set a proper working level 3~4 is ok


Set working time 6~12 min


5.3 Ion+/-



Press choose  to adjust the working time.

Press  to decrease the energy, and press  to increase the energy.

Press  to select the model. three models A,B,C can be selected.

Press  to start or pause the operation.

Operation guide:



Choose the ion working mode: Mode A positive ion
Mode B negative ion
Mode C lifting



Set a proper working level 3~4 is ok


Set working time 6~12 min

5.4 Ice cold/hot handle



Press choose  to adjust the working time.

Press  to decrease the temperature, and press  to increase the temperature.

Press  to start or pause the operation.

Operation guide:



Set a proper temperature for treatment , it has 6 level : -10°C , -5°C , 0°C , 36°C , 38°C , 40°C

Set working time at 6~12min

Hot temperature for open pores, promote blood circulation, better absorption of nutrients

Cold temperature for calm down skin, shrink pores, help lock in moisture

5.5 Spray gun

Please fill some toner or other liquid nutrient to the spray bottle
Press the spray gun start key, it is the red button ,on the back of the machine



Please use your finger to press like below picture, toner is atomized



Spray down from the forehead, please keep the bottle from face at a proper distance
After using, please turn off the machine
Do a hydrating mask after the treatment

6.Illustration of Operation




1.Neck

Move in small circular motions from the inside of the neck to the outside, then lift from the neck upward to the cheek, neck treatment area can be divided into three lines(up,medium,down areas) to operate
Attention: strength should be top heavier and down lighter, operation electrodes should tight skin, repeat 8-10 times




2.Chin

Make operation head move in small circular motions by pressing tightly against the double chin and then pull back to the back of the ear, repeat 5-8 times




Make operation head move in small circular motions by clinging to the chin and circle out to the ear beads. Then pulled down to the chin with light force,repeat 5-8 times



3.Mouth

Make operation head move in small circular motion by clinging to the corner of mouth, and circle out to the ear beads. Then pulled down to the chin with light force,repeat 5-8 times




4.Cheek

Make operation head move in small circular motion by closing to the light bone of nose, circle out to the ear beads. Then pulled down to the chin with light force, repeat 5-8 times



5.Forehead

Make operation head move in small circular motion by closing to the light bone of nose, circle out to the ear beads. Then pulled down to the chin with light force, repeat 5-8 times



6.Eyes' around

Dredge the meridian and collaterals of eyes, repeat 5-8 times




Move in small circular motion in the lower eye then lift to the temple, repeat 5-8 times




Move in small circular motion from the corner of the eye to the Jingming points to the eyebrow, repeat 5-8 times

Belly



1. Move slowly from lower waist to the end of groin, repeat 3 times
2. Move slowly from lower waist to pot belly, along the gap of ribs, repeat 3 times
3. Move slowly from lower waist to belly button line, repeat 3 times




4. Separately move in small circular motions in thickness fat of body, each one in 10 circles
5. Improve anti-cellulite and the growth of collagen, tighten skin



6. From small to large by clockwise direction, circling around the navel, promote intestinal peristalsis and improve constipation

Waist



1. Knead fat back and forth in each side waist(decompose the cellulite effectively and contour waist) repeat 10 times



2. Move 10 circles in renal region
3. Separately move in small circular motions in thickness fat of buttock, each one in 10 circles
4. Move and write "8" in two sides of waist (repeat 10 times)




Lie on your side:

5. Move from the thickness fat of buttock to groin, repeat 8-10 times
6. Move slowly from the thickness fat of buttock near to axillary fossa, repeat 8-10 times

Foreleg



1. Move one time slowly from knee to groin, line by line



2. Move in horizontal from inside of leg to middle, down to up, slower the moving, repeat 2 times
3. Move in horizontal from outside of leg to middle, down to up, slower the moving, repeat 2 times

Caution: Lower the energy 1-3 Level when treating inside of leg



4. Separately move in small circular motions in thickness fat of body, each one in 10 circles (Improve anti-cellulite and the growth of collagen, tight skin)



5. In inside of leg, move along the muscle texture from knee to crotch, repeat back and forth movement in 4 times
6. In outside of leg, move along the muscle texture from knee to crotch, repeat back and forth movement in 4 times



7. Move from knee to groin, repeat 8 times

Backleg



1. Move slowly from the bottom of crus to popliteal fossa, repeat 3 times



2. Separately move in small circular motions in thickness fat of body, each one in 10 circles




3. Push to the up of popliteal fossa from the root of the calf, repeat 5 times

Backleg




1. Push to the bottom of buttocks from up of popliteal fossa slowly line by line, repeat once time



2. Move in horizontal from inside of thigh to middle, down to up, slower the moving, repeat 2 times
 3. Move in horizontal from outside of thigh to middle, down to up, slower the moving, repeat 2 times
- Caution: Lower the energy 1-3 Level when treating inside of thigh



4. Separately move in small circular motions in thickness fat of body, each one in 10 circles (Improve anti-cellulite and the growth of collagen, tight skin)



5. In inside of leg, move along the muscle texture from knee to bottom of buttocks, repeat back and forth movement in 4 times
 6. In outside of leg, move along the muscle texture from knee to bottom of buttocks, repeat back and forth movement in 4 times
- Push to the above of eight-liao point from up of knee fast line by line, repeat 8 times


Back




1. From the Da Zhui point down to governor vessel, eight-liao poin, urinary bladde and up to Da Zhui point then slide to the Da Zhui point from the governor vessel .
Caution: Repeat 3 times movement on urinary bladde of each side after governor vessel



2. From the Da Zhui point down to governor vessel, eight-liao poin, urinary bladde and up to Da Zhui point then slide to the Da Zhui point from the governor vessel .
Caution: Repeat 3 times movement on urinary bladde of each side after governor vessel



3. Move in small circle motion on dazhui point then go to daban point massage on circle motion, repeat 2 times from left to right
Caution: Move 5 times on seams of the scapula)



4. After playing five circle in the cardiopulmonary area, five circles should be played on the cardiopulmonary area


Back



5. Move in small circle motion from down to up, up to down, repeat 6 times, one side finished, the other side start




6. Move and write "8" in two sides of waist,back, repeat 6 times



7. Move from above of buttock to underarm until two sides finished, repeat 3 times


Outside Arm



1. Move from outside of arm to underarm, down to up, slower the moving, repeat 2 times



2. Separately move in small circular motions in thickness fat of body, each one in 10 circles (Improve anti-cellulite and the growth of collagen, tight skin)




3. Move from top of outside of arm to elbow then up to underarm from two sides then down,repeat 6 times

Outside Arm




4. From up of elbow to underarm which near to back



5. Move from inside of elbow to underarm, down to up, fast the moving, repeat 2 times



6. Separately move in small circular motions in thickness fat of body, each one in 10 circles (Improve anti-cellulite and the growth of collagen, tight skin)



7. Move from inside of elbow to underarm, down to up, fast the moving, repeat 2 times