

9in1 Multifunction Dermabrasion Machine
Model:SR-VS908



Catalog

Preface	-----	3
Functions	-----	3
Specifications	-----	3
Operating Instruction	-----	3
Diamond Dermabrasion	-----	3
Ultrasonic	-----	5
Skin Scrubber	-----	6
Hot & Cold treatment	-----	6
Rotary Brush Function	-----	7
Galvanic Function	-----	7
Vac/Spary Function	-----	8
High Frequency Function	-----	8
Safety Warning	-----	8
Caution	-----	9
Trouble Shooting	-----	9

Preface

The Microdermabrasion was designed by Italian Florence's Mattioli at first, until now had over 20 years' history. This kind of technology only supplies to the dermatology expert and doctor at first uses, but more and more SPA House & Beauty Saloon also gradually uses, it was called the best method of medical beauty. This method successfully has developed in Europe for several years, cures patient thousands of examples, and has obtained the extremely satisfactory effect.

Now we want to introduce the Diamond Dermabrasion, it is a creative improvement of the Micro-crystal Dermabrasion. The Diamond Dermabrasion provided a non-surgical skin refinish procedure, by using sterile diamond heads to abrade or rub off the top skin layer, then Vacuuming out the particles along with any dirt and dead skin back up. This procedure removes skin debris, imperfections, blemishes, wrinkles and unwanted pigmentation on the skin. Moderate usage of Diamond Dermabrasion, accompany with skin products that penetrate through top skin layer that reach the dermis level, helps to replenish natural nutrients, restore cell activity, and promote skin health. The results after Micro Dermabrasion treatment is usually a healthy, glowing, beautiful skin.

Functions

1. Diamond Dermabrasion
2. Ultrasonic
3. Skin Scrubber
4. Hot & Cold
5. Rotary Brush
6. High Frequency
7. Vacuum
8. Spray
9. Galvanic

Specifications

Voltage:

110V/60Hz □ 220V/50Hz □

Power: 150W

Operating Instruction

Diamond Dermabrasion

1. Push in the movable ring of copper nozzle at the hand piece, and pull out

the 6x4 size vacuum hose. Then push in the black plastic ring on the vacuum hose.

2. Put a drop of cream or milk lotion on the vacuum hose nozzle, push in the movable ring and insert in the vacuum

hose to the limit.

3. Insert a cotton swab into the Diamond Head, then tight up the Diamond Head on the hand piece. Remove the cotton swab after every treatment.

4. Turn on the power switch.

5. Adjust the Vacuum intensity by Vacuum Regulator according to the following reference table, the intensity displays on Vacuum gauge and check with the client for comfort.

Recommended Vacuum intensity on treatment area

Treatment Area	Normal Skin (inHg)	Sensitive Skin (inHg)
Face	12-22	8-12
Forehead	12-22	10-12
Chest	15-20	15-18
Nose	10-15	10-14
Neck	18-18	10-15
Near Eye Area	8-12	8-10

Treatment Schedule Table


Treatment	Frequency (days)	Total Visits
Face Peeling	6-10	3-5
Fine Lines	7-13	6-8
Neck	12-14	3-5
Acne Scars	7-12	8-10
Hyperpigmentation	7-12	8-10

6. Use the incline area of Hand piece to contact with the skin surface. Proceed with a slow sliding motion. Do not place the hand piece at one spot for too long, doing so will cause harm to the skin.

7. Adjust the crystal flow rate by the Vacuum Regulator slowly from minimum to maximum. For commonly use, do not set the flow rate to maximum unless it is of necessity. This is for better comfort intensity.

8. Gently stretch the skin with thumb and finger and glide the Hand-piece 2 to 3 times over the area being treated. Then cross again in different directions to prevent streaking lines as in the reference picture. Brush strokes or polishing motions may be used in different sites.

DIAGRAM OF SUGGESTED MOVEMENTS AND AREAS FOR DERMABRASION



Ultrasonic function

I. Preparation

- (1). Before start the ultrasonic, cleanse the area to be worked. After cleansed, apply suitable lotion and gel to the skin.
- (2). Press Power On/Off for Ultrasonic to start the function.
- (3). Output intensity display on Intensity window, Press Intensity “up” and “down” button to set output intensity, output intensity classify into 9 levels, from minimum 0 to maximum 1.8 Watts/cm. Each level is 0.2 Watts/cm. For Facial work, low to mid intensity is recommended. For Body work, mid to maximum intensity is recommended.

II. Face Treatment

- (1). Connect Face Electrode (big one) into Face Electrode output.
- (2). Press Intensity “up” and “down” button to set output intensity to your need.
- (3). Press Pulse/Continuing mode to select Pulse or Continuing wave, Pulse/Continuing wave indicator will light on.
- (4). Apply the Electrode working on the face.

III. Eye Treatment

- (1). Connect Eye Electrode (small one) into Electrode output. (Should take out the Face Electrode first)
- (2). Press Intensity “up” and “down” button to set output intensity to your need.
- (3). Press Pulse/Continuing mode to select Pulse or continuing wave, Pulse/Continuing wave indicator will light on.
- (4). Apply the Electrode working around the Eye area.

***** Important:**

1. Clean skin with clear water.
2. Never stop the electrode on the skin in one spot.
3. Always keeps the electrode moving across the skin.
4. Never point the electrode towards chest area.
5. Do not immerse the handle or the machine into water. Doing so could cause serious injury or death.
6. Store in a dry cool place.

Skin Scrubber

Functions and Characters:

Deep clear and clean up aged horniness.

Reduce melanin to get a whiten skin.

In-depth eliminates the wrinkle, strengthens the skin elasticity.

Improve blood circulation and lymph metabolism.

Strengthen the skin cell reproductive property and the permeability.

Accelerate the absorption of nutrition.

Digital control, with double function integration.

No effect and negative for skin.

Skin Scrubber Cautions:

Please use this machine together with clear water, toilet gel or astringent. This may get better effect.

Massage, wrinkles with the transformed waves or the continuous wave.

Please keep the probe and mainframe clean and sterilize before and after using.

It is not allowed to use this machine on area around eyes and wound place.

Children, pregnant women and people with mental disease are not allowed to use it.

The working time can be adjusted according to the request.

It is normal if there is a sound of “zhizhizhi” while the machine works, as its working frequency is something like sound frequency.

If the machine is broken and doesn't work, please contact the manufacturer and agent.

Hot & Cold treatment

1. Connect the hot & cold plug to the machine; please note the socket is in the

“Hot & Cold” Area.

2. Start the machine: To turn on the machine, press the “on” button. And then start the hot & cold treatment function.
3. Please wait 3 minutes, and then you can apply the treatment. For the machine need time to warming up or cool down the treatment probe.
4. Normally one treatment program is about 5 ~10 minutes. Please note that do not stay in one treatment place over time when applying the treatment to avoid any hurt caused by cool head.
5. After treatment, please turn off the machine and then pull out the plug, make clean the probe.

***Please clean the treatment probe before treatment and after treatment.**

Rotary Brush Function

1. Make sure that the power switch is in the “Off” position.
2. Connect the Brush handle wire to Brush Handle.
3. Connect the other end of wire to the machine.
4. Simply pick the desired accessory to be used for the treatment and insert it into the head of the gun making sure that the little pins line up to the slots on the head, allowing the piece to lock into place. Note, that if the accessory is not locked into place properly, it will not spin when any pressure is applied during use.
5. Press the “on/off” button, and the machine attachment on the handle should begin to rotate.
6. To control the speed of the attachment rotation, adjust the knob on the left side of the face plate “Min/Max”. Turning this knob will all the user to control the speed.
7. On the handle, there is a trigger which can be used to stop the attachment from rotating and also to reverse the direction of rotation. Make sure to STOP The rotation before reversing the direction to prevent damage to the handle..

Galvanic Function

8. Connect the Negative & positive (Red) wire into the machine using the plug and fasten the securing ring to hold the plug in place.
9. Connect the Electrode handle to the end of the red wire.
10. Connect Black wire to the machine and connect the straight indifferent electrode to the end of the Black wire for the client to hold.
11. Chose the electrode that best desired treatment and insert it into the Electrode handle and tighten the secure head to prevent dropping the electrode.
12. Select the current you wish to use, “-“or “+”.Do not set the intensity too high to start with.(Notice: The machine has auto protection to shut off the galvanic

function when the output current reached excessive intensity to the treatment area. To restart The function, turn the intensity knob left to the off Position, then turn right to reactivate the function.)

13. To select the desired current, simply press the button on the front panel as drawn on the machine.

Vac/Spary Function

14. Start the machine: To turn on the machine, press the Vac/Spray power switch.

15. Connect the hoses to the machine, black for suction and white for spray.

16. Determine the choice of function and use the accessories accordingly, Glass ventouses for suction and bottles for spray.

17. Depending on the function being used, disconnect the opposite hose for best results.

18. To increase vacuum or spray, turn the intensity control to the right and to decrease the vacuum or spray, turn the control to the left.

19. When the machine is not in use make sure the power switch is off.

High Frequency Function

20 Connect the High Frequency handle's wire to the machine (it only fits one way) making sure to fasten it into place the securing ring.

21. Simply insert the desired glass electrode (metal end first) into the hole at the end of the handle.

22. Using the power switch, turn the machine to the "on" position and the glass electrode should light up.

23. To control the intensity of the electrode, use the knob located at the center of the face plate labeled "Min/Max". Turning this knob allows the user to adjust the current to match the desired need.

24. When the machine is not in use, make sure the power switch is in the "Off" position.

Safety Warning

Important Notice of Caution before using this machine:

Never operate this equipment on a person has or suffers from the following medical conditions described below. (The salon or operator may consider having a client sign a waiver to the effect that he/she is free of any poor health conditions.)

1. Heart Disease
2. Heart Murmur or Irregular Heart Rate
3. Has a Pacemaker or other type of heart rate balancer.
4. Pregnancy
5. Epilepsy

6. Alcoholism
7. Ulcers
8. Disease of Stomach, Intestine, Liver, or Pancreas.
9. Breast Cancer.
10. Undergoing treatment for any form of Cancer.
11. Had major surgery within the past year.
12. Asthma
13. Cold, Flu, or Upset Stomach
14. Any serious illness or disease not listed above.

Do not use on these areas:

1. Sensitive Skin
2. Skin Rash
3. Damaged, Inflamed, or Infected Skin.
4. Surgical scar less than 12 months old.
5. Skin Swollen.
6. The probes can not point to Eye ball.
7. Genitals.

Caution

1. Always begin treatment conservatively and observe the effects by the touch and appearance of the skin. A vacuum power that is too high or too big may cause an inflammatory response.
2. Edema and erythema can occur following treatment, but will resolve in few hours.
3. Overly aggressive techniques can cause mild flaking over the next 3-4 days. Epidermis that is treated too aggressively can lead to bleeding in the superficial dermis.
4. Diamond Dermabrasion should not be used on sensitive areas such as eyelids. The eyelids should be closed at all times during the procedure and may be covered with damp cotton pads.
5. The following areas are considered inappropriate for treatment: Skin tumor, moles, birthmark, angioma, lips, eyes, herpes, skin cancer, and liver spots.
6. Do not modify or attempt to repair the machine, contact the technician for maintenance.

Trouble Shooting

If the machine does not function properly, please follow these procedures before requesting maintenance.

1. Low suction power/ low vacuum pressure:

Please check again the vacuum hose should be connected well to the socket and then hand pieces. And then turn the Vacuum Regulator to the maximum,

and then start the Diamond Dermabrasion, and then use finger to block the hole of the hand pieces. Normally the vacuum Gauge can reach the 24 inch Hg.

**** Please check the O ring in the Diamond pens & O ring in the socket!**

If you cannot settle the problem finally, please contact your distributor. Don't open the machine by yourself without the technician guidance.

2. No response when turn on the power.

Please exchange the fuse. And connect the power cord correctly. If cannot settle the problem, please contact to your distributor.

3. Supersonic hand piece cannot work properly.

Please connect again the plug. If cannot settle the problem, please contact to your distributor.

Hot & Cold treatment cannot use.

Please connect again the plug. If cannot settle the problem, please contact to your distributor.