



## **Important Notice:**

1. Please read the user manual carefully before any operations.
2. Do not put any liquid onto or into the machine.
3. Avoid operations on pregnant women, or those people with epilepsy, pacemakers or any serious psychological or medical condition, or infections on the area treated.
4. Please disassemble or assemble the equipment by professional personnel. Please check if there is any connection inside the machine dropped when some function does not work.
5. Clean and sterilize the spare parts after each operation.

## ✧ **Working Principle:**

Equipped with 10 groups of pads with BIO-Electricstimulation and Infrared Heating,which work on the special points of body at the same time so that to achieve the effects of Loosing Weight,Breast Enhancement,Muscle Tightening and Removing the internal toxin body.

## ✧ **Product Features:**

- 1.Advanced controller with micro-computer system and digital display.
- 2.Ten pairs of silicon gel electrode pads with ultra strong ELECTRICAL CONDUCTIVITY and INFRARED.
- 3.Ten groups of electrode pulse outputs , independent working or Simultaneous working.
- 4.Five different pulse waveforms,adjustable massage intensity, auto circulation.
- 5.Attaching the electrode pads by elastic bandages, without shedding.

## ✧ **Treatment Applications and Effects:**

1. Reduce cellulite on Arms,Abdnomen,Buttocks,Thighs and Calfs for losing weight purposes.
2. Enhance the blood circulation of shoulder,neck,back,hip and waist and vertebrae,removing body fatigue and stress.
3. Stimulate tensional tight muscle to do contracting exercises , lift and firm slack muscle and muscle mass for slimming and shaping body line.
4. Effectively treating cervical and lumber cervical and lumbar vertebrae strains due to work ;
5. Relieving arthritis,arthralgia,rheumatism.
6. Unclogging mammary gland,developing strong and beautiful breasts.
7. Faradic stimulation to prevent breast cancer and relative diseases.

## ✧ **Applicable Symptoms:**

- ✓ Physiological obesity and dietary obesity
- ✓ Muscle sag
- ✓ Flat and saggy breast
- ✓ Cervical spondylosis
- ✓ Periarthritis of Shoulder
- ✓ Soft tissues、 nerve injury
- ✓ Chronic lumbar muscle strain
- ✓ Pain in waist and lower extremities
- ✓ Sciatica pain、 neurasthenia
- ✓ Rheumatic arthritis
- ✓ Shoulder,hip and waist,back and legs recovering therapy



## ✧ **Technical Parameters:**



**EMS BODY SLIMMING EQUIPMENT****Manual Instructions**

Rated Voltage	AC110V±10%,60Hz±1Hz
Rated Power	<=40W
Working Current	<=50MA
Pulse Output	10 channels
Maximum Treatment Time	90 minutes
Infrared Release	9~14μm
Heating Temperature	0~80°C,adjustable
Pulse Waveforms	5 modes
Pulse Frequency	30Hz/s~800Hz/s
Pulse Power	0~270 angles,adjustable
Pulse Speed	3 levels,optional
Shink Time	0.1~5s
Loose Time	0.1~5s
Machine Dimension	40cm(L)x27cm(W)x13cm(H)
Packing Size	47cm(L)x35cm(W)x23cm(H)
Net.Weight	7.5KG
Gross.Weight	8.5KG

☆ **Panel Introduction:**



1.START/STOP button	2,3.START/STOP indicator light
4.Pulse SPEED button	5,6,7. Speed Indicator Light(3 levels)
8.ERROR warning light	9.TREATMENT TIMER
10,11.Time UP/DOWN	12.Pulse Waveform SELECT button
13,14,15,16.F1~F4 waveforms	17.Cycle Mode (from F1 to F4)
18.Cycle Indicator Light	19~28.Pulse Power Dials
29~38.Pulse Outles	39.Temperature Dial

## ✧ Operation Instructions:

1. Connect Power Wire to the machine.
2. Connect electrode lines to electrode pads, then connect the lines to the machine.
3. Put the electrode pads on the treating area of body symmetrical and wrap the pads with bandages (**Warning: Never put the electrode pads on the heart area**)
4. Keep the ten groups **Pulse Power Dials** off before turning on the machine.
  - **If one power dial is on, the machine would send out alarm voice with error lights flashing**
  - **Once meeting with the alarm, please turn off all the power dials.**
5. Press **Power Switch** to start the machine.
6. Set the **Treatment Timer** to the desired time using the ▲ and ▼ arrows.
7. Press **SELECT** to choose the desired pulse waveform (F1, F2, F3, F4).
8. Press **AUTO** to enter the circulatory pulse waveform (F1~F4).
9. Press **SPEED** to select the desired pulse speed.
10. Press **START** to begin the treatment.
11. Control the pulse intensity of each channel by adjusting the **Pulse Power Dial**.
  - **The pulse intensity should be adjusted from weak to strong.**
12. Control the infrared temperature by adjusting the **Temperature Dial**.
  - **The temperature should be adjusted from low to high.**
13. Press **STOP** to pause the treatment.
14. After the operation, turn off the machine and then take down the electrode pads.

## ✧ **Pulse Waveforms Introductions**

### ✓ **Program 1-Thin Wave:**

This waveform is best for enhancing and building the breast. It can activate the growth of Mammary Gland Cells. To lift and firm the pectoralis major, pectoralis minor muscle and serratus anterior by its electric stimulation so that to build a strong breast support. It would have an obvious effects at promoting the growth of breast, improving the slack and sagged breast.

### ✓ **Program 2-Thick Wave:**

This waveform can improve blood circulation, promote metabolism, firm and tighten skin.

### ✓ **Program 3-Dilatational Wave:**

Tightening and relaxing muscle to make it do aerobic exercise in order to decompose and consume fat and redundant energy inside body, prevent the accumulation of fat, then decompose the formative solid fat, keeping it at the free state.

### ✓ **Program 4-Ups and Downs Wave:**

Forcing the fat to decompose heat, moisture and compressing the plump fat cells, which is further useful for the elimination and consumption of metabolic rubbish, toxin, free fat.

### ✓ **Program 5-Cycle Waveform P1~P4**

## ✧ **Recommended Operations:**

- ① Apply some gel or cream on the electro pads and fix them up with the bandage on the treating area of body;
- ② Attach the pads on the body as followings:

● **Picture 1: Neck and Shoulder**

Massage the neck and shoulder, stimulate the blood circulation so that to reduce fatigue and relieve stress.



● **Picture 2: Neck and Shoulder and Arms**

Massage the neck, shoulder, stimulate the blood circulation so that to reduce fatigue and relieve stress; Removing the cellulite of arms and tightening the saggy muscle so that to slim arms.



● **Picture 3: Breast**

Unclogging mammary gland, stimulating the development of breasts and adjusting the secretion of Hormone.




● **Picture 4: Belly**

Removing the cellulites and lifting the saggy muscle



● **Picture 5: Back**

Exercise the back muscle to promote blood circulation and alleviate fatigue and reduce ache and pain; Remove the redundant cellulite and tighten the back muscle.



● **Picture 6: Hips**

Eliminate the cellulites and lift up the saggy muscle.



● **Picture 7: Legs**



Remove the cellulites and lift up the saggy muscle



✧ **Packing List:**

Item Name	Unit	Quantity	Remark
Mainframe	set	1	
Accessories Holder	Set	1	
Square Electrode Pads	pair	10	
Electrode Lines	piece	10	
Elastic Bandage	set	1	
Fuse and Power Wire	set	1	
English Manual Instructions	Copy	1	