

User Manual

MS-11Y6



Preface

To our valued users,

Thanks for choosing our latest high frequency beauty equipment. This is the latest skin cleansing and acne removing product with integration technology, which applies high frequency technology. MS-11Y6 is a professional beauty equipment, which needs to be operated by people with professional training. Any improper use of this instrument may bring adverse consequence to human body. Therefore, we advise all the people to read this manual thoroughly and strictly follow the instructions specified hereby before operation.

We believe our quality product will offer you good returns, and our perfect after-sales service will make you worry-free.

Thank you and best regards.

Table of Contents

Preface
Part I
Brief Introduction
Advantages
Skin Cleansing & Acne Removing
Head Therapy
Part II
Detailed Operations
Technical Parameters
Contraindications
Precautions for Use
Troubleshooting & Solutions
FAQs
Packing List
Operational Diagram

Part I

Brief Introduction

At present equipment with high frequency technology is relatively popular and extensively used in beauty field. It can effectively replace some traditional acne removing programs in beauty salons, and can be operated easily, conveniently and swiftly, and solves acne problems for people pursuing beauty, and works safely and effectively, and requires no skin breaking, no injection, no medication, no operation, and without side-effects, and is operated externally during the whole process and has an immediate effect.

Advantages

1. It's has operating heads with multiple shapes, which can be changed according to different demands.
2. Electric current with high frequency can release ozone, which can diminish inflammation and do sterilization for skin.
3. It suits acne skin, and can effectively repair acne marks and holes while removing acne.
4. Specially-designed comb-like tube head helps lessen dandruff and improve hair quality.
5. Comfortable, painless, and non-invasive during the process.
6. There is no repairing period, you can leave immediately, and it will not affect your normal work and life after treatment.
7. No consumption, low cost, but with quick returns.
8. Unevenness, bleeding, and swelling will not occur after treatment.

Principles

High Frequency

To increase the activity of skin stratum basal, improve skin texture and skin's defense capability against bacteria with the use of high frequency. The formation of ozone, which has the effect of sterilization and anti-bacteria. It limits sebum's secretion, exsiccated and heals pustule infection. The expansion of skin pores will release impurity, improve blood circulation of lymph and vein, enhance skin texture, and keep the balance of water and grease.

Skin Cleansing & Acne Removing

Effects

1. Diminish inflammation and do sterilization for deep skin.
2. Remove acne, fade acne marks, and lighten and rejuvenate skin.
3. Improve coarse and large pores.
4. Remove acne, and repair dark acne skin.
5. Do sterilization with ozone, and remove acne and kill mites.
6. Accelerate blood circulation, and increase skin metabolism.
7. Lighten skin color and refine pores.

Indications

1. People with exuberant grease secretion and large pores.

2. People with acne skin, and who have acne over and over again.
3. People with extensive acne on face.
4. People who have dark acne after got acne.
5. People with clogged pores and accumulated mites.

Contraindications

1. People who have severely sensitive skin, or in allergic period.
2. People who have prosthesis in the nose, or with extensive trauma on face.
3. Pregnant woman, or people with severe heart disease, severe Three Highs (hypertension, hyperlipidemia and hyperglycemia), or serious diseases, etc.
4. People with skin diseases and infectious diseases.
5. People who just had mesotherapy, micro-needling, or with trauma on face.
6. People who have unrealistic illusion about the effects.

Notes (after treatment)

1. Drink more water, and protect yourself from the sun.
2. Do not use irritating products that contain scrub, exfoliation or AHA (Alpha Hydroxyl Acid).
3. Avoid staying up late, smoking and drinking alcohol. Eat more fruits and vegetables, and light food.
4. Apply moisturizing mask one time each day within three days.
5. Wash your face with warm water within three days.
6. Don't eat and drink too much, and avoid eating spicy food.

Head Therapy

Effects

1. Improve Qi and blood, and accelerate head's qi and blood circulation.
2. Accelerate hair growth.
3. Lessen dandruff, prevent hair loss or greasy hair.
4. Improve dry and broken hair.
5. Improve hair quality, and prevent gray hair.
6. Improve insomnia and dreaminess, and dizziness and headache.

Indications

1. People with dizziness and headache, and insomnia and dreaminess.
2. People who always feel scalp itchy.
3. People who always have hair loss or broken hair.
4. People with poor hair quality and splitting hair.
5. People with dry and oily hair.
6. People with relatively less hair.

Contraindications

1. People with abnormally sensitive skin.
2. People with festered and infected scalp, and with contact skin infectious diseases.
3. People in pregnancy, with severe heart disease, with serious Three Highs (hypertension, hyperlipidemia and hyperglycemia), serious diseases, etc.

4. People who are overly aging.
5. People who had head operation, or in surgical recovery.
6. People who have unrealistic illusion about the effects.

Notes (after head therapy)

1. It's prohibited to stand up immediately.
2. Don't wash the face and head within 6 hours.
3. Avoid catching a cold and having electric fan and AC blowing exactly at the person.

Part II

1. Detailed Operations



Have the handle wire and power cord installed, and then press the power switch.

Detailed Operations for Transparent Tube with High Frequency

Select the proper operating head



and install it on the



handle (as picture).

Twist the part that connect handle's end part with connecting wire, and then adjust the handle to operate.

The Functions for Different Transparent Tube Heads

Comb-shaped Tube Head: used for head, suitable for people with insomnia and dreaminess, and wanting to improve hair quality.

Curved-shape Tube Head: used for sterilization and acne removal of nasal fossa.

Tongue-shaped Tube Head: used for face, diminishing inflammation and removing acne on small area.

Mushroom-shaped Tube Head: used for face, diminishing inflammation and removing acne on big area.

2. Technical Parameters

Output Frequency: 7W

Output Voltage: DC 15V

Power Supply: 110-240V AC power

Frequency: 50HZ

3. Contraindications

People with the following situations should use the instrument with caution. Please consult a doctor or professional before using this instrument. The details are as follow:

- (1) Women in pregnancy, or lactation;
- (2) People with heart disease, or with heart pacemaker;
- (3) People with unhealed surgical wound, or in surgical recovery;
- (4) People with epilepsy, severe diabetes or hyperthyroidism;
- (5) People with malignant tumors, hemophiliac or severe bleeding;
- (6) People with skin diseases or infectious diseases should use with caution;
- (7) People with metal implants inside the body, such as stent, pacemaker, etc., or who are allergic to metal.

4. Precautions for Use

- (1) Plug with ground pin must be used, and ensure power socket is already factually grounded before using the instrument.
- (2) To ensure voltage of the instrument is adaptive. If the voltage of local power supply is unstable, we suggest users add voltage regulator with matching power between the

mains supply and the instrument.

- (3) To guarantee the curative effect and normal service life of the instrument, please uniformly use specified parts provided or suggested by original manufacturer.
- (4) The instrument can't be placed in damp places or near water sources, and also can't be exposed to sunlight directly.
- (5) Don't place the instrument near a strong heat source since this may affect its service life and normal use.
- (6) Please remove all the metal objects from the body before treatment so as to avoid unexpected situations which may affect the curative effect.
- (7) Please restrain from using the instrument aiming at eyes, thyroid, parathyroid, testicles, abdomen of pregnant woman, pacemaker, and etc.
- (8) People who are suffering from illness should use with caution unless he gets permission from a doctor.
- (9) Please turn off the power switch of the instrument if no one uses it, and ensure the main power is off after someone used it and before he left, thus the safety of the instrument can be guaranteed.
- (10) Using the instrument and training the operators in strictly accordance with the instructions specified in the manual.
- (11) Start with the lowest energy level when just operated, and increase the level gradually after the client adjusting himself to it.
- (12) It's prohibited to use the instrument on the rim of eyes, ears, the periphery of parts that have skin diseases, and parts that have plastic surgery.
- (13) Clean the instrument with normal saline after operation so as to ensure its cleanliness and hygiene, as well as extending its service life.

5. Troubleshooting & Solutions

- (1) The instrument cannot be started, and the lights on its back is not working?
 - A. Make sure the power cord is plugged into a good power socket.
 - B. Check whether the fuse tube on the back of the instrument is loose or burnt out.
- (2) The handle is not working?
 - A. Check whether the start button is on.
 - B. If the handle is broken, then you have to contact the manufacturer for repair or reorder.
- (3) The probe is without output?
 - A. Check whether the handle is correctly installed. If no, install it correctly and then operate it.
 - B. If the operating head is damaged, then replace a new one.
 - C. Check whether the power cord is connected, and whether the switch or the switch on the handle is on.

6. FAQs

- (1) Q: Is there any dependency for this treatment?
A: There is no dependency. It helps skin improve in accordance with skin problems. It helps skin absorption, solves problems like acne and coarse skin. If you quit in later stage, the skin will only back to previous conditions of secretion and senility, and it will not cause dependency.

(2) Q: Can I put on makeup after the treatment?

A: You can wear makeup without questions. The skin will become rejuvenated and plump. It's more appropriate to put on makeup at this moment. It's a kind of care that without broken skin or wound, and it's safe, healthy, fast and effective.

(3) Q: How long does it take to see the effects?

A: The causes for acne vary from people to people. There are internal causes and external causes. The management time for acne is related with how serious the acne is, the self-repairing capability of skin, regular management and skin care habits in daily life. Acnes removing should be done step by step.

(4) Q: Do acne will come out again after treatment?

A: There are many factors causing acne, and it has something to do with life style, dietary habits, and skin care habits. If you don't want to have acne you should manage from the following aspects, namely, remove excess grease from hair follicle, keep hair follicle unobstructed, and keep the balance of water and oil.

7. Packing List

Main Machine x1

Comb-shaped Tube Head x1

Curved-shape Tube Head x1


Tongue-shaped Tube Head x1

Mushroom-shaped Tube Head x1

High Frequency Handle x1

Power Cord x1

8. Operational Diagram

Parameters Setting	Product	Technique	Diagram
Skin Cleansing & Rejuvenating: 45 Minutes/ Once a Week			
Adjust handle in accordance with individual demands as long as have treatment time controlled.	Makeup Remover + Facial Cleanser + Toner + Facial Mask + Instrument	<ol style="list-style-type: none"> 1. Remove makeup, 2 minutes. 2. Clean whole face with facial cleanser, 3 minutes. 3. Apply toner on whole face, 1 minute. 4. Apply educe liquid and use hot steam (apply the liquid on parts with large pores or with severe clogging. Cheeks can be applied with toner, and face with 	<p>Technique 9</p>  <p>Technique 7</p>

sensitive skin can be applied with toner and used with cold steam. Hot steam and educe liquid is not advised.), 5 to 8 minutes.

5. Getting rid of ripe acne with acne pin.

6. **Techniques for Transparent Tube with High Frequency:**

Select tongue-shaped tube head, and press the switch on the handle to start.

7. Starting with chin, moving upwards and circlewise, and moving one line after another. Repeat 3 times.

8. Move circlewise from forehead towards hairline. Repeat 3 times. (between the eyebrows can be exercised with more treatment. Repeat 1 to 2 times.)

9. Starting with chin, lifting upwards, and lifting one line after another. Repeat 3 times.

10. Lift from forehead towards hairline. Repeat 3 times. (between the eyebrows can be exercised with more treatment. Repeat 1 to 2 times.)

11. Treat the other side.

12. Choose curved-shape tube head, and press the switch on the handle to start.

13. Techniques: Treat nasal fossa with stamp method.

14. Select mushroom-shaped tube head, and treat whole face with stamp method. Repeat 3 to 5 times.

15. Apply toner, 1 minute.

16. Apply facial mask, and wait for 15 minutes.

17. Remove the mask and clean the face, 2 minutes.

18. Apply toner, essence, facial cream to prevent against the sun.

19. Treatment done.



Technique 8



Course of treatment recommended:

Ten treatments is a course. After one treatment, pores will become very clean, and lots of acne will disappear. After one course, skin color will turn bright little by little, pores will be improved, acne will lessen, and acne marks will become light. After three courses, acne, by and large, will not grow again, and acne marks will become invisible as long as keep diet under control and develop good life style and skin care habits.

