



MS-2183

# Instructions

## Warning

Dear users:

Thanks for choosing our company's product, in order to fully make use of the products, we sincerely suggest follow:

- ◆◆ Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆ Please follow the instruction to fit and operate the machine correctly.
- ◆◆ Please do not remove and change any accessories of the machine.
- ◆◆ It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆ It forbids others to change and maintain the machine except the professional people.
- ◆◆ Please contact us if there is any trouble during operation.

This product does not belong to medical, please don't use for disease treatment!

## System Overview



## Package Listing



large



middle



small

Install cotton filter:



Twist open the IN connector, place cotton filter inside, then turn clockwise to twist the cover back on.

Note:


1. Make sure the cotton filter is placed inside the IN connector or massage oil may be sucked into the machine and cause damage.
2. Replace cotton filter after each treatment to keep it clean.

## How To Test Suction

1. Connect with power cord.

2. Click  to turn on the device.

3. Twist the INTENSITY button  clockwise to the maximum intensity.

4. Place your finger on the IN hole  to see if there's suction.



5. Insert one end of the Y-shaped tube into the IN hole, and then



connect two other ends with two vacuum cups. Click start to turn on the device. Make sure the two cups are in full contact with the skin and then see if there's



suction.

Note: Do not start with high intensity for beginners. Start low and then increase slowly after adaptation.

## How To Connect Tubes?



1. This Y-shaped tube is to be connected with



9 pairs of vacuum cups for scraping, cupping, massage, relaxation and detoxification.

Large vacuum cups can be used on large treatment areas, such as thighs, buttocks, waist, abdomen, arms, etc. Small vacuum cups can be used on smaller treatment areas, such as

arms, calves, face, neck, etc.



2. Connect one end of the Y-shaped tube with the IN hole.




3. Connect two other ends with two vacuum cups.



Turn on the device and adjust intensity.

Note: Make sure to connect Y-shaped ends with two vacuum cups; make sure the cups are in full contact with the skin to avoid possible leakage.



4. When using only one vacuum cup, tie a knot on one of the tubes so



that it doesn't leak air.



5. Y-shaped ends are to be connected with



breast vacuum cups for breast enlargement, buttock augmentation, massage and body shaping.



6. Connect one end of the Y-shaped tube with the IN hole.



7. Connect two other ends with two breast vacuum cups. Turn on the device and adjust intensity.




two breast vacuum cups.

Note: Make sure to connect Y-shaped ends with two breast vacuum cups; make sure the cups are in full contact with the skin to avoid possible leakage.



8. When using only one breast vacuum cup,

tie a knot on one of the



tubes so that it doesn't leak air.



9. This hose can be used to connect



grease cups for lymphatic drainage, meridian dredging, and massage.

Large grease cups can be used on large treatment areas, such as thighs, buttocks, waist, abdomen, arms, etc. Small grease cups can be used on smaller treatment areas, such as arms, calves, face, neck, etc.




10. Connect one end of the hose with

the IN hole, and another end



with grease cup. When using one grease cup only,



tie a knot



on one tube, so it doesn't leak air.

## How To Adjust Parameters ?



1. Power Button



2. INTENSITY +/-

Twist clockwise to increase intensity and counterclockwise to decrease intensity



3. VACUUM METER

The pointer will swing left and right as suction changes



4. Mode Switching

Turn left or right to change mode. The indicator light turns on when specific mode is selected.



## Specification:

Model: MS-2183

Voltage: 220V/110V

Input Power:

Power: 12W

Vacuum Suction: <-80kPa

Material: ABS


Package size: 40\*35\*28cm(15.7\*13.8\*11 inch)

G.W.: 5.4KG

## Functions:

### 1 Breast Lifting

Step 1. Massage breast with massage oil or cream for 2-3 minutes firstly



Step 2. Choose breast cup according to body condition of customer, connect two cups with Small Y type hose.



Step 3. Make customer half-sleep to stretch out breast, cover breast with plastic cloth, don't apply on people with sunken breast.

Step 4. Turn INTENSITY knob to adjust suction of cups and chose the BREAST LIFTING or TREATMENT mode

Step 5. Put both cups on breasts and make sure cup touch skin closely and no air leak. Turn INTENSITY button from left to right till the breast lifting, then adjust INTENSITY knob to get more suction until get the best comfortable feeling of customer.

Step 6. After 20 minutes, turn INTENSITY knob to left and the cups will loose, then apply lifting cream/oil on breasts and massage breasts for 2-3 minutes.

## 2 Butt Lifting

Step 1. Massage butt with massage oil or cream for 2-3 minutes firstly.

Step2. Choose butt cup according to body condition of customer, connect two cups with Small Y type hose.

Step3. Make customer Lying on the bed with face down

Step4. Turn INTENSITY knob to adjust suction of cups, chose the BREAST LIFTING or TREATMENT mode

Step5. Put cups on butts and make sure cup touch skin closely and no air leak. Turn INTENSITY button from left to right till the butts lifting, then adjust INTENSITY knob to get more suction until get the best comfortable feeling of customer.

Step6. After 20 minutes, turn INTENSITY knob to left and the cups will loose, then apply lifting cream/oil on butt and massage breasts for 2-3 minutes.

### 3 Rhythm Massage Treatment (Fingertip Moving)

#### 3.1 For Shoulder, back, arm, waist, thigh, leg

Step1. Massage back with massage oil or cream for 2-3 minutes firstly.

Step2. Choose appropriate cup size, connect the cups with silicon tube



Step3. Turn INTENSITY knob to adjust suction of cup, chose BODY RHYTHM-1 mode

Step4. Put ONE CUP on back and make sure cup touch skin closely and no air leak. Turn INTENSITY button from left to right till the skin lifting, then adjust INTENSITY knob to get more suction until get the best comfortable feeling of customer.

Step5. Move cup slowly from neck to shoulder to back to waist to butt, back and forth for about 3-4 times, turn INTENSITY

knob to left and the cup will loose, then apply lifting cream/oil on skin and massage arm for 2-3 minutes.



### 3.2 For face, chin, neck

- Step 1. Massage arm with massage oil or cream for 2-3 minutes firstly.
- Step 2. Choose appropriate cup size, connect the cup with silicon tube.
- Step3. Turn INTENSITY knob to adjust suction of cup, chose FACE RHTHM-1/FACE RHTHM-2 mode
- Step 4. Put a fingertip cup on arm and make sure cup touch skin closely and no air leak. Turn INTENSITY button from left to right till the skin lifting, then adjust INTENSITY knob to get more suction until get the best comfortable feeling of customer.



Step 5. Move cup slowly from arm to underarm, back and forth for about 10-15minutes, turn INTENSITY knob to left and the cup will loose, then apply lifting cream/oil on skin and massage arm for 2-3 minutes.

#### 4 Scrub and cupping Treatment

##### 4.1 For Shoulder, back, arm, waist, thigh, leg

Step 1. Massage back with massage oil or cream for 2-3 minutes firstly.

Step 2. Choose appropriate cup size, connect the cups and big Y type hose.



Step 3. Turn INTENSITY knob to adjust suction of cups, chose BODY RHYTHM-1 or BODY RHYTHM-1 mode

Step 4. Put both cups on back and make sure cup touch skin closely and no air leak. Turn INTENSITY button get the best comfortable feeling suction and then move cups on other treatment area after skin loose automatically, repeat operation until finish back treatment





Step 5. Apply lifting cream/oil on skin and massage back for 2-3 minutes after treatment

#### 4.2 For Buttock Shape

- Step 1. Massage buttock with massage oil or cream for 2-3 minutes firstly.
- Step 2. Choose appropriate cup size, connect cup and suitable type hose
- Step 3. Turn INTENSITY knob to adjust suction of cups, chose body RHTHM-1/body RHTHM-2 mode
- Step 4. Pull up cup from the base to the thigh and pull from the sides of the hips to the hips, both for 3-5 times, and make sure cup touch skin closely and no air leak. Turn INTENSITY button get the best comfortable feeling suction. Then Select the right cup, secure the hips, and adjust the right amount of energy for 5-10 minutes.


#### 4.3 For Abdomen

- Step 1. Massage Abdomen with massage oil or cream for 2-3 minutes firstly.
- Step 2. Choose appropriate cup size, connect cup and suitable type hose
- Step 3. Turn INTENSITY knob to adjust suction of cups, chose body RHTHM-1/body RHTHM-2 mode
- Step 4. Use cup push along six meridians from up to down to the groin for 3-5 times, and make sure cup touch skin closely and no air leak. Turn

INTENSITY button get the best comfortable feeling suction. Then Select the right cup use in the abdomen, such as acupuncture points as xiayuan, tianshu, daheng, guan yuan, waist side, suction and waiting for 5-10 minutes (depending on the customer).

#### 4.4 For face, chin, neck

- Step 1. Massage face with massage oil or cream for 2-3 minutes firstly.
- Step 2. Choose appropriate cup size, connect both cups and big Y type hose
- Step 3. Turn INTENSITY knob to adjust suction of cups, chose FACE RHTHM-1/FACE RHTHM-2 mode
- Step 4. Put both cups on face and make sure cup touch skin closely and no air leak. Turn INTENSITY button get the best comfortable feeling suction and then move cups on other treatment area after skin loose automatically, repeat operation until finish back treatment



Step 5. Apply lifting cream/oil on skin after treatment

### Recommend Treatment Time

- Scrub Treatment  
4 treatments/course, one treatment/month



Note: shower/bath/swimming is forbid after treatment,12h interval will be okay

- Rhythm Massage Treatment  
4 treatments/course, one treatment/week
- Breast & Butt Lifting  
4 treatments/course, one treatment/week

## Attention

1. Keep the machine level to avoid shaking and collision.
2. Turn off the power and pull out plug when washing the machine.
3. People with serious heart disease, diabetes mellitus and hepatitis are not suitable for treatment with the machine.
4. Women during menstruation and pregnancy are not suitable for treatment.
5. Scrub time should be controlled within 30 minutes with appropriate strength that customers feel comfortable.
6. Don't apply the machine half hour before or after meals, people with empty stomach, fatigue; insufficient sleeping and weak health are not suitable for treatment.
7. People with and blood productivity and spontaneous hemorrhage are not suitable for treatment, including hemophilia, purpuric disease and leukemia.
8. People with serious skin allergic symptoms or acariasis are not suitable for treatment.
9. People with malignant skin tumor or partial skin cancer, fracture, varicosity, partial skin elasticity loss
10. are not suitable for treatment. People in drunk, thirsty or hungry situation should pay special attention.
11. Wash filter regularly to ensure the machine effect and lifetime.
12. When pairing with essential oil, don't apply too much essential oil or cosmetics to avoid damage the machine.
13. After application of the machine, turn the air pressure to the lowest point.
14. Don't apply lymphatic detoxifying treatment when you are annoyed with lymphadenitis.
15. Patients with myocardial infarction, phlebitis depravation and intravenous tumor are not suitable for treatment.
16. Don't apply the treatment within three weeks after surgery for hyperthyroidism, pregnant women.
17. Patient with serious heart disease, diabetes mellitus and hepatitis are not suitable for scrub treatment.
18. Women during menstruation are not suitable for scrub treatment. Pregnant women can have light scrub treatment to promote blood circulation.
19. The scrub time shall not be too long, about half hour is ok with appropriate strength that customers feel comfortable.
20. Don't apply the treatment half hour before or after meals, people with empty

stomach, fatigue; insufficient sleeping and weak health are not suitable for scrub treatment.

21. After scrub treatment, have a glass of warm water after test for a while, and then have warm bath to enhance the effect.
22. You would feel a little bit aching one or two days after scrub treatment that is normal and would automatically disappear.

## Guarantee

### User file

|                                |  |     |                |     |  |
|--------------------------------|--|-----|----------------|-----|--|
| User name                      |  | Sex |                | Age |  |
| Address                        |  |     |                |     |  |
| Telephone                      |  |     |                | Zip |  |
| Purchasing Address             |  |     |                |     |  |
| Retail Price                   |  |     | Order ID       |     |  |
| Commodity name & specification |  |     | Model NO.      |     |  |
| Dealer Signature               |  |     | User Signature |     |  |

Dear Customer:

Thanks for using our Vacuum Cupping Massage Machine. We really hope you could share your purchase experience to other customer if you are 100% satisfied with the machine. Really appreciate your support, time and effort.

If you have any question about the machine, we will prove the following warranty policy:

1. 2 months upon date of purchase, if there is function problem, we should offer complete or partial return/replacement.
2. 1 year upon date of purchase, we will guarantee to assist in maintenance, including providing accessories.
3. All of damage by man-made, we will not guarantee