

---

# USER MANUAL

## MS-22Y7E



---

## Preface

Dear user:

Welcome to use our latest product MS-22Y7E, which replaces all the items in the traditional beauty salon and can be done from head to toe. It not only replaces manual services, but it also has magical effects. MS-22Y7E is a professional beauty equipment that requires professional who is trained to operate. If the instrument is used incorrectly, it may cause adverse consequences to the human body. Therefore, we hope that before operating the MS-22Y7E, please read this manual in detail and follow the instructions in the manual to operate the instrument.

We believe that our quality products will bring you good income, and our perfect after-sales service will make you worry-free.

Thank you!

---

# Catalog

Preface

Part I

Background

Functions

Part II

Detailed Operation

Technical Parameters

Taboos Crowds

Precautions

Trouble Shooting

Packing List

Operation Diagram

---

# Part I

## I. Background

It replaces the lymphatic detoxification program of traditional beauty salons. It can be done from head to toe, facial lift, eye bags, dark circles, lymphatic drainage of the neck, neck and shoulders, breasts, breasts, slimming, hips, abdomen Shape, stomach cold, palace cold, slimming and so on. And there is also a magical effect, which replaces the traditional artificial service, which has the functions of traditional scraping, massage, cupping, massage, and walking. Safe, effective, simple, convenient, no injection, no medication, no surgery, no side effects, full-time external use and immediate results. The negative pressure health instrument can actively exert dynamic negative pressure on the body to give the body effective deep stretching, massage movement, deep dredge, improve the obstruction phenomenon and improve the sub-health state of the body.

### I. Principle

Through vacuum suction and negative pressure probe massage, the body is deeply massaged, and the health and acupoints are used to regulate the excretion of toxins from the body.

The upward aerobic exercise of the skin and fat enhances the function of the lymphatic system and the circulation of blood, and the effect of cupping with negative pressure suction is obvious.

Large size slider is suitable for large area scraping of back, thigh and abdomen. Medium size is suitable for arms, abdomen, shoulder and neck. The small size is suitable for face scraping and lifting, neck lymphatic drainage and other small areas.

### II. Effects

1. Chest, solve chest nodules, hyperplasia, block.
2. Buttocks, low temperature, fat, loose.
3. Scraping, detoxification, cupping, detoxification, phlegm and dampness.
4. Regulate endocrine and enhance physical fitness.
5. Improve immunity and delay aging.
6. Meridian dredge, detoxification.
7. Lymphatic detoxification, scraping.
8. Improve skin tone and promote blood circulation.
9. Accelerate metabolism, soften tissue and reduce fat accumulation.

---

10. Lose weight, slim down, and keep fit.

### III. Treatments area

The instrument can be operated in a wide range of parts, such as face, chest, buttocks, back, abdomen, legs, arms, sucking, detoxification, weight loss and so on.

The scraping is introduced from the inside to the body surface, and finally is emitted from the body surface. The stimulating effect of sucking promotes the peripheral blood flow, which in turn makes the blood flow of the whole body also benefit. Sucking is also a kind of physical therapy like acupuncture. It is a kind of therapy that uses negative pressure to make it adsorb to the skin and cause congestion. It can slowly discharge the moisture and toxins in the lesion, and promote local blood circulation to relieve pain. The purpose of restoring function can be to chill and damp, dredge meridians, qi and activating blood, reduce swelling and pain, detoxify and diarrhea, adjust the balance of yin and yang of the human body, contact fatigue and enhance physical fitness.

## II. Functions

### Facial Care

#### I. Main functions

1. Improve wrinkles and fine lines.
2. Tighten facial skin tissue.
3. Improve allergic skin immunity.
4. Improve skin pigmentation.
5. Accelerate the discharge of aged products.
6. Cervical lymphatic drainage.

#### II. Application

1. A person with a lack of luster on the face, a rough skin and a blackhead.
2. People who drink regularly and often suffer from insomnia.
3. The face is prone to long spots and long acne.
4. The contour of the face is not clear, and the meridians are blocked.
5. Double chin, sagging skin.

---

### III. Taboo crowds

1. Facial allergies or wounds.
2. Severe sensitive muscles.
3. Pregnancy, recovery period.
4. The face has a pustule.
5. Just doing a complete surgery.

### IV. Note:

It is forbidden to expose to the sun, to hydrate more, to keep the skin moist and refreshed, and to maintain the skin regularly. It can be combined with professional skin cleansing and hydrating programs to promote more effective absorption of nutrients.

## III. Back Lymphatic Detoxification

### I. Mainfunction

1. Use negative pressure to smother it on the skin to cause congestion. This treatment method can remove cold and dampness, dredge meridians, eliminate stasis, qi and blood circulation, reduce swelling and pain, remove toxic and diarrhea, adjust human body yin and yang balance, relieve fatigue Enhance physical function and thus achieve the purpose of strengthening the body and eliminating evil and curing the disease. Therefore, many diseases can be treated with cupping therapy.
2. It can also promote blood circulation in the body, eliminate toxins from the body, relieve fatigue, effectively improve dizziness, fatigue, help regulate physical and mental fatigue, accelerate metabolism and waste excretion, and reduce toxicity to skin and internal organs.

### II. Application

- 1.The person who is nervous, anxious, and prone to temper.
2. People with back pain, muscle aches, insomnia, etc.
3. Working in the office, sedentary, easy to sleep and dreamers.

### III. Taboo crowds

1. Three high populations and heart patients.
2. Pregnancy period, lactation period.

- 
3. The surgical wound is healing or recovering from surgery.
  4. Patients with epilepsy and severe diabetes and hyperthyroidism.
  5. Patients with malignant tumors and patients with hemophilia or severe bleeding.
  6. Patients with skin diseases and infectious diseases.
  7. Patients with severe gynecological diseases.
  8. The population during the treatment of gynecological diseases.
  9. Patients with lymphoma or signs of cancer.
  10. People who are too weak.
  11. Drunk, full, fasting, thirst, overworked.

#### IV. Note

1. Do not take a bath until 4-6 hours (or not taking a bath the same day).
2. Avoid blowing and getting cold.
3. Drink more hot water to help detox metabolism.
4. Cupping time is not as long as possible, according to the actual situation.

## IV. Systemic Lymphatic Drainage

### I. Working principles

The lymphatic system is a circulatory system, and the immune system. It is also the body's defense station and recycling garbage station. It prevents external bacterial infection and recovers garbage and toxins from the body. Lymph is also the body's largest detoxification system.

Corresponding to the scavenger in the human body, it is the main medium for collecting toxins from various organs of the human body. If the lymphatic system slows down, it will cause toxins to accumulate in the body, and people will easily become edema and gain weight. The principle of weight loss of lymphatic drainage massage is mainly through the massage and stimulation of lymph glands, dredging the lymphatic circulation system, promoting the excretion of toxins in the body, accelerating fat burning, and achieving the purpose of losing weight. Lymphatic detoxification also has a very important role in allowing our blood vessels to be well cleaned, lymphatic drainage, lymphatic vessels and veins. When the blood volume is too much, the lymphatic vessels can act as a fluid to reduce the burden on the blood vessels. effect.

As the age of the human lymphatic system increases, the speed of circulation may decrease, and the effect of detoxification will also decrease. The appropriate massage can effectively accelerate the circulation rate of the lymphatic system, so that the toxin can be excreted as soon as possible.

---

## II. Main functions

1. Accelerate lymphatic flow.
2. Accelerate the process of natural detoxification, help the body to remove waste through the lymphatic system.
3. Stimulate lymphoid activity.
4. Increase lymphocyte productivity.
5. Improve body metabolism.
6. Promote the function of the immune system.
7. Scraping, sucking, draining moisture and toxins.
8. Dredge meridians, health detoxification, and enhance the body's resistance.

## III. Application

1. People with lack of luster on the face, skin becomes very rough.
2. People who drink often and often suffer from insomnia.
3. People with waist and abdomen obesity can quickly improve edema and consume fat through lymphatic detoxification.
4. people who have bad breath and are often constipated.
5. Poor mental state, poor resistance, easy to catch cold and sick people.

## IV. Taboo crowds

1. Three high population and heart disease patients.
2. Pregnancy period, lactation period.
3. The surgical wound is healing or recovery period.
4. Patients with epilepsy and severe diabetes and hyperthyroidism.
5. Patients with malignant tumors and patients with hemophilia or severe bleeding.
6. Patients with skin diseases and infectious diseases.
7. Patients with severe gynecological diseases.
8. The population during the treatment of gynecological diseases.
9. Have lymphatic cancer patients or signs of cancer.
10. People who are too weak.
11. People who are drunk, full, fasting, thirsty, overworked.

## V. Note

1. Do not bathing until 4-6 hours after the completion.
2. Avoid blowing cold.
3. Avoid eating spicy and greasy things.
4. Drink more hot water to help metabolism.



- 
5. Avoid staying up late, smoking.
  6. Cupping is not the longer the better, depending on the actual situation.

## **V. Breast Therapy**

### **I. Main function**

1. The effect of negative pressure on the human body is attracted by the special tool in the breast through the negative pressure inside the instrument, which produces different negative pressure and suction in the breast, so that the breast is relatively lifted, so that the local capillaries shrink and congest, stimulate the organs, and enhance Cell viability, promote blood circulation and improve body resistance.
2. During the operation, the negative growth pressure and the disappearance of negative pressure cause the local pores to continuously open and close, promote skin respiration, increase oxygen absorption, and accelerate waste removal.
3. It has the functions of warming meridians, replenishing qi and activating blood circulation, promoting blood circulation and removing blood stasis, dissolving phlegm and dispelling phlegm.
4. Evacuate the nodules, relieve breast hyperplasia, and prevent breast disease.

### **II. Application**

1. The breast has nodules, chest pain.
2. Breast size is asymmetrical.
3. There is chest hyperplasia, temper.
4. The chest is meridian, the areola is dark.
5. The axillary lymph is blocked, there is a side milk.
6. The chest is not full, losing flexibility.

### **III. Taboo crowds**

1. Menstrual period, pregnancy, lactation.
2. Patients with diseases.
3. People who have had chest surgery.
4. People with skin problems or chest allergies.
5. Severe hyperplasia, patients with tumors and chest diseases.
6. Heart disease, high blood pressure patients.
7. People who are too weak.

---

## IV. Note

1. Drink warm water and keep warm.
2. Wear stereotypes and comfortable underwear, do not press hard on the breast  
Bathing after 3.4-6 hours.
4. Avoid blowing and getting cold.
5. The initial use of 1-3 times is not too strong, so as not to cause local congestion.
6. Breast asymmetry should operate on the small side first on the other side of the operation.
7. Chest stiffness after breast augmentation surgery, unsuitable for too much intensity.
8. Choose the peak of hormone secretion (the fourth day after menstrual period) for better breast treatment.

## VI. Buttocks Therapy

### I. Main function

1. The effect of negative pressure on the human body is attracted by the special tool on the buttocks through the negative pressure inside the instrument, which produces different negative pressure and suction on the buttocks, so that the buttocks are relatively lifted, so that the local capillaries shrink and congest, stimulate the organs and strengthen Cell viability, promote blood circulation and improve body resistance.
2. Mechanical negative pressure stimulation, the origin of the central nervous system through the reflex pathway, regulating the nerve activity tends to balance.
3. During the operation, the inspiratory contraction, the negative growth pressure and the disappearance of the negative pressure cause the local pores to continue to open and close, promote skin respiration, increase the absorption of oxygen, and accelerate the removal of waste.
4. Increase the elasticity of the buttocks and dredge the meridians.

### II. Application

1. The hips are slack, flat, and the hips are not good.
2. The hip depression affects the beauty.
3. Need more perfect and healthy hips and body.
4. Fat accumulation in the buttocks, loose.
5. People who are often sedentary and have low hips.
6. Dysmenorrhea, people with inflammation.

---

### III. Taboo crowds

1. Menstrual period, pregnancy, lactation.
2. People with diseases such as heart disease and high blood pressure.
3. People with serious diseases, such as acute infectious diseases, skin diseases, etc.
4. Frail elderly, over-aged people.
5. Have had surgery, or are recovering from surgery.
6. People who are too weak.
7. People whose hips are too old.

### IV. Note

1. Do not Bathing until 4-6 hours.
2. When you start using 1-3 times, the strength is not too large, so as to avoid local congestion.
3. You can wear hips and underwear to prevent slack or sagging.
4. Avoid wearing miniskirts and shorts.

## VII. Abdominal Cupping Treatment

### I. Main function

1. Abdominal cupping is the use of external force to open the pores so that the waste in the body can be excreted through the pores, and also has a certain renewal effect on some aging cells. By local conditioning, local blood and body circulation can be accelerated to accelerate metabolism and achieve weight loss.
2. Through these principles, cupping can reduce the fatty stomach of the stomach. Cupping can also have a good therapeutic effect on colds. It can regulate internal organs, balance yin and yang, and improve autoimmune function.

### II. Application

1. People who want to lose weight and build strong muscles.
2. Constipation or indigestion.
3. Abdominal cold, palace cold.
4. Sedentary, with fat, fat accumulation.

### III. Taboo crowds

1. Menstrual period, pregnancy, lactation.

- 
2. People with heart disease, high blood pressure and other diseases.
  3. People with serious diseases such as acute infectious diseases, skin diseases, etc.
  4. Frail elderly, over-aged people.
  5. People who have undergone surgery or are recovering.
  6. People who are too weak.
  7. People who are drunk, full, fasting, thirsty, overworked.

#### IV. Note

1. Do not Bathing until 4-6 hours.
2. The initial strength of 1-3 times is not too large, so as not to adapt or cause local congestion.
3. Keep your abdomen warm and avoid blowing cold.
4. Avoid wearing the navel.
5. It is not easy to operate after hunger and fullness, and it is operated after half an hour after a meal.
6. Avoid eating cold and spicy objects, drink warm water.
7. Cupping is not as long as possible, depending on the situation.

---

## Part II


### I. Detailed Operation

Connect the power line and turn on the power button, it will light up. Connect the tube to the holes.





Turn on or turn off the machine



Button For Adjust Vacuum Energy

## VIII. Technical Parameters

Package size: 25\*25\*20cm

Voltage: 220V-240V / 50Hz; 100V-120V / 60Hz

Power: 50W

Gross weight: 3.5 kg

## IX. Taboo Crowds

1. Patients with the following symptoms should be cautious when using the device.
2. Please consult the doctor or professional for details in the following situations.
3. Pregnant or breastfeeding woman.
4. Heart disease or pacemaker.
5. Patients whose surgical wounds have not healed.
6. Patients with epilepsy, severe diabetes and hyperthyroidism.
7. Patients with malignant tumors, patients with hemophilia or patients with major bleeding.
8. Patients with skin diseases and infectious diseases should use caution.

## X. Precautions

1. Keep the machine level and avoid vibration and collision.
2. When washing the machine, turn off the power and unplug it.
3. People with severe heart disease, diabetes and hepatitis should not be treated with this machine.
4. Women during menstruation and pregnancy should not be treated.
5. The operating time should be controlled within 30 minutes and make the customer feel comfortable.
6. Do not use this machine before or after half an hour after meals, fasting, fatigue, fatigue, lack of sleep and physical weakness are not suitable for treatment.
7. Patients with hematological and spontaneous bleeding such as hemophilia, purpura, and leukemia should not be treated.
8. People with severe skin allergies or rickets should not be treated.

- 
9. Malignant skin tumors or partial skin ulcers, fractures, varicose veins, and partial skin loss are not suitable for treatment.
  10. Special attention should be paid to drunkenness, thirst or hunger.
  11. Clean the instrument head regularly to ensure the machine's effectiveness and service life.
  12. When using with essential oils, do not use too much essential oil or cosmetics to avoid damaging the machine.
  13. After using the unit, adjust the air pressure to the lowest level.
  14. Do not use lymphatic drainage therapy when you are suffering from lymphangitis.
  15. Patients with myocardial infarction, phlebitis, and venous tumors should not be treated.
  16. Hyperthyroidism, pregnant women should not use this product within three weeks after surgery.
  17. Patients with severe heart disease, diabetes, and hepatitis should not be scrubbed.
  18. Women are not suitable for treatment during menstruation. Pregnant women can use light scrubs to promote blood circulation.
  19. The operating time should not be too long, about half an hour, and the strength is moderate, so that customers feel comfortable.
  21. One or two days after the operation, you will feel a little pain. This is normal and will disappear automatically afterwards.
  22. Remember that the liposuction head and the soft hose must be inserted into the corresponding suction holes.
  23. In order to use the instrument better, be sure to install the filter before use.
  24. To extend the life of your instrument, please ensure that it is kept within 30 minutes of each use.
  25. At the beginning of the operation, the suction is adjusted from small to large to avoid discomfort.

## **XI. Trouble Shooting**

1. The instrument cannot be started, and the button light on the back of the instrument does not light.
  - A. Make sure the power cord is connected to an active power outlet.
  - B. The fuse tube on the back of the instrument is loose or burnt.
2. The instrument has no suction or suction is too small.
  - A. Please turn off the instrument and check the filter element of the instrument. It may be necessary to replace the filter.
  - B. Please check the rubber ring of the handle plug, which may be caused by the rubber ring wear.
  - C. Please check if the oil filter cup outside the filter is tight and the rubber ring is worn. Air leaks at this location may cause a difference in air pressure.

D. If the above method cannot be processed, please contact the instrument dealer for assistance.

3. The instrument can be started, but the information displayed is wrong.

A. Unplug the power cord from the back of the instrument and wait for about 1 minute before inserting and restarting the instrument.

B. If the above method cannot be processed, please contact the instrument dealer for assistance.

## XII. Packing List

1 x main machine

1 x tubes


1 x large scraping and detoxification operation head

1 x medium scraping and detoxification operation head



1 x small scraping and detoxification operation head

1 x power line

## XIII. Operation Diagram

| Parameter Adjustment  | product  | Technique   | Diagram   |
|---|--|---|---|
| <b>Face (Operating time is 45 minutes); once a week</b>       |  |   |   |
| Suction is slowly adjusted from weak to strong<br>Energy: 1-5 | cleansing oil<br>+facial cleanser<br>+<br>massagecream<br>+ instrument | <ol style="list-style-type: none"> <li>1. Facial makeup removal, cleaning 5 minutes.</li> <li>2. Full-face clean-up of skin water.</li> <li>3. Hands on the whole face acupoints, Chengjiang, Dicang, Daying, chiach'e, hsia kuan, Yingxiang, cheek crevice, temple, Ching Ming hole, Cuanzhu points, Yuyao, Sizhukong, apple of the eye seam, Chengqi, ssu pai, ) 3 times.</li> <li>4. Hands alternately lift the three lines of the face to the underarms, 3 times.</li> <li>5. Operate the other side.</li> <li>6. Instrument operation: in the jaw line combined with the hand pull to the</li> </ol> |  |










|  |  |   |   |
|--|--|---|---|
|  |  | <p>Armpit discharge<br/>, 3 times.</p> <p>7. Lifting the line up on the face, 3 times.</p> <p>8. 8. Lifting the hairline and the sides on the forehead, 3 times.</p> <p>9. Take a massage at the forehead, 3 times.</p> <p>10. Operate the other side.</p> <p>11. Choose a suitable .<br/>transparent can in the nose to absorb and put a massage, you can suck blackheads, 3-5 times.</p> <p>12. Wash your face.</p> <p>13. End.</p> |  <hr/>  |
|--|--|---|---|

**Treatment advice**  
 10 times for a course of treatment, once finished rosy in the back, blood circulation is accelerated, after a course of treatment, the skin begins to metabolize toxins, the acne muscles are improved, the blackheads are slowly reduced, and the contours of the skin are enhanced under three courses of treatment, and the resistance is enhanced. The skin is healthy and shiny.

Suggested price: 698/time 4980/ten times

**Back(Operating time is 60 minutes); once a week**

|   |   |  |   |
|---|---|--|---|
| <p>Suction is slowly adjusted from weak to strong<br/>Energy: 1-5</p> | <p>Lymphatic detoxification equipment<br/>+ cupping<br/>+ essential oil</p> | <p>1. Manipulation: hands licking the oil to the waist to point to the wind and wind, 3 times.</p> <p>2. Push the Du Meridian and the bladder through both hands, 3 times.</p> <p>3. Hands alternately push the Du Meridian and Bladder, 3 times.</p> <p>4. The palms of both hands are alternately pushed from the Dazhui to push the Du Meridian and the Bladder Meridian and pull it from the waist side to discharge through the armpits, first left and then right, 3 times.</p> <p>5. Push your shoulders with your thumb and sew your shoulders to your underarms. After 3 passes, drain from your armpits to your fingers.</p> <p>6. The palms of both hands alternately</p> |  <hr/>  <hr/>  |
|---|---|--|---|

|  |  |  |  |
|--|--|--|--|
|  |  | <p>start from the Dazhui. The ribs are pushed along the line and pushed to the waist. The hands are discharged from the waist to the armpit. 3 times.</p> <p>7.The hands of the tiger's mouth alternately from the hairline through the large plate to the finger discharge, 3 times.</p> <p>8.Operate the other side,</p> <p>9.Instrument operation: lymphatic detoxification instrument is pushed from the vertebra to the caudal vertebra, 3 times,</p> <p>10.The instrument pushes the bladder to the caudal vertebra (first left and right) 3 times.</p> <p>11.The instrument pushes the shoulder to squat to the armpit, 3 times.</p> <p>12.The instrument pushes the line to the waist in the direction of the ribs, 3 times.</p> <p>13.The instrument starts to push the large plate through the neck, 3 times.</p> <p>14.On the other side of the instrument operation,</p> <p>15.Going the tank: If the body is heavy, you can go to the tank, the bladder and the governor go back and forth 3 times, the shoulder blade is back and forth 3 times, and the whole back goes 3 times in the direction of the rib (depending on the customer's stress).</p> <p>16.Choose the right cupping head and go back and forth in the neck and slab position, 3-5 times.</p> <p>17.Fixed cans: Dazhui Shoulder wells Tianzong points Feishu points Bagua area Arms 5-10 minutes (according to customer's stress).</p> <p>18.Take the can and wipe it with a towel.</p> <p>19. End.</p> |  <hr/>  <hr/>  <hr/>  |
|--|--|--|--|

Treatment advice

10 times for a course of treatment, once finished, the body out of the body, expelling body moisture and

wind evil, a treatment to the back becomes easy, relieve the symptoms of shoulder and neck pain, after three courses of wealth disappeared, blood circulation is normal, sleep quality Good, physical enhancement.



Suggested price: 698/time 4980/ten times








**Body(Operating time is 120 minutes); once a week**








Suction is slowly adjusted from weak to strong Energy: 1-5








equipment + essential oil


1. face + neck remover, clean.
2. Full face rubbing lotion.
3. Hands on the whole face acupoints Chengjiang, Dicang, Daying, chiach'e, hsia kuan, Yingxiang, cheek crevice, temple, Ching Ming hole, Cuanzhu points, Yuyao, Sizhukong, apple of the eye seam, Chengqi, ssu pai, ) 3 times.
4. Hands alternately lift the three lines of the face to the underarms, 3 times,
5. Operate the other side.
6. Neck: Put your hands on the chest and rub the oil to the back of the neck and wind up the wind, 3 times.
7. First left and right: hands alternately from the lower jaw lymph through the back of the ear to the armpit, 3 times
8. Unilateral tiger mouth ringing push the neck side three lines to the armpit, 3 times.
9. One-sided finger kneeling hot the neck side and push the lymph of the neck to the underarm, 3 times.
10. Comfort your hands to the point of Feng Chi Feng Fu, 3 times.
11. The other side is the same as above.
12. Select the small head of the instrument on the face from the chin and lift it up one line and lift it through the ear to the armpit, 1-2 times (according to the customer's condition).
13. Select the appropriate instrument to slide up and on both sides of the forehead, 1-2 times (according to










|  |  |   |  |
|--|--|---|--|
|  |  | <p>customer's situation).</p> <p>14. Lymphatic detoxification instrument unilaterally pulls from the mandibular lymph with the hand alternately through the back of the ear to the armpit, 3 times.</p> <p>15. The instrument is arranged in a circle on the neck side and circled to the armpit, 3 times.</p> <p>16. Sliding the inner clavicle to detox, 3-5 times.</p> <p>17. Select the appropriate small head to start from the ear to the underarm, 2-3 times (according to customer situation).</p> <p>18. The instrument operates on the other side as above.</p> <p>19. Underarm manipulation, first left and then right: Use both hands to save oil under the armpit, 3 times.</p> <p>20. One hand palm root knead the pole spring point, 3-5 times.</p> <p>21. The hands of the tiger's mouth alternately push the gallbladder through the armpit to the finger discharge, 3 times.</p> <p>22. Both hands alternately push the inner clavicle lymph to the finger discharge, 3 times.</p> <p>23. The hands of the tiger's mouth alternately push the arm three meridians (pulmonary pericardium through the heart) 3 times.</p> <p>24. Comfort the entire arm to the finger discharge, 3 times.</p> <p>25. Kneeling finger rubbing arm.</p> <p>26. Palm roots press the pole springs, 3-5 times.</p> <p>27. Overall comfort to finger discharge, 1 time.</p> <p>28. Instrument operation, first left and then right.</p> <p>29. Push the gallbladder from the chest to the armpit and appease 3 times by hand.</p> |        |
|--|--|---|--|


|  |  |  |
|--|--|--|
|  | <p>30. Push three meridians from the elbow to the armpit 3 times.</p> <p>31. Operate the other side, <i>ibid</i>.</p> <p>32. From the liver and gallbladder to the armpit, 3 times.</p> <p>33. Select the appropriate instrument to walk the cans in the arm position, 1-2 times (according to customer situation).</p> <p>34. Choose the right cupping head inside and outside the arm, liver and gallbladder, set the can, 5-10 minutes (according to customer situation).</p> <p>35. Inguinal lymphatic drainage method:<br/>first left and right: hands 3-5 times in the whole leg.</p> <p>36 points, press the door, 3-5 times.</p> <p>37. Hands and three mouths alternately push the four meridians (spleen and stomach through the liver and gallbladder) 3 times.</p> <p>38. The whole appeasement 3 times</p> <p>39. Both hands are transferred to the groin 3-5 times</p> <p>40. The thumb of both hands alternately dials the groin 3-5 times.</p> <p>41. Two-handed twist twists and turns 3 times.</p> <p>42. Hands rubbing the hot thigh roots.</p> <p>43. The instrument pushes four meridians to the groin, 3 times.</p> <p>44. Circle back and forth in the groin position, 3 times.</p> <p>45. The instrument and the hand alternately push the groin, 3 times.</p> <p>46. Push the line from the middle of the thigh to the groin while appeasing, 3 times.</p> <p>47. Operate the other side (<i>ditto</i>).</p> <p>48. Walk the cans at the thigh position (one line next to the line) 1-2 times (according to customer's stress).</p> <p>49. Set the cans in the legs, thighs, Zusanli, Sanyinjiao, about 10 minutes.</p> |  <hr/>  <hr/>  <hr/>  <hr/>  <hr/>  <hr/>  |
|--|--|--|

|  |   |  |
|--|---|--|
|  | <p>50. Remove the can and end.</p> <p>51. The customer is kneeling, left leg: hands from the calf oil to the thigh bag back to the heel 3 times.</p> <p>52. The palms of both hands alternately push the leg pack from bottom to top and return to the heel for 3 times.</p> <p>53. The hands of the tiger mouth alternately push the bladder tendon from the bottom to the top - kidney tendons - liver tendons - gallstones each 3 times.</p> <p>54. Push the axillary part 3 times with both hands.</p> <p>55. Hands alternately twisted and twisted 3 times from bottom to top.</p> <p>56. Appease 3 times.</p> <p>57. Pushing the bladder ribs from the bottom to the top with both hands - kidney ribs - liver ribs - ribs 3 times each.</p> <p>58. End of appeasement.</p> <p>59. Instruments: first calf, thigh, detoxification instrument from the bottom to the armpit, push the bladder through the kidney - kidney - liver - gallbladder 2-3 times.</p> <p>60. The instrument pushes four meridians from the armpit to the root of the thigh, 2-3 times.</p> <p>61. The instrument extends four meridians from the calf and pushes the detox to the root of the thigh, 3 times.</p> <p>62. Choose the right cupping head and extend the four meridians in the leg for 1-2 times (according to the customer's condition).</p> <p>63. Set the cans in the legs, the center of the vulture (armpit), the thigh, the Chengfu (middle of the thigh), the Chengshan (middle calf), the Yongquan (foot) for about 10 minutes (according to the customer) Happening).</p> <p>64. End of cans.</p> |        |
|--|---|--|

|   |                                  |  |  |
|---|----------------------------------|--|--|
|   |                                  |  |  |
| <p>Treatment advice</p> <p>10 times for a course of treatment, once finished, the whole body is smooth, the lymph of the whole body begins to circulate, the metabolism of the body is accelerated after a course of treatment, the nodules of the lymphatic block begin to disappear slowly, the quality of sleep and gynecological problems are improved after three courses of treatment, and the menstruation is normal. The size of the arm thighs becomes smaller. The axillary lymph node bulge basically disappeared.</p> <p>Suggested price: 980/time 6980/ten times</p> |                                  |  |  |
| <p><b>Breast(Operating time is 60 minutes); once a week</b></p>   |                                  |  |  |
| <p>Suction is slowly adjusted from</p>  | <p>equipment + essential oil</p> | <p>1.Standing on the head of the bed, start with the hands of the sputum in the middle of the sputum to the underarm</p> |  |

|                                       |  |  |  |
|---------------------------------------|--|--|--|
| <p>weak to strong<br/>Energy: 1-5</p> |  | <p>bag and lift the suspensory ligament 3 times.</p> <p>2.Acupuncture points with both hands: Shanzhong, Rugen, Dabao, Yingchuang, Zhongfu, Yunmen, 3 times.</p> <p>3. Appease action: start half a circle from the squat to the squat to the lifting ligament 3 times.</p> <p>4. The palms of both hands overlap. Starting from the middle of the hole, you can draw 8 words on both sides of the chest. 3 times.</p> <p>5. The palms of both hands are pulled from the position of the secondary milk to the suspensory ligament. First left and right 10 times.</p> <p>6. Appease action 3 times.</p> <p>7. Standing on the side of the customer, first left and right: both hands palms to raise the chest 3 times.</p> <p>8. Two-handed tiger mouth pushes the breast duct around the chest: (one round of the chest) 3 times</p> <p>9. Thumbs up to clear the chest nodules (first alternately in the circle) 3 times</p> <p>10. Both hands palms appease the entire chest.</p> <p>11. Operate the other side (same as 6.7.8.9.).</p> <p>12. Choose the right cupping head and gently lift the can on the chest, 1-2 times.</p> <p>13. Set the cans in the chest, the warehouse corner (below the collarbone), the Zhongfu point (near the armpit position) for about 5-10 minutes (according to the customer's condition).</p> <p>14. The nipple retracted can be selected in the appropriate position of the nipple, about 3-8 minutes (according to customer situation).</p> <p>16. Remove the instrument cup and Wipe clean.</p> |  <hr/>  <hr/>  <hr/>  <hr/>  <hr/>  <hr/>  |
|---------------------------------------|--|--|--|




|  |  |          |   |
|--|--|----------|---|
|  |  | 17. End. |  |
|--|--|----------|---|





**Treatment advice**

10 times for a course of treatment, once the chest is done, the chest has a lifting effect, a treatment has a firming effect on the chest, the breast and the expansion and sagging will slowly tighten and concentrate, and the three courses will improve the nodules and hyperplasia. The chest is healthy, away from the disease, the chest is tight and full, and the effect will be more obvious.

Suggested price: 698/time 4980/ten times

**Buttocks(Operating time is 60 minutes); once a week**

|   |                           |   |   |
|---|---------------------------|---|---|
| Suction is slowly adjusted from weak to strong<br>Energy: 1-5 | equipment + essential oil | <ol style="list-style-type: none"> <li>1. Technique: Stand sideways with both hands and squeegee from the gossip area to the waist and lift it up from the waist along the buttocks (ie, appease action) 3 times.</li> <li>2. Push the gossip area with both hands and thumb 3 times.</li> <li>3. Appease the movement 3 times after the point: Shenshu point -Baliao- Changqiang - Huantiao - Cheng Fu point 3 times.</li> <li>4. Appease action 3 times.</li> <li>5. First left and right, hands extending from bottom to top of thighs-bladder tendon-kidney meridian-liver meridian-bile channel pushed straight to the veins 3 times each.</li> <li>6. The hands are extended from the bottom of the thigh from bottom to top - the bladder - the kidney - the liver - the gallbladder is divided to the veins 3 times.</li> </ol> |  |
|---|---------------------------|---|---|


|  |  |   |  |
|--|--|---|--|
|  |  | <p>7. The palms of both hands overlap and push up from the roots of the thighs to the veins (lift) 3 times.</p> <p>8. Hands from the bottom of the buttocks from bottom to top, a line of straight to the highest point of the buttocks (plastic) back and forth 3 times.</p> <p>9. Repeat action 7.</p> <p>10. Appease action.</p> <p>11. Operation on the other side: Same as 5.6.7.8.9.10.</p> <p>12. Hands and hot gossip area.</p> <p>13. Choose the right cupping device to carry the can: pull up from the root of the thigh, 3-5 times.</p> <p>14. Lifting 3-5 times from the hips to the hips.</p> <p>15. Fixed cans: Select a suitable cup, fix the buttocks, adjust the appropriate energy, about 5-10 minutes.</p> <p>16. Take off the cup.</p> <p>17. End.</p> |     |
|--|--|---|--|








**Treatment advice**







10 times for a course of treatment, after the completion of the hips, the buttocks have the effect of lifting and unblocking. After a course of treatment, the excess fat in the buttocks begins to disappear, the surrounding fat begins to concentrate, and the hip curve is formed after three courses of treatment, which increases the elasticity of the buttocks and improves the gynecological palace cold. problem.

Suggested price: 698/time 4980/ten times

**Abdomen(Operating time is 60 minutes); once a week**

|   |                                  |  |   |
|---|----------------------------------|--|---|
| <p>Suction is slowly adjusted from weak to strong<br/>Energy: 1-5</p> | <p>Equipment + essential oil</p> | <p>1. Abdominal ringing oyster sauce 3 times.</p> <p>2. Hands back and forth 3 times.</p> <p>3. Hands alternately knead belly fat by chiropractic for 3 times, soothing excessive.</p> <p>4. Hands alternately pull the waist with 16 veins on both sides.</p> |  |
|---|----------------------------------|--|---|

|  |  |  |  |
|--|--|--|--|
|  |  | <p>5. Both hands are drawn to the waist and then lifted from the waist side to the navel 3 times.</p> <p>6. Hands overlap with a small circle clockwise 3 times intestines.</p> <p>7. The whole hands are circled and appease.</p> <p>8. Acupuncture points: Shangyu Zhongyu Xia Shen Shen Qi Hai Guan Yuan Intermediate Tianshu Daheng Bone. 2 times.</p> <p>9. Hands and circles to appease 3 times.</p> <p>10. Hands alternately push the meridians from top to bottom (Renmai - side by side 2 inch kidney - side by side 4 inch stomach - side by side 6 inch spleen - side by side 8 inch liver - both sides of the gallbladder) 3 times.</p> <p>11. Push your thumb straight from the middle to the navel and slide it to the waist and lift it up to the groin, 3 times.</p> <p>12. Hands and circles to appease 3 times.</p> <p>13. Appease the platoon to the end of the groin.</p> <p>14. Combine the hand to pull the ulnar vein to the groin, first left and right, 6-8 times.</p> <p>15. The instrument pushes six meridians from top to bottom to the groin, 3 times.</p> <p>16. Choose the right tank in the abdomen, squat, Tianshu, Daheng, Guanyuan, waist side, about 5-10 minutes (according to customer situation).</p> <p>17. Take the can end.</p> |        |
|--|--|--|--|

|  |  |  |  |
|--|--|--|--|
|  |  |  |       |
|--|--|--|--|

**Treatment advice**

10 times for a course of treatment, after a single operation, the abdomen is detoxified. After a course of treatment, the excess fat in the abdomen is slowly reduced, and the digestion and absorption capacity is improved. The three courses of treatment improve the abdominal coldness, improve immunity and strengthen the body.

Suggested price: 698/time 4980/ten times

---