

# User Manual

MS-5481



# Preface

Dear Users,

We're pleased to present to you our latest 5-in-1 comprehensive beauty machine MS-5481, which includes Face RF, Tripolar RF, Body RF, Vacuum Bipolar RF and 40K technologies for body slimming and skin firming. It is a professional beauty machine, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore, we advise any personnel to read this manual thoroughly and follow the instructions strictly before operation. We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

# Table of Contents

|                                |  |
|--------------------------------|--|
| Preface                        |  |
| Part I                         |  |
| Introduction                   |  |
| Advantages                     |  |
| Working Principles             |  |
| Face Firming and Anti-ageing   |  |
| Eye Anti-wrinkle               |  |
| Neck Care                      |  |
| Body Management                |  |
| (1). Waist&Abdomen Shaping     |  |
| (2). Arm Shaping               |  |
| (3). Back Shaping              |  |
| (4). Hip Shaping               |  |
| (5). Leg Shaping               |  |
| (6). Breast Shaping            |  |
| Part II                        |  |
| 1. Detailed Operation          |  |
| 2. Technical Parameters        |  |
| 3. Dos and Don'ts              |  |
| 4. Troubleshooting & Solutions |  |
| 5. FAQs                        |  |
| 6. Packing List                |  |
| 7. Simplified Treatment Steps  |  |
| 8. Operational Diagrams        |  |

# Part I

## Introduction

Face RF, Tripolar RF, Body RF, Vacuum Bipolar RF and 40K are currently among the most popular approaches for skin care, skin firming, body shaping and anti-ageing. So, this machine is an effective replacement of all the facial and body care programs in beauty salons. Moreover, it is convenient and easy to operate. It solves both skin and body management problems for beauty-seekers. This equipment requires no injection, no medication, and no operation, for which it has no side effects. It is used externally during the whole process and has an immediate effect. RF heats deeper skin layers, stimulates collagen to regenerate and reorganize in dermis, and promotes blood circulation in skin bottom, which can achieve the results of skin firming and anti-aging. Ultrasounds with specific frequencies target those parts in which there are stubborn fat to bring customers amazing results.

## Advantages

1. This 5-in-1 beauty machine is used for facial and body care.
2. The ultrasonic wave can powerfully burn and remove fat, which surpasses liposuction to solve the problems for obese people.
3. There are several hand pieces for treatment. You can use corresponding one based on different needs and different treatment areas.
4. Single hand piece has various functions. The RF has more powerful and even energy release.
5. The perfect combination of ultrasound and radio frequency solves major problems in losing weight. It helps with weight loss and meanwhile firms skin, increases skin elasticity, promotes metabolism, and ultimately strengthen people's constitutions.
6. It adopts the advanced technology of ultrasonic wave fat-burning.
7. It is painless and non-invasive treatment process. Once the treatment is over, customer can go back home and there is no downtime, which means it will not affect customer's life and work.
8. No consumption, low cost and quick returns.
9. Wider treatment range for body and face.
10. Unevenness, bleeding, and swelling will not appear after treatment.
11. Various handpieces have the energy release of red light, which can dissolve fat and massage body, and accelerate metabolism and kill bacteria at the same time. The effect is remarkable with a higher comfort level.

## Working Principles

### Face RF

RF wave can directly penetrates the skin and uses the resistance formed by the skin to produce energy, which raises the temperature of the bottom layer of the skin. Using the principle of dermis collagen protein tightening and stimulating collagen protein regeneration, it has two therapeutic effects: timely skin firming and lasting collagen protein regeneration.

After treatment, an obvious sense of contraction can be felt immediately at the treatment area, and the effect of skin firming and lifting can be seen about a month later, and the effect will become more and more obvious. The effect of RF treatment can last for about 18 months to 3 years depending on the area of treatment and individual maintenance methods, which, however, may be longer depending on individual maintenance methods. The latest research points out that better results can be achieved after repeated treatment.

### Biological Effects

The waves vibrating millions of times per second can penetrate the epidermis and then directly act on the dermis to heat up collagen cells in a deep layer, and stimulate the growth of collagen fibers. In this way, it can make the skin firmer and more elastic. In the process of treatment, the RF waves can pass through the epidermis and act on the deep dermis to produce heat energy rapidly. When collagen tissue is heated to 45 °C ~ 60 °C, it will contract immediately and then the collagen will be stimulated to regenerate. At the same time, biological heat can effectively accelerate the blood flow and make free fatty acid (ffa) released to promote the dissolution of superficial fat to achieve the purpose of tightening and lifting saggy part of the body.



### Tripolar RF

RF wave can directly penetrates the skin and uses the resistance formed by the skin to produce energy, which raises the temperature of the bottom layer of the skin. Using the principle of dermis collagen protein tightening and stimulating collagen protein regeneration, it has two therapeutic effects: timely skin firming and lasting collagen protein regeneration.

After treatment, an obvious sense of contraction can be felt immediately at the treatment area, and the effect of skin firming and lifting can be seen about a month later, and the effect will become more and more obvious. The effect of RF treatment can last for about 18 months to 3 years depending on the area of treatment and individual maintenance methods, which, however, may be longer depending on individual maintenance methods. The latest research points out that better results can be achieved after repeated treatment.

### Biological Effects

The waves vibrating millions of times per second can penetrate the epidermis and then directly act on the dermis to heat up collagen cells in a deep layer, and stimulate the growth of collagen fibers. In this way, it can make the skin firmer and more elastic. In the process of treatment, the RF waves can pass through the epidermis and act on the deep dermis to produce heat energy rapidly. When collagen tissue is heated to 45 °C ~ 60 °C, it will contract immediately and then the collagen will be stimulated to regenerate. At the same time, biological heat can effectively accelerate the blood flow and make free fatty acid (ffa) released to promote the dissolution of superficial fat to achieve the purpose of tightening and lifting saggy part of the body.



### Body RF

Based on the RF thermal effect, a specific radio frequency wave is produced at a specific depth of the skin, which acts directly on the dermis. In this way, it heats the tissue to promote the catabolism of subcutaneous fat, and stimulate the regeneration and reorganization of collagen

and elastic fibers at the same time. With the continuous regeneration and reorganization of collagen in the dermis, the saggy skin will be lifted and firmed after several courses of treatment are taken.

### Biological Effects

The waves vibrating millions of times per second can penetrate the epidermis and then directly act on the dermis to heat up collagen cells in a deep layer, and stimulate the growth of collagen fibers. In this way, it can make the skin firmer and more elastic. In the process of treatment, the RF waves can pass through the epidermis and act on the deep dermis to produce heat energy rapidly. When collagen tissue is heated to  $45\text{ }^{\circ}\text{C} \sim 60\text{ }^{\circ}\text{C}$ , it will contract immediately and then the collagen will be stimulated to regenerate. At the same time, biological heat can effectively accelerate the blood flow and make free fatty acid (ffa) released to promote the dissolution of superficial fat to achieve the purpose of tightening and lifting saggy part of the body.




### 40K

The principle of ultrasonic wave is used to make human adipose cells take a strong impact and produce friction movement among them after strong sound waves enter into the human body, which can effectively consume heat and water of cells and shrink adipose cells. In addition, when the sound wave vibrates, it can produce a strong impact on cells. And then the adipose cells will burst instantly, and they will be reduced so as to achieve the effect of fat removal.

**Ultrasonic Cavitation:** Tens of thousands of tiny bubbles, namely cavitation bubbles, are produced by vibration of liquid. These bubbles grow in the negative pressure region formed by the longitudinal transmission of ultrasound, which will close rapidly in the positive pressure region to be compressed and stretched under alternating positive and negative pressures. When the bubbles burst, there will be huge instantaneous pressure of up to tens of MPa to hundreds of MPa, and strong vibration and noise will be produced.

**Advantages:** It only aims for low-density adipose tissue for its specific frequency, and protects high-density tissue such as vascular and nerve tissue.



### Vacuum Bipolar RF

Massaging skin and muscles with a special vacuum suction handle can effectively improve the humoral fluidity of human cells, and increase the movement of cells to improve skin elasticity.

Meanwhile, it can accelerate the blood circulation of microvessels, expel toxins from the body through the normal circulation of the lymphatic system, and reduce the probability of occurrence of adverse conditions such as spots, pigmentation and congestion in various parts of the body. The kneading sensation generated by vacuum suction can increase the tissue activity of skin and muscle, and therefore help to reduce tough tissue and increase the elasticity of skin tissue, so that both body slimming and shaping can be realized at the same time. In addition, the movement of vacuum cup can also stimulate sympathetic nervous system to lower skin sensitivity. The suction and release process of air pressure can not only improve capillary system but also promote the flow between deeper veins and lymphoid gland network to strengthen blood vessels and alleviate the condition of varicose veins.

**Advantages:** it is superior to the previous vacuum working mode which is single and simple. Using different modes for different parts can be more efficient and more effective for slimming and shaping body, and it will never harm human body. Uniquely designed RF tips integrate RF and vacuum as a whole, in which RF and vacuum can exist independently and also are complementary to each other. Compared to the ordinary one-polar RF, it can ideally shape body in a more effective, faster and evener way

## **Vacuum Physical Effects**

### **1. Skin layer**

The suction and release of air pressure:

- (1) It can improve the fluidity of cells, thus increasing the movement of cells, and can treat diseases related to blood stasis and blood stagnation.
- (2) It can drain excess water of lymph glands and veins.

Effects: redundant water is discharged from tissue fibers.

### **2. Vascular layer**

The suction and release of air pressure:

- (1) Improve the blood circulation of microvessels for the capillary system, and also promote the flow between the deeper veins and lymphatic networks.
- (2) Remove harmful substance in the body.

Effects: Strengthen blood vessels, eliminate toxins and relieve the condition of varicose veins.

### **3. Fibrous layer**

The suction and release of air pressure:

It can deliver two kinds of effects in different tissues (including skin, muscle and so on).

- (1) Repair cells and increase their activity.
- (2) Restore and improve skin elasticity.
- (3) Stimulate the production of bone collagen and improve skin plumpness.
- (4) Improve the oxygen absorption of skin and increase the consumption of carbon dioxide.

Effects: Break the hard fibers of cellulite, and make it more elastic so as to achieve the effect of body-shaping.

### **4. Nerve layer**


The suction and release of air pressure:

It stimulates the surface and deep layer of sympathetic nervous system.

- (1) Lower the sensitivity of skin.
- (2) Restore skin elasticity and help resist tissue fibrosis.

Effects: Repair cells and lower skin sensitivity.





## Face Firming&Anti-ageing

### 1. Handpiece

#### Face RF

### 2. Effects

1. Firm skin and smooth out fine lines
2. Relive saggy condition of skin
3. Build face contour and make face have stereoscopic appearance
4. Moisturize skin and promote skin's absorption
5. Accelerate metabolism and blood circulation
6. Stimulate collagen regeneration and slow down skin ageing

### 3. Applicable Range

1. Those whose faces are dark and dull.
2. Those whose skin is flabby and sagging.
3. Those who have fine lines, nasolabial folds , wrinkles at eye corners.
4. Those whose face contours are not clear.
5. Those who have rough skin and large pores.
6. Those who work for a long time in the office and are exposed to ultraviolet rays.
7. Those who have sagging skin, edema or puffiness after giving birth

### 4. Inapplicable Range:

1. Those who have just undergone plastic surgery with prostheses and metal materials, etc., implanted at the operation area.
2. Those who have used injection products recently, such as hyaluronic acid, skin booster, and wrinkle-removal products.
3. Those who have hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
4. Those who are in allergic period and have severely sensitive skin and metal allergies.
5. Those who have lesions on the skin.
6. Those who are of advanced years.
7. Those who are in menstruation, pregnancy, lactation, surgical recovery.
8. Those who have skin diseases and infectious diseases.
9. Those who have unrealistic expectations about the effect.

### 5. Do's and Don'ts after operation

1. Wash face with warm water within three days
2. Stay hydrated and use suncream
3. Do not use irritating products such as AHA, scrubbing and exfoliating products within one



week

4. Do not stay up late, avoid drinking alcohol and smoking cigarettes and eat more fruits, vegetables and foods with less grease.
5. Apply face mask at least three times within a week

## **Eye Anti-wrinkle**

### **1. Handpiece**

#### **Tripolar RF**

### **2. Effects**

1. Relieve eye fatigue, dark circles, eye bags and edema.
2. Reduce lines around eyes and crow's feet.
3. Reduce dark spots and accelerate blood circulation.
4. Accelerate metabolism and prevent pigmentation.
5. Supplement eyes with nutrition and moisturize skin.
6. Firm and soften skin, lift the corners of the eyes.
7. Accelerate blood circulation around eyes and help skin absorb nutrients efficiently.

### **3. Applicable Range**

1. Those with wrinkles, fine lines, bags and dark circles around the eyes.
2. Those with dry skin, dry lines and dynamic wrinkles.
3. Those with eyes prone to fatigue and dryness.
4. Those who frequently face computer screens and mobile phones.
5. Those who often stay up late and have dark circles.
6. Those who often stay in a dry or hot environment.

### **4. Inapplicable Range**

1. Those who had just undergone plastic surgery.
2. Those who have hypertension, heart disease, diabetes, severe hyperthyroidism and malignant tumors.
3. Those who are having an allergic reaction and have severe sensitive skin or are allergic to metals.
4. Those who have skin lesions.
5. Those who are of advanced years.
6. Those who are pregnant or convalescents.
7. Those who have skin diseases and contagious diseases.
8. Those who are sensitive to electric currents

### **5. Do's and Don'ts after Operation**

1. Avoid excessive exposure to the sun. Suncream is encouraged.
2. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
3. Avoid facial washing with overheated water, taking a sauna, going for hot springs or other violent exercise within 7 days after operation.
4. Moisturizers and eye masks are recommended to use and it will be preferable to use eye masks 3 times or more a week because the treated area is relatively dry.

# Neck Care

## 1. Handpiece

Face RF

## 2. Effects

1. Reduce the flabbiness, roughness and dullness of the neck skin.
2. Firm skin and increase skin elasticity.
3. Remove double chin.
4. Accelerate lymph detox and improve facial skin.
5. Prevent neck and lymph diseases.

## 3. Applicable Range

1. Those with rough and dull skin on the neck.
2. Those with lymphatic blocking
3. Those with sagging and inelastic neck skin.
4. Those who often bow their heads
5. Those who are not satisfied with their skin colour.

## 4. Inapplicable Range

1. Those who have just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin.
4. Those with skin lesions.
5. Those who are of advanced years.
6. Those who are pregnant or convalescents.
7. Those with skin disease and other contagious diseases.

## 5. Do's and Don'ts after Operation

1. Avoid sunburn. Keep the neck warm.
2. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
3. Avoid facial washing with overheated water, taking a sauna, go for hot springs or violent exercise within 7 days after operation.
4. Moisturizers and eye masks are recommended to use and it will be preferable to use eye masks 3 times or more a week because the treated area is relatively dry. Essence or neck cream is encouraged.
5. Avoid bowing head for a long time.

# Body Management

## (1). Waist&Abdomen Shaping

Handpieces: Body RF+40K+Vacuum Bipolar RF

## Effects

1. Relieve cold hands and cold feet, cold womb or cold body of women.
2. Tighten the skin on the waist and abdomen.
3. Reduce lumbar and abdominal fat.
4. Tighten skin, relieve stretch marks, obesity lines and increase skin elasticity.
5. Accelerate metabolism, relieve constipation and improve intestinal peristalsis.

## Applicable Range

1. Those with cold feet, cold hands, cold womb and cold body.
2. Those with lumbar and abdominal fat or who have sagging skin after birth.
3. Those who often sit for a long time, or have unsatisfactory waistlines.
4. Those with stretch marks.
5. Those with constipation or obstruction of abdominal meridians.

## Inapplicable Range

1. Those with "3 Hs"(hypertension, hyperglycemia, hyperlipidemia) and who have heart disease.
2. Those in pregnancy, menstruation or lactation.
3. Those whose incisions are healing or who are on the road to recovery.
4. Those with epilepsy and severe diabetes and hyperthyroidism.
5. Those with malignant tumors, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those with severe gynecological diseases.
8. Those whose gynecological diseases are being treated.

## Do's and don'ts after Operation

1. Keep abdomen warm. Avoid getting a cold.
2. Avoid eating and drinking too much. No alcohol, no spicy and greasy food. Staying up late should be avoided. Drink warm water.
3. Take a bath after 4-6 hours.
4. Sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.
5. Before going to bed at night, you can rub your abdomen clockwise with your hands, so as to help with metabolism for weight loss.

## (2). Arm Shaping

### Handpieces: Body RF+40K+Vacuum Bipolar RF

## Effects

1. Stimulate collagen growth in the bottom of the skin to shape and firm it.
2. Relieve the condition of saggy skin.
3. Relieve the condition of flabby arms and thick arms.
4. Tighten skin.
5. Accelerate blood circulation and dredge meridians and collaterals.

## Applicable Range

1. Those with thick arms and want to look better in clothes.
2. Those with flabby arms.

3. Those with sagging arm skin.
4. Those whose arms are prone to pain and numbness.
5. Those who often hold baby or children

### **Inapplicable Range**

1. Those who have just had plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
3. Those who are during allergic period and has severely sensitive skin.
4. Those with skin trauma or cut.
5. Those who are of advanced years.
6. Those who are pregnant and who are recovering from surgery.
7. Those with skin diseases and infectious diseases.

### **Do's and Don'ts after Operation**

1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours.
3. Drink plenty of warm water to stay hydrated.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise within one week after operation.

### **(3) . Back Shaping**

Hip is located in the middle of human body, and are the key hub of qi and blood circulation of meridians. It is the main switch of six meridians and also the bridge connecting the qi and blood circulation of upper and lower parts of human body, which plays an important role in building S-shaped curve and promote female attractiveness.

### **Handpieces: Body RF+Vacuum Bipolar RF**

#### **Effects**

1. Relieve shoulder and back soreness, and improve the lump on back of the neck.
2. Dredge channels and collaterals, and improve channels and collaterals blocking.
3. Increase blood circulation and metabolism.
4. Improve blood supply to the head and sleep.
5. Regulate the functions of viscera and strengthen the body.
6. Firm skin and prevent skin from sagging.
7. Reduce excess fat on the back and shape back.

#### **Applicable Range**

1. Those with sore shoulders and backs and stiff necks.
2. Those who suffer from insomnia, dream a lot or have a worse memory.
3. Those who are prone to fatigue, drowsiness and poor circulation of qi and blood.
4. Those with thick back and want to look better in clothes.
5. Those who have a lump on the back of the neck.

#### **Inapplicable Range**

1. Those with metal implants in the body, such as stents and pacemakers, are allergic to metal.
2. Those are during pregnancy, menstruation and lactation.

1. Those whose surgical wound is healing or convalescence
2. Those who has "three Hs" , heart disease, epilepsy and severe diabetes and hyperthyroidism.
3. Those with malignant tumor and hemophilia or severe bleeding.
4. Those with skin diseases and infectious diseases.
5. Those who are too weak.
6. Those who are drunk, thirsty, overworked, on a full or empty stomach.

## **Do's and Don'ts after Operation**

1. Keep warm, avoid getting a cold and drink plenty of hot water.
2. Take shower in 4-6 hours after operation.
3. Avoid staying up late, drinking and overeating.
4. Avoid raw, cold and spicy foods and get enough sleep.
5. Avoid shoulder-baring and backless clothes

## **(4). Hip Shaping**

### **Handpieces: Body RF+Vacuum&RF**

#### **Effects**

1. Improve blood circulation and speed up metabolism.
2. Improve menstrual pain, irregular menstruation, abnormal leucorrhea and other female diseases.
3. Improve sleep quality, improve female sexual function and tighten the vagina.
4. Maintain the normal function of ovary, stimulate the secretion of glands, and make couple's relationship closer.
5. Make complexion ruddy, fade color spots and return to young state.
6. Shape hips, improve sagging and outward expansion of hips, tighten skin and increase elasticity.

#### **Applicable Range**

1. Those with saggy hips and fat accumulation.
2. Those with stretch marks induced by obesity and obesity.
3. Those whose butt shape is not good-looking, flat and soft outward expansion.
4. Those with cold and cool hips with low hip temperature.
5. Those with cold womb, dysmenorrhea, irregular menstruation, gynecological inflammation and other problems.
6. Those with decreased estrogen levels and poor sex lives.

#### **Inapplicable Range**

1. Those who are in menstruation, pregnancy, lactation and operation recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc
3. Those with skin diseases, patients with infectious diseases and skin sensitive period.
4. Those with wounds during recovery from surgery.
5. Those with allergic and severely sensitive skin.
6. Those who have just had liposuction.
7. Those who are of advanced years.

## **Do's and Don'ts after Operation**

1. Keep hips warm and avoid wearing miniskirts and shorts.

2. Take a shower in 4-6 hours after operation.
3. Drink plenty of hot water to avoid getting a cold.
4. Avoid staying up late, drinking and overeating.
5. Avoid eating raw, cold and spicy foods and get enough sleep.
6. Sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.

## **(5). Leg Shaping**

### **Handpieces: Body RF+40K+Vacuum Bipolar RF**

#### **Effects**

1. Tighten skin and prevent sagging.
2. Stimulate collagen regeneration and improve fat lines.
3. Increase leg circulation and detoxification, metabolism.
4. Promote blood circulation and remove blood stasis, dredge channels and collaterals, and prevent varicose veins.
5. Tighten excessive protruding fat on legs and get rid of fat on thick thighs.

#### **Applicable Range**

1. Those with poor circulation of lower limbs, edema and obese people.
2. Those who have low immunity and feel discomfort and pain all over the body are prone to colds.
3. Those with constipation and has coarse, flabby skin.
4. Those whose channels and collaterals of the legs are blocked, and the proportion of the legs is poor.

#### **Inapplicable Range**

1. Those who are in menstruation, pregnancy or lactation.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc
3. Those with skin diseases or infectious diseases.
4. Those who have incisions on their bodies.
5. Patients with severe varicose veins and tumors.
6. Those who are having allergic reactions or who have severely sensitive skin.
7. People who have just had liposuction.
8. Those who are of advanced years.
9. Pregnant women or who are on the road to recovery.

#### **Do's and Don'ts after Operation**

1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours later.
3. Drink more warm water to keep hydrated and speed up metabolism.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
6. Wear pants after operation. Avoid wearing miniskirts and shorts.

## **(6). Breast Shaping**

## Handpieces: Body RF+40K+Vacuum&RF

### Effects

1. Correct breast shape and relieve accessory breasts
2. Relieve nodules on breasts, slight hyperplasia and distending pain in breast
3. Relieve outward expansion of breasts
4. Relieve irregular menstruation, speckles on the face and the lack of elasticity of breasts.
5. Relive breast atrophy, breast sagging and obstruction of lactiferous ducts

### Applicable Range

1. Those who have accessory breast and unsatisfactory breast shape
2. Those who have nodules on breasts, slight hyperplasia and distending pain in breast during the menstruation
3. Those who have free fat, saggy breast and outward expansion of breasts
4. Those who have low immunity
5. Those who have irregular menstruation, speckles on the face and breast inelasticity
6. Those who think that their breasts are underdeveloped, or who have postpartum breast atrophy, breast relaxation, blockage of lactiferous ducts and so on.

### Inapplicable Range

1. Patients with cardiovascular and cerebrovascular diseases, diabetes, abnormal blood coagulation and impaired function of important organs.
2. Those who are taking anticoagulants, vasodilators, corticosteroids and other drugs or who have been taking them for a long time.
3. Those with a focus of infection on the skin of the chest.
4. Patients with severe hyperplasia, fibroma and cyst of the chest
5. Those who are in pregnancy and lactation

### Do's and Don'ts after Operation


1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours later.
3. Drink more warm water to keep hydrated and speed up metabolism.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
6. Wear shaping and comfortable bra and don't squeeze your chest.



# Part II

## 1. Detailed Operation

After checking all are connected properly, the following interface will appear when power switch is pressed



### 1. Function Selection



time display and adjusting buttons



button for decreasing time



button for increasing time



exit button




stop button



start button

### 2. Detailed Operation of 40K

Choose  to go to the following interface





energy level and adjusting button



button for lowering energy level



button for raising energy level



mode selection



A: constant working mode, the handpiece releases energy constantly after start button is pressed. Suitable for those with thick fat accumulation.

B: discontinuous working mode, the handpiece releases energy discontinuously after start button is pressed. Suitable for those who have unwanted fat on some areas and want to shape body.

### 3. Detailed operation of Vacuum Bipolar RF



Choose **Vacuum+Bipolar RF** and then go to the following interface



energy level and adjusting buttons



button for lowering energy level



button for raising energy level



suction and release time and adjusting buttons



button for increasing suction or release time




button for decreasing suction or release time




Knob for adjusting suction level of the Vacuum Bipolar RF handpiece.

Turning it clockwise is to increase suction and turning it anticlockwise is to decrease suction. It is recommended that you should not set the suction at too high level at the beginning in case there is discomfort.

#### 4. Detailed Operation of Tripolar RF



Choose to go the following interface



energy level and adjusting buttons



button for lowering energy level




button for raising energy level

#### 5. Detailed Operation of Body RF



Choose to go to the following interface



energy level and adjusting buttons





button for lowering energy level



button for raising energy level

## 6. Detailed operation of Face RF

Choose [Face RF](#) to go the following interface



energy level and adjusting buttons



button for lowering energy level



button for raising energy level

## 2. Technical Parameters

Power supply Input:100V-240V

Power: 380W

Unoisetion 40K Cavitation Probe

Frequency: 40KHz

Power: 30W

Multipolar RF Head for Body

Frequency: 1MHz

Power: 70W

Bipolar RF with Vacuum

Frequency: 1MHz

Power: 70W

Vacuum: <-80kPa

Pressure: >250kPa

Air flow: >10L/minute

Noise level: <70dB (30cm away)

Quadrupole RF Head for Face and Body

Frequency: 1MHz

Power: 60W

Tripolar RF Head for Face and Eyes' Around

Frequency: 1MHz

Power: 20W

Dimensions: 36x40x27cm

Net Weight: 5.4kgs

Gross Weight: 6.35kgs

Input Power: 100V-240V

Power Consumption: Max 150W

Display: Touch Screen

### 3. Dos and Dont's

1. The plug with grounding pin must be used before the instrument is used and the power socket of the instrument is properly grounded.
2. Ensure that the voltage of the instrument is suitable. If the voltage of the local power supply is not stable, we suggest that the user add a voltage regulator with matching power between the mains and the instrument.
3. In order to ensure the treatment effect and normal service life of the instrument, please use the specified parts provided or suggested by the original manufacturer.
4. The instrument should not be placed in damp places or near water sources, and should not be directly exposed to sunlight.
5. Do not place the instrument near a strong heat source, as this may affect the life and normal use of the instrument.
6. Before treatment, please remove all metal objects from the body first, so as to avoid unexpected conditions and influence the curative effect.
7. Please restrain from using the device on eyes, thyroid, parathyroid, testicles, abdomen of pregnant women, and pacemakers.
8. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
9. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use so as to ensure the safety of using the electrical products.

10. Set the energy at lowest level at the beginning and slowly raise it after the customer adapts to it.
11. The handle can be used alone. Better results will be delivered if one handle works in conjunction with other handles.
12. Make sure the handle is in full contact with skin during operation in case skin is heated unevenly.
13. The treatment area should be kept moisturized during the operation and do not perform the treatment on dry skin.
14. When the treatment receiver is taking other weight-loss drugs, it is recommended to stop taking them for 1 to 2 months before treatment. Extension of the course of treatment is recommended if the treatment receiver want to use this machine immediately.
15. 40K handle should not be used on head, chest, breast, heart and back.
16. Clean the instrument with normal saline after operation to ensure its cleanliness to prolong its service life.

## 4. Troubleshooting & Solutions

**1. The instrument cannot be started, and the button light on the back of the instrument does not work?**

- A. Make sure the power cord is connected to the socket with normal power supply.
- B. Whether the fuse on the back of the instrument is loose or blown.

**2. The machine can be started but there is no or weakened RF energy output?**

- A: Please check whether the handpiece is connected to the machine properly.
- B: The wires in the handpiece are loose or burnt out.

**3. No suction or very little suction?**

- A. Please check whether the handpiece is properly connected to the machine.
- B. Please check whether the container of the filter is screwed tight and rubber rings are worn down. This case can occur if the area is not airtight.
- C. If methods above can not solve such problem, please contact machine distributor for help.

**4. I can start the instrument, but the screen shows erroneous message?**

- A. Please unplug the power cord behind the instrument and wait for about 1 minute before restarting the instrument.
- B. If the above method cannot solve the problem, please contact the instrument dealer for assistance.

## 5. FAQs

**1. Q: How long can I start seeing results of RF treatment?**

A: Normally, you can see the treatment results that day when you take the operation or within one week. Skin collagen tissue is heated and then contracts, and you can obviously feel the sensation of skin being tightened. Because RF works to stimulate subcutaneous tissue to promote collagen regeneration, so the more you take the treatment, the more noticeable the effects will be.

**2. Q: Is RF harmful to skin?**

A: RF treatment is non-surgical operation. It stimulates collagen regeneration at the skin bottom and boosts metabolism, and therefore has no damage to the skin. Local redness and hotness are

the normal conditions of accelerated blood circulation, which will be relieved after a moment. So there's no need to worry about that.

**3. Q: How long does the body-firming treatment take?**

A: One treatment session is 60 minutes. We will combine professional techniques and instruments to achieve remarkable results.

**4. Q: What functions does this instrument have?**

A: It has the functions of dissolving fat and shaping body, smoothing&rejuvenating skin, soothing&repairing skin, firming &lifting skin and anti-aging. It can be operated on every part of the body. In terms of facial care, it can form face contour, alleviate wrinkles and saggy skin. In terms of body care, it can reduce regional fatness and build S-shaped curve, and promote metabolism and detoxification of the whole body. Meanwhile, it can also enhance the functions of viscera and solve sub-health problems.

**5. Q: Which one is better in terms of weight loss, liposuction or this machine?**

A: The principle of liposuction is to suck out the excess fat in one part of the body through negative pressure suction so as to achieve the purpose of rapid weight loss. It's a surgical treatment requiring anesthesia and so there is recovery period, which has risks despite its immediate results. However, our instrument is risk-free without side effects. It can not only help you reduce fat but also stimulate collagen regeneration to lift and tighten skin. It can also build perfect curves by lifting skin, and accentuate women's sexual charm.

**6. Q: Will I regain weight after operation?**

A: Ultrasound works to dissolve fat to reduce weight, and users are not prone to regain their weight. Because what's dissolved through ultrasonic treatment is fat, not water. Fat needs a longer period of time to accumulate, so normally people will not regain weight.

**7. Q: Do I need to be on a diet?**

A: Pay attention to what you eat. Because radio frequency and ultrasonic treatment work to dissolve fat and accelerate metabolism. It's best to avoid eating spicy, oily, or fried food in case they affect regular metabolism. Proper workout for sweating helps with fat metabolism process. In such way, the results will be more noticeable.

**8. Q: Does it have side effects on the body?**

A: RF beauty treatment is non-invasive and is one of the most effective and the safest approach to remove wrinkles and shape body. Generally speaking, side effects won't appear. A small number of people may experience transient redness or swelling that will disappear after a few hours. Those with dry skin may experience atrophy after treatment, which will disappear after three days because collagen regeneration will make skin plump up though RF heat can cause moisture loss of skin. Overall, there are no side effects on people's health.

**9. Q: Can ultrasound treatment be performed on every part of the body?**

A: Ultrasound treatment works to rupture fat cells, by which fat is released from cells and then metabolized through lymphatic system. So its power is very high. Heart is sensitive to high frequency sound waves since it is a hollow organ. After being vibrated by sound waves, the heart will reflect the waves due to the different transmission of sound waves by the muscle tissue and the blood in the heart. This force can pull valves away from the heart muscle. If you use the instrument on eyes, it can lead to retinal detachment, so the ultrasound treatment should not be performed on the areas of eyes and heart. (not on the areas of waist, back and chest)



**10. Q: Does ultrasound treatment have side effects on human body?**

A: Ultrasonic beauty treatment is non-surgical and non-invasive program without anaesthesia. It relies on ultrasound (mechanical sound wave) - dissolve fat through cavitation - ultrasound focusing. And the results are fat smash - fat dissolving - skin firming and body sculpting respectively. It only targets low-density fat cells but leaves high-density tissue like blood vessels intact, and therefore has no side effects on human body. Slight tinnitus will probably occur during the treatment, which is within normal range. So there's no need to worry about that.

**11. Q: Why does tinnitus occur?**

A: It's because ultrasound comes with a high frequency vibration, which could go up to 20kHz. It penetrates into 20mm-deep skin fat layers through ultrasound focusing. Highly intensified ultrasound can cause high-speed friction among fat cells, which can heat and dissolve them. And during the process, sound may be produced so that customer will have the sensation of tinnitus.

**12. Q: What is collagen?**

A: Collagen is a kind of biological macro molecule substance, a kind of white, opaque and non-branched fibrous protein. It can supplement the nutrients needed by all layers of the skin, enhance the collagen activity in the skin, retain moisture, nourish skin, slow down ageing, relieve facial relaxation, provide nutrients for the hair and so on. Collagen is a nutrient that human body must replenish to delay ageing. With the growth of age, collagen will gradually be lost. For women, collagen will begin to decrease since the age of 20. At the age of 25, collagen loss will reach its peak. At the age of 40, the collagen content is less than half of that one has at age of 18. It is the loss of collagen and water that causes the wrinkles on the face of the elderly. The breakage of collagen fibers and elastic fibers will cause skin tissue to be oxidized, atrophied and collapsed. The skin will be dry, wrinkled, flabby and inelastic. Therefore, in order to delay ageing, there must be plenty of collagen.

**13. Q: How does vacuum detoxify?**

A: Vacuum can cause capillary congestion, stimulate cell to increase vitality. The suction and release of pressure during operation can make local pores continue to open and close, which can promote skin breathing, increase the amount of skin oxygen, and speed up waste removal. So it achieves the effects of promoting qi and blood circulation, dispelling cold and damp, dredging meridians and collaterals, clearing toxins, dissipating heat and so on.

## **6. Packing List**

- Multipolar RF Head for Body x 1
- 40KHz Cavitation Head x 1
- Vacuum Bipolar RF Head x1
- Three Pole RF Head for Face and eyes' around x 1
- Quadrupole RF Head for Face and Body x1
- Power cord x 1

## **7. Simplified Treatment Steps**

### **Facial skin firming&anti-aging/eye wrinkle removal/neck care**

Remove makeup--clean the target area--apply toner--massage techniques--Face RF/Tripolar RF--apply face mask, eye mask or neck mask--apply water, essence, face cream(or eye cream and neck cream) and suncream

Needed skin products: RF cream/massage cream

### Waist&abdomen shaping/arm shaping/leg shaping

Massage techniques--Body RF—40K—Vacuum Bipolar RF—wipe the target area clean

Needed skin products: gel, RF cream/essential oil

### Back shaping/butt shaping

Massage techniques—Body RF—Vacuum Bipolar RF—wipe the target area clean




Needed skin products:RF cream/essential oil



### Breast-shaping

Massage techniques—Vacuum Bipolar RF—wipe the target area clean

Needed skin products:RF cream/essential oil

## 8. Operational Diagrams






| Parameter Adjustment                                                            | Product                                                                                          | Techniques                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | Diagrams                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Face Firming and Anti-ageing: 60 min once a week</b></p>                  |                                                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                                                                                                                                                                                                                                                                                                             |
| <p><b>Face RF</b><br/>Recommended energy level: 3~7<br/>Time: 10-15 minutes</p> | <p>Makeup Remover+Face Cleanser+Cold&amp;Hot Steam+Massage Cream+Base Mask+Essence+Face Mask</p> | <ol style="list-style-type: none"> <li>1. Remove makeup and clean face, 5 minutes.</li> <li>2. Apply water mask and hot steam, 10 minutes</li> <li>3. Clean the face, 2 minutes.</li> <li>4. Apply toner , 1 minute</li> <li>5. Apply massage cream evenly to face and soothe face, 3 times.</li> <li>6. Apply digital pressing technique on the following acupoints (RN-24, DU-26, ST-4, ST-6, LI-20, BL-1, BL-2, EM3, TE-23, the temple, GB-1, ST-1, ST-2), 3 times.</li> <li>7. Soothe face, 3 times</li> <li>8. Alternately apply lifting technique from chin to ear lobe, mouth corner to Ermen(SJ21), and nose wing to temple, lift eye corner from lower eyelid to temple with both hands, 3 times.</li> <li>9. Repeat the technique on the other side.</li> </ol> | <p>Technique5, 7</p>  <hr/> <p>Technique6</p>  <hr/> <p>Technique8, 17</p>  <hr/> <p>Technique10, 12, 19</p> |

|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                                  |
|--|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | <p>10. Apply lifting technique towards hairline on forehead, 3 times</p> <p>11. Apply plucking technique with fingers alternately to lift face.</p> <p>12. Apply plucking technique with one hand towards hairline on forehead.</p> <p>13. One side, apply plucking technique with both hands alternately to lift face, 2-3 minutes</p> <p>14. Repeat the technique on the other side.</p> <p>15. Move in a zigzag motion on forehead with middle and ring fingers, slide to the front side of ears and lift back and forth for 3 times and then slide to the back of ears.</p> <p>16. Clean face, 2 minutes.</p> <p>17. Apply essence evenly to face, 1 minute.</p> <p>18. <b>Face RF</b> operation: lift from chin to earlobe, mouth corner to ear center and nose wing to temple line by line, 3 times.</p> <p>19. Repeat the techniques above on the other side.</p> <p>20. Move the handpiece towards hairline on forehead, 3 times.</p> <p>21. Clean face, 2 minutes.</p> <p>22. Apply face mask, 15 minutes.</p> <p>23. Clean face, 2 minutes.</p> <p>24. Apply toner, essence, face cream and suncream.</p> <p>25. The end.</p> |  <hr/> <p>Technique15</p>  |
|--|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**Treatment Suggestions**

One course of treatment consists of ten treatment sessions. After one treatment session, skin will be firmed and rejuvenated as the nutrients it needs are absorbed. After one course of treatment, skin will be lifted with glow. Two courses will make skin smooth and plump, leaving face contour clearer. Three courses will promote skin metabolism and detox processes, help fade pigmentation and skin darkness, slow down skin ageing and prevent skin from sagging, and stimulate collagen regeneration to leave skin smooth and firm with glow.







## Eye Anti-wrinkle: 25 min 2-3 times a week

|                                                                                     |                                                                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Tripolar RF</b><br/>Recommended energy level: 3~7<br/>Time: 10-15 minutes</p> | <p>Make-up remover+face cleanser+essence/eye cream+eye mask+MS-5481</p> | <ol style="list-style-type: none"> <li>1. Remove makeup and clean face, 5 minutes</li> <li>2. Apply toner, 1 minute</li> <li>3. Apply essence(eye cream) evenly to eye area in circular motion, 3 minutes.</li> <li>4. <b>Tripolar RF</b> operation: adjust energy level, mode and time: 10 minutes.</li> <li>5. Move the instrument to lift skin from lower eyelid to eye corner, 3-6 times.</li> <li>6. Move the instrument to lift skin from lower eyelid to temple, 3-6 times.</li> <li>7. Move the instrument to lift skin from lower eyelid to temple in small circles, 3-6 times.</li> <li>8. Move the instrument to lift skin from lower eyelid to temple, 3-6 times.</li> <li>9. Lift from brow bone to hairline.</li> <li>10. Repeat the operation on the other side.</li> <li>11. Wash eye area clean, 1 minute</li> <li>12. Apply eye mask, 15 minutes</li> <li>13. Remove mask and wash the area clean, 2 minutes</li> <li>14. Apply eye essence and eye cream.</li> <li>15. The end.</li> </ol> | <p>Technique3, 8</p>  <p>Technique5, 11</p>  <p>Technique10</p>  <p>Technique12</p>  <p>Technique14</p>  |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### Treatment Suggestions

It is recommended to take this treatment 2-3 times per week. After one treatment, eye area is lifted and firmed and blood circulation is accelerated. After one month, fine wrinkles and black circles are reduced and skin colour is brightened. After three months, skin in eye area is totally firmed and rejuvenated with glossiness. If the customers keep using this machine, eye aging can be relieved and prevented.

## Neck Care: 60 min once a week





|                                                                                   |                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Face RF</b><br/> Recommended energy level: 3~7<br/> Time: 10-15 minutes</p> | <p>Makeup Remover+Face Cleanser+Massage Cream+Essence+Neck Mask</p> | <ol style="list-style-type: none"> <li>1. Remove makeup and clean neck, 5 min.</li> <li>2. Apply toner, 1 min.</li> <li>3. Neck massage, rub oil into chest with both hands towards the back of the neck and press Gb-20 and Du-16 acupoint, 3 times.</li> <li>4. Alternately lift the lower jaw and double chin area to pass the back of ear to the armpit, 3 times.</li> <li>5. Push downwards with thumb and index finger splaying, tracing 3 meridians from the neck side to the armpit, 3 times.</li> <li>6. Push downwards with 4 fingers, tracing 3 meridians from neck side to the armpit, 3 times.</li> <li>7. Push downwards with bending fingers till the neck turns hot and slide to the armpit, 3 times.</li> <li>8. Push downwards with thumb and index finger splaying from external collarbone to armpit, 3-5 times.</li> <li>9. Do the same on the other side.</li> <li>10. The end</li> <li>11. <b>RF</b> Operation, use one hand and the handpiece to lift double chin area and then pass the back of ear to armpit, 3 times.</li> <li>12. Move the handpiece in circles along three lines on one side of neck to armpit, 3 times.</li> <li>13. Slide the handpiece from inner and outer sides of clavicle to the armpit, 3-5 times.</li> <li>14. Move the handpiece in circles on the whole neck area, 3 times.</li> <li>15. Do the same on the other side.</li> </ol> | <p>Technique3</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5, 6, 7</p>  <hr/> <p>Technique11</p>  <hr/> <p>Technique12</p>  <hr/> <p>Technique14</p>  |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|  |  |                                                                                                                                                                                                                                                                          |  |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|  |  | <p>16. Wipe neck clean with warm towel, 2 minutes</p> <p>17. Apply neck mask, 15 minutes</p> <p>18. Wipe neck clean with warm towel, and put the warm towel on neck for 5 minutes.</p> <p>19. Apply toner, neck essence, neck cream and suncream</p> <p>20. The end.</p> |  |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|

**Treatment Suggestions**



A course of treatment consists of ten treatment sessions. After one treatment, skin becomes rejuvenated and smooth. After a full course, neck skin becomes smooth, and the neck lymphatic circulation is accelerated. After 2 courses of treatment, skin will be tightened, neck lines will be reduced, and the the condition of double chin will be relieved. After 3 courses of treatment, skin will be firmed with glow, lymphatic detox will be accelerated and some problems will be resolved such as acne and skin dullness. Moreover, skin aging will be prevented and skin will remain in a young state.

**Waist&Abdomen Shaping: 60 min once a week**

|                                                                                                                                                                                                                                                                                       |                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
| <p><b>Body RF</b><br/>Recommended energy level: 3~7<br/>Time: 10-15 minutes</p> <p><b>40K</b><br/>Energy level:3~7<br/>Time: 15-20 min<br/>Mode: A/B</p> <p><b>Vacuum Bipolar RF</b><br/>Energy level:3-7<br/>Time: 10-15 minutes<br/>Suction time and release time can be freely</p> | <p>Massage Cream(Essential Oil)+Gel+MS-5481</p> | <ol style="list-style-type: none"> <li>1. Rub oil into abdomen with hands moving in circular motion, 3 times.</li> <li>2. Rub abdomen back and forth with both hands, 3-5 times.</li> <li>3. Knead abdomen with both hands using chiropractic technique, 3 times.</li> <li>4. Lift Meridian BV(Belt Vessel) on both sides of waist with both hands alternately, 16 times.</li> <li>5. Move hands in 8-shaped motion with right hand and left hand moving clockwise and anti-clockwise respectively on waist, 3 times.</li> <li>6. Overlap hands and message the intestinal canal clockwise, 3 times.</li> <li>7. Move hands in circular motion to soothe the treatment area, 3 times.</li> </ol> | <p>Technique1、 7、 10</p>  |
|                                                                                                                                                                                                                                                                                       |                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | <p>Technique2</p>         |
|                                                                                                                                                                                                                                                                                       |                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | <p>Technique3</p>         |
|                                                                                                                                                                                                                                                                                       |                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | <p>Technique4</p>         |

|                                                                                              |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|----------------------------------------------------------------------------------------------|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>adjusted.<br/>When release time is zero, the handpiece is at Constant Suction status.</p> |  | <p>8. Apply digital pressing technique on the following acupoints: RN13, RN12, RN10, RN8, RN6, RN4, RN3, ST25, SP15. 2 times.</p> <p>9. Starting from RN3 point, push to belly navel, slide to waist and then lift upwards to groin with both thumbs, 3 times.</p> <p>10. Soothe the treatment area with hands and then slide to groin.</p> <p>11. <b>Body RF</b> operation: move the handpiece slowly along ascending colon, transverse colon and descending colon to rectum,3-5 times.</p> <p>12. Move the handpiece in rhombus-shaped motion on abdomen, 3-5 times</p> <p>13. Move the handpiece up and down on waist side, 3-5 times.</p> <p>14. Use the handpiece to lift Meridian BV(Belt Vessel) on both sides of waist 16 times.</p> <p>15. Use one hand and the handpiece to lift Meridian BV(Belt Vessel) on both sides of waist 16 times.</p> <p>16. <b>40K</b> operation: move the handpiece slowly along ascending colon, transverse colon and descending colon to rectum,3-5 times.</p> <p>17. Move the handpiece in rhombus-shaped motion on abdomen, 3-5 times</p> <p>18. Move the handpiece up and down on waist side, 3-5 times.</p> <p>19. Use the handpiece to lift Meridian BV(Belt Vessel) on both sides of waist 16 times.</p> <p>20. Use one hand and the handpiece to lift Meridian BV(Belt Vessel) on both sides of waist 16 times.</p> <p>21. <b>Vacuum Bipolar RF</b> operation: starting from one</p> | <div data-bbox="1182 197 1318 405"> </div> <p>Technique5</p> <div data-bbox="1182 461 1318 669"> </div> <p>Technique6</p> <div data-bbox="1182 725 1318 934"> </div> <p>Technique8</p> <div data-bbox="1182 990 1318 1198"> </div> <p>Technique9</p> <div data-bbox="1182 1254 1318 1462"> </div> <p>Technique11、 16</p> <div data-bbox="1182 1518 1318 1727"> </div> <p>Technique12、 17</p> <div data-bbox="1182 1805 1318 2036"> </div> <p>Technique13、 18、 24</p> |
|----------------------------------------------------------------------------------------------|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|










|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                            |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | <p>side of the waist, lift the handpiece to abdomen and groin, 3 times.</p> <p>22. Repeat the previous operation on the other side, 3 times.</p> <p>23. Suction&amp;release mode, move the handpiece slowly on abdomen, 3 times.</p> <p>24. Constant suction mode, move the handpiece up and down to shape waist, 3-5 times.</p> <p>25. Use one hand and the handpiece to lift Meridian BV(Belt Vessel) on both sides of waist 16 times.</p> <p>26. Wipe the area clean with towel, the end.</p> |  <p>Technique14、 15、 19、 20、 25</p>  |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment session, metabolism will be accelerated. After one course, fat will be slowly metabolized, body will start to slim down. Constipation and other issues will be relieved at the same time. After two courses, the effects will gradually become obvious, skin will be tightened and lifted, and skin collagen will be increased. Stretch marks and flabby skin are also alleviated. Three courses will help consolidate the effects. Unwanted fat will gradually disappear, and waist and abdomen curves will start to take shape.

**Arm Shaping: 70 min once a week**

|                                                                                                                                                                                                                               |                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                     |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Body RF</b><br/>Recommended energy level: 3~7<br/>Time: 10-15 minutes</p> <p><b>40K</b><br/>Energy level:3~7<br/>Time: 15-20 min<br/>Mode: A/B</p> <p><b>Vacuum Bipolar RF</b><br/>Energy level:3-7<br/>Time: 10-15</p> | <p>Massage Cream(Essential Oil)+Gel+MS-5481</p> | <ol style="list-style-type: none"> <li>1. Left-right order: lay the arm flat, rub oil into skin from lower arm to the entire arm and slide to the fingers(soothing technique), 3 times.</li> <li>2. Push the entire arm with both palms, 3 times.</li> <li>3. Push hands upwards tracing three meridians: Large Intestine Meridian(LI) - Triple Energizer Meridian(TE) - Small Intestine Meridian(SI) to armpits with thumbs and index fingers splaying, 3 times.</li> <li>4. Soothe the area, 3 times</li> <li>5. Rub three meridians on arms with bending fingers</li> </ol> | <p>Technique1</p>  <p>Technique2, 3, 4, 5, 6</p>  <p>Technique7</p>  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|                                                                                                                                                                                                      |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>minutes<br/>Suction<br/>time and<br/>release time<br/>can be<br/>freely<br/>adjusted.<br/>When<br/>release time<br/>is zero, the<br/>handpiece<br/>is at<br/>Constant<br/>Suction<br/>status.</p> |  | <p>back and forth to warm them up, 3 times.<br/>6. Soothe the area, 3 times<br/>7. Lay inner arm upwards, and push hands tracing three yin meridians on inner arm: Lung Meridian(LU) - Pericardium Meridian(PC) - Heart Meridian(HT) to armpit with thumbs and index fingers splaying, 3 times respectively.<br/>8. Rub three meridians back and forth with bending fingers to warm them up, 3 times.<br/>9. Soothe the area and slide to fingers.<br/>10. Wipe arms clean with towel, and then apply gel evenly to the entire arm.<br/>11. Repeat above techniques on the other arm.<br/>12. <b>Body RF</b> operation: lay customer's arm flat, push from fat part of lower arm along three meridians to armpit, 5-8 times.<br/>13. Move the handpiece in spiral curve on upper arm towards armpit, 5-8 times.<br/>14. Repeat the operation on flabby area, 3 times.<br/>15. Move the handpiece along three meridians to armpit, 3 times.<br/>16. Let inner side of arm face upwards and push three meridians on upper arm to armpit, 3 times.<br/>17. Move the handpiece in spiral curve on fat area of upper arm towards armpit, 3 times.<br/>18. Move the handpiece back and forth on the fat area of upper arm, 5-8 times.<br/>19. Repeat the techniques above on the other arm.<br/>20. <b>40K</b> operation: lay</p> | <p>Technique8, 16, 18, 24, 26</p>  <p>Technique9</p>  <p>Technique12, 15, 20, 23</p>  <p>Technique13, 14, 21, 22, 32</p>  |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|







|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |  |
|--|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|  |  | <p>customer's arm flat, push from fat part of lower arm along three meridians to armpit, 5-8 times.</p> <p>21. Move the handpiece in spiral curve on upper arm towards armpit, 5-8 times.</p> <p>22. Repeat the operation on flabby area, 3 times.</p> <p>23. Move the handpiece along three meridians to armpit, 3 times.</p> <p>24. Let inner side of arm face upwards and push three meridians on upper arm to armpit, 3 times.</p> <p>25. Move the handpiece in spiral curve on fat area of upper arm towards armpit, 3 times.</p> <p>26. Move the handpiece back and forth on the fat area of upper arm, 5-8 times.</p> <p>27. Repeat the techniques above on the other arm.</p> <p>28. <b>Vacuum Bipolar RF</b> operation.</p> <p>29. let customer lay her arm flat, move the handpiece back and forth on upper arm, 5-8 times.</p> <p>30. Move the handpiece in spiral curve on upper arm, 5-8 times.</p> <p>31. Let the inner side of arm face upwards, move the handpiece on the inner side of upper arm up and down, 5-8 times.</p> <p>32. Move the handpiece in small spiral curve on the the inner side of upper arm, 5-8 times.</p> <p>33. Wipe the area clean with towel, the end.</p> |  |
|--|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|







**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment session, it will feel warm and tight, and the arms will be relaxed. After one course of treatment, excessive fat on arms will be reduced and the skin will appear to be tight and plump. Two courses of treatment will help shape arms, dredge meridians, and enhance

constitution to strengthen the effects. 3 courses will help consolidate effects and prevent arms from regaining fat.

## Back Shaping: 60 min once a week

|                                                                                                                                                                                                                                                                                |                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Body RF</b><br/>Energy level: 3~7<br/>Time: 15-20 min</p> <p><b>Vacuum Bipolar RF</b><br/>Energy level:3-7<br/>Time: 15-20 minutes<br/>Suction time and release time can be freely adjusted. When release time is zero, the handpiece is at Constant Suction status.</p> | <p>Massage Cream(Essential Oil)+Gel+MS-5481</p> | <ol style="list-style-type: none"> <li>1. Massage techniques</li> <li>2. Apply oil to back and press Gb-20 and Du-16 acupoints.</li> <li>3. Pluck Dabanjing(the area connecting neck and shoulder) (from hairline), 3-5 times</li> <li>4. Move thumb outwards from Bladder Meridian(BL) to the sacral region and then to Gb-20 and Du-16 points.</li> <li>5. Starting from neck to caudal vertebra, move hands in S-shape, 3 times.</li> <li>6. Push Bladder Meridian(BL) with thumb in left-right order, 3 times.</li> <li>7. Push Bladder Meridian(BL) to the sacral region with both thumbs at the same time, 3 times.</li> <li>8. Push Bladder Meridian(BL) in three bending fingers, 3 times.</li> <li>9. Push along the medial border of scapula with both hands alternately in left-right order, 3-6 times.</li> <li>10. Push medial border of scapula with both hands horizontally, 3-6 times.</li> <li>11. Soothe the back with both hands and push Gb-20 and Du-16 points, 3 times.</li> <li>12. Overlap both thumbs and push SI-11 point and slide to the arm, 3 times.</li> <li>13. Do the same on the other side,3 times</li> <li>14. Rub Meridian GV and Bladder Meridian(BL) with hands to warm them up.</li> <li>15. The end</li> <li>16. <b>Body RF</b> handpiece</li> </ol> | <p>Technique2, 11</p>  <hr/> <p>Technique 3, 18, 25</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5</p>  <hr/> <p>Technique6, 7, 8, 14, 16, 23</p>  <hr/> <p>Technique 9, 19, 26</p>  <hr/> <p>Technique10</p> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|




|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | <p>operation, move the handpiece along Meridian GV and Bladder Meridian(BL) respectively from neck to the sacral region, 3-5 times.</p> <p>17. Move the handpiece in circles on Du-14 point and sacral region, 3-5 times respectively.</p> <p>18. From neck, move the handpiece on Dabanjing(the area connecting neck and shoulder) back and forth, 3-5 times.</p> <p>19. Move the handpiece to lift the medial border of scapula back and forth in left-right order, 3-5 times.</p> <p>20. Move the handpiece in 8-shaped motion downwards to sacral region, 3 times</p> <p>21. Starting from the left side, move the handpiece in annular motion upwards to armpit, 3 times.</p> <p>22. Apply lifting technique by moving the handpiece from one side of the waist upwards to armpit, 3-5 times.</p> <p>23. <b>Vacuum Bipolar RF</b> handpiece operation, move the handpiece along Meridian GV and Bladder Meridian(BL) respectively from neck to the sacral region, 3-5 times.</p> <p>24. Move the handpiece in circles on Du-14 point and sacral region, 3-5 times respectively.</p> <p>25. From neck, move the handpiece on Dabanjing(the area connecting neck and shoulder) back and forth, 3-5 times.</p> <p>26. Move the handpiece to lift the medial border of scapula back and forth in</p> |  <hr/> <p>Technique12, 13</p>  <hr/> <p>Technique17, 24</p>  <hr/> <p>Technique20, 27</p>  <hr/> <p>Technique21, 28</p>  <hr/> <p>Technique22, 29</p>  <hr/> |
|--|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|





|  |  |                                                                                                                                                                                                                                                                                                                                                                                                      |  |
|--|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|  |  | <p>left-right order, 3-5 times.</p> <p>27. Move the handpiece in 8-shaped motion downwards to sacral region, 3 times</p> <p>28. Starting from the left side, move the handpiece in annular motion upwards to armpit, 3 times.</p> <p>29. Apply lifting technique by moving the handpiece from one side of the waist upwards to armpit, 3-5 times.</p> <p>30. Wipe back clean with towel, the end</p> |  |
|--|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|

**Treatment suggestions**

A course of treatment consists of 10 treatment sessions. After one treatment session, back will be significantly relaxed, neck stiffness will also be alleviated, and skin will be lifted. A course will help relax shoulder, reduce back fat, relieve the lump on the nape of the back, and shape the back line. After two courses, the lump will appear to be smaller, meridians will be dredged, sub-health problems will be solved, and sleep quality will be improved. 3 courses will help with back thinning, back line refinement, strengthening visceral function, improving digestion, detoxification and metabolism(one course of treatment is to reduce fat, two courses are to strengthen effects, and three courses help consolidate effect and prevent the reappearance of above condition.)

**Hip Shaping: 60 min once a week**

|                                                                                                                                                                                                                                                                                   |                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                                            |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Body RF</b><br/>Energy level:3~7<br/>Time: 15-20 min</p> <p><b>Vacuum Bipolar RF</b><br/>Energy level:3-7<br/>Time: 15-20 minutes<br/>Suction time and release time can be freely adjusted.<br/>When release time is zero, the handpiece is at Constant Suction status.</p> | <p>Massage<br/>Cream(Essential Oil)+Gel+MS-5481</p> | <ol style="list-style-type: none"> <li>1. Massage techniques.</li> <li>2. Standing on the side, rub oil into skin by moving both hands clockwise and anticlockwise respectively from sacral region on hips , 3 times (this is soothing technique.)</li> <li>3. Apply pushing technique on the sacral region with two thumbs, 3 times.</li> <li>4. Soothe the treatment area for 3 times and then apply digital pressing technique on the following points: Bl-23, Baliao points(BL31-34),</li> </ol> | <p>Technique2, 5, 11</p>  <hr/> <p>Technique3</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique6, 14, 17</p> |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | <p>Du-1, Gb-30, BL-36, 3 times.</p> <p>5. Soothe the treatment area, 3 times.</p> <p>6. Left-right order, apply pushing technique from the top of thigh along the four meridians -- Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) to Meridian BV respectively with both hands, 3 times each.</p> <p>7. Alternately apply pushing technique with both hands along Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) from the top of thigh to Meridian BV respectively, 3 times each.</p> <p>8.Overlap both palms and then apply pushing technique from the top of thigh to Meridian BV (lifting), 3 to 5 times.</p> <p>9. Apply pushing technique from both sides of hips to the top of hips(shaping) back and forth, 3 times.</p> <p>10. Repeat technique 7</p> <p>11. Soothe the treatment area.</p> <p>12. Repeat the techniques above on the other side.</p> <p>13. The end.</p> <p>14. Body RF operation: lift from the top of thigh to Meridian BV line by line, 3 times.</p> |  <p>Technique7, 10</p>  <p>Technique9, 15, 18</p>  <p>Technique16</p>  |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|











|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |  |
|--|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|  |  | <p>15. Lift upwards from both sides of the buttocks to the highest point of the buttocks 3 times.</p> <p>16. Move the handpiece in small spiral curve on hip, 3-5 times</p> <p>17. Lift from the top of thigh to Meridian BV line by line, 3 times.</p> <p>18. Lift upwards from both sides of the buttocks to the highest point of the buttocks 3 times.</p> <p>19. Do the same on the other side.</p> <p>20. The operation of Vacuum Bipolar RF is the same with Body RF.</p> <p>21. Wipe the area clean with towel, the end.</p> |  |
|--|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|

**Treatment Suggestions**

One course of treatment consists of 10 treatment sessions. After one treatment, buttocks will experience the sensations of being heated and lifted. After one course of treatment, buttocks start to show obvious curve lines, excessive fat begin to disappear slowly, and skin appear to be tight. After 3 courses, the conditions of cold body and uterus, other gynecological problems such as irregular menstruation are relieved, which can help consolidate the effect to make females more attractive.

**Leg Shaping: 100 min once a week**

|                                                                                                                                                                      |                                                 |                                                                                                                                                                                                                                                                                                                                   |                                                                                                                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------|
| <p><b>Body RF</b><br/>Energy level:3~7<br/>Time: 10-15 min</p> <p><b>40K</b><br/>Energy level:3~7<br/>Time: 15-20 min<br/>Mode: A/B</p> <p><b>Vacuum Bipolar</b></p> | <p>Massage Cream(Essential Oil)+Gel+MS-5481</p> | <p>1. Customer lies on his front, back side of the leg, left-right order: rub oil into the skin from calf to thigh and then move back to heel, 3 times</p> <p>2. Alternate the heels of both palms to push the entire leg upwards and move back to heel, 3 times</p> <p>3. Move hands upwards to push four meridians: Bladder</p> | <p>Technique1、 2、 3、 4、 7</p>  |
|                                                                                                                                                                      |                                                 |                                                                                                                                                                                                                                                                                                                                   | <p>Technique5</p>              |

|                                                                                                                                                                                                         |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>RF</b><br/> Energy level:3-7<br/> Time: 10-15 minutes<br/> Suction time and release time can be freely adjusted.<br/> When release time is zero, the handpiece is at Constant Suction status.</p> |  | <p>Meridian(BL) - Kidney<br/> Meridian(KI) - Liver<br/> Meridian(LV) - Gallbladder<br/> Meridian(GB) with thumb and index finger splaying.</p> <p>4. Push popliteal fossa with both hands alternately, 3 times.<br/> 5. Move both hands upwards alternately in undulating motion, 3 times.<br/> 6. Soothe the area, 3 times<br/> 7. Push four meridians upwards with bending fingers of both hands, 3 times.<br/> 8. Soothe the area, 3 times<br/> 9. Repeat the operation on the other side.<br/> 10. Wipe the treatment area clean with hot towel, 2 minutes.<br/> 11. Apply gel evenly to the treatment area, 2 minutes.<br/> 12. <b>Body RF</b> operation: push Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) upwards to popliteal fossa respectively, 3 times.<br/> 13. Move in small spiral curve on calf to dissolve fat, 3 times<br/> 14. Move the instrument upwards to popliteal fossa and soothe the treatment area with hand, 3 times.<br/> 15. Push 4 meridians</p> | <p>Technique12、 18</p>  <p>Technique13、 19</p>  <p>Technique14、 18、 24</p>  <p>Technique15、 17、 21、 23、 25</p>  <p>Technique16、 22、 26</p>  <p>Technique27</p>  <p>Technique28、 29、 30、 31、 32</p> |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | <p>from popliteal fossa to the top of thigh respectively, 3 times.</p> <p>16. Move the handpiece in small spiral curve from popliteal fossa to the top of thigh to dissolve fat, 3 times.</p> <p>17. Push meridians from popliteal fossa to the top of thigh, 3 times.</p> <p>18. <b>40K</b> operation: push Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) upwards to popliteal fossa respectively, 3 times.</p> <p>19. Move in small spiral curve on calf to dissolve fat, 3 times</p> <p>20. Move the instrument upwards to popliteal fossa and soothe the treatment area with hand, 3 times.</p> <p>21. Push 4 meridians from popliteal fossa to the top of thigh respectively, 3 times.</p> <p>22. Move the handpiece in small spiral curve from popliteal fossa to the top of thigh to dissolve fat, 3 times.</p> <p>23. Push meridians from popliteal fossa to the top of thigh, 3 times</p> <p>24. <b>Vacuum&amp;Bipolar RF</b> operation: push Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder</p> | <div data-bbox="1209 203 1342 405" data-label="Image"> </div> <p data-bbox="1129 421 1362 495">Technique33、<br/>35、 40、 42、 47</p> <div data-bbox="1209 510 1342 712" data-label="Image"> </div> <p data-bbox="1129 728 1318 801">Technique36、<br/>43、 48</p> <div data-bbox="1209 817 1342 1019" data-label="Image"> </div> <p data-bbox="1129 1034 1318 1108">Technique37、<br/>44、 49</p> <div data-bbox="1209 1124 1342 1326" data-label="Image"> </div> <p data-bbox="1129 1344 1318 1417">Technique38、<br/>45、 50</p> <div data-bbox="1209 1433 1342 1635" data-label="Image"> </div> |
|--|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |  |
|--|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|  |  | <p>Meridian(GB) upwards to popliteal fossa respectively, 3 times.</p> <p>25. Push 4 meridians from popliteal fossa to the top of thigh respectively, 3 times.</p> <p>26. Move the handpiece in small spiral curve from popliteal fossa to the top of thigh to dissolve fat, 3 times.</p> <p>27. Use the handpiece and one hand to lift from both sides of thigh to the middle line by line, 3-5 times.</p> <p>28. The front of leg, customer lies on her back: rub oil into skin from foot to the top of thigh (soothing technique), 3 times.</p> <p>29. Alternate heels of both palms to push leg towards the top of thigh, 3 times</p> <p>30. Push four meridians: Spleen Meridian (SP) - Stomach Meridian (ST) - Liver Meridian(LV) - Gallbladder Meridian(GB) to the top of thigh with thumb and index finger splaying.</p> <p>31. Push four meridians with bending fingers of both hands, 3 times.</p> <p>32. <b>Body RF</b> operation: lift along the four meridians from calf to knee, 3 times(if there is not too much unwanted fat on calf, perform the operation from thigh)</p> <p>33. Move in annular</p> |  |
|--|--|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|


|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |  |
|--|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|  |  | <p>motion from the two sides of calf to knee, 3 times(if there is not too much unwanted fat on calf, perform the operation from thigh)</p> <p>34. Lift from knee to the top of thigh line by line, 3 times</p> <p>35. Move in small circles on thigh, 3 times.</p> <p>36. Lift from the two sides of thigh to the middle with one hand and the handpiece, 3 times</p> <p>37. Lift from knee to the top of thigh in annular motion, 3 times.</p> <p>38. Lift from knee to the top of thigh line by line, 3 times.</p> <p>39. <b>40K</b> operation: lift along the four meridians from calf to knee, 3 times(if there is not too much unwanted fat on calf, perform the operation from thigh)</p> <p>40. Move in annular motion from the two sides of calf to knee, 3 times(if there is not too much unwanted fat on calf, perform the operation from thigh)</p> <p>41. Lift from knee to the top of thigh line by line, 3 times</p> <p>42. Move in small circles on thigh, 3 times.</p> <p>43. Lift from the two sides of thigh to the middle with one hand and the handpiece, 3 times</p> <p>44. Lift from knee to</p> |  |
|--|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|








|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |  |
|--|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
|  |  | <p>the top of thigh in annular motion, 3 times.</p> <p>45. Lift from knee to the top of thigh line by line, 3 times.</p> <p>46. <b>Vacuum&amp;Bipolar RF</b> operation: lift from knee to the top of thigh line by line, 3 times.</p> <p>47. Move in small circles on thigh, 3 times.</p> <p>48. Lift from the two sides of thigh to the middle with one hand and the handpiece, 3 times</p> <p>49. Lift from knee to the top of thigh in annular motion, 3 times.</p> <p>50. Lift by moving the handpiece from knee to the top of thigh line by line, 3 times.</p> <p>51. Do the same on the other side.</p> <p>52. Wipe legs clean with towel, the end.</p> |  |
|--|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|


**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment, leg will be relaxed, and circulation will be accelerated. After one course of treatment, leg will become thin, there will be contraction of skin collagen due to heat effect, and tightening effect can obviously be felt on the skin. After 2 courses of treatment, there will be obvious results with the slimness and firmness of legs promoted. 3 courses will help consolidate the effects. This treatment can stimulate dermis to promote constant collagen regeneration. So the treatment results will be more obvious.

**Breast Shaping: 60 min once a week**

|                                                                  |                                             |                                                                                                                                                                    |                                                                                                         |
|------------------------------------------------------------------|---------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| <p><b>Body RF:</b><br/>Energy level: 3~7<br/>Time: 15-20 min</p> | <p>Massage Cream(Essential Oil)+MS-5481</p> | <p>1. Standing at the bedside, use both hands to rub oil into skin from Shanzhong(RN17) to armpit and then lift the suspensory ligament (massage) for 3 times.</p> | <p>Technique1</p>  |
|                                                                  |                                             |                                                                                                                                                                    | <p>Technique2</p>                                                                                       |

|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | <p>2. Use both thumbs to press acupoints: Shanzhong(RN17), Rugen(ST18), Dabao(SP21), Yingchuang(ST16), Zhongfu(LU1), Yunmen(LU2) 3 times</p> <p>3. Massage the area, 3 times</p> <p>4. Overlap two palms and move from Shanzhong(RN17) in a 8-shaped motion back and forth, 3 times</p> <p>5. Alternate hands to push from accessory breast towards suspensory ligament, in left-right order</p> <p>6. Massage the area, 3 times</p> <p>7. Sit beside the customer, use both palms to massage and lift breasts in a left-right order, 3-5 times</p> <p>8. Alternately push the lactiferous ducts (around the breast) with thumbs and index fingers splaying 3-5 times</p> <p>9. Alternately move in circles with two thumbs to dredge nodules on the chest, 3 times</p> <p>10. Massage and lift breast with both hands</p> <p>11. Repeat the techniques on the other breast.</p> <p>12. <b>Body RF</b> handpiece operation: left-right order, lift towards nipple around breast in conjunction with soothing technique of the other hand, 5-8</p> |  <p>Technique4</p>  <p>Technique5</p>  <p>Technique7, 13</p>  <p>Technique8, 12</p>  <p>Technique9, 14</p>  <p>Technique15</p>  <p>Technique16</p> |
|--|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|  |  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                                                                                                                                                     |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | <p>times</p> <p>13. Lift in circles on the breast towards collarbone in conjunction with soothing technique of the other hand, 5-8 times</p> <p>14. Move in small circles to dredge the area where nodules exist, 3-5 times</p> <p>15. Treatment receiver lies on her side and lift up her arm, move the handpiece in circular motion on accessory breast to burn fat, 5-8 times.</p> <p>16. Push the handpiece from the accessory breast in the armpit to breast (shaping and narrowing the accessory breast) 5-8 times</p> <p>17. Repeat the techniques on the other side.</p> <p>18. Wipe breasts clean, the end.</p> |  <p>The diagram shows a handpiece with a circular tip being applied to a breast. Three red arrows point upwards from the handpiece towards the breast, indicating the direction of the treatment or the area being targeted.</p> |
|--|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**Treatment Suggestions**

One course of treatment consists of ten treatment sessions. After one treatment, there will be sensation of heat on the breast, blood circulation will be promoted and there will be some degree of lifting effect. After one course, the effect will be more obvious and the nodules will be relieved. After three courses, the skin will be firmed and the breast will be well shaped with increased elasticity. The functions of endocrine system will be improved to make females more attractive.