

# User Manual

MS-97D1



## Preface

Dear Users,

We're pleased to present to you our latest 8-in-1 comprehensive beauty machine for skin and body care, which includes face RF handle, eye RF handle, body RF handle, vacuum RF handle, 40K handle, BIO handle, hot&cold handle and ultrasound handle. It is a professional beauty machine, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore, we advise any personnel to read this manual thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

# Table of Contents

Preface	
Part I	
Introduction	
Advantages	
Working Principles	
Face Firming and Anti-ageing	
Face Relaxing and Repairing	
Face Firming and Whitening	
Eye Anti-ageing	
Neck Care	
Body Management	
(1). Working Principles	
(2). Waist&Abdomen Shaping	
(3). Arm Shaping	
(4). Back Shaping	
(5). Hip Shaping	
(6). Leg Shaping	
(7). Hip Shaping	
Part II	
1. Detailed Operations	
2. Technical Parameters	
3. User Contraindications	
4. Dos and Don'ts	
5. Troubleshooting & Solutions	
6. FAQs	
7. Packing List	
8. Operational Diagrams	

# Part I

## Brief Introduction

Multi-polar RF handle, vacuum&RF handle, 40K ultrasound fat-burning handle, hot&cold handle, BIO handle and ultrasound handle are among the most popular tools for skin firming, body shaping and anti-ageing. So, this machine is an effective replacement of all the facial and body care programs in beauty salons. Moreover, it is convenient and easy to operate. It solves both skin and body management problems for beauty-seekers. This equipment requires no injection, no medication, and no operation, for which it has no side effects. It is used externally during the whole process and has an immediate effect. RF heats deeper skin layers, stimulates collagen to regenerate and reorganize in dermis, and promotes blood circulation in hypodermis, which can achieve the results of skin firming and anti-aging. Ultrasounds with specific frequencies target those parts in which there are stubborn fat to bring customers amazing results.

## Advantages

1. 8-in-1 beauty machine that works on facial and body care.
2. The ultrasonic wave can powerfully burn and remove fat, which surpasses liposuction to solve the problems for obese people.
3. Various working heads can be chosen to be used in different parts according to different needs.
4. One working head has various functions. The RF has more powerful and even energy release.
5. The perfect combination of ultrasound and radio frequency solves major problems in losing weight. It helps with weight loss and meanwhile firms skin, increases skin elasticity, promotes metabolism, and ultimately strengthen people's constitutions.
6. It adopts the advanced technology of ultrasonic wave fat-burning.
7. It is painless and non-invasive treatment process. Once the treatment is over, customer can go back home and there is no recovery period, which will not affect customer's life and work.
8. No consumption, low cost and quick returns.
9. Wider treatment range for body and face.
10. Unevenness, bleeding, and swelling will not appear after treatment.
11. Various handles have the energy release of red LED light, which can dissolve fat and massage body, and accelerate metabolism and sterilize skin at the same time. The effect is remarkable with a higher comfort level.

## Face Firming and Anti-ageing

### 1. Working Principles

#### Face RF


RF wave can directly penetrates the skin and uses the resistance formed by the skin to produce energy, which raises the temperature of the bottom layer of the skin. Using the principle of dermis collagen protein tightening and stimulating collagen protein regeneration, it has two therapeutic effects: timely skin firming and lasting collagen protein regeneration.

After treatment, an obvious sense of contraction can be felt immediately at the treatment area, and the effect of skin firming and lifting can be seen about a month later, and the effect will become more and more obvious. The effect of RF treatment can last for about 18 months to 3

years depending on the area of treatment and individual maintenance methods, which, however, may be longer depending on individual maintenance methods. The latest research points out that better results can be achieved after repeated treatment.

### Biological Effects

The waves vibrating millions of times per second can penetrate the epidermis and then directly act on the dermis to heat up collagen cells in a deep layer, and stimulate the growth of collagen fibers. In this way, it can make the skin firmer and more elastic. In the process of moving the handle, the RF waves have passed through the epidermis and acted on the deep dermis to produce safe biothermal energy rapidly. When collagen tissue is heated to  $45^{\circ}\text{C} \sim 60^{\circ}\text{C}$  in the natural friction, it will contract immediately and stimulate the continuous proliferation of collagen. At the same time, biological heat can effectively accelerate the blood flow in adipocytes and made free fatty acid (ffa) released to promote the dissolution of superficial fat to achieve the purpose of tightening and lifting saggy part of the body.





### BIO

It stimulates the energy contained in the ATPs of cells to restore the normal operation and function of cells after its imitating human biological current enters into muscle cells through skin to provide energy for muscle and body. It can improve cell permeability to make it 30%-40% higher than usual. It can help the body better absorb nutrition, water, oxygen etc., and eliminate wastes, toxins and retained liquids as soon as possible. It can promote the synthesis of muscle fibers to make it 68% higher than usual, keep the skin moisturized, soft, firm and shiny, and increase collagen synthesis by 73% to reduce wrinkles and make the skin softer, elastic and shiny.

## Skin facial lifting

Microcurrent therapy -acne treatment & skin firming

Microcurrent-safe and painless facial helps erase fine lines and wrinkles, while firming your skin and defining your features. It improves muscle tone, reduces puffiness, increases cellular activity, and tightens pores. the overall result is a healthier and younger looking skin, no matter your age.

## Ultrasound

It makes high-frequency ultrasonic wave of 28,000 HZ/ second penetrate into the subcutaneous tissue at a depth of 4~6mm. Ultrasonic wave has high energy and strong permeability, which can drive facial cells to vibrate accordingly. It can not only drive pore grease out, but also change the volume of cells to improve blood circulation, which can stimulate cell regeneration and make the skin shiny and elastic. Secondly, because of the heat effect caused by its high-frequency vibration, it also has the function of quickly introducing nutrition into the skin, which can make the skin better absorb the nutrients of the skin care products.

### 2. Effects:

1. Activate cells and increase skin's elasticity
2. Firm skin and smooth out fine lines
3. Relieve saggy condition of skin
4. Build face contour and make face have stereoscopic appearance
5. Shrink pores and plump up cheeks
6. Moisturize skin and promote skin's absorption
7. Accelerate metabolism and blood circulation
8. Stimulate collagen regeneration and slow down skin ageing

### 3. Applicable range

1. Those with saggy skin
2. Those with fine line on face, nasolabial folds and dynamic wrinkles
3. Those with unsatisfactory face contour and double chin
4. Those with rough skin and large pores
5. Those with saggy skin and edema after birth
6. Those with low skin metabolism and poor absorption

### 4. Inapplicable Range:

1. Those who have just undergone plastic surgery with prostheses and metal materials, etc., implanted at the operation area.
2. Those who have experience in injection products recently, such as hyaluronic acid, or those who have mesotherapy, wrinkle removal, or plastic surgery.
3. Those who have hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
4. Those who are in allergic period and have severely sensitive skin and metal allergies.
5. Those who have trauma or wound on the skin.

6. Those who are of advanced years.
7. Those who are in menstruation, pregnancy, lactation, surgical recovery.
8. Those who have skin diseases and infectious diseases.
9. Those who have unrealistic illusions about the effect.
10. Those who are sensitive to electric currents.

## 5. Do's and Don'ts after Operation

1. Wash face with warm water within three days
2. Stay hydrated and use suncream
3. Do not use irritating products such as AHA, scrubbing and exfoliating products within one week
4. Do not stay up late, avoid drinking alcohol and smoking cigarettes and eat more fruits, vegetables and foods with less grease.
5. Apply face mask at least three times within a week

# Face Soothing and Repairing

## 1. Working Principles

### Ultrasound

It makes high-frequency ultrasonic wave of 28,000 HZ/ second penetrate into the subcutaneous tissue at a depth of 4~6mm. Ultrasonic wave has high energy and strong permeability, which can drive facial cells to vibrate accordingly. It can not only drive pore grease out, but also change the volume of cells to improve blood circulation, which can stimulate cell regeneration and make the skin shiny and elastic. Secondly, because of the heat effect caused by its high-frequency vibration, it also has the function of quickly introducing nutrition into the skin, which can make the skin better absorb the nutrients of the skin care products.

### Cold&Hot Handle

It uses the principle of thermal expansion and cold contraction. Our skin reduces the metabolic rate at 5 degrees Celsius, slows down the movement of cells and slows down the sensitivity of nerve endings, which can inhibit pain and swelling of skin for allergic reactions. Due to thermal expansion and cold contraction, the temperature of 5 degrees Celsius causes the blood vessels to contract, resulting in pore contraction to restrain thick pores and firm the skin. Therefore, it can smooth facial wrinkles and crow's feet.

Physiological effects of Cold Hammer:


1. It can constrict blood vessels, resulting in tightening of pores, which can greatly improve the size of pores.
2. It can strengthen myofibril, make the skin firm and elastic, and smooth the wrinkles and crow's feet.
3. It can constrict blood vessels, and reduce varicose veins and red veins effectively. It is easier to purify and detoxify venous blood in terms of harmful substances, and it has a significant effect on reduction and elimination of spots.
4. It can reduce the metabolic rate, inhibit and slow down the cell activity, reduce the sensitivity of the peripheral nerve to play a calming and soothing effect, and relieve skin swelling and pain, especially for allergic skin.

Physiological Effects of Hot Hammer:

1. 40°C makes people feel warm and comfortable, which will be able to relax muscles, keep the balance of vegetative nervous system and relieve pain and soreness. It can promote blood circulation, increase the oxygen content in the blood, and open pores and extend blood vessels to make skin absorb nutrients quickly, which is a milder and more comfortable way for skin absorption.

3. It can increase the metabolic rate of the skin, increase the heat energy content of the tissue, and make skin naturally eliminate the waste and toxins, which has a certain effect on the elimination of bags under the eyes and dark circles.

The alternating use of cold and heat hammer has multiplying effect on skincare, and can activate cells, and improve and purify the skin, which is a good instrument for health and beauty treatment.



## 2. Effects

1. Reduce acne and relieve pimples
2. Diminish inflammation, kill bacteria and soothe skin
3. Remove dirt in pores and promote skin absorption
4. Regulate grease secretion, shrink pores and smooth skin
5. Promote skin metabolism
6. Relieve redness, swelling, acne and allergic reactions on face.

## 3. Applicable Range

1. Those with large pores and rough skin
2. Those with pimples and closed acne
3. Those with poor skin absorption and low metabolism
4. Those with excessive grease secretion
5. Those who are prone to pimples and acne.
6. Those who have sensitive skin, acne-prone skin or whose skin is vulnerable to redness, swelling, and allergies.

## 4. Inapplicable Range

1. Those who have had invasive surgery and the wounds have not healed
2. Those with skin diseases and infectious diseases
3. Those who have used injection products recently, such as hyaluronic acid, skin booster, and something that are used for reducing wrinkles or who have plastic surgery.
4. Those who have sensitive skin or are having allergic reactions
5. Those who are sensitive to current.





## Ultrasound

It makes high-frequency ultrasonic wave of 28,000 HZ/ second penetrate into the subcutaneous tissue at a depth of 4~6mm. Ultrasonic wave has high energy and strong permeability, which can drive facial cells to vibrate accordingly. It can not only drive pore grease out, but also change the volume of cells to improve blood circulation, which can stimulate cell regeneration and make the skin shiny and elastic. Secondly, because of the heat effect caused by its high-frequency vibration, it also has the function of quickly introducing nutrition into the skin, which can make the skin better absorb the nutrients of the skin care products.

## Cold&Hot Handle

It uses the principle of thermal expansion and cold contraction. Our skin reduces the metabolic rate at 5 degrees Celsius, slows down the movement of cells and slows down the sensitivity of nerve endings, which can inhibit pain and swelling of skin for allergic reactions. Due to thermal expansion and cold contraction, the temperature of 5 degrees Celsius causes the blood vessels to contract, resulting in pore contraction to restrain thick pores and firm the skin. Therefore, it can smooth facial wrinkles and crow's feet.

Physiological effects of Cold Hammer:

1. It can constrict blood vessels, resulting in tightening of pores, which can greatly improve the size of pores.
2. It can strengthen myofibril, make the skin firm and elastic, and smooth the wrinkles and crow's feet.
3. It can constrict blood vessels, and reduce varicose veins and red veins effectively. It is easier to purify and detoxify venous blood in terms of harmful substances, and it has a significant effect on reduction and elimination of spots.
4. It can reduce the metabolic rate, inhibit and slow down the cell activity, reduce the sensitivity of the peripheral nerve to play a calming and soothing effect, and relieve skin swelling and pain, especially for allergic skin.


Physiological Effects of Hot Hammer:

1. 40°C makes people feel warm and comfortable, which will be able to relax muscles, keep the balance of vegetative nervous system and relieve pain and soreness. It can promote blood circulation, increase the oxygen content in the blood, and open pores and extend blood vessels to make skin absorb nutrients quickly, which is a milder and more comfortable way for skin absorption.
3. It can increase the metabolic rate of the skin, increase the heat energy content of the tissue, and make skin naturally eliminate the waste and toxins, which has a certain effect on the elimination of bags under the eyes and dark circles.

The alternating use of cold and heat hammer has multiplying effect on skincare, and can activate cells, and improve and purify the skin, which is a good instrument for health and beauty treatment.

## Skin tightening

Use the cold therapy hammer to shrink pore after skin moisturizing and nutrition importation, this cooling therapy can promote skin tightening and enhance acne relieving and soothing



## 2. Effects

1. Firm skin and smooth out fine lines
2. Relive saggy condition of skin
3. Shrink pores and plump up cheeks
4. Moisturize skin and promote skin's absorption
5. Accelerate metabolism and blood circulation
6. Stimulate collagen regeneration and slow down skin ageing
7. Improve skin colour and make it have more clarity.

## 3. Applicable Range:

1. Those with saggy skin
2. Those with fine line on face, nasolabial folds and dynamic wrinkles
3. Those with dull skin and uneven skin colour
4. Those with sallow skin
5. Those with rough skin and large pores
6. Those with saggy skin and edema after birth
7. Those with low skin metabolism and poor absorption

## 4. Inapplicable Range:

1. Those who have just undergone plastic surgery with prostheses and metal materials, etc., implanted at the operation area.
2. Those who have experience in injection products recently, such as hyaluronic acid, or those who have mesotherapy, wrinkle removal, or plastic surgery.
3. Those who have hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
4. Those who are in allergic period and have severely sensitive skin and metal allergies.
5. Those who have trauma or wounds on the skin.
6. Those who are of advanced years.
7. Those who are in menstruation, pregnancy, lactation, surgical recovery.
8. Those who have skin diseases and infectious diseases.
9. Those who have unrealistic illusions about the effect.

## 5. Do's and Don'ts after operation

1. Wash face with warm water within three days
2. Stay hydrated and use suncream

3. Do not use irritating products such as AHA, scrubbing and exfoliating products within one week
4. Do not stay up late, avoid drinking alcohol and smoking cigarettes and eat more fruits, vegetables and foods with less grease.
5. Apply face mask at least three times within a week

## Eye Anti-wrinkle

### 1. Working Principles

#### RF

In the process of treating the eye, the RF instrument uses the radio frequency wave to make the bottom layer of the eye skin quickly heat up and promote the collagen growth of the skin, so as to achieve the effect of firming the skin around the eye, removing the bags under the eyes, removing dark circles, and relieving bloodshot eyes lines and reducing fine lines.




#### BIO

It stimulates the energy contained in the ATPs of cells to restore the normal operation and function of cells after its imitating human biological current enters into muscle cells through skin to provide energy for muscle and body. It can improve cell permeability to make it 30%-40% higher than usual. It can help the body better absorb nutrition, water, oxygenated blood, etc., and eliminate wastes, toxins and retained liquids as soon as possible. It can promote the synthesis of muscle fibers to make it 68% higher than usual, keep the skin moisturized, soft, firm and shiny, and increase collagen synthesis by 73% to reduce wrinkles and make the skin softer, elastic and shiny.

## Skin facial lifting

### Microcurrent therapy -acne treatment & skin firming

Microcurrent-safe and painless facial helps erase fine lines and wrinkles, while firming your skin and defining your features. It improves muscle tone, reduces puffiness, increases cellular activity, and tightens pores. The overall result is a healthier and younger looking skin, no matter your age.



## 2. Effects

1. Relieve eye fatigue, dark circles, eye bags and edema.
2. Reduce lines around eyes and crow's feet.
3. Reduce dark spots and accelerate blood circulation.
4. Accelerate metabolism and prevent pigmentation.
5. Supplement eyes with nutrition and moisturize skin.
6. Firm and soften skin, lift the corners of the eyes.
7. Accelerate blood circulation around eyes and help skin absorb nutrients efficiently.

## 3. Applicable Range

1. Those with wrinkles, fine lines, bags and dark circles around the eyes.
2. Those with dry skin, dry lines and dynamic wrinkles.
3. Those with eyes prone to fatigue and dryness.
4. Those who frequently face computer screens and mobile phones.
5. Those who often stay up late and have dark circles.
6. Those who often stay in a dry or hot environment.

## 4. Inapplicable Range

1. Those who had just undergone plastic surgery.
2. Those who have hypertension, heart disease, diabetes, severe hyperthyroidism and malignant tumors.
3. Those who are having an allergic reaction and have severe sensitive skin or are allergic to metals.
4. Those who have skin trauma or wound.
5. Those who are of advanced years.
6. Those who are pregnant or convalescents.
7. Those who have skin diseases and contagious diseases.
8. Those who are sensitive to electric currents

## 5. Do's and Don'ts after Operation

1. Avoid excessive exposure to the sun. Sun protection is encouraged.
2. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.

3. Avoid facial washing with overheated water, taking a sauna, going for hot springs or other violent exercise within 7 days after operation.
4. Moisturizers and eye masks are encouraged to use. It will be preferable to use eye masks 3 times or more a week.

## Neck Care

### 1. Working Principles

#### RF

RF waves directly penetrate the skin, making use of the resistance formed by the skin to produce energy to raise the temperature of the bottom layer of the skin, heat the tissue quickly and continuously and promote the growth of collagen, which achieves the effect of firming, lifting and smoothing skin.



It makes high-frequency ultrasonic wave of 28000 HZ/ second penetrate into the subcutaneous tissue at a depth of 4~6mm. Ultrasonic wave has high energy and strong permeability, which can drive facial cells to vibrate accordingly. It can not only drive pore grease out, but also change the volume of cells to improve blood circulation, which can stimulate cell regeneration and make the skin shiny and elastic. Secondly, because of the heat effect caused by its high-frequency vibration, it also has the function of quickly introducing nutrition into the skin, which can make the skin better absorb the nutrients of the skin care products.

### 2. Effects

1. Reduce the flabbiness, roughness and dullness of the neck skin.
2. Firm skin and increase skin elasticity.
3. Remove double chin.
4. Accelerate lymph detox and improve facial skin.
5. Prevent neck and lymph diseases.

### 3. Applicable Range

1. Those with rough and dull skin on the neck.
2. Those with lymphatic blocking
3. Those with sagging and inelastic neck skin.
4. Those who often bow their heads
5. Those who are not satisfied with their skin colour.

## 4. Inapplicable Range

1. Those who have just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin.
4. Those with skin trauma or wound.
5. Those who are of advanced years.
6. Those who are pregnant or convalescents.
7. Those with skin disease and other contagious diseases.

## 5. Do's and Don'ts after Operation


1. Avoid sunburn. Keep the neck warm.
2. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
3. Avoid facial washing with overheated water, taking a sauna, go for hot springs or violent exercise within 7 days after operation.
4. Moisturizing and mask for the neck is encouraged, Use masks at least 3 times a week. Essence or neck cream is encouraged.
5. Avoid bowing head for a long time.

# Body Management

## (1) Working Principles

### Body RF

RF waves directly penetrate the skin, making use of the resistance formed by the skin to produce energy to raise the temperature of the bottom layer of the skin, heat the tissue quickly and continuously and promote the growth of collagen, which achieves the effect of firming, lifting and smoothing skin. The perfect combination of 40K and RF can not only help people lose weight, but also increase skin elasticity and firm skin.



### 40K

The principle of ultrasonic wave is used to make human adipose cells take a strong impact and produce friction movement among adipose cells after strong sound waves enter into the human body, which can effectively consume heat and water of cells and shrink adipose cells. In addition, when the sound wave vibrates, it can produce a strong impact on cells. And then the adipose

cells will burst instantly, and they will be reduced so as to achieve the effect of fat removal.


**Ultrasound Cavitation:** Tens of thousands of tiny bubbles, namely cavitation bubbles, are produced by vibration of liquid. These bubbles grow in the negative pressure region formed by the longitudinal transmission of ultrasound, which will close rapidly in the positive pressure region to be compressed and stretched under alternating positive and negative pressures. When the bubbles burst, there will be huge instantaneous pressure of up to tens of MPa to hundreds of MPa, and strong vibration and noise will be produced.

**Advantages:** It only aims for low-density adipose tissue for its specific frequency, and protects high-density tissue such as vascular and nerve tissue.

## 02

### Dissolve fat cells- ultrasound unoisetion cavitation 2.0

Unoisetion cavitation 2.0 vibration: Ultrasound waves penetrate into the fat layer of body to produce cavitation effect, effectively making fat cells burst and emulsify into glyceride



**Experiment:** Put a piece of pork ointment on the fat of the unoisetion cavitation 2.0 probe, it was dissolved gradually by ultrasonic vibration.



## Vacuum&RF

Massaging skin and muscles with a special vacuum suction handle can effectively improve the humoral fluidity of human cells, and increase the movement of cells to improve skin elasticity. Meanwhile, it can accelerate the blood circulation of microvessels, expel toxins from the body through the normal circulation of the lymphatic system, and reduce the probability of occurrence of adverse conditions such as spots, pigmentation and congestion in various parts of the body. The kneading sensation generated by vacuum suction can increase the tissue activity of skin and muscle, and therefore help to reduce tough tissue and increase the elasticity of skin tissue, so that both body slimming and shaping can be realized at the same time. In addition, the movement of vacuum cup can also stimulate sympathetic nervous system to lower skin sensitivity. The suction and release process of air pressure can not only improve capillary system but also promote the flow between deeper veins and lymphoid gland network to strengthen blood vessels and alleviate the condition of varicose veins.

**Advantages:** it is superior to the previous vacuum working mode which is single and simple. Using different modes for different parts can be more efficient and more effective for slimming and shaping body, and it will never harm human body. Uniquely designed RF tips integrate RF and vacuum as a whole, in which RF and vacuum can exist independently and also are complementary to each other. Compared to the ordinary one-polar RF, it can ideally shape body in a more effective, faster and even way



## Vacuum Physical Effects

### 1. Skin layer

The suction and release of air pressure:

- (1) It can improve the fluidity of cells, thus increasing the movement of cells, and can treat diseases related to blood stasis and blood stagnation.
- (2) It can drain excess water of lymph glands and veins.

Effects: redundant water is discharged from tissue fibers.

### 2. Vascular layer

The suction and release of air pressure:

- (1) Improve the blood circulation of microvessels for the capillary system, and also promote the flow between the deeper veins and lymphatic networks.
- (2) Remove harmful substance in the body.

Effects: Strengthen blood vessels, eliminate toxins and relieve the condition of varicose veins.

### 3. Fibrous layer

The suction and release of air pressure:

It can deliver two kinds of effects in different tissues (including skin, muscle and so on).

- (1) Repair cells and increase their activity.
- (2) Restore and improve skin elasticity.
- (3) Stimulate the production of bone collagen and improve skin plumpness.
- (4) Improve the oxygen absorption of skin and increase the consumption of carbon dioxide.

Effects: Break the hard fibers of cellulite, and make it more elastic so as to achieve the effect of body-shaping.

### 4. Nerve layer

The suction and release of air pressure:

It stimulates the surface and deep layer of sympathetic nervous system.

- (1) Lower the sensitivity of skin.
- (2) Restore skin elasticity and help resist tissue fibrosis.

Effects: Repair cells and lower skin sensitivity.



## (2) Waist&Abdomen Shaping Effects

1. Relieve cold hands and cold feet, cold womb or cold body of women.

2. Tighten the skin on the waist and abdomen.
3. Reduce lumbar and abdominal fat.
4. Tighten skin, relieve stretch marks, obesity lines and increase skin elasticity.
5. Accelerate metabolism, relieve constipation and improve intestinal peristalsis.

### **Applicable Range**

1. Those with cold feet, cold hands, cold womb and cold body.
2. Those with lumbar and abdominal fat or who have sagging skin after birth.
3. Those who often sit for a long time, or have unsatisfactory waistlines.
4. Those with stretch marks.
5. Those with constipation or obstruction of abdominal meridians.

### **Inapplicable Range**

1. Those with "3 Hs"(hypertension, hyperglycemia, hyperlipidemia) and who have heart disease.
2. Those in pregnancy, menstruation or lactation.
3. Those whose incisions are healing or who are on the road to recovery.
4. Those with epilepsy and severe diabetes and hyperthyroidism.
5. Those with malignant tumors, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those with severe gynecological diseases.
8. Those whose gynecological diseases are being treated.

### **Do's and don'ts after Operation**

1. Keep abdomen warm. Avoid getting a cold.
2. Avoid eating and drinking too much. No alcohol, no spicy and greasy food. Staying up late should be avoided. Drink warm water.
3. Take a bath after 4-6 hours.
4. Sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.
5. Before going to bed at night, you can rub your abdomen clockwise with your hands, so as to help with metabolism for weight loss.

## **(2) Arm Shaping**

### **Effects**

1. Stimulate collagen growth in the bottom of the skin to shape and firm it.
2. Improve sagging skin.
3. Improve flabby arms and thick arms.
4. Soothe and firm sagging skin.
5. Tighten skin.
6. Accelerate blood circulation and dredge meridians and collaterals.

### **Applicable Range**

1. Those with thick arms and want to look better in clothes.
2. Those with flabby arms.
3. Those with sagging arm skin.
4. Those whose arms are prone to pain and numbness.
5. Those who often hold baby or children

### **Inapplicable Range**

1. Those who have just had plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
3. Those who are during allergic period and has severely sensitive skin.
4. Those with skin trauma or cut.
5. Those who are of advanced years.
6. Those who are pregnant and who are recovering from surgery.
7. Those with skin diseases and infectious diseases.

### **Do's and Don'ts after Operation**

1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours.
3. Drink plenty of warm water to stay hydrated.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise within one week after operation.

### **(3) Back Shaping**

#### **Effects**

1. Relieve shoulder and back soreness, and improve the lump on back of the neck.
2. Dredge channels and collaterals, and improve channels and collaterals blocking.
3. Increase blood circulation and metabolism.
4. Improve blood supply to the head and sleep.
5. Regulate the functions of viscera and strengthen the body.
6. Firm skin and prevent skin from sagging.
7. Reduce excess fat on the back and shape back.

#### **Applicable Range**

1. Those with sore shoulders and backs and stiff necks.
2. Those who suffer from insomnia, dream a lot or have a worse memory.
3. Those who are prone to fatigue, drowsiness and poor circulation of qi and blood.
4. Those with thick back and want to look better in clothes.
5. Those who have a lump on the back of the neck.

#### **Inapplicable Range**

1. Those with metal implants in the body, such as stents and pacemakers, are allergic to metal.
2. Those are during pregnancy, menstruation and lactation.
1. Those whose surgical wound is healing or convalescence
2. Those who has "three Hs", heart disease, epilepsy and severe diabetes and hyperthyroidism.
3. Those with malignant tumor and hemophilia or severe bleeding.
4. Those with skin diseases and infectious diseases.
5. Those who are too weak.
6. Those who are drunk, thirsty, overworked, on a full or empty stomach.

### **Do's and Don'ts after Operation**

1. Keep warm, avoid getting a cold and drink plenty of hot water.
2. Take shower in 4-6 hours after operation.
3. Avoid staying up late, drinking and overeating.
4. Avoid raw, cold and spicy foods and get enough sleep.
5. Avoid shoulder-baring and backless clothes

## **(4) Hip Shaping**

Hip is located in the middle of human body, and are the key hub of qi and blood circulation of meridians. It is the main switch of six meridians and also the bridge connecting the qi and blood circulation of upper and lower parts of human body, which plays an important role in building S-shaped curve and promote female attractiveness.

### **Effects**

1. Improve blood circulation and speed up metabolism.
2. Improve menstrual pain, irregular menstruation, abnormal leucorrhea and other female diseases.
3. Improve sleep quality, improve female sexual function and tighten the vagina.
4. Activate the warm nest function, stimulate the secretion of glands, and increase the feelings between husband and wife.
5. Make complexion ruddy, fade color spots and return to young state.
6. Shape hips, improve sagging and outward expansion of hips, tighten skin and increase elasticity.

### **Applicable Range**

1. Those with saggy hips and fat accumulation.
2. Those with stretch marks induced by obesity and obesity.
3. Those whose butt shape is not good-looking, flat and soft outward expansion.
4. Those with cold and cool hips with low hip temperature.
5. Those with cold womb, dysmenorrhea, irregular menstruation, gynecological inflammation and other problems.
6. Those with decreased estrogen levels and poor sex lives.

### **Inapplicable Range**

1. Those who are in menstruation, pregnancy, lactation and operation recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc
3. Those with skin diseases, patients with infectious diseases and skin sensitive period.
4. Those with wounds during recovery from surgery.
5. Those with allergic and severely sensitive skin.
6. Those who have just had liposuction.
7. Those who are of advanced years.

### **Do's and Don'ts after Operation**

1. Keep hips warm and avoid wearing miniskirts and shorts.
2. Take a shower in 4-6 hours after operation.
3. Drink plenty of hot water to avoid getting a cold.
4. Avoid staying up late, drinking and overeating.
5. Avoid eating raw, cold and spicy foods and get enough sleep.
6. Sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.

## **(5) Leg Shaping**

### **Effects**

1. Tighten skin and prevent sagging.
2. Stimulate collagen regeneration and improve fat lines.
3. Increase leg circulation and detoxification, metabolism.

4. Promote blood circulation and remove blood stasis, dredge channels and collaterals, and prevent varicose veins.
5. Tighten excessive protruding fat on legs and get rid of fat on thick thighs.

### **Applicable Range**

1. Those with poor circulation of lower limbs, edema and obese people.
2. Those who have low immunity and feel discomfort and pain all over the body are prone to colds.
3. Those with constipation and has coarse, flabby skin.
4. Those whose channels and collaterals of the legs are blocked, and the proportion of the legs is poor.

### **Inapplicable Range**

1. Those who are in menstruation, pregnancy or lactation.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc
3. Those with skin diseases or infectious diseases.
4. Those who have incisions on their bodies.
5. Patients with severe varicose veins and tumors.
6. Those who are having allergic reactions or who have severely sensitive skin.
7. People who have just had liposuction.
8. Those who are of advanced years.
9. Pregnant women or who are on the road to recovery.

### **Do's and Don'ts after Operation**

1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours later.
3. Drink more warm water to keep hydrated and speed up metabolism.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
6. Wear pants after operation. Avoid wearing miniskirts and shorts.

## **(7) Breast Shaping**

### **Effects**

1. Correct breast shape and relieve accessory breasts
2. Relieve nodules on breasts, slight hyperplasia and distending pain in breast
3. Relieve outward expansion of breasts
4. Relieve irregular menstruation, speckles on the face and the lack of elasticity of breasts.
5. Relive breast atrophy, breast sagging and obstruction of lactiferous ducts

### **Applicable Range**

1. Those who have accessory breast and unsatisfactory breast shape
2. Those who have nodules on breasts, slight hyperplasia and distending pain in breast during the menstruation
3. Those who have free fat, saggy breast and outward expansion of breasts
4. Those who have low immunity
5. Those who have irregular menstruation, speckles on the face and breast inelasticity

6. Those who think that their breasts are underdeveloped, or who have postpartum breast atrophy, breast relaxation, blockage of lactiferous ducts and so on.

### **Inapplicable Range**

1. Patients with cardiovascular and cerebrovascular diseases, diabetes, abnormal blood coagulation and impaired function of important organs.
2. Those who are taking anticoagulants, vasodilators, corticosteroids and other drugs or who have been taking them for a long time.
3. Those with a focus of infection on the skin of the chest.
4. Patients with severe hyperplasia, fibroma and cyst of the chest
5. Those who are in pregnancy and lactation


### **Do's and Don'ts after Operation**

1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours later.
3. Drink more warm water to keep hydrated and speed up metabolism.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
6. Wear shaping and comfortable bra and don't squeeze your chest.


# Part II

## 1. Detailed Operation


After checking that the instrument is connected properly, you will see the following starting interface when the power switch is pressed.




### ① Function Selection



time display and adjustment




button for decreasing time



button for increasing time



energy level display




buttons for adjusting energy level



button for lowering energy level




button for raising energy level



home button



reset button



start/stop button





knob for adjusting suction level(turning clockwise is to increase suction, and turning anticlockwise is to decrease suction.

## ② Detailed Operation of 40K




choose **40K UNOISETION** to enter the following interface



time display and adjustment



energy level display



buttons for adjusting energy level





two modes




M1 is the mode in which handle releases energy constantly(constant working mode)





M1 is the mode in which handle releases energy discontinuously(discontinuous working mode)

### ③ Detailed Operation of Eye RF




choose to enter the following interface




time display and adjustment



energy level display




buttons for adjusting energy level



two modes



NOR is fixed mode, the red light is on and RF handle directly reaches preset energy level.




PRO is smart mode, the red light is on and RF handle slowly reaches preset energy level after it touches skin.

#### ④ Detailed Operation of Face RF




choose to enter the following interface



time display and adjustment



energy level display




buttons for adjusting energy level



two modes



NOR is fixed mode, the red light is on and RF handle directly reaches preset energy level.





PRO is smart mode, the red light is on and RF handle slowly reaches preset energy level after it touches skin.

### ⑤ Detailed Operation of Body RF




choose BODY RF to enter the following interface



time display and adjustment



energy level display



buttons for adjusting energy level



two modes



NOR is fixed mode, the red light is on and RF handle directly reaches preset energy level.




RF handle slowly reaches preset energy level after it touches skin.

### ⑥ Detailed Operation of Vacuum&RF



choose **BODY VACUUM&RF** to enter the following interface



time display and adjustment




energy level display




buttons for adjusting energy level




two modes



NOR is fixed mode, the red light is on and RF handle directly reaches preset energy level.




RF handle slowly reaches preset energy level after it touches skin.



suction time display and adjustment(suction time is longer than

release time)



release time display and adjustment("0" indicates the machine is


in Constant Suction mode)




Suction&Release mode includes 5 specific modes



M1 is Constant Suction mode, in which there is constant suction in the Vacuum&RF handle(suitable for those who need to shape body and cupping therapy for detox)



M2 is Suction&Release mode with lower frequency of release(suitable for those who need massage and lower-intensity-than-M1 cupping therapy for detox)



M3 is Suction&Release mode with higher frequency of release than M2(suitable for those who need massage and lower-intensity-than- M2 cupping therapy for detox)



M4 is Suction&Release mode with higher frequency of release than M3(suitable for those who need to massage body and accelerate blood circulation)





M5 is the Suction&Release mode with highest frequency of release(suitable for those who want massage for relaxation)


### ⑦ Detailed Operation of Face Ultrasound




choose FACE ULTRASOUND to enter the following interface



time display and adjustment



energy level display





buttons for adjusting energy level

## ⑧ Detailed Operation of Cold&Hot Handle




choose **HOT&COLD HANDLE** to enter the following interface



time display and adjustment




energy level display



buttons for adjusting energy level



hot and cold handle functions



Hot handle(can be used to open up pores and introduce nutrients into skin





Cold handle(can be used to calm and soothe skin, remove acne and repair skin, contract pores and so on)


## ⑨ Detailed Operation of BIO




choose **BIO** to enter the following interface



time display and adjustment



energy level display



buttons for adjusting energy level



mode display and selection



buttons for mode selection



M1 with low-frequency of electric current release, suitable for the first-time users



M2 with higher frequency of electric current release than M1, suitable for dull skin.



M3 with higher frequency of electric current release than M2, suitable for the skin which is rough and lacks nutrients



M4 is constant working mode, suitable for dry, rough and dull skin.

## 2. Technical Parameters

40K Unoisetion Cavitation 2.0 Handle:

Frequency:40KHz

Power:15W

Vacuum& RF Handle for body

Frequeuncy:1MHz

Power:50W

Vacuum: <-80kPa

Pressure: >250kPa

Air flow: >10L/minute

Noise level: <70dB (30cm away)

6-Polar for body

Frequency:1MHz

Power:30W

Quadrupole RF Handle for face

Frequency:1MHz

Power:30W

Bipolar RF Handle for the eye area

Frequency:1MHz

Power:10W

Hot & Cold probe for body and face

Cold: Temperature:0~5 degree centigrade

Hot: 42 degree centigrade

BIO microcurrent Probe:

Power supply output:24V

Frequency:1.5K

Power:1W

## 3. User Contraindications

People described as below should be cautious about using this instrument. Please consult a doctor or a professional before using the instrument. The details are as follows:

1. Those who are in pregnancy or lactation.
2. Those with heart disease or who are equipped with pacemaker.
3. Those whose wounds have not healed and who are recovering from the operation.
4. Those with epilepsy, severe diabetes and hyperthyroidism.
5. Those with malignant tumor, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those who have metals in their bodies such as stents, pacemakers or who are allergic to metals.
8. Those who are sensitive to electric currents.

## 4. Dos and Dont's

1. The plug with grounding pin must be used before the instrument is used and the power



socket of the instrument is properly grounded.

2. Ensure that the voltage of the instrument is suitable. If the voltage of the local power supply is not stable, we suggest that the user add a voltage regulator with matching power between the mains and the instrument.
3. In order to ensure the treatment effect and normal service life of the instrument, please use the specified parts provided or suggested by the original manufacturer.
4. The instrument should not be placed in damp places or near water sources, and should not be directly exposed to sunlight.
5. Do not place the instrument near a strong heat source, as this may affect the life and normal use of the instrument.
6. Before treatment, please remove all metal objects from the body first, so as to avoid unexpected conditions and influence the curative effect.
7. Please restrain from using the device on eyes, thyroid, parathyroid, testicles, abdomen of pregnant women, and pacemakers.
8. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
9. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use, so as to ensure the safety of using the electrical products.
10. Set the energy at lowest level at the beginning and slowly raise it after the customer adapts to it.
11. Clean the instrument with normal saline after operation to ensure its cleanliness to prolong its service life.
12. The handle can be used alone. Better results will be delivered if one handle works in conjunction with other handles.
13. Make sure the handle is in full contact with skin during operation in case skin is heated unevenly.
14. The treatment area should be kept moisturized during the operation and do not perform the treatment on dry skin.
15. When the treatment receiver is taking other weight-loss drugs, it is recommended to stop taking them for 1 to 2 months before treatment. Extension of the course of treatment is recommended if the treatment receiver want to use this machine immediately.
16. 40K handle should not be used on head, chest, breast, heart and back.

## **5. Troubleshooting & Solutions**

### **1. The instrument cannot be started, and the button light on the back of the instrument does not work?**

- A. Make sure the power cord is connected to the socket with normal power supply.
- B. Whether the fuse on the back of the instrument is loose or blown.

### **2. No RF output from the instrument?**

- A. Please check whether the cord connects the handle to the machine properly
- B. Please check whether the treatment area has been cleaned. Grease or essential oil products may cause poor contact between the handle and human body to cause this problem.

### **3. Weakened RF output?**

- A. Please check whether the cord connects the handle to the machine properly
- B. Please check whether non-conductive grease and other substances are on the treatment handle, which can cause poor contact between the handle and skin to weaken the output.

C. Please check whether the handle used is the one specified by the instrument.

#### **4. No suction or very little suction?**

A. Please turn off the instrument and check whether the filter used in the instrument should be replaced

B. Please check the rubber seal in the handle, because the rubber seal may be worn to cause air leakage.

C. Please check whether the cup containing the filter is screwed tight and whether the rubber seal is worn, for air leakage may lead to poor air pressure.

D. If the methods above cannot solve this problem, please contact the instrument dealer for assistance.

#### **5. I can start the instrument, but the screen shows erroneous message?**

A. Please unplug the power cord behind the instrument and wait for about 1 minute before restarting the instrument.

B. If the above method cannot solve the problem, please contact the instrument dealer for assistance.

## **6. FAQs**

### **1. Q: How long can I start seeing results of RF treatment?**

A: Normally, you can see the treatment results that day when you take the operation or within one week. Skin collagen tissue is heated and then contracts, and you can obviously feel the sensation of skin being tightened. Because RF works to stimulate subcutaneous tissue to promote collagen regeneration, so the more you take the treatment, the more noticeable the effects will be.

### **2. Q: Is RF harmful to skin?**

A: RF treatment is non-surgical operation. It stimulates collagen regeneration at the skin bottom and boosts metabolism, and therefore has no damage to the skin. Local redness and hotness are the normal conditions of accelerated blood circulation, which will be relieved after a moment. So there's no need to worry about that.

### **3. Q: How long does the body-firming treatment take?**

A: One treatment session is 60 minutes. We will combine professional techniques and instruments to achieve remarkable results.

### **4. Q: What functions does this instrument have?**

A: It has the functions of dissolving fat and shaping body, smoothing&rejuvenating skin, relaxing and repairing skin, firming &lifting skin and anti-aging. It can be operated on every part of the body. In terms of facial care, it can form face contour, alleviate wrinkles and saggy skin. In terms of body care, it can reduce regional fatness and build S-shaped curve, and promote metabolism and detoxification of the whole body. Meanwhile, it can also enhance the functions of viscera and solve sub-health problems.

### **5. Q: Which one is better in terms of weight loss, liposuction or this machine?**

A: The principle of liposuction is to suck out the excess fat in one part of the body through negative pressure suction so as to achieve the purpose of rapid weight loss. It's a surgical treatment requiring anesthesia and so there is recovery period, which has risks despite its immediate results. However, our instrument is risk-free without side effects. It can not only help

you reduce fat but also stimulate collagen regeneration to lift and tighten skin. It can also build perfect curves by lifting skin, and accentuate women's sexual charm.

**6. Q: Will I regain weight after operation?**

A: Ultrasound works to dissolve fat to reduce weight, and users are not prone to regain their weight. Because what's dissolved through ultrasonic treatment is fat, not water. Fat needs a longer period of time to accumulate, so normally people will not regain weight.

**7. Q: Do I need to be on a diet?**

A: Pay attention to what you eat. Because radio frequency and ultrasonic treatment work to dissolve fat and accelerate metabolism. It's best to avoid eating spicy, oily, or fried food in case they affect regular metabolism. Proper workout for sweating helps with fat metabolism process. In such way, the results will be more noticeable.

**8. Q: Does it have side effects on the body?**

A: RF beauty treatment is non-invasive and is one of the most effective and the safest approach to remove wrinkles and shape body. Generally speaking, side effects won't appear. A small number of people may experience transient redness or swelling that will disappear after a few hours. Those with dry skin may experience atrophy after treatment, which will disappear after three days because collagen regeneration will make skin plump up though RF heat can cause moisture loss of skin. Overall, there are no side effects on people's health.

**9. Q: Can ultrasound treatment be performed on every part of the body?**

A: Ultrasound treatment works to rupture fat cells, by which fat is released from cells and then metabolized through lymphatic system. So its power is very high. Heart is sensitive to high frequency sound waves since it is a hollow organ. After being vibrated by sound waves, the heart will reflect the waves due to the different transmission of sound waves by the muscle tissue and the blood in the heart. This force can pull valves away from the heart muscle. If you use the instrument on eyes, it can lead to retinal detachment, so the ultrasound treatment should not be performed on the areas of eyes and heart. (not on the areas of waist, back and chest)

**10. Q: Does ultrasound treatment have side effects on human body?**

A: Ultrasonic beauty treatment is non-surgical and non-invasive program without anaesthesia. It relies on ultrasound (mechanical sound wave) - dissolve fat through cavitation - ultrasound focusing. And the results are fat smash - fat dissolving - skin firming and body sculpting respectively. It only targets low-density fat cells but leaves high-density tissue like blood vessels intact, and therefore has no side effects on human body. Slight tinnitus will probably occur during the treatment, which is within normal range. So there's no need to worry about that.

**11. Q: Why does tinnitus occur?**

A: It's because ultrasound comes with a high frequency vibration, which could go up to 20kHz. It penetrates into 20mm-deep skin fat layers through ultrasound focusing. Highly intensified ultrasound can cause high-speed friction among fat cells, which can heat and dissolve them. And during the process, sound may be produced so that customer will have the sensation of tinnitus.

**12. Q: What is collagen?**

A: Collagen is a kind of biological macro molecule substance, a kind of white, opaque and non-branched fibrous protein. It can supplement the nutrients needed by all layers of the skin, enhance the collagen activity in the skin, retain moisture, nourish skin, slow down ageing, relieve facial relaxation, provide nutrients for the hair and so on. Collagen is a nutrient that human body must replenish to delay ageing. With the growth of age, collagen will gradually be

lost. For women, collagen will begin to decrease since the age of 20. At the age of 25, collagen loss will reach its peak. At the age of 40, the collagen content is less than half of that one has at age of 18. It is the loss of collagen and water that causes the wrinkles on the face of the elderly. The breakage of collagen fibers and elastic fibers will cause skin tissue to be oxidized, atrophied and collapsed. The skin will be dry, wrinkled, flabby and inelastic. Therefore, in order to delay ageing, there must be plenty of collagen.


**13. Q: How does vacuum detoxify?**

A: Vacuum can cause capillary congestion, stimulate cell to increase vitality. The suction and release of pressure during operation can make local pores continue to open and close, which can promote skin breathing, increase the amount of skin oxygen, and speed up waste removal. So it achieves the effects of promoting qi and blood circulation, dispelling cold and damp, dredging meridians and collaterals, clearing toxins, dissipating heat and so on.

**7. Packing List**

- Machine Body x1
- Fittings Holder x3
- Fuse x2
- Filter x1
- Power Cord x1
- 2-polar RF Handle x1
- 4-polar RF Handle x1
- 6-polar RF Handle x1
- Vacuum&RF Handle x1
- BIO Handle x1
- 40K Handle x1
- Ultrasound Handle x1
- Hot&Cold Handle x1

**8. Operational Diagrams**

Parameter Adjustment	Product	Techniques	Diagrams
<b>Face Firming and Anti-ageing: 60 minutes once a week</b>			
<b>Face RF</b> Recommend energy level:30-70% Time: 10-15 minutes Mode: NOR/PRO	Makeup Remover+Face Cleanser+Cold&Hot Steam+Massage Cream+Base Mask+Essence+Face Mask	<ol style="list-style-type: none"> <li>1. Remove makeup and clean face, 5 minutes.</li> <li>2. Apply water mask and hot steam, 10 minutes</li> <li>3. Clean the face, 2 minutes.</li> <li>4. Apply toner , 1 minute</li> <li>5. Apply massage cream evenly to face and soothe face, 3 times.</li> </ol>	Technique5, 7  Technique6





<p><b>BIO</b>  Recommend ed energy level:30-70%  Time: 10-15 minutes  Mode:  M1\M2\M3\M4</p> <p><b>Ultrasound</b>  Recommend ed energy level:30-70%  Time: 10-15 minutes</p>		<p>6. Apply digital pressing technique on the following acupoints (RN-24, DU-26, ST-4, ST-6, LI-20, BL-1, BL-2, EM3, TE-23, the temple, GB-1, ST-1, ST-2), 3 times.</p> <p>7. Soothe face, 3 times</p> <p>8. Alternately apply lifting technique from chin to ear lobe, mouth corner to Ermen(SJ21), and nose wing to temple, lift eye corner from lower eyelid to temple with both hands, 3 times.</p> <p>9. Repeat the technique on the other side.</p> <p>10. Apply lifting technique towards hairline on forehead, 3 times</p> <p>11. Apply plucking technique with fingers alternately to lift face.</p> <p>12. Apply plucking technique with one hand towards hairline on forehead.</p> <p>13. One side, apply plucking technique with both hands alternately to lift face, 2-3 minutes</p> <p>14. Repeat the technique on the other side.</p> <p>15. Move in a zigzag motion on forehead with middle and ring fingers, slide to the front side of ears and lift back and forth for 3 times and then slide to the back of ears.</p> <p>16. Wash face clean, 2 minutes</p> <p>17. Apply essence evenly to face, 1 minute</p> <p>18. <b>Face RF</b> operation, apply lifting technique from jaw to earlobe, mouth corner to ear center, and nose wing to temple by moving the handle line by line, 3 times.</p> <p>19. Repeat the technique on the other side.</p> <p>20. Apply lifting technique by moving the handle towards</p>	<div data-bbox="1161 197 1283 394" data-label="Image"> </div> <p data-bbox="1082 412 1244 483">Technique 8, 18, 23, 24</p> <div data-bbox="1161 501 1283 698" data-label="Image"> </div> <p data-bbox="1082 712 1244 784">Technique 10, 12, 20, 25</p> <div data-bbox="1161 801 1283 999" data-label="Image"> </div> <p data-bbox="1082 1012 1244 1039">Technique15</p> <div data-bbox="1161 1057 1283 1254" data-label="Image"> </div> <p data-bbox="1082 1267 1244 1294">Technique27</p> <div data-bbox="1161 1312 1283 1509" data-label="Image"> </div> <p data-bbox="1082 1523 1244 1550">Technique29</p> <div data-bbox="1161 1568 1283 1765" data-label="Image"> </div>
--	--	--	---

		<p>hairline on forehead, 3 times.</p> <p>21. Wash face clean, 2 minutes</p> <p>22. Apply essence evenly to face, 1 minute</p> <p>23. <b>BIO</b> handle operation: lift line by line from jaw to earlobe, corner of mouth to Ermen(SJ21), and nose wing to temple, 3 times.</p> <p>24. Use one hand and the handle to lift from double chin area to the back of ear, jaw to Ermen(SJ21), mouth corner to temple, nose wing to temple, 3 times</p> <p>25. Lift towards hairlines on forehead, 3 times.</p> <p>26. Do the same on the other side.</p> <p>27. <b>Ultrasound</b> handle operation: lift line by line from jaw to the back of ear, mouth corner to Ermen(SJ21), and nose wing to temple by moving the handle in circular motion, 3 times</p> <p>28. Do the same on the other side.</p> <p>29. Move the handle in circular motion horizontally on forehead, 3 times.</p> <p>30. Apply face mask, 15 minutes.</p> <p>31. Wash face clean, 2 minutes.</p> <p>32. Apply toner, essence, face cream and suncream.</p> <p>33. The end.</p>	
--	--	--	--

**Treatment Suggestions**

One course of treatment consists of ten treatment sessions. After one treatment session, skin will be firmed and rejuvenated as the nutrients it needs are absorbed. After one course of treatment, skin will lifted with glow. Two courses will make skin smooth and plump, leaving face contour clearer. Three courses will promote skin metabolism and detox processes, help fade pigmentation and skin darkness, slow down skin ageing and prevent skin from sagging, and stimulate collagen regeneration to leave skin smooth and firm with glow.

## Face Relaxing and Repairing: 60 minutes once a week






<p><b>Ultrasound</b> Recommended energy level:30-70% Time: 10-15 minutes</p> <p><b>Cold handle</b> Recommended energy level:30-70% Time: 10-15 minutes</p>	<p>Makeup Remover+ Face Cleanser+ Soothing and Repairing Gel+Cold&amp;Hot Hammer+Soothing and Repairing essence+Soothing and Repairing Mask+MS-97D1</p>	<ol style="list-style-type: none"> <li>1. Remove makeup and clean face, 5 minutes</li> <li>2. Apply toner, 1 minute</li> <li>3. Apply gel for soothing and repairing skin to face in conjunction with cold steam, 5-8 minutes</li> <li>4. Remove mask and clean face, 2 minutes</li> <li>5. Apply some essence for soothing and repairing skin to face, 1 minute</li> <li>6. <b>Ultrasound</b> operation: move the handle upwards in circles from jaw to earlobe, mouth corner to ear center, and nose wing to temple, 3 times</li> <li>7. Move the handle in circular motion towards hairline on forehead, 3 times.</li> <li>8. Do the same on the other side.</li> <li>9. Move the handle in small circles on nose area, 3-5 times</li> <li>10. Wash face clean, 2 minutes</li> <li>11. Apply some essence for soothing and repairing skin to face, 1 minute</li> <li>12. <b>Cold handle</b> operation:move the handle upwards in circles from jaw to earlobe, mouth corner to ear center, and nose wing to temple, 3 times</li> <li>13. Do the same on the other side.</li> <li>14. Move the handle in circular motion horizontally on forehead, 3 times.</li> <li>15. Apply face mask, 15 minutes.</li> <li>16. . Wash face clean, 2 minutes.</li> </ol>	<p>Technique6</p>  <p>Technique8</p>  <p>Technique12</p>  <p>Technique14</p> 
--	---	--	--

		17. Apply toner, essence, face cream and suncream. 18. The end.	
--	--	--	--

**Treatment Suggestions**

One course of treatment consists of ten treatment sessions. After one treatment session, skin will be smoothed and hydrated. After one course, skin will be return to normal and less prone to allergic reactions caused by seasonal factors, cosmetic products and dryness of skin. After two courses, skin will be supple and velvety. There will be ruddy cheeks thanks to promoted blood circulation. After three courses, effects will be consolidated.

**Face Firming and Whitening: 60 minutes once a week**

<p><b>Face RF</b> Recommend ed energy level:30-70% Time: 10-15 minutes Mode: NOR/PRO <b>Ultrasound</b> Recommend ed energy level:30-70% Time: 10-15 minutes</p> <p><b>Cold&amp;Hot handle</b> Recommend ed energy level:30-70% Time: 10-15 minutes</p>	<p>Makeup Remover+ Face Cleanser+ Cold&amp;Hot Hammer+Massage Cream+ Base Mask+Whitening Essence+ Anti-wrinkle Essence+Facial Mask</p>	<ol style="list-style-type: none"> <li>1. Remove makeup and clean face, 5 min.</li> <li>2. Apply base mask and use hot steam, 10 min.</li> <li>3. Clean the face, 2 min.</li> <li>4. Apply toner,1 minute</li> <li>5. Apply massage cream evenly on face and soothe face, 3 times.</li> <li>6. Press the following acupoints (RN-24, DU-26, ST-4, ST-6, LI-20, BL-1, BL-2, EM3, TE-23, GB-1, the Temple, ST-1, ST-2), 3 times.</li> <li>7. Soothe the whole face, 3 times.</li> <li>8. Alternate both hands to lift in 3 lines from chin to ear lobe, corner of mouth to ear gate, nose wing to temple, lower eyelid to temple, lift corners of eyes, 3 times.</li> <li>9. Do the same to the other side.</li> <li>10. Lift towards hairlines on forehead, 3 times.</li> <li>11. Alternate hands in lifting the face as if plucking the string, 3-5 times.</li> <li>12. Alternate fingers in plucking forehead towards hairline with one hand, 3-5 times</li> <li>13. One side,pluck and lift face by alternating hands,2-3 minutes</li> <li>14. Do the same to the other</li> </ol>	<p>Technique5, 7</p>  <hr/> <p>Technique6</p>  <hr/> <p>Technique 8, 18, 19</p>  <hr/> <p>Technique10, 21</p>  <hr/> <p>Technique15</p>  <hr/> <p>Technique24, 29</p>
--	--	--	--






		<p>side.</p> <p>15. Move in a zigzag shape on forehead with middle and ring fingers, slide to the front side of ears and lift back and forth for 3 times and pass behind ears.</p> <p>16. Wash face clean, 2 minutes</p> <p>17. Apply anti-wrinkle evenly to face, 1 minute.</p> <p>18. <b>Face RF</b> handle operation: lift line by line from jaw to earlobe, mouth corner to ear center, and nose wing to temple, 3 times.</p> <p>19. Use one hand and the handle to lift from double chin area to the back of ear, jaw to Ermen(SJ21), mouth corner to temple, nose wing to temple, 3 times</p> <p>20. Do the same on the other side.</p> <p>21. Lift towards hairline on forehead, 3 times.</p> <p>22. Wash face clean, 2 minutes</p> <p>23. Apply whitening essence evenly to face, 1 minute.</p> <p>24. <b>Ultrasound</b> handle operation: move the handle upwards in circular motion from jaw to the back of ear, mouth corner to Ermen(SJ21), and nose wing to temple, 3 times</p> <p>25. Move the handle in circular motion towards hairline on forehead, 3 times.</p> <p>26. Do the same to the other side.</p> <p>27. Wash face clean, 2 minutes</p> <p>28. Apply whitening essence evenly to face, 1 minute.</p> <p>29. <b>Cold&amp;Hot handle</b> operation: move the handle upwards in circular motion</p>	<div data-bbox="1161 203 1283 398" data-label="Image"> </div> <hr/> <p data-bbox="1082 416 1299 450">Technique 25, 31</p> <div data-bbox="1161 465 1283 660" data-label="Image"> </div>
--	--	--	---



		<p>from jaw to the back of ear, mouth corner to Ermen(SJ21), and nose wing to temple, 3 times</p> <p>30. Do the same to the other side.</p> <p>31. Move the handle in circular motion horizontally on forehead, 3 times.</p> <p>32. Apply face mask, 15 minutes.</p> <p>33. Wash face clean, 2 minutes.</p> <p>34. Apply toner, essence, face cream and suncream.</p> <p>35. The end.</p>	
--	--	---	--

**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one session, skin will be firmed and hydrated with necessary nutrients absorbed. After a full course, skin will be lifted with glow. Two courses help increase skin smoothness and plumpness and make face contour clearer. Three courses help promote skin metabolism and detoxification, fade pigmented spots and brighten dull skin, prevent skin sagging and aging and make the collagen at the bottom layer of skin regenerate to make skin be rejuvenated, smoothed and firmed.

**Eye Anti-wrinkle: 30 minutes 2-3 times a week**


<p><b>Eye RF</b> Recommended energy level:30-70% Time: 10-15 minutes Mode: NOR/PRO</p> <p><b>BIO</b> Recommended energy level:30-70% Time: 10-15 minutes Mode: M1\M2\M3\M4</p>	<p>Makeup Remover+Face Cleanser+Essence /Eye Cream+MS-97D1+Eye Mask</p>	<p>1. Remove makeup and clean face, 5 minutes</p> <p>2. Apply toner, 1 minute</p> <p>3. Apply essence(eye cream) evenly to eye area in circular motion, 3 times.</p> <p>4. <b>Eye RF handle</b> operation:adjust energy level, mode and time: 10 minutes.</p> <p>5. Move the instrument to lift skin from lower eyelid to eye corner, 3-6 times.</p> <p>6. Move the instrument to lift skin from lower eyelid to temple, 3-6 times.</p> <p>7. Move the instrument to lift skin from lower eyelid to temple in small circles, 3-6 times.</p> <p>8. Move the instrument to</p>	<p>Technique3, 12</p> 
			<p>Technique5, 14</p> 
			<p>Technique6, 8, 15</p> 






		<p>lift skin from lower eyelid to temple, 3-6 times.</p> <p>9. Lift from brow bone to hairline.</p> <p>10. Repeat the operation on the other side.</p> <p>11. Wash eye area clean, 1 minute</p> <p>12. Apply essence evenly to eye area in circular motion, 1 minute</p> <p>13. <b>BIO</b> handle operation:adjust energy level, mode and time: 10 minutes.</p> <p>14. Move the instrument to lift skin from lower eyelid to eye corner, 3-6 times.</p> <p>15. Move the instrument to lift skin from lower eyelid to temple, 3-6 times.</p> <p>16. Lift from brow bone to hairline.</p> <p>17. Repeat the operation on the other side.</p> <p>18. Wash eye area clean, 1 minute</p> <p>19. Apply eye mask, 15 minutes</p> <p>20. Remove mask and wash the area clean, 2 minutes</p> <p>21. Apply eye essence and cream.</p> <p>22. The end</p>	<p>Technique7</p>  <hr/> <p>Technique9, 16</p> 
--	--	--	---

**Treatment Suggestions**

It is recommended to take this treatment 2-3 times per week. After one treatment, eye area is lifted and firmed and blood circulation is accelerated. After one month, fine wrinkles and black circles are reduced and skin colour is brightened. After three months, skin in eye area is totally firmed and rejuvenated with glossiness. If the customers keep using this instrument, eye aging can be relieved and prevented.

**Neck Care: 60 minutes once a week**

<p><b>Face RF</b>          Recommended energy level:30-70%          Time: 10-15 minutes</p>	<p>Makeup Remover+Face Cleanser+Massage Cream+Essence+ Neck Mask</p>	<p>1. Remove makeup and clean neck, 5 min.</p> <p>2. Apply toner, 1 min.</p> <p>3. Neck massage, rub oil into chest with both hands towards the back of the neck</p>	<p>Technique3</p> 
---	--	--	---



<p>Mode: NOR/PRO</p> <p><b>Ultrasound</b> Recommend ed energy level:30-70% Time: 10-15 minutes</p>		<p>and press Gb-20 and Du-16 acupoint, 3 times.</p> <p>4. Alternately lift the lower jaw and double chin area to pass the back of ear to the armpit, 3 times.</p> <p>5. Push downwards with thumb and index finger splaying, tracing 3 meridians from the neck side to the armpit, 3 times.</p> <p>6. Push downwards with 4 fingers, tracing 3 meridians from neck side to the armpit, 3 times.</p> <p>7. Push downwards with bending fingers till the neck turns hot and slide to the armpit, 3 times.</p> <p>8. Push downwards with thumb and index finger splaying from external collarbone to armpit, 3-5 times.</p> <p>9. Do the same on the other side.</p> <p>10. The end</p> <p><b>RF</b> Operation, use one hand and the handle to lift double chin area and then pass the back of ear to armpit, 3 times.</p> <p>12. Move the handle in circles along three lines on one side of neck to armpit, 3 times.</p> <p>13. Slide the handle from inner and outer sides of clavicle to the armpit, 3-5 times.</p> <p>14. Move the handle in circles on the whole neck area, 3 times.</p> <p>15. Do the same on the other side.</p> <p>16. Wash neck clean, 1 minute.</p> <p>17. Apply essence evenly to neck, 1 minute</p> <p><b>Ultrasound</b> handle operation, use one hand and</p>	<p>Technique4</p>  <hr/> <p>Technique5, 6, 7</p>  <hr/> <p>Technique11, 18</p>  <hr/> <p>Technique12, 19</p>  <hr/> <p>Technique14, 21</p> 
--	--	---	---









		<p>the handle to lift double chin area and then pass the back of ear to the armpit, 3 times.</p> <p>19. Move the handle in circles along three lines on one side of neck to armpit, 3 times.</p> <p>20. Slide the handle from inner and outer sides of clavicle to the armpit, 3-5 times.</p> <p>21. Move the handle in circles on the whole neck area, 3 times.</p> <p>22. Do the same on the other side.</p> <p>23. Wash neck clean, 1 minute</p> <p>24. Apply neck mask, 15 minutes</p> <p>25. Wash neck clean, 1 minute</p> <p>26. Apply toner, neck essence, neck cream and suncream</p> <p>27. The end.</p>	
--	--	---	--



**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment, skin becomes rejuvenated and smooth. After a full course, neck skin becomes smooth, and the neck lymphatic circulation is accelerated. After 2 courses of treatment, skin will be tightened, neck lines will be reduced, and the the condition of double chin will be relieved. After 3 courses of treatment, skin will be firmed with glow, lymphatic detox will be accelerated and some problems will be resolved such as acne and skin dullness. Moreover, skin aging will be prevented and skin will remain in a young state.

**Waist&Abdomen Shaping: 60 minutes once a week**

<p><b>40K</b>          Recommend energy level:30-70%          Time:15-20 minutes          Mode:          M1:constant working mode          M2: discontinuou</p>	<p>Massage Cream(Essential Oil)+Gel+MS-97D1</p>	<ol style="list-style-type: none"> <li>1. Massage techniques</li> <li>2. Rub oil into abdomen with hands moving in circular motion, 3 times.</li> <li>3. Rub abdomen back and forth with both hands, 3-5 times.</li> <li>4. Knead abdomen with both hands using chiropractic technique, 3 times.</li> <li>5. Lift Meridian BV(Belt Vessel) on both sides of</li> </ol>	<p>Technique2, 8, 11</p>  <hr/> <p>Technique3, 19</p> 
---	---	--	--




<p>s working mode</p> <p><b>Vacuum&amp;RF:</b>  <b>Energy level:</b>  30-70%  Time: 15-20 minutes  Mode:  NOR/PRO  Suction and release time can be freely adjusted.  When release time is zero, it is in Constant Suction mode</p>		<p>waist with both hands alternately, 16 times.</p> <p>6. Move hands in 8-shaped motion with right hand and left hand moving clockwise and anti-clockwise respectively on waist, 3 times.</p> <p>7. Overlap hands and message the intestinal canal clockwise, 3 times.</p> <p>8. Move hands in circular motion to soothe the treatment area, 3 times.</p> <p>9. Apply digital pressing technique on the following acupoints: RN13, RN12, RN10, RN8, RN6, RN4, RN3, ST25, SP15. 2 times.</p> <p>10. Starting from RN3 point, push to belly navel, slide to waist and then lift upwards to groin with both thumbs, 3 times.</p> <p>11. Soothe the treatment area with hands and then slide to groin.</p> <p>12. The end</p> <p>13. <b>40K</b> handle operation: starting from one side of the waist, lift the handle to abdomen and groin, 3 times.</p> <p>14. Repeat the previous operation on the other side, 3 times.</p> <p>15. Move the handle on abdomen in small circular motion, 3 times.</p> <p>16. Move the handle on abdomen in big circular motion, 3 times.</p> <p><b>17. Vacuum&amp;RF</b> handle operation, Constant Suction mode, one side, lift from waist side to abdomen line by line, 3 times</p> <p>18. Do the same on the other side</p> <p>19. Lift Meridian BV(Belt Vessel) horizontally back and forth, 3 times</p>	<p>Technique4</p>  <hr/> <p>Technique5</p>  <hr/> <p>Technique6</p>  <hr/> <p>Technique7</p>  <hr/> <p>Technique9</p>  <hr/> <p>Technique10</p>  <hr/> <p>Technique 13, 17, 22</p>  <hr/> <p>Technique15</p> 
--	--	--	--

		<p>20. Move the cup up and down to shape waist, 3 times</p> <p>21. Move cup in big circles on abdomen, 3 times</p> <p>22. One side, lift from waist side to abdomen line by line, 3 times</p> <p>23. Wipe the are clean with towel, the end.</p>	<p>Technique16, 21</p>  <hr/> <p>Technique20</p> 
--	--	--	---

**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment session, metabolism will be accelerated. After one course, fat will be slowly metabolized, body will start to slim down. Constipation and other issues will be relieved at the same time. After two courses, the effects will gradually become obvious, skin will be tightened and lifted, and skin collagen will be increased. Stretch marks and flabby skin are also alleviated. Three courses will help consolidate the effects. Unwanted fat will gradually disappear, and waist and abdomen curves will start to take shape.

**Arm Shaping: 60 minutes once a week**

<p><b>40K</b> Energy level:30-70% Time:15-20 minutes Mode: M1:constant working mode M2: discontinuous working mode</p> <p><b>Body RF:</b> Energy level:30-70% Time: 15-20 minutes Mode:NOR/P RO</p>	<p>Massage Cream(Essential Oil)+Gel+MS-97D1</p>	<ol style="list-style-type: none"> <li>1. Massage techniques</li> <li>2. Left-right order: lay the arm flat, rub oil into skin from lower arm to the entire arm and slide to the fingers(soothing technique), 3 times.</li> <li>3. Push the entire arm with both palms, 3 times.</li> <li>4. Push hands upwards tracing three meridians: Large Intestine Meridian(LI) - Triple Energizer Meridian(TE) - Small Intestine Meridian(SI) to armpits with thumbs and index fingers splaying, 3 times.</li> <li>5. Soothe the area, 3 times</li> <li>6. Rub three meridians on arms with bending fingers back and forth to warm them up, 3 times.</li> </ol>	<p>Technique2, 3, 4, 5, 6, 7</p>  <hr/> <p>Technique8, 9</p>  <hr/> <p>Technique10</p>  <hr/> <p>Technique13, 16, 21, 24</p>
---	---	--	---

7. Soothe the area, 3 times

8. Lay inner arm upwards, and push hands tracing three yin meridians on inner arm: Lung Meridian(LU) - Pericardium Meridian(PC) - Heart Meridian(HT) to armpit with thumbs and index fingers splaying, 3 times respectively.

9. Rub three meridians back and forth with bending fingers to warm them up, 3 times.

10. Soothe the area and slide to fingers

11. Repeat above techniques on the other arm.

12. Wipe arms clean with towel, and then apply gel evenly to the entire arm.

13. **40K** operation: lay customer's arm flat, push three meridians from fat part of lower arm to armpit, 3 times.

14. Move the handle in annular motion along three meridians to armpit, 3 times.

15. Redouble the operation on flabby arms, 3 times.


16. Move the handle along three meridians to armpit, 3 times.

17. Lay inner arm upwards and push three meridians on upper arm to armpit, 3 times.

18. Move the handle in circles along three meridians on upper arm to armpit, 3 times.

19. Move the handle along three meridians from upper arm to armpit.

20. Repeat those



Technique14, 15, 22, 23

Technique17, 25



Technique19, 27













		<p>techniques on the other arm.</p> <p>21. <b>Body RF</b> operation: lay customer's arm flat, push three meridians from fat part of lower arm to armpit, 3 times.</p> <p>22. Move the handle in annular motion along three meridians to armpit, 3 times.</p> <p>23. Redouble the operation on flabby arms, 3 times.</p> <p>24. Move the handle along three meridians to armpit, 3 times.</p> <p>25. Lay inner arm upwards and push three meridians on upper arm to armpit, 3 times.</p> <p>26. Move the handle in circles along three meridians on upper arm to armpit, 3 times.</p> <p>27. Move the handle along three meridians from upper arm to armpit.</p> <p>28. Repeat those techniques on the other arm.</p> <p>29. Wipe arms clean with towel, the end.</p>	
--	--	--	--





**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment session, it will feel warm and tight, and the arms will be relaxed. After one course of treatment, excessive fat on arms will be reduced and the skin will appear to be tight and plump. Two courses of treatment will help shape arms, dredge meridians, and enhance constitution to strengthen the effects. 3 courses will help consolidate effects and prevent arms from regaining fat.

**Back Shaping: 60 minutes once a week**






<p><b>Body RF:</b> Energy level:30-70% Time: 15-20 minutes Mode:NOR/P</p>	<p>Massage Cream(Essential Oil)+MS-97D1</p>	<p>1. Massage techniques 2. Apply oil to back and press Gb-20 and Du-16 acupoints. 3. Pluck Dabanjing(the area connecting neck and shoulder) (from hairline), 3-5</p>	<p>Technique2, 11</p>
---	---	---	-----------------------



<p>RO</p> <p><b>Vacuum&amp;RF:</b>  <b>Energy level:</b>  30-70%  Time: 15-20 minutes  Mode:  NOR/PRO  Suction&amp;Release mode:  M1/M2/M3/M4/M5  Suction and release time can be freely adjusted.  When release time is zero, it is in Constant Suction mode</p>		<p>times</p> <p>4. Move thumb outwards from Bladder Meridian(BL) to the sacral region and then to Gb-20 and Du-16 points.</p> <p>5. Starting from neck to caudal vertebra, move hands in S-shape, 3 times.</p> <p>6. Push Bladder Meridian(BL) with thumb in left-right order, 3 times.</p> <p>7. Push Bladder Meridian(BL) to the sacral region with both thumbs at the same time, 3 times.</p> <p>8. Push Bladder Meridian(BL) in three bending fingers, 3 times.</p> <p>9. Push along the medial border of scapula with both hands alternately in left-right order, 3-6 times.</p> <p>10. Push medial border of scapula with both hands horizontally, 3-6 times.</p> <p>11. Soothe the back with both hands and push Gb-20 and Du-16 points, 3 times.</p> <p>12. Overlap both thumbs and push SI-11 point and slide to the arm, 3 times.</p> <p>13. Do the same on the other side, 3 times</p> <p>14. Rub Meridian GV and Bladder Meridian(BL) with hands to warm them up.</p> <p>15. The end</p> <p>16. <b>Body RF</b> handle operation, move the handle along Meridian GV and Bladder Meridian(BL) respectively from neck to the sacral region, 3-5 times.</p> <p>17. Move the handle in circles on Du-14 point and sacral region, 3-5 times respectively.</p> <p>18. From neck, move the handle on Dabanjing(the area connecting neck and shoulder) back and forth, 3-5</p>	 <p>Technique 3, 18, 25</p>  <p>Technique 4</p>  <p>Technique 5</p>  <p>Technique 6, 7, 8, 14, 16, 23</p>  <p>Technique 9, 19, 26</p>  <p>Technique 10</p>  <p>Technique 12, 13</p> 
---	--	---	---

		<p>times.</p> <p>19. Move the handle to lift the medial border of scapula back and forth in left-right order, 3-5 times.</p> <p>20. Move the handle in 8-shaped motion downwards to sacral region, 3 times</p> <p>21. Starting from the left side, move the handle in annular motion upwards to armpit, 3 times.</p> <p>22. Apply lifting technique by moving the handle from one side of the waist upwards to armpit, 3-5 times.</p> <p>23. <b>Vacuum&amp;RF</b> handle operation, move the handle along Meridian GV and Bladder Meridian(BL) respectively from neck to the sacral region, 3-5 times.</p> <p>24. Move the handle in circles on Du-14 point and sacral region, 3-5 times respectively.</p> <p>25. From neck, move the handle on Dabanjing(the area connecting neck and shoulder) back and forth, 3-5 times.</p> <p>26. Move the handle to lift the medial border of scapula back and forth in left-right order, 3-5 times.</p> <p>27. Move the handle in 8-shaped motion downwards to sacral region, 3 times</p> <p>28. Starting from the left side, move the handle in annular motion upwards to armpit, 3 times.</p> <p>29. Apply lifting technique by moving the handle from one side of the waist upwards to armpit, 3-5 times.</p> <p>30. Wipe back clean with towel, the end</p>	<p>Technique17, 24</p>  <hr/> <p>Technique20, 27</p>  <hr/> <p>Technique21, 28</p>  <hr/> <p>Technique22, 29</p> 
<p><b>Treatment suggestions</b>  A course of treatment consists of 10 treatment sessions. After one treatment session,</p>			

back will be significantly relaxed, neck stiffness will also be alleviated, and skin will be lifted. A course will help relax shoulder, reduce back fat, relieve the lump on the nape of the back, and shape the back line. After two courses, the lump will appear to be smaller, meridians will be dredged, sub-health problems will be solved, and sleep quality will be improved. 3 courses will help with back thinning, back line refinement, strengthening visceral function, improving digestion, detoxification and metabolism(one course of treatment is to reduce fat, two courses are to strengthen effects, and three courses help with effect consolidation and prevent relapse.)

## Hip Shaping: 60 minutes once a week

<p><b>Body RF:</b> Energy level:30-70% Time: 15-20 minutes Mode:NOR/P RO</p> <p><b>Vacuum&amp;RF:</b> <b>Energy level:</b> 30-70% Time: 15-20 minutes Mode: NOR/PRO Suction&amp;Release mode: M1/M2/M3/M4/M5 Suction and release time can be freely adjusted. When release time is zero, it is in Constant Suction mode</p>	<p>Massage Cream(Essential Oil)+MS-97D1</p>	<ol style="list-style-type: none"> <li>1. Standing on the side, rub oil into skin by moving both hands clockwise and anticlockwise respectively from sacral region on hips , 3 times (this is soothing technique.)</li> <li>2. Apply pushing technique on the sacral region with two thumbs, 3 times.</li> <li>3. Soothe the treatment area for 3 times and then apply digital pressing technique on the following points: Bl-23, Baliao points(BL31-34), Du-1, Gb-30, BL-36, 3 times.</li> <li>4. Soothe the treatment area, 3 times.</li> <li>5. Left-right order, apply pushing technique from the top of thigh along the four meridians -- Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) to Meridian BV respectively with both hands, 3 times each.</li> <li>6. Alternately apply pushing technique with both hands along Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) from the top of thigh to Meridian BV respectively, 3 times each.</li> <li>7. Overlap both palms and then apply pushing technique from the top of</li> </ol>	<p>Technique1, 4</p>  <p>Technique2</p>  <p>Technique3</p>  <p>Technique 6, 12, 15</p>  <p>Technique7</p>  <p>Technique 8, 13, 16</p>
---	---	---	---








		<p>thigh to Meridian BV (lifting), 3 to 5 times.</p> <p>8. Apply pushing technique from both sides of hips to the top of hips(shaping) back and forth, 3 times.</p> <p>9. Repeat technique 7</p> <p>10. Soothe the treatment area.</p> <p>11. Repeat the techniques above on the other side.</p> <p>12. <b>Body RF</b> handle operation, lift from the top of thigh to Meridian BV line by line, 3 times.</p> <p>13. Lift upwards from both sides of the buttocks to the highest point of the buttocks 3 times.</p> <p>14. Move in small circles on hip, 3-5 times</p> <p>15. Lift from the top of thigh to Meridian BV line by line, 3 times.</p> <p>16. Lift upwards from both sides of the buttocks to the highest point of the buttocks, 3 times.</p> <p>17. Do the same on the other side.</p> <p>18. <b>Vacuum RF</b> operation, the techniques are the same with Body RF handle operation</p> <p>19. Wipe hips clean with towel, the end.</p>	 <p>Technique14</p> 
--	--	--	--






**Treatment Suggestions**

One course of treatment consists of 10 treatment sessions. After one treatment, buttocks will experience the sensations of being heated and lifted. After one course of treatment, buttocks start to show obvious curve lines, excessive fat begin to disappear slowly, and skin appear to be tight. After 3 courses, the conditions of cold body and uterus, other gynecological problems such as irregular menstruation are relieved, which can help consolidate the effect to make females more attractive.

**Leg Shaping: 60 minutes once a week**

40K Energy	Massage Cream(Essential	1. Massage techniques 2. Customer lies on his front,	Technique2, 3, 4, 5, 7, 8, 9
------------	-------------------------	---	------------------------------

<p>level:30-70%  Time:15-20 minutes  Mode:  M1:constant working mode  M2: discontinuous working mode</p> <p><b>Vacuum&amp;RF:</b>  <b>Energy level:</b> 30-70%  Time: 15-20 minutes  Mode: NOR/PRO  Suction&amp;Release mode: M1/M2/M3/M4/M5  Suction and release time can be freely adjusted.  When release time is zero, it is in Constant Suction mode</p>	<p>Oil)+RF  Cream+Gel+MS-97D1</p>	<p>back side of the leg, left-right order: rub oil into the skin from calf to thigh and then move back to heel, 3 times</p> <p>3. Alternate the heels of both palms to push the entire leg upwards and move back to heel, 3 times</p> <p>4. Move hands upwards to push four main meridians: Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) with thumbs and index fingers splaying.</p> <p>5. Push popliteal fossa with both hands alternately, 3 times.</p> <p>6. Move both hands upwards alternately in an undulating motion, 3 times.</p> <p>7. Soothe the area, 3 times</p> <p>8. Push four main meridians upwards with bending fingers of both hands, 3 times.</p> <p>9. Soothe the area, 3 times</p> <p>10. Repeat the operation on the other side.</p> <p>11. Wipe leg clean with towel and apply gel evenly.</p> <p>12. <b>40K</b> operation: push Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) upwards to popliteal fossa respectively, 3 times.</p> <p>13. Move in small circles on calf to dissolve fat, 3 times</p> <p>14. Move the handle upwards to popliteal fossa and soothe the treatment area with hand, 3 times.</p> <p>15. Push 4 meridians from popliteal fossa to the top of thigh respectively, 3 times.</p> <p>16. Move in small circles from popliteal fossa to the top of thigh to dissolve fat, 3 times.</p>	 <p>Technique6</p>  <p>Technique 12, 14, 21</p>  <p>Technique 13, 22, 23</p>  <p>Technique 15, 17, 24</p>  <p>Technique 16, 18, 25</p>  <p>Technique26</p>  <p>Technique27</p>
---	---------------------------------------	--	--

		<p>17. Push meridians from popliteal fossa to the top of thigh, 3 times</p> <p>18. Move the handle in small circles on thigh to dissolve fat, 3-5 times</p> <p>19. Repeat the techniques above on the other leg.</p> <p>20. Wipe leg clean with towel and apply RF cream evenly.</p> <p>21. <b>Vacuum&amp;RF handle</b> operation: lift from heel to popliteal fossa line by line with one hand and the cup, 3 times</p> <p>22. Move cup in small circles on calf, 3 times</p> <p>23. Move the cup on calf from top to bottom back and forth,3-5 times.</p> <p>24. Lift from popliteal fossa to the top of thigh, 3 times.</p> <p>25. Move cup in small circles on thigh, 3 times.</p> <p>26. Lift from two sides of leg to the middle, 3 times.</p> <p>27. Lift upwards from two sides of thigh in annular motion</p> <p>28. Do the same on the other side.</p> <p>29. Massage the front side of legs: rub oil into skin from feet to thigh root(soothe the treatment area), 3 times</p> <p>30. Push four meridians: Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) to the top of thigh with thumb and index finger splaying.</p> <p>31. Push four meridians with bending fingers of both hands,3 times</p> <p>32. The end</p> <p>33. Do the same on the other side.</p> <p>34. Wipe legs clean with towel and apply gel evenly</p>	 <p>Technique 29, 30, 31, 37, 40, 43</p>  <p>Technique39, 46</p>  <p>Technique38, 44</p>  <p>Technique45</p> 
--	--	--	---

		<p>35. <b>40K</b> handle operation: lift the handle towards knee from calf while tracing 4 meridians, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>36. Move the handle on two sides of calf in annular motion, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>37. Lift the handle from knee to the top of thigh line by line, 3 times.</p> <p>38. Move the handle in small circles on thigh, 3 times.</p> <p>39. Lift the handle from knee to the top of thigh in annular motion, 3 times.</p> <p>40. Lift the handle from knee to the top of thigh line by line, 3 times.</p> <p>41. Wipe leg clean with towel and apply RF cream evenly.</p> <p>42. operation: lift the cup towards knee from calf while tracing 4 meridians, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>43. Lift from knee to the top of thigh,3 times</p> <p>44. Move in small circles on thigh, 3 times.</p> <p>45. Lift from two sides of thighs to the middle, 3-5 times.</p> <p>46. Lift upwards on two sides of thigh in annular motion.</p> <p>47. Do the same on the other thigh.</p> <p>48. Wipe legs clean with towel, the end.</p>	
--	--	--	--








**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment, leg will be relaxed, and circulation will be accelerated. After one course of treatment, leg will become thin, there will be contraction of skin collagen due to heat effect, and tightening effect can obviously be felt on the skin. After 2 courses of treatment, there will be obvious results with the slimness and firmness of legs promoted. 3 courses will help consolidate the effects. This treatment can stimulate dermis to promote constant



collagen regeneration. So the treatment results will be more obvious.

## Breast Shaping: 60 minutes once a week

<p><b>Body RF:</b> Energy level:30-70% Time: 15-20 minutes Mode:NOR/P RO</p>	<p>Massage Cream(Essential Oil)+MS-97D1</p>	<ol style="list-style-type: none"> <li>1. Standing at the bedside, use both hands to rub oil into skin from Shanzhong(RN17) to armpit and then lift the suspensory ligament (massage) for 3 times.</li> <li>2. Use both thumbs to press acupoints: Shanzhong(RN17), Rugen(ST18), Dabao(SP21), Yingchuang(ST16), Zhongfu(LU1), Yunmen(LU2) 3 times</li> <li>3. Massage the area, 3 times</li> <li>4. Overlap two palms and move from Shanzhong(RN17) in a 8-shaped motion back and forth, 3 times</li> <li>5. Alternate hands to push from accessory breast towards suspensory ligament, in left-right order</li> <li>6. Massage the area, 3 times</li> <li>7. Sit beside the customer, use both palms to massage and lift breasts in a left-right order, 3-5 times</li> <li>8. Alternately push the lactiferous ducts (around the breast) with thumbs and index fingers splaying 3-5 times</li> <li>9. Alternately move in circles with two thumbs to dredge nodules on the chest, 3 times</li> <li>10. Massage and lift breast with both hands</li> <li>11. Repeat the techniques on the other breast.</li> <li><b>12. Body RF</b> handle operation: left-right order,</li> </ol>	<p>Technique1</p> 
			<p>Technique2</p> 
			<p>Technique4</p> 
			<p>Technique5</p> 
			<p>Technique7, 13</p> 
			<p>Technique8, 12</p> 
			<p>Technique9, 14</p> 

		<p>lift towards nipple around breast in conjunction with soothing technique of the other hand, 5-8 times</p> <p>13. Lift in circles on the breast towards collarbone in conjunction with soothing technique of the other hand, 5-8 times</p> <p>14. Move in small circles to dredge the area where nodules exist, 3-5 times</p> <p>15. Treatment receiver lies on her side and lift up her arm, move the handle in circular motion on accessory breast to burn fat, 5-8 times.</p> <p>16. Push the handle from the accessory breast in the armpit to breast(shaping and narrowing the accessory breast) 5-8 times</p> <p>17. Repeat the techniques on the other side.</p> <p>18. Wipe breasts clean, the end.</p>	<div data-bbox="1141 212 1308 392" data-label="Image"> </div> <p data-bbox="1082 414 1252 459">Technique15</p> <div data-bbox="1141 504 1308 627" data-label="Image"> </div> <p data-bbox="1082 683 1252 728">Technique16</p> <div data-bbox="1141 772 1308 896" data-label="Image"> </div>
--	--	---	---

**Treatment Suggestions**

One course of treatment consists of ten treatment sessions. After one treatment, there will be sensation of heat on the breast, blood circulation will be promoted and there will be some degree of lifting effect. After one course, the effect will be more obvious and the nodules will be relieved. After three courses, the skin will be firmed and the breast will be well shaped with increased elasticity. The functions of endocrine system will be improved to make females more attractive.