

User Manual

MS-98B3



Preface

Dear users:

We're pleased to present to you our 8-in-1 beauty machine featured with weight loss, figure slimming, anti-aging, repairing, lifting, tightening, body shaping and so on. It is a multifunctional beauty machine, targeted at both face and body, mainly focusing on lifting and tightening of face, removing wrinkles, and body shaping. We aim at a safe and, at the same time, effective beauty result by using high-end technologies. RF and ultrasonic beauty machines are anti-aging machines for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore we advise any personnel read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

Table of Contents

Preface
Part I
Introduction
Advantages
Working Principles & Treatments
Facial Anti-aging
Facial Tightening and Anti-Wrinkle
Anti-aging around Eyes
Anti-aging around Neck
Body-shaping: Waist & Abdomen
Body-shaping: Breasts
Body-shaping
Arms
Back
Buttocks
Legs
Body Sculpting & Slimming
Overall Facial Care
Part II
1.Detailed Operations
2.Technical Parameters
3. Safety Precautions
4. Dos and Don'ts
5. Troubleshooting & Solutions
6.FAQs
7.Packing List
8.Operational Diagrams

Part I

Brief Introduction

RF and ultrasonic beauty machines are currently the most popular beauty machines that work on tightening, body-shaping and anti-aging. It's an effective replacement of all the facial and body care programs in beauty salons. Moreover, they are convenient and easy to operate. It solves both skin and body management issues for beauty-lovers. This equipment requires no injection, no medication, no operation, thus having no side-effects. It is used externally during the whole process and has an immediate effect. RF machines heats deeper skin layers, stimulates collagen cells to reform in dermis, increases blood circulation in hypodermis, thus achieving skin tightening and anti-aging. Coupled with cold hammer, which soothes skin and shrinks pores, this equipment takes care of your skin in all aspects.

Advantages

1. 8-in-1 multifunctional beauty machine that works on both of your facial and body care.
2. Ultrasonic fat burning replaces and transcends liposuction, bringing good news for treating obesity.
3. Various work heads can be replaced in different parts according to different requirements.
4. The perfect combination of ultrasound and RF solves the dilemma of losing weight. It helps with weight loss while at the same time maintains firm skin, increases skin elasticity, metabolism, and ultimately better health.
5. Adopting the cutting-edge technology of ultrasonic fat burning.
6. Suitable for all skin types. It also improves and soothes skin.
7. Painless, non-invasive during the process. Skipping any recovery periods without affecting any work and life plans.
8. No consumption, low cost and quick returns.
9. Wider treatment range and faster and more visible effects.
10. Unevenness, bleed, swelling and stasis will not appear after treatment.
11. Ultrasound, vacuum and radio frequency all have the energy release of red LED respectively. They can dissolve fat and massage while at the same time accelerate metabolism, and sterilize skin. The effect is remarkable with a higher comfort level.
12. BIO microcurrent can increase cell function, cell vitality and skin rejuvenation.

Facial Anti-aging

Principles

RF

RF works through the circuits formed by electrodes in pair. The optimum frequency of radio frequency biothermal effect is 3 MHz, and the maximum speed of extreme change is 8 million charged particles in tissues per second. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy. During the process, RF promotes the metabolism and proliferation of skin collagen through high electromagnetic wave change, therefore heating up dermis layer of the skin, thus achieving effects like tightening, lifting, wrinkle removal, anti-aging and smooth skin.

Biological Effect

Radio frequency can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, making the skin more compact and elastic.

Every time energy comes into the skin, it feels a temporary sense of heat. This means that collagen has been stimulated to produce fever, resulting in tightening of the skin. While operating sliding, the radio frequency wave emitted by the therapeutic head of the super-frequency wave system, under the cold protection of the epidermis, has acted on the deep dermis through the epidermis, rapidly heating water molecules to produce safe Bio-thermal energy. When the natural friction of collagen tissue is heated up to 45°C~60°C, it will shrink immediately and stimulate the collagen to grow continuously. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, and release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

Face Ultrasound


Ultrasonic facial beauty instrument is based on the theory of high frequency vibration. It reaches deep into the subcutaneous tissues, 4-6 mm to be exact, with the strong vibration wave of 2.8MHz. With its high energy and permeability, it vibrates tissues in every layers of your skin, exporting oil and dirt from pores. Also, it changes cells volumes so as to improve blood circulation and stimulates proliferation of cells, making it lustrous and elastic. Moreover, thanks to its high frequency vibration induced warm effect, it helps better absorption of nutrients into the skin at a higher speed.

Cold Hammer

Cold hammer uses the principle of air conditioning. Special wafer ceramics are used to absorb energy and refrigerate at -5°C in an instant. The method of physical heat accumulation produced by wafer energy absorption makes one side produce heat up

to 40 degrees. The fan discharges excess heat energy and achieves the effect of constant temperature and freezing alternately. It also promotes metabolism, alleviates allergic reaction, shrinks pore, calms skin, whitens and tenders skin.

Biological effects: During the cooling process, tissue temperature cools down, blood vessels and veins contract, metabolic rate reduces, cell activity inhibits, peripheral nerve sensitivity reduces, thus having a sedative and soothing effect, counteracting or reducing the symptoms of allergic reaction, such as redness, swelling, and pain, therefore achieving certain effects, such as calming the skin, reducing redness and swelling, shrinking pores and so on.



Effects

1. Deeply channel nutrients into skin and promote absorption.
2. Improve blood circulation and help metabolism.
3. Remove dirt and oil secretion from pores.
4. Shrink pores, calm skin, reduce redness and swelling.
5. Stimulate collagen regeneration and improve facial relaxation.
6. Improve facial fine lines and nourish skin.
7. Improve uneven, dull and sallow skin tone.

Applicable Range

1. Those with facial dullness, skin relaxation and sagging.
2. Those with dry, dehydrated and malnourished skin.
3. Those with uneven, dull and sallow skin tone.
4. Those with rough skin and large pores.
5. Those who are under long-term exposure of UV radiation at workplace.
6. Those with skin relaxation and sagging after child delivery.
7. Those with poor absorption and poor metabolism.

Inapplicable Range

1. Those who have just undergone plastic surgery which have implanted prostheses, metal materials and so on.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.
4. Those who have recently received certain injections, such as hyaluronic acid, intradermal injection, or plastic surgery.
5. Those with skin trauma or wound.
6. Those in menstruation, pregnancy, lactation, surgical recovery.
7. Those with skin diseases and infectious diseases.

Notes after Operation

1. Do not wash your face with overheated water within three days (lukewarm or cold water is recommended).
2. Enhance moisturizing and sunscreen.
3. Do not soak in hot springs, steam saunas or exercise strenuously, etc. within three days.
4. Use facial masks at least three times a week.
5. Within the first week after operation, it is better not to use alcohol, AHA or other exfoliating products.
6. Spicy, greasy food should be avoided. Staying up late, smoking and drinking are prohibited. Eat more vegetables, fruits, and less greasy food.
7. Avoid food that will cause '3 Hs'(high blood pressure, high cholesterol and high blood sugar), less greasy food is advised.

Facial Tightening and Anti-Wrinkle

Principles

Face RF

Radio frequency can alternate electrical polarities of biological tissues on the treated area millions of times per second and change directions in the same frequency concurrently. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy, which stimulates the dermis to produce more collagen and fill up the gap of the lost collagens, therefore plumping up the skin and revive elasticity.

During multipolar RF treatment, polarities alternate consecutively, which produces a higher intensity energy release, heats up the bottom layers of skin. Through fast and consecutive heating of tissues, collagen regeneration is boosted, thus having a more

visible result in a shorter period of time. It also comes with a wider and more intensifies treatment scope.

Biological Effect

Radio frequency can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, making the skin more compact and elastic.

Every time energy comes into the skin, it feels a temporary sense of heat. This means that collagen has been stimulated to produce fever, resulting in tightening of the skin. While operating sliding, the radio frequency wave emitted by the therapeutic head of the super-frequency wave system, under the cold protection of the epidermis, has acted on the deep dermis through the epidermis, rapidly heating water molecules to produce safe Bio-thermal energy. When the natural friction of collagen tissue is heated up to 45°C~60°C, it will shrink immediately and stimulate the collagen to grow continuously. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, and release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

BIO

The combination of micro-current and inductive current has a strong cosmetic effect. The function of microcurrent is to increase the synthesis of ATP (adenosine triphosphate), which is five times higher than usual. It can provide muscle and body energy, make you look energetic; Improve the permeability of cells, 30% - 40% higher than usual, help the human body better absorb nutrients and water, oxygen-containing blood, etc., and remove waste and toxins as soon as possible. Promote the synthesis of muscle fibers, 68% higher than usual, keep the skin moist, soft, compact and lustrous; Increase the level of collagen synthesis to 73%, reduce wrinkles, make the skin more soft, elastic and lustrous. The function of induced current is to balance cell's electric potential, increase cell ATP, accelerate blood circulation, improve collagen synthesis and increase skin elasticity.

BIO Cell Activation Instrument, which imitates human biological current, enters muscle cells through skin, stimulates ATP contained in cells, and restores normal operation and function of cells. It collects the latest technology to shape facial contour, thin face, reduce double chin and wrinkles and shrink pores for customers in the shortest time possible. BIO makes facial muscles tight and delicate, and rosy cheeks will appear naturally . It is more effective on loose skin, so as to achieve the goal of shaping a beautiful face.

Effects

1. Activate cells and increase skin elasticity.

2. Improve loose, soft and inelastic skin.
3. Lift the cheeks to make the face more three-dimensional.
4. Improve facial fine lines, nasionobial folds, wrinkles and dynamic wrinkles.
5. Improve facial contour and tighten skin.
6. Stimulate collagen proliferation and prevent skin aging.

Applicable Range

1. Those with sagging, soft and inelastic skin.
2. Those with indistinct facial contour.
3. Those with dull, coarse skin and large pores.
4. Those with facial fine lines, nasionobial folds, wrinkles and dynamic wrinkles.
5. Those who are under long-term exposure of UV radiation at workplace.
6. Those with skin relaxation, sagging and swelling after child delivery.

Inapplicable Range

1. Those who have just undergone plastic surgery which have implanted prostheses, metal materials and so on.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.
4. Those with skin trauma or wound.
5. Those in menstruation, pregnancy, lactation, surgical recovery.
6. Those with skin diseases and infectious diseases.
7. Those who are unrealistic about the effects.

Notes after Operation

1. Do not wash your face with overheated water within three days (lukewarm or cold water can be used).
2. Strengthen moisturizing and sunscreen.
3. Do not soak in hot springs, steam saunas or exercise strenuously, etc. within three days.
4. Use facial masks at least three times a week.
5. Spicy, greasy food should be avoided. Staying up late, smoking and drinking are prohibited. More vegetables, fruits, and less greasy food.
6. Avoid food that will cause 'three highs'(high blood pressure, high cholesterol and high blood sugar), less greasy food is advised.

Anti-aging around Eyes

Principles

Eye RF

RF works through the circuits formed by electrodes in pair. Radio frequency waves penetrate skin directly and produce thermal energy created by strong vibration from tissue resistance (1 MHz). During the treatment, RF device stimulates new collagen production by heating up hypodermis around eyes, which helps to tighten the skin, remove dark circles and bags, treat couperose and fine lines.

Biological Effects

RF can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, make the skin more compact and elastic.

Every time energy comes into the skin, it feels a temporary sense of heat. This means that collagen has been stimulated to produce fever, resulting in tightening of the skin. While operating sliding, the radio frequency wave emitted by the therapeutic head of the super-frequency wave system, under the cold protection of the epidermis, has acted on the deep dermis through the epidermis, rapidly heating water molecules to produce safe Bio-thermal energy. When the natural friction of collagen tissue is heated up to 45°C~60°C, it will shrink immediately and stimulate the collagen to continue. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

BIO

The combination of micro-current and inductive current has a strong cosmetic effect. The function of microcurrent is to increase the synthesis of ATP (adenosine triphosphate), which is five times higher than usual. It can provide muscle and body energy, make you look energetic; Improve the permeability of cells, 30% - 40% higher than usual, help the human body better absorb nutrients and water, oxygen-containing blood, etc., and remove waste and toxins as soon as possible. Promote the synthesis of muscle fibers, 68% higher than usual, keep the skin moist, soft, compact and lustrous; Increase the level of collagen synthesis to 73%, reduce wrinkles, make the skin more soft, elastic and lustrous. The function of induced current is to balance cell's electric potential, increase cell ATP, accelerate blood circulation, improve collagen synthesis and increase skin elasticity.

BIO Cell Activation Instrument, which imitates human biological current, enters muscle cells through skin, stimulates ATP contained in cells, and restores normal operation and function of cells. It collects the latest technology to shape facial contour, thin face, reduce double chin and wrinkles and shrink pores for customers in the

shortest time possible. BIO makes facial muscles tight and delicate, and rosy cheeks will appear naturally . It is more effective on loose skin, so as to achieve the goal of shaping a beautiful face.

Effects

1. Relieve eye fatigue, dark circles, eye bags and edema.
2. Reduce lines around eyes and crow's feet.
3. Reduce dark spots and accelerate blood circulation.
4. Accelerate metabolism and prevent hyperpigmentation.
5. Supplement nutrition to the eyes and moisturize skin.
6. Tighten and soften skin, lift the corners of the eyes.
7. Accelerate blood circulation around eyes and help skin absorb nutrients efficiently.

Applicable Range

1. Those with wrinkles, fine lines, bags and dark circles around the eyes.
2. Those with dry skin, dry lines and dynamic wrinkles.
3. Those with eyes prone to fatigue and dryness.
4. Those frequently face computer screens and mobile phones.
5. Those who often stay up late and have dark circles.
6. Those who often stay in a dry or hot environment.

Inapplicable Range

1. Those who had just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.
4. Those with skin trauma or wound.
5. Those who are overaging.
6. Those who are pregnant or convalescents.
7. Those with skin disease and contagious diseases.

Notices after Operation

1. Avoid excessive exposure to the sun. Sun protection is encouraged.
2. After operation, the treated area is relatively dry, moisturizing and sun protection is encouraged.
3. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
4. Avoid facial washing with overheated water, sauna steaming, hot springs or other strenuous exercise within 7 days after operation.

5. Moisturizers and eye masks are encouraged to use, preferably 3 times or more on eye masks a week.

Anti-aging on Neck

Principles

Face RF

RF instrument can achieve the effect of tightening skin, lifting skin and shaping face. The optimum frequency of radio frequency biothermal effect is 3 MHz, and the maximum speed of extreme change is 8 million charged particles in tissues per second. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy.

RF wave directly penetrates the skin and generates energy by the resistance action of the skin, which makes the temperature of the skin bottom rise, rapidly and continuously heats the tissue and promotes the growth of collagen in the skin. The principle is tightening collagen and stimulating the proliferation of collagen in the dermis and then is utilized to enhance the tightness and durability in time. After 2-6 months, collagen will gradually proliferate and recombine, making sagging or relaxed skin tight and elevated.

Biological Effect

Radio frequency can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, making the skin more compact and elastic.

Every time energy comes into the skin, it feels a temporary sense of heat. This means that collagen has been stimulated to produce fever, resulting in tightening of the skin. While operating sliding, the radio frequency wave emitted by the therapeutic head of the super-frequency wave system, under the cold protection of the epidermis, has acted on the deep dermis through the epidermis, rapidly heating water molecules to produce safe Bio-thermal energy. When the natural friction of collagen tissue is heated up to 45°C~60°C, it will shrink immediately and stimulate the collagen to continue to grow. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

BIO

BIO Cell Activation Instrument, which imitates human biological current, enters muscle cells through skin, stimulates ATP contained in cells, and restores normal

operation and function of cells. It collects the latest technology to shape facial contour, thin face, reduce double chin and wrinkles and shrink pores for customers in the shortest time possible. BIO makes facial muscles tight and delicate, and rosy cheeks will appear naturally . It is more effective on loose skin, so as to achieve the goal of shaping a beautiful face.

Effects

1. Reduce fine lines and wrinkles on the neck.
2. Improve the flabbiness, roughness and dullness of the neck skin.
3. Tighten skin and increase skin elasticity.
4. Remove double chin.
5. Accelerate lymphatic detoxification and improve facial skin.
6. Prevent cervical and lymphatic diseases.

Applicable Range

1. Those with fine lines and wrinkles on the neck.
2. Those with loose and inelastic neck skin.
3. Those with dark skin.
4. Those who often bow their heads.
5. Those who are dissatisfied with their neck skin color.

Inapplicable Range

1. Those who have just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin.
4. Those with skin trauma or wound.
5. Those who are overaging.
6. Those who are pregnant or convalescents.
7. Those with skin disease and other contagious diseases.

Notes after Operation

1. Sunscreen is encouraged . Keep the neck warm.
2. After operation, the treated area is relatively dry, moisturizing and sun protection is encouraged.
3. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
4. Avoid facial washing with overheated water, sauna steaming, hot springs or strenuous exercise within 7 days after operation.
5. Moisturizing and mask for the neck is encouraged, at least 3 times masks a week.

Essence or neck cream is encouraged.

6. Avoid bowing head for a long time.

Body-shaping: Waist & Abdomen


Principles

40K

Using the sound wave with a frequency of 4 MHz emitted by the strong sound wave head gathering strong sound waves, human adipocytes can produce a strong impact and friction movement between adipocytes after entering the human body, which can effectively consume heat, water and shrink adipocytes. In addition, when sound waves vibrate, it can make adipocytes smaller. Strong crack will come along, cells will burst instantly, adipocytes reduce, so as to achieve the effect of fat removal.

Cavitation Principle of Ultrasound: Tens of thousands of tiny bubbles, namely cavitation bubbles, are produced by vibration of liquid. These bubbles grow in the negative pressure region formed by the longitudinal propagation of ultrasound, and close rapidly in the positive pressure region, thus being compressed and stretched under alternating positive and negative pressures. The bubbles will be compressed until burst, which will generate huge instantaneous pressure, generally up to tens of MPa to hundreds of MPa, and produce strong vibration and noise.

Advantages: It only aims for low-density adipose tissue in a specific frequency, and protects high-density tissue such as vascular and nerve tissue.



Vacuum & RF

Massage skin and muscles with the special negative pressure suction head, which can effectively improve body fluid flow and increase cell activity to achieve effects like improving skin elasticity, accelerating blood circulation of capillary blood vessels,

clearing out toxins through lymphatic drainage, and improving and lowering the occurrence of dark spots, hyperpigmentation and blood stasis. The kneading action generated by vacuum suction can increase the tissue activity of skin and muscle, thus helping to reduce the cellular tissue and increase the elasticity of skin tissue, so that both slimming and shaping can be completed at the same time. In addition, vacuum negative pressure movement can also stimulate the surface and deep sympathetic nervous system, improve skin sensitivity. The suction and release process of air pressure can not only improve the capillary system but also improve the flow between deeper veins and lymph network, strengthen blood vessels, and improve varicose veins.

Advantages: breaking through the previous single and simple working mode of negative air pressure, using different modes for different parts can be more efficient and more effective, slimming, shaping, and will never harm human body. Unique RF function head design, RF and negative pressure are relatively independent but complementary metabolic system as a whole, RF has negative pressure, negative pressure has RF, each other as a whole. Compared to the ordinary single-stage RF can be more effective, faster and more uniform perfect remodeling body.

Vacuum Physical Effects

1. Skin layer

The suction and release of air pressure:

(1) It can improve the fluidity of cells, thus increasing the movement of cells, and can treat diseases related to blood stasis and blood stagnation.

(2) Draining lymph glands and veins.

Effects: redundant water is discharged from tissue fibers.

2. Vascular layer

The suction and release of air pressure:

(1) Improve the blood circulation of microvessels, not only for the capillary system, but also the flow between the deeper veins and lymphatic networks.

(2) Remove excess toxins in the body.

Effects: Strengthen blood vessels, eliminate toxins and relieve varicose veins.

3. Fibrous layer

The suction and release of air pressure:

(1) Repair cell tissue, increase its activity. Lift and revive skin elasticity.


(2) Stimulate the production of bone collagen and improve skin plumpness.

(3) Improve the oxygen support capacity of skin and increase the consumption of carbon dioxide.

Effects: Break the hard fibers of cellulite, make it more elastic, so as to achieve the effect of body-shaping.

Expel Waste and Tighten Skin-Vacuum Sextupolar 3D SMART RF

Quadrupolar Vacuum&Sextupolar 3D SMART RF radio frequency combines the mechanism of advanced Sextupolar 3D SMART RF energy and vacuum therapy, can promote tissue metabolism, repel the cellulite, which make fat granule in different depths and shocking make fusion energy, so that the combustion of fatty acids at the same time broken out of a pot, not only the burning consumption of fatty acids but also more directly emptying so rapidly reduce the size of granule.



Effects

1. Relieve cold hands and cold feet, cold womb or cold body of women.
2. Tighten the skin on the waist and abdomen.
3. Reduce lumbar and abdominal fat.
4. Tighten skin, relieve stretch marks, obesity lines and increase skin elasticity.
5. Accelerate metabolism, relieve constipation and encourage intestinal peristalsis.

Applicable Range

1. Those with cold hands and cold feet and cold womb.
2. Those with lumbar and abdominal fat or who have sagging skin after birth.
3. Those sitting for a long time, or with unideal waistlines.
4. Those with striae due to obesity and pregnancy.
5. Those with constipation or obstruction of abdominal meridians.

Inapplicable Range

1. Those with "3 Hs"(hypertension, hyperglycemia, hyperlipidemia) and who have heart disease.
2. Those in pregnancy, menstruation or lactation.
3. Those who have just undergone surgical wounds or convalescents.
4. Those with epilepsy and severe diabetes and hyperthyroidism.
5. Those with malignant tumors, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those with severe gynecological diseases.
8. Those whose gynecological diseases are being treated.

Notes after Operation

1. Avoid crop top or anything exposing waist and abdomen as well as cold and windy environment.
2. Avoid binge eating or drinking. No alcohol, no spicy and greasy food. Staying up late should be avoided. Drink warm water.
3. Keep abdomen warm. Take a bath after 4-6 hours.
4. Steam sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.
5. Before going to bed at night, you can rub your abdomen clockwise with your hands, so as to better the effect of weight loss and metabolism.

Body-shaping: Breasts

Principle

Body RF

RF heat effect produces specific radio frequency waves in a specific depth of the subcutaneous skin, which directly acts on the dermis through the skin epidermis, thus heating the tissue, promoting the catabolism of subcutaneous fat, and stimulating the proliferation and reorganization of collagen and elastic fibers. With the continuous regeneration and rearrangement of collagen in the dermis, it can be achieved. Let the flabby skin feel the effect of lifting upward and compacting immediately after the treatment. At the same time, massage can promote breast blood circulation, accelerate blood circulation, dredge congested blood clots and silt, and free fat to the chest under the condition of vibration. At the same time, more nutrients are absorbed into the breasts, which not only prevents atypical lobular hyperplasia and breast cancer, but also enables you to have natural, healthy, beautiful and strong breasts, and makes beautiful female friends more confident.

Effects

1. Adjust chest type and reduce accessory breast.
2. Relieve thoracic nodules, slight hyperplasia, alleviate menstrual breast pain.
3. Reduce chest expansion.
4. Alleviate irregular menstruation and reduce facial spots or skin without elasticity.
5. Alleviate postpartum breast atrophy, breast relaxation, breast duct blockage.

Applicable Range

1. Those with bad breast shape and accessory breasts.
2. Those with nodules in the chest, slight hyperplasia and painful breast during menstruation.
3. Those with extended chest fat extended.

4. Those with low immunity.
5. Those with irregular menstruation and facial spots or whose skin is inelastic.
6. Those with improperly developed mammary glands.

Inapplicable Range

1. Those with cardiovascular and cerebrovascular diseases, diabetes mellitus, abnormal coagulation function and impaired vital organ function.
2. Those who are long-term or ongoing users of anticoagulants, vasodilators, corticosteroids and other drugs.
3. Those with infectious lesions on the chest skin.
4. Patients with severe thoracic proliferative or fibroma or cyst.
5. Those in pregnancy and lactation.

Notes after Operation

1. Keep warm.
2. Wear styled and comfortable underwear. Don't press your chest hard.
3. Avoid getting cold. Take a bath after 4-6 hours.

Body-shaping

Body RF

RF instrument can achieve the effect of tightening skin, lifting skin and shaping face. The optimum frequency of radio frequency biothermal effect is 3 MHz, and the maximum speed of extreme change is 8 million charged particles in tissues per second. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy.

RF wave directly penetrates the skin and generates energy by the resistance action of the skin, which makes the temperature of the skin bottom rise, rapidly and continuously heats the tissue and promotes the growth of collagen in the skin. The principle is tightening collagen and stimulating the proliferation of collagen in the dermis and then is utilized to enhance the tightness and durability in time. After 2-6 months, collagen will gradually proliferate and recombine, making sagging or relaxed skin tight and elevated.

Biological Effect

Radio frequency can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, making the skin more compact and elastic.

Every time energy comes into the skin, it feels a temporary sense of heat. This means that collagen has been stimulated to produce fever, resulting in tightening of the skin. While operating sliding, the radio frequency wave emitted by the therapeutic head of the super-frequency wave system, under the cold protection of the epidermis, has acted on the deep dermis through the epidermis, rapidly heating water molecules to produce safe Bio-thermal energy. When the natural friction of collagen tissue is heated up to $45^{\circ}\text{C}\sim 60^{\circ}\text{C}$, it will shrink immediately and stimulate the collagen to continue to grow. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

Body sculpture & cellulite treatment


Body 8-Polar **3D SMART RF** & Unoisetion cavitation 2.0 & Vacuum
Sextupolar **3D SMART RF** + Laser Photon

01

Soften Fat Cells-Body 8-Polar **3D SMART RF**


Using a high power of 1MHz 8-Polar 3D SMART RF radio frequency which is the same power for professional use, with a patent pending technology that generates heat deep under the skin, promoting metabolism by its hyper-thermal effect, this procedure enhances blood/lymph flow and is used for the purpose of fat reduction.

1




Before

2



In progress

3



After


40K

Using the sound wave with a frequency of 4 MHz emitted by the strong sound wave head gathering strong sound waves, human adipocytes can produce a strong impact and friction movement between adipocytes after entering the human body, which can effectively consume heat, water and shrink adipocytes. In addition, when sound waves vibrate, it can make adipocytes smaller. Strong crack will come along, cells will burst instantly, adipocytes reduce, so as to achieve the effect of fat removal.

Cavitation Principle of Ultrasound: Tens of thousands of tiny bubbles, namely cavitation bubbles, are produced by vibration of liquid. These bubbles grow in the negative pressure region formed by the longitudinal propagation of ultrasound, and close rapidly in the positive pressure region, thus being compressed and stretched under alternating positive and negative pressures. The bubbles will be compressed

until burst, which will generate huge instantaneous pressure, generally up to tens of MPa to hundreds of MPa, and produce strong vibration and noise.

Advantages: It only aims for low-density adipose tissue in a specific frequency, and protects high-density tissue such as vascular and nerve tissue.



Vacuum Slimming

Massage skin and muscles with the special negative pressure suction head, which can effectively improve body fluid flow and increase cell activity to achieve effects like improving skin elasticity, accelerating blood circulation of capillary blood vessels, clearing out toxins through lymphatic drainage, and improving and lowering the occurrence of dark spots, hyperpigmentation and blood stasis. The kneading action generated by vacuum suction can increase the tissue activity of skin and muscle, thus helping to reduce the cellular tissue and increase the elasticity of skin tissue, so that both slimming and shaping can be completed at the same time. In addition, vacuum negative pressure movement can also stimulate the surface and deep sympathetic nervous system, improve skin sensitivity. The suction and release process of air pressure can not only improve the capillary system but also improve the flow between deeper veins and lymph network, strengthen blood vessels, and improve varicose veins.

Advantages: Breaking through the previous single and simple working mode of negative air pressure, using different modes for different parts can be more efficient and more effective, slimming, shaping, and will never harm human body. Unique RF function head design, RF and negative pressure are relatively independent but complementary metabolic system as a whole, RF has negative pressure, negative pressure has RF, each other as a whole. Compared to the ordinary single-stage RF can be more effective, faster and more uniform perfect remodeling body.

Vacuum Physical Effects

1. Skin layer

The suction and release of air pressure:

(1) It can improve the fluidity of cells, thus increasing the movement of cells, and can treat diseases related to blood stasis and blood stagnation.

(2) Draining lymph glands and veins.

Effects: redundant water is discharged from tissue fibers.

2. Vascular layer

The suction and release of air pressure:

(1) Improve the blood circulation of microvessels, not only for the capillary system, but also the flow between the deeper veins and lymphatic networks.

(2) Remove excess toxins in the body.

Effects: Strengthen blood vessels, eliminate toxins and relieve varicose veins.

3. Fibrous layer

The suction and release of air pressure:

(1) Repair cell tissue, increase its activity.

(2) Lift and revive skin elasticity.

(3) Stimulate the production of bone collagen and improve skin plumpness.

(4) Improve the oxygen support capacity of skin and increase the consumption of carbon dioxide.

Effects: Break the hard fibers of cellulite, make it more elastic, so as to achieve the effect of body-shaping.

4. Nerve layer


The suction and release of air pressure:

(1) Stimulates the surface and deep sympathetic nervous system.

(2) Improves skin sensitivity.

(3) Repair skin elasticity and resist tissue fibrosis.

Effects: Repair and improve skin sensitivity.



Body-shaping: Arms

Effects

1. Stimulate collagen hyperplasia in the bottom of the skin to shape and firm it.
2. Improve loose skin.
3. Improve flabby arms and thick arms.
4. Soothe and firm loose skin.
5. Dissolve the fat, and say goodbye to flabby arms.
6. Accelerate blood circulation and dredge channels and collaterals.

Applicable Range

1. Those with thick arms and want to look better in clothes.
2. Those with jiggly and flabby arms.
3. Those with loose arm skin.
4. Those whose arms are prone to pain and numbness.

Inapplicable Range

1. Those who have just had plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc.
3. Those who are during allergic period and has severely sensitive skin.
4. Those with skin trauma or cut.
5. Those who are overly old.
6. Those who are pregnant and who are recovering from surgery.
7. Those with skin diseases and infectious diseases.

Notes after Operation

1. Keep warm after operation. Do not eat cold food.
2. 4-6 hours later in the shower.
3. Drink plenty of warm water to stay hydrated.
4. Refuse to overeat or stay up late.
5. Avoid sauna, hot spring or strenuous exercise within one week after operation.

Body-shaping: Back

Effects

1. Relieve shoulder and back soreness, and improve the lump on back of the neck.
2. Dredge channels and collaterals, and improve channels and collaterals blocking.
3. Increase blood circulation and metabolism.
4. Improve blood supply to the head and sleep.
5. Regulate the functions of viscera and strengthen the body.
6. Firm skin and prevent sagging.
7. Improve excess fat on the back and shape the back.

Applicable Range

1. Those with sore shoulders and backs and stiff necks.
2. Those with insomnia, dreaminess and fading memory.
3. Those who are prone to fatigue, drowsiness and poor circulation of qi and blood.
4. Those with thick back and want to look better in clothes.
5. Those who has a lump on the back of the neck.

Inapplicable Range

1. Those with metal implants in the body, such as stents and pacemakers, are allergic to metal.
2. Those are during pregnancy, menstruation and lactation.
3. Those whose surgical wound is healing or convalescence
4. Those who has "three Hs" , heart disease, epilepsy and severe diabetes and hyperthyroidism.
5. Those with malignant tumor and hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those who are too weak.
8. Those who are drunk, full, empty, thirsty and overworked.

Notes after Operation

1. Keep warm, avoid cold air and drink plenty of hot water.
2. Shower 4-6 hours after operation.
3. Avoid staying up late, drinking and overeating.
4. Avoid raw, cold and spicy foods and get enough sleep.
5. Avoid shoulder-baring and backless clothing.

Body-shaping: Buttocks

Buttocks are located in the middle of human body, and are the key hub of meridian qi and blood operation. They are the main switch of six meridian and also the bridge connecting the upper jiao qi and lower jiao qi blood operation of human body. They are the important factor of the body management whose S curve increases the feminine charm.

Effects

1. Improve blood circulation and speed up metabolism.
2. Improve menstrual pain, irregular menstruation, abnormal leucorrhea and other female diseases.
3. Improve sleep quality, improve female sexual function and tighten the vagina.
4. Activate the warm nest function, stimulate the secretion of glands, and increase the feelings between husband and wife.
5. Make complexion ruddy, fade color spots and return to young state.
6. Shape hips, improve sagging and outward expansion of hips, tighten skin and increase elasticity.

Applicable Range

1. Those with saggy hips and fat accumulation.
2. Those with stretch marks induced by obesity and obesity.
3. Those whose butt shape is not good-looking, flat and soft outward expansion.
4. Those with cold and cool hips with low hip temperature.
5. Those with dysmenorrhea, dysmenorrhea, irregular menstruation, gynecological inflammation and other problems.
6. Those with decreased estrogen levels and poor sex lives.

Inapplicable Range

1. Those who are in menstruation, pregnancy, lactation and operation recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc
3. Those with skin diseases, patients with infectious diseases and skin sensitive period.
4. Those with wounds during recovery from surgery.
5. Those with allergic and severely sensitive skin.
6. Those who have just had liposuction.
7. Those who are overly old.

Notes after Operation

1. Keep hips warm and avoid wearing miniskirts and shorts.
2. Shower 4-6 hours after operation.
3. Drink plenty of hot water to avoid cold air.
4. Avoid staying up late, drinking and overeating.
5. Avoid eating raw, cold and spicy foods and get enough sleep.
6. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.

Body-shaping: Legs

Effects

1. Tighten skin and prevent sagging.
2. Stimulate collagen regeneration and improve fat lines.
3. Increase leg circulation and detoxification, metabolism.
4. Promote blood circulation and remove blood stasis, dredge channels and collaterals, and prevent varicose veins.
5. Tighten excessive protruding fat on legs and rid thick thighs.

Applicable Range

1. Those with poor circulation of lower limbs, edema and obese people.
2. Those who have low immunity and feel discomfort and pain all over the body are prone to colds.
3. Those with constipation and has coarse, flabby skin.
4. Those whose channels and collaterals of the legs are blocked, and the proportion of the legs is poor.

Inapplicable Range

1. Those who are in menstruation, pregnancy, lactation and operation recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc
3. Those with skin diseases, patients with infectious diseases and skin sensitive period.
4. Those with wounds during recovery from surgery.
5. Patients with severe varicose veins and tumors.
6. Those with allergic and severely sensitive skin.
7. People who have just had liposuction.
8. Those who are overly old.
9. Those who are pregnant or recovering from surgery.

Notes after Operation

1. Keep warm after operation. Do not eat cold food.
2. 4-6 hours later in the shower.
3. Drink more warm water to replenish water and speed up metabolism.
4. Refuse to overeat or stay up late.
5. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
6. Wear pants when you're done. Avoid miniskirts and shorts.

Body Sculpting & Slimming

Effects

1. Relieve shoulder and back soreness, and improve the lump on back of the neck.
2. Dredge channels and collaterals, and improve channels and collaterals blocking.
3. Accelerate blood circulation and dredge channels and collaterals.
4. Promote blood circulation and remove blood stasis, dredge channels and collaterals, and prevent varicose veins.
5. Regulate the functions of viscera and strengthen the body.
6. Improve cold hands, feet, uterus and body of female.
7. Stimulate collagen regeneration, firming and shaping.
8. Improve loose and soft waist and abdomen skin.
9. Improve the jiggly and thick arm.
10. Improve unwanted fat on the back, waist and abdomen.
11. Tighten skin and prevent sagging.
12. Improve stretch marks, fat marks, and increase skin elasticity.
13. Accelerate metabolism, improve constipation and intestinal peristalsis.

Applicable Range

1. Those with cold hands, feet, womb and body.
2. Those who are sedentary and don't look good at the waist.
3. Those with unwanted fat on waist and abdomen, prominent belly, loose skin after delivery.
4. Those with fat marks, stretch marks.
5. Those with constipation, abdominal meridian blocking.
6. Those with big arms, thick back, poor leg shape ratio and want to look better in clothes.
7. Those with jiggly and flabby arms.
8. Those with loose skin on arms and legs.
9. Those whose hands and legs are easily sore and numb.
10. Those with sore shoulders and backs and stiff necks.
11. Those with insomnia and dream, memory decline of the population.
12. Those who are prone to fatigue, drowsiness and poor circulation of qi and blood.

13. Those with thick back and unattractive clothes.
14. Those with poor circulation of lower limbs, edema and obese people.
15. Those who have low immunity and feel discomfort and pain all over the body are susceptible to colds.

Inapplicable Range

1. Those who are during pregnancy, menstruation and lactation.
2. Those with hypertension, heart disease, diabetes, thyroid, epilepsy and severe hyperthyroidism.
3. Those with malignant tumor and hemophilia or severe bleeding.
4. Those who have just had liposuction.
5. Those with severe varicose veins and tumors.
6. Those whose surgical wound is healing or recovering.
7. Those with skin diseases and infectious diseases.
8. Those with severe gynecological diseases.
9. Those with gynecological diseases being treated.
10. Those with allergic and severely sensitive skin.
11. Those with skin trauma or cut.
12. Those who are overly old.
13. Those with metal implants in the body, such as stents and pacemakers, are allergic to metal.
14. Those who are too weak.
15. Those who over-drink, full, empty, thirsty and overworked.

Notes after Operation

1. Avoid wearing clothes that show navel, shoulder or back, miniskirts or shorts, and avoid getting cold from the wind.
2. Avoid overeating, staying up late, drinking and eating raw, cold, spicy and greasy food.
3. Keep warm after operation, and do not bathe until 4-6 hours apart.
4. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
5. Before going to bed at night, you can knead your abdomen clockwise with your hands to make the metabolic effect of weight loss better.
6. Drink plenty of hot water to hydrate and speed up metabolism.

Overall Facial Care

Effects

1. Deeply introduce nutrition and boost absorption.
2. Shrink pores, calm skin, reduce swelling and redness.
3. Activate cells and increase skin elasticity.

4. Improve soft, sagging and inelastic skin.
5. Lift cheeks to make them look more three dimensional.
6. Improve fine lines, canthus lines, crow's feet, dynamic wrinkles and fine lines on neck.
7. Improve facial outline, lift and tighten skin.
8. Stimulate collagen regeneration and prevent skin aging.
9. Relieve eye fatigue, dark circles, bags and edema.
10. Fade pigmented spots and accelerate blood circulation.
11. Accelerate metabolism and prevent hyperpigmentation.
12. Improve sagging, coarse and dull neck skin.
13. Reduce double chin.
14. Accelerate lymphatic detoxification and improve facial skin quality.
15. Prevent neck and lymphatic diseases.

Applicable Range

1. Those with sagging, soft and inelastic skin.
2. Those with dry, dehydrated and malnourished skin.
3. Those with uneven, dull and sallow skin tone.
4. Those with rough skin and large pores.
5. Those with skin relaxation, sagging and swelling after child delivery.
6. Those with poor absorption and poor metabolism.
7. Those with indistinct facial contour.
8. Those with wrinkles, fine lines, dynamic wrinkles and nasionobial folds.
9. Those who are under long-term exposure of UV radiation at workplace.
10. Those with wrinkles, fine lines, bags under the eyes, and dark circles under the eyes.
11. Those with dry eyes and eyestrain.
12. Those who face computer screens and cell phones for long hours.
13. Those who often stay up late with dark circles under their eyes.
14. Those who always be in a dry or hot environment.
15. Those with fine lines and wrinkles on necks.
16. Those who often bow their heads.

Inapplicable Range

1. Those who have just undergone plastic surgery which have implanted prostheses, metal materials and so on.
2. Those who recently experienced certain injections, such as hyaluronic acid, intradermal injection, wrinkle removal, or plastic surgery.
3. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
4. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.

5. Those with skin trauma or wound.
6. Those who are overaging.
7. Those in menstruation, pregnancy, lactation, surgical recovery.
8. Those with skin diseases and infectious diseases.
9. Those who are unrealistic about the effects.

Notes after Operation

1. Do not wash your face and neck with overheated water (with either warm or cold) within 7 days.
2. Strengthen hydrating and moisturizing, avoid overexposure to sunlight, and always apply sunscreen.
3. It is recommended to apply facial mask + eye mask + neck mask at least 3 times a week.
4. Within 1-3 days after operation, it is best not to use products such as alcohol, AHAs or exfoliating scrub and so on.
5. Avoid sauna, hot spring or strenuous exercise within 7 days after operation.
6. Avoid spicy and greasy food, staying up late, smoking and drinking; eat more fruits and vegetables, less greasy food.
7. Avoid food that will lead to "3Hs"(hypertension, hyperglycemia, hyperlipidemia), and mainly consume less greasy food.
8. Keep your neck warm.
9. Apply neck serum or cream.
10. Avoid head-bowing for long hours.





Part II

1. Detailed Operation

After checking that the instrument is connected properly, the following starting interface will appear when the power supply is turned on.



1.1 Function Selection



1.2 Detailed Operation: 40K



Select to go to the interface below



40K Time




Time -




Time +



40K Energy



Energy +/-



40K Consecutive Mode (Constant working mode is suitable for those with thick layers of fat and want to lose weight.)



40K Intermittent Mode (Suitable for those who wish to consolidate the effect and lose weight while enjoy the process.)

1.3 Detailed Operation: Body Vacuum&RF




Select to go to the interface below



Vacuum&RF Time



Vacuum&RF Energy



Energy +/-



RF Mode Selection (NOR is the set mode. The red LED is on and energy is at the set level and stay at the set level after turning on the device. PRO is the smart mode. The temperature goes up slowly to the set energy and red LED is on after the probe comes in full contact with the skin.)



Vacuum Suction Time (Select suction time subjecting to your comfort level. It is advised to start from the minimal number and add up slowly.)



Vacuum&RF Suction Time -



Vacuum&RF Suction Time +



Vacuum Release Time (Release time refers to the interval period between suction and release. 0 means direct suction.)



Vacuum&RF Release Time -



Vacuum&RF Release Time +



Vacuum&RF Mode



M1: Direct Suction Mode (Suitable for those who opt for detoxification and scraping)



M2: Slow Intermittent mode (Suitable for those who opt for detoxification through massage and slight petechiae)



M3: Fast Intermittent Mode (Suitable for those who opt for detoxification through massage and slight petechiae)



M4: Faster Intermittent Mode (Suitable for those who opt for detoxification through massage and accelerate the circulation)



M5: Fastest Intermittent Mode (Suitable for those who opt for relaxation through massage)

1.4 Detailed Operation: Eye RF




Select to go to the interface below




RF Time




RF Energy




Energy +/-



RF Mode Selection (NOR is the set mode. The red LED is on and energy is at the set level and stay at the set level after turning on the device. PRO is the smart mode. The temperature goes up slowly to the set energy and red LED is on after the probe comes in full contact with the skin.)



Red Light Display (The icon turns gray when red light's off)




Click to start

1.5 Detailed Operation: Face Ultrasound




Select to go to the interface below




Ultrasound Time



Ultrasound Energy



Ultrasound Energy +/-



Click to start

1.6 Detailed Operation: Face RF




Select to go to the interface below




RF Time




RF Energy




Energy +/-



RF Mode Selection (NOR is the set mode. The red LED is on and energy is at the set level and stay at the set level after turning on the device. PRO is the smart mode. The temperature goes up slowly to the set energy and red LED is on after the probe comes in full contact with the skin.)



Red Light Display (The icon turns gray when red light's off)



Click to start

1.7 Detailed Operation: Cold&Hot Hammer




Select to go to the interface below




Cold & Hot Hammer Time



Cold & Hot Hammer Energy



Cold & Hot Hammer Energy +/-




Modes Selection of Hot & Cold Hammer



Hot Mode(open up pores, heat up skin and lead-in nutrition.)



Cold Mode(shrink pores, calm skin and repair sensitive skin.)



Click to start

1.8 Detailed Operation: Body RF




Select to go to the interface below



RF Time




RF Energy




Energy +/-



RF Mode Selection (NOR is the set mode. The red LED is on and energy is at the set level and stay at the set level after turning on the device. PRO is the smart mode. The temperature goes up slowly to the set energy and red LED is on after the probe comes in full contact with the skin.)



Red Light Display (The icon turns gray when red light's off)




Click to start

1.9 Detailed Operation: BIO




Select to go to the interface below



BIO Time



BIO Energy



BIO Energy +/-



BIO Mode



BIO Energy +



BIO Energy -

2. Technical Parameters

Power supply Input:100V-240V

Power :190W

40K Unoisetion Cavitation 2.0 Handle

Frequency:40KHz

Power:20W

3-Polar RF Handle for eye around with Photon

Frequency:1MHz

Power:30W

Quadrupole RF Handle for face with Photon

Frequency:1MHz

Power:50W

Sextupolar RF Handle for body with Photon

Frequency:1MHz

Power:50W

Vacuum&Bipolar RF Handle for body with Photon

Frequency:1MHz

Power:10W

Vacuum: <-80kPa

Pressure: >250kPa

Air flow: >10L/minute

Noise level: <70dB (30cm away)

face ultrasound handle

Frequency: 1KHz

Power:20W

hot&cold handle

Temperature:0~5 degree centigrade

Bio face lift probe with Photon

Frequency:1MHz

Power:30W

3. Safety Precautions

Patients with the following symptoms should be cautious about using this instrument. Please consult a doctor or a professional before using the instrument. The details are as follows:

1. Those who are pregnant or lactating.
2. Those with heart disease or who are equipped with pacemaker.
3. Those whose wounds have not healed and who are recovering from the operation.
4. Those with epilepsy, severe diabetes and hyperthyroidism.
5. Those with malignant tumor, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Use with caution for those who are allergic to electric currents.

4. Dos and Dont's

1. Use a grounding pin and make sure the power socket of the device is properly connected before each operation.
2. Ensure that the voltage of the device is appropriate. If the voltage of the local power supply is unstable, we recommend that users increase to the matching power voltage.
3. In order to ensure the therapeutic effect and normal service life expectancy of the device, please use the specified accessories provided or recommended by the original manufacturer.
4. The device should not be placed in wet environment or near water, nor should they be directly exposed to sunlight.
5. Do not place the device near a strong heat source, because this may affect the life and normal use of the device.
6. Before treatment, please remove all metal objects from the body first, so as to avoid unexpected conditions and influence the curative effect.
7. Please restrain from using the device on eyes, thyroid, parathyroid, testicles, abdomen of pregnant women, and pacemakers.
8. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
9. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use, so as to ensure the safety of the electrical products.
10. Use device or train device operators in strict accordance with instructions in the user manual.
11. When taking other weight-loss drugs, it is recommended to stop taking them for 1 to 2 months before treatment. Extension of the course of treatment is recommended if you wish to lose weight immediately.
12. Don't use the instrument with an empty stomach. After a full meal, wait at least 1 hour before a course of weight loss treatment.
13. When operating, the instrument should fully touch the skin to avoid uneven heat.
14. Don't stay at one place for over 5 seconds to avoid scald.
15. Brief customers what they might feel before BIO and 40K treatment to avoid unexpected situations.
16. Start from the lowest energy level and slowly add up.
17. When using this device, the operating parts must be kept moist and dry skin should be avoided.
18. Clean the instrument with normal saline after operation to ensure its cleanliness and hygiene and prolong its service life.

5. Troubleshooting & Solutions

1. The instrument cannot be started, but the key light on the back of the instrument does not work?

- A. Make sure the power cord is connected to a valid power socket.
- B. Whether the fuse tube on the back of the instrument is loose or burnt out.

2. No RF output of the instrument?

- A. Please check whether the plug connecting the instrument handle and the body is tightly connected.
- B. Please check whether the treatment area has been cleaned. Grease or oil essential oil products may cause poor contact between the treatment handle of the instrument and human body, resulting in no rf output.

3. Weakened RF output?

- A. Please check whether the plug connecting the instrument handle and the body is tightly connected.
- B. Please check whether non-conductive grease and other substances are stained on the treatment handle, which will cause poor contact and weaken the output.
- C. Please check whether the product used is the adaptive product specified by the instrument.

4. No suction or very little suction?

- A. Please turn off the instrument and check the filter element specially used for the instrument, which may need to be replaced.
- B. Please check the rubber ring of the handle plug, because the rubber ring may be worn and cause air leakage.
- C. Please check whether the oil filter cup outside the filter is tightened and whether the rubber ring is worn. It is possible that air leakage at this position may lead to poor air pressure.
- D. If the above methods cannot address the issue, please contact the instrument dealer for assistance.

5. I can start the instrument, but the screen show's error message?

- A. Please unplug the power plug behind the instrument and wait for about 1 minute before restarting the instrument.
- B. If the above methods cannot be handled, please contact the instrument dealer for assistance.

6.FAQs

1. Q: How long can I start seeing results of RF treatment?

A: Under normal circumstances within that day or a week. Skin collagen tissue is heated and produces contractions, and you can obviously feel skin tightening. Because RF works by stimulating subcutaneous tissue and bring about sustainable collagen regeneration, so the more you do it, the more visible effects you shall receive.

2. Q: Is RF harmful to skin?

A: RF treatment falls under the category of non-surgery operation. It stimulates collagen regeneration at the skin bottom and boosts metabolism, therefore has no damage to the skin. Partial redness and fever is the normal phenomenon of accelerated blood circulation, and will be resolved by oneself after a moment, so there's no need to worry.

3.Q: How long does the body firming treatment take?

A: One operation time is 60 minutes. We will combine professional techniques and instruments to achieve remarkable results.

4.Q: What are all the functions of this instrument?

A: Dissolve fat, body shaping, smooth skin, firming, lifting, and anti-aging. It can be operated all over the body. In terms of facial care, it can contour facial outline, alleviate wrinkles and saggy skin. In terms of body care, it can reduce regional fatness and build S curve, promote the metabolism and detoxification of the whole body. Meanwhile, it can also enhance the functions of viscera and regulate the sub-health of the body.

5.Q:Which one is better in terms of weight loss, liposuction or this machine?

A: The principle of liposuction is to suck out the excess fat in one part of the body through the method of negative pressure suction, so as to achieve the purpose of rapid weight loss on parts of the body. It achieves fast results, requires anesthesia, recovery period. It falls under the category of surgery operation and comes with risks. However, our machine is risk-free and has no side effects. While losing weight, it can stimulate collagen regeneration to lift and tighten skin. By lifting, it can also shape perfect curves and accentuate women's natural charm.

6.Q: Will I experience rebound after operation?

A: RF works by dissolving fat to reduce weight, and are not prone to rebound. Because what's dissolved through RF treatment is fat, not water. The formation of fat needs a longer period of time to accumulate, so it is not prone to rebound.

7.Q: Do I need to be on a diet?

A: Pay attention to what you eat. Because radio frequency and ultrasonic treatment end with blasting adipose and accelerated metabolism. It's best to avoid to eating spicy, oily, or fried food, to avoid affecting regular metabolism. Also proper workout and sweating help with the process, in a way weight loss will be more visible as well.

8.Q: Does it have any side effects on the body?

A: RF beauty treatment is non-invasive and so far the most safe and effective method to remove wrinkles and shape body. Generally speaking side effects won't appear. A small number of people may experience transient redness or swelling that will disappear after a few hours. Those with dry skin may experience atrophy after the

initial treatment. The skin will lose moisture due to RF heat. However, the skin will become plump at the beginning of collagen regeneration and these symptoms will disappear after three days. There are no side effects on health.

9.Q: Can RF instrument be used for chest shaping?

A: As the growth of age and female physiology period, sub-health state, gravity, stimulation during lactation, a large number of nutrition loss... All the above will cause collagen loss, decreased blood flow, damaged elastic tissue, flabby ligament gland, and prolapsed breasts. Through RF energy conduction, collagen protease can be promoted to repair elastic tissue and ligament glands, thus achieving the effect of chest lifting and shaping.

10. Q: Can ultrasound be operated all over the body?

A: Ultrasonic cavitation ruptures fat cells and causes mechanical disruption of fat cell membrane. Liquefied fat cells then cleared away by metabolic pathways. Our heart is a hollow organ, and sensitive to high frequency sound waves. Back and forth reflex occur due to the fact that muscle tissue and blood in the heart do not conduct sound waves in the same way. The reflex force can detach cardiac valves from cardiac muscles. If directed at the eye, it can cause retinal detachment. In a word, avoid eye and heart area when using ultrasonic beauty device. (It is also advised not to use it on waist, back and chest.)

11.Q: Does ultrasound have side effects on human body?

A: Ultrasonic beauty machine requires no surgery, no opening up and no anesthesia, therefore is indeed non-invasive. It relies on ultrasound (mechanical sound wave) - blast fat through cavitation - ultrasound focusing. And the results are fat smash - fat dissolving - skin firming and body sculpting respectively. It only targets at low density fat cells while leaving high density tissues like blood vessels intact, therefore having no side effects on human body. Slight tinnitus might appear to happen during the treatment, which is within normal range so there's no need to worry.

12.Q: Why does tinnitus occur?

A: It's because ultrasound comes with a high frequency vibration, which could go up to 20kHz. It also goes deep into skin fat layers, which could go down to 20mm below epidermis. Highly intensifies ultrasound cause high speed friction among fat cells, resulting in them bursting and heating up, and eventually emulsified. And during the process, ultrasound can cause microvibration to the tissue that make us feel like as "tinnitus".

13.Q: What is collagen?

A: Collagen is a kind of biological macromolecule substance, a kind of white, opaque and non-branched fibrous protein. It can supplement the nutrition needed by all layers of the skin, enhance the collagen activity in the skin, lock moisture, nourish the skin, delay aging, beautify, relieve facial relaxation, provide nutrients for the hair and

so on. Collagen is a nutrient that human body must replenish to delay aging. With the growth of age, collagen will gradually be lost. At the age of 20, women have begun to age, collagen gradually decrease. At the age of 25, collagen loss reaches its peak. At the age of 40, the content is less than half of which one has at age 18. It is the loss of collagen and water that causes the wrinkles on the face of the elderly. The breakage of collagen fibers and elastic mesh will cause skin tissue to be oxidized, atrophied and collapsed. The skin will be dry, wrinkled, flabby and inelastic etc,. Therefore, in order to delay aging, collagen must be supplemented.

14. Q: Why do I need breast maintenance?

A: The lymph nodes of bosom rank the top among other parts of the body, which also lead to toxin accumulation very easily, and urbanite life rhythm is very fast now, working pressure is very big, which can bring about the hyperplasia of different level. The saying goes that bosom is the cradle of children, the garden of the woman, the symbol of female identity! Breast maintenance can make your figure more beautiful, add scores for your whole temperament, thus turn head rate increases 10 times certainly! Good figure and temperament will also ensure husband's loyalty!
The evolution of breast cancer: breast milk residual, secretions, toxins -- nodules and blockages -- lobular hyperplasia, ductal hyperplasia, breast hyperplasia -- duct blockage and adhesion -- fibroma, breast cyst -- breast cancer.

15.Q: Why do I need hip maintenance?

A: Because it can help the body lymphatic system discharge poison, improve the gynaecological disease, and avoid harm brought by buttock blockage. Hip impassability surely bring about gynaecological issues. What's inside buttocks is pelvic cavity, bowel. It connects daimai, lumbar vertebra, sciatic nerve from above, uterus ovary, uterine adenexa in front, anus, vagina, inguinal lymphatic from below.
If suffering from squeeze of shangjiao and xiaojiao meridian channels, hip would be most vulnerable to cold, wet and blood stasis. When hip is affected by cold, it can cause meridian contraction of hips, dysmenorrhea, irregular menstruation, dark menstruation color, clot, deterred blood flow and meridian channels. It also affects fertility.


16.Q: How does negative pressure detoxify?

A: Negative pressure can congest capillary, stimulate cell to increase vitality. The suction and release of pressure during operation can make local pores continue to open and close, therefore promote skin breathing, increase the amount of skin oxygen, speed up the elimination of waste. So it achieves the effects of promoting qi and blood circulation, dispelling coldness and dampness, dredging channels and collaterals, clearing toxins, heat and so on.


7. Packing List

- 1 x Main Machine
- 1 x 40K Unoisetion Cavitation 2.0 Handle with Photon
- 1 x 3-Polar RF Handle for eye around with Photon
- 1 x Quadrupole RF Handle for face with Photon
- 1 x Sextupolar RF Handle for body with Photon
- 1 x Vacuum&Bipolar RF Handle for body with Photon
- 1 x Face ultrasound handle
- 1 x Hot&cold handle
- 1 x Bio face lift probe with Photon
- 2 x Holder
- 1x Power Cord

8. Operational Diagrams

Parameter Adjustment	Product	Techniques	Diagrams
Facial Anti-aging: 60 min, once/week			
RF Energy: 30%-80%	Makeup remover+	1. Remove makeup and cleanse face, 5 min.	Technique6
Mode: PRO	Face wash+	2. Apply toner, 1 min.	
Ultrasound Energy: 30%-80%	Cleaning base mask+	3. Apply cleaning base mask+hot steam, 10 min.	
	Cold&Hot steam+	4. Wash face clean, 3 min.	
	Essence+	5. Apply massage cream evenly on face and caress face, 3 times.	
	Device+	6. Push the following acupoints (Ren-24, Du-26,	Technique11,15
	Facial mask		



<p>Mode: Consecutive Cold&Hot Hammer Energy: 30%-80%</p> <p>Mode: Cold</p>		<p>St-4, St-6, LI-20, BI-1, BI-2, EM3, TE-23, Gb-1, the Temple, St-1, St-2), 3 times.</p> <p>7. Caress the whole face, 3 times.</p> <p>8. Use both hands to lift face as if plucking the string, 1 min.</p> <p>9. Alternate hands in lifting the face as if plucking the string on one side, 1 min.</p> <p>10. Do the same on the other side.</p> <p>11. Alternate both hands to lift in 3 lines from chin to ear lobe, corner of mouth to ear gate, nose wing to temple, lower eyelid to temple, lift canthus, 3 times.</p> <p>12. Lift towards hairlines on forehead, 3 times.</p> <p>13. Do the same on the other side.</p> <p>14. RF Operation: Move the probe from jaw to earlobe, mouth corner to ear gate, nose wing to temple in circles, 3 times.</p> <p>15. Brace the skin with the probe and another free hand to keep the treatment surface tight and lifted; move from jaw to earlobe, mouth corner to ear gate, nose wing to temple, lower eyelid to temple to lift the treatment area, 3 times.</p> <p>16. Combined with hands to lift towards hairline on forehead, 3 times.</p> <p>17. Do the same on the other side.</p> <p>18. Clean face, 3 min.</p> <p>19. Apply essence evenly on</p>	<div data-bbox="1145 228 1289 452" data-label="Image"> </div> <p data-bbox="1056 488 1264 519">Technique12,16</p> <div data-bbox="1145 564 1289 788" data-label="Image"> </div> <p data-bbox="1056 824 1225 855">Technique14</p> <div data-bbox="1145 878 1289 1102" data-label="Image"> </div> <p data-bbox="1056 1115 1225 1146">Technique20</p> <div data-bbox="1145 1182 1289 1406" data-label="Image"> </div> <p data-bbox="1056 1442 1225 1473">Technique21</p> <div data-bbox="1145 1518 1289 1742" data-label="Image"> </div> <p data-bbox="1056 1778 1225 1809">Technique23</p>
--	--	---	--





		<p>the face, 1 min.</p> <p>20. Ultrasonic Operation: Move the probe from jaw to earlobe, mouth corner to ear gate, nose wing to temple in circles, 3 times.</p> <p>21. Move the probe horizontally in circles on forehead, 3 times.</p> <p>22. Do the same on the other side.</p> <p>23. Apply hydrating mask+cold hammer to cover the whole face in circles, 3 times.</p> <p>24. Remove the mask and wash face clean, 3 min.</p> <p>25. Apply toner, essence, cream and sunscreen.</p>	
--	--	---	---

Suggested Treatment:

10 treatments account for one full course. After one time, skin feels heated and appears to be clean and glossy. After a full course, skin is translucent, smooth, fine and full nourished with increase collagen regeneration. Two courses help increase skin elasticity, increase skin metabolism, detoxification, fade pigmented spots and brighten dull skin. Three courses help consolidate the effect and prevent skin sagging and laxity.

Facial Tightening and Anti-wrinkle: 60 min, once/week

<p>RF Energy: 30%-80%</p> <p>Mode: PRO</p> <p>BIO Energy: 30%-80%</p> <p>Mode: Consecutive</p>	<p>Makeup remover+ Face wash+ Cleansing facial mask+ Essence+ Device+ Facial mask</p>	<ol style="list-style-type: none"> 1. Remove makeup and cleanse face, 5 min. 2. Apply toner, 1 min. 3. Apply cleaning base mask+hot steam, 10 min. 4. Wash face clean, 3 min. 5. Apply massage cream evenly on face and caress face, 3 times. 6. Push the following acupoints (Ren-24, Du-26, St-4, St-6, LI-20, BI-1, BI-2, EM3, TE-23, Gb-1, the Temple, St-1, St-2), 3 times. 7. Caress the whole face, 3 	<p>Technique6</p>  <p>Technique11</p>  <p>Technique12</p>
--	---	---	---



	<p>times.</p> <p>8. Use both hands to lift face as if plucking the string, 1 min.</p> <p>9. Alternate hands in lifting the face as if plucking the string on one side, 1 min.</p> <p>10. Do the same on the other side.</p> <p>11. Alternate both hands to lift in 3 lines from chin to ear lobe, corner of mouth to ear gate, nose wing to temple, lower eyelid to temple, lift canthus, 3 times.</p> <p>12. Lift towards hairlines on forehead, 3 times.</p> <p>13. Do the same on the other side.</p> <p>14. RF Operation: Move the probe from jaw to earlobe, mouth corner to ear gate, nose wing to temple in circles, 3 times.</p> <p>15. Brace the skin with the probe and another free hand to keep the treatment surface tight and lifted; move from jaw to earlobe, mouth corner to ear gate, nose wing to temple, lower eyelid to temple to lift the treatment area, 3 times.</p> <p>16. Combined with hands to lift towards hairline on forehead, 3 times.</p> <p>17. Do the same on the other side.</p> <p>18. Clean face, 3 min.</p> <p>19. Apply essence evenly on the face, 1 min.</p> <p>20. BIO Operation: Move the probe from jaw to earlobe, mouth corner to ear gate,</p>	 <hr/> <p>Technique14,20</p>  <hr/> <p>Technique15,21</p>  <hr/> <p>Technique16,22</p>  <hr/>
--	--	--







		<p>nose wing to temple in circles, 3 times.</p> <p>21. Move the probe horizontally in circles on forehead, 3 times.</p> <p>22. Do the same on the other side.</p> <p>23. Apply hydrating mask+cold hammer to cover the whole face in circles, 3 times.</p> <p>24. Remove the mask and wash face clean, 3 min.</p> <p>25. Apply toner, essence, cream and sunscreen.</p>	
--	--	---	--

Suggested Treatment:

10 treatments account for one full course. After one time, the skin will be tightened. After a course of treatment, the skin will be plump, firm and ruddy. 2 courses help increase skin elasticity, accelerate metabolism and detoxification, lighten dark spots, brighten skin, and lower the rate of skin sensitivity. 3 courses help stabilize skin state, prevent skin aging and laxity, and help define facial outlines.

Anti-aging around Eyes: 60min, once/week

<p>RF Energy: 30%-80%</p> <p>Mode: PRO</p> <p>BIO Energy: 30%-80%</p> <p>Mode: Intermittent</p>	<p>Makeup remover+ Face wash+ Toner+ Eye essence+ Device+ Eye mask</p>	<ol style="list-style-type: none"> 1. Remove makeup and clean face, 5 min. 2. Apply toner, 1 min. 3. Apply essence evenly around eyes and gently massage to rub it in, 1 min. 4. Push the following acupoints(BI-1, BI-2, EM3, TE-23 , the Temple, Gb-1, St-1) using middle and ring fingers, 3 times. 5. Move hands in circular motion on eyes and push temples, 3 times. 6. Move middle and ring fingers horizontally below lower eyelids and slide to back of the ear, 3 times. 7. Lift middle and ring fingers 	<p>Technique3,5,10</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique6</p>
---	--	---	--






		<p>from inner eye corner to upper eyelids and slide to back of the ear, 3 times.</p> <p>8. Lift eye corners with “scissor” hand gesture, 3-5 times.</p> <p>9. Repeat the previous operation on the other side, 3 times.</p> <p>10. Move hands in circular motion on eyes and push temples, 3-5 times.</p> <p>11. RF operation: move the probe beneath lower eyelid in small circles and move all the way till the temple, 3-5 times.</p> <p>12. Combined with hand, move the probe from lower eyelid to temple to lift the treatment area, 3-5 times.</p> <p>13. Lift from the beginning of the eyebrows to the end, 3-5 times.</p> <p>14. Combined with hand, move the probe to lift the corner of the eyes towards the hairline, 3-5 times.</p> <p>15. Repeat the same technique with BIO operation.</p> <p>16. Do the same on the other side.</p> <p>17. Apply eye mask, 15 min.</p> <p>18. Wash eyes and face clean, 2 min.</p> <p>19. Apply toner, facial essence, eye essence, cream and sunscreen.</p>	 <p>Technique7</p>  <p>Technique8</p>  <p>Technique11</p>  <p>Technique12,14</p>  <p>Technique13</p> 
--	--	---	--

Suggested Treatment:

10 treatments account for one full course. Once done, the eyes are lifted and tightened, increasing blood circulation to the eyes. After one course, it lightens fine

lines, dark circles under the eyes and brightens the skin. After 3 courses, the skin around the eyes is firm, tender and shiny, improving and preventing the aging of the eyes.

Anti-aging on Neck: 60 min, once/week





<p>RF Energy: 30%-80%</p> <p>Mode: PRO</p> <p>BIO Energy: 30%-80%</p> <p>Mode: Intermittent</p>	<p>Makeup remover+ Face wash+ Toner+ Massage cream+ Essence+ Device+ Neck mask</p>	<ol style="list-style-type: none"> 1. Remove makeup and cleanse face, 5 min. 2. Apply toner, 1 min. 3. Neck massage: apply and rub oil on the chest with both hands to the back of the neck and press Gb-20 and Du-16 acupoint, 3 times. 4. Alternately lift the lower jaw and double chin to pass behind the ear and under the armpit, 3 times. 5. Push downwards with hand spread and hukou pointing down, tracing 3 meridian channels from the neck side to the armpit, 3 times. 6. Push downwards with 4 fingers, tracing 3 meridian channels from neck side to the armpit, 3 times. 7. Push downwards with kneeling fingers till the neck turns hot and slide to the armpit, 3 times. 8. Push downwards with hand spread and hukou pointing down from external collarbone to armpit, 3-5 times. 9. Do the same on the other side. 10. Wash neck clean. 2 min. 11. Apply essence evenly on neck, 1 min. 12. RF operation: coupled 	<p>Technique3</p>  <hr/> <p>Technique4,12</p>  <hr/> <p>Technique5,6,7,13</p>  <hr/> <p>Technique15</p>  <hr/> <p>Technique16</p> 
---	--	---	--

		<p>with hand massage, lift double chin to pass behind the ear and under the armpit, 3 times.</p> <p>13. Move the device in circles on one side of neck to pass under armpit, 3 times.</p> <p>14. Slide the device from internal and external clavicle to the armpit, 3-5 times.</p> <p>15. Move the device in circles on the whole neck area, 3 times.</p> <p>16. Coupled with hand, lift neck in lines, 3-5 times.</p> <p>17. Do the same on the other side.</p> <p>18. Repeat the same techniques with BIO operation.</p> <p>19. Apply neck mask, 15 min.</p> <p>20. Wash neck clean, 2 min.</p> <p>21. Apply essence and neck cream.</p>	
--	--	---	--

Suggested Treatments:

10 treatments account for one full course. After one treatment, skin becomes tender, smooth. After a full course, neck skin becomes smooth, speed up the neck lymphatic circulation. After 2 courses it tightens skin, fade neck lines, improve the double chin. After 3 course, skin is fine, firm and glowing, which also promote lymph drainage, improve facial dullness, some problems such as acne, and delay the aging of the skin, make skin restore young state.

Body-shaping on Waist&Abdomen: 60 min, once/week

<p>Vacuum&RF Energy: 30%-80%</p> <p>Suction Time: 0.3-1.0</p> <p>Release Time: 0-0.5</p> <p>Mode: M1</p> <p>40K Energy: 30%-80%</p> <p>Mode: Consecutive</p>	<p>Essential oil+ Gel+ Device+ Towel</p>	<ol style="list-style-type: none">1. Apply and rub oil on abdomen with hands, 3 times.2. Rub stomach back and forth with both hands, 3-5 times.3. Knead abdomen with both hands using chiropractic techniques, 3 times.4. Lift Meridian BV(Belt Vessel) on both sides of waist with both hands alternately, 16 times.5. Move hands as if writing an "8" number on waist, 3 times.6. Overlap hands and message the intestinal canal clockwise, 3 times.7. Move hands in circular motion and caress the treatment area, 3 times.8. Push on the following acupoints: Ren-13, Ren-12, Ren-10, Ren-8, Ren-6, Ren-4, Ren-3, St-25, SP-15, Ren-2, 2 times.9. Starting from Ren-3 point, push both thumbs to belly navel, slide them to waist and the lift upwards to groin, 3 times.10. Caress the treatment area with hands and slide to	<p>Technique1,7,10</p>  <hr/> <p>Technique2</p>  <hr/> <p>Technique3</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5</p>
--	--	---	---

groin.

11. Apply gel evenly on abdomen, 1 min.

12. 40K operation: starting from one side of the waist, lift the device to abdomen and groin, 3 times.

13. Repeat the previous operation on the other side, 3 times.

14. Move the device on abdomen in small circular motions, 3 times.

15. Move the device on abdomen in big circular motions, 3 times.

16. Move the probe to shape the waist as if writing "8", 3-5 times.

17. Move the probe up and down to shape the waist, 3-5 times.

18. Combined with hand, lift Meridian BV on both sides of the waist, 16 times.

19. Wipe clean abdomen with hot towel, 2 min.


20. Apply essential oil evenly on the abdomen.

21. Repeat the same technique with vacuum&RF operation.

22. Clean with fresh towel.



Technique6




Technique8



Technique9







Technique12



Technique14








Technique15






			 <hr/> <p>Technique16</p>  <hr/> <p>Technique17</p>  <hr/> <p>Technique18</p>  <hr/>
--	--	--	---

Suggested Treatment:

Ten treatments account for one full course of treatment. Once done, abdomen is heating up, metabolism is accelerated. After one course, fat slowly reduce through metabolism, body starts to slim down. Constipation and other issues are addressed at the same time. After 2 courses the effect gradually becomes obvious, skin is tightened and lifted, skin collagen increased. Stretch marks and flabby skin are also improved. 3 courses help reinforce the effect. Redundant adipose gradually disappear, lumbar abdomen curve start to take in shape, and slim waist starts to show.

Body-shaping on Breasts: 60 min, once/week






<p>RF Energy: 30%-80%</p> <p>Mode: PRO</p>	<p>Massage cream(essen tial oil)+ Device</p>	<ol style="list-style-type: none"> 1. Stand at the head of the bed, both hands from Ren-17, zhong point began to rub oil to armpit up at the same time lift the suspension ligament (this is to caress treatment area) 3 times. 2. Both hands thumb point: Ren-17, St-18, SP-21, St-16, LU-1, LU-2, 3 times. 3. Caress the treatment area, 3 times. 4. Overlap both palms from Ren-17 point to chest both sides back and forth as if writing an "8" number, 3 times. 5. Lift the accessory breast to the suspensory ligament alternately with both hands, first left, then right, 10 times. 6. Caress the treatment area, 3 times. 7. Sit beside the customer, first left, then right, with both hands to soothe and lift the chest, 3-5 times. 8. Alternately push the breast duct around the chest (one circle of the chest), 3-5 times. 9. Use the thumb to dredge the nodules in the chest, 3 times. 10. Lift the whole chest with hands, 3-5 times. 11. The technique on the other side is the same as above. 12. RF operation: first left, then right, from the bottom 	<p>Technique1,3,6,10</p>  <hr/> <p>Technique2</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5</p>  <hr/> <p>Technique7</p>  <hr/> <p>Technique8,12</p>
--	--	--	---



		<p>to the top along the direction of the chest to the nipple, combined with hand pacification, 5-8 times.</p> <p>13. Draw half a circle along the chest to lift and pull from the bottom to the collarbone, combining with hand comfort, 5-8 times.</p> <p>14. Make small circles to clear the nodule site for 3-5 times.</p> <p>15. The customer lies on his side and stands up with her arms. The device burns fat in circles at the position of auxiliary breast, 5-8 times.</p> <p>16. The instrument is pushed from the position of underarm auxiliary breast to the chest (shaping and receiving auxiliary breast), 5-8 times</p> <p>17. The technique on the other side is the same as above.</p> <p>18. Clean with hot towel.</p>	 <hr/> <p>Technique9,14</p>  <hr/> <p>Technique13</p>  <hr/> <p>Technique15</p>  <hr/> <p>Technique16</p> 
--	--	---	---

Suggested Treatment:

10 times account for a course of treatment. After a full course, breasts will experience heat, blood circulation is boosted, and a certain improvement effect will be seen. After a course of treatment the effect is more obvious, nodules are improved. 3 courses help tighten the skin and consolidate shaping effect, and increase the elasticity of the breasts. Regulating endocrine makes women more attractive.

Body-Shaping on Arms: 60 min, once/week





<p>RF Energy: 30%-80%</p> <p>Mode: PRO</p> <p>40K Energy: 30%-80%</p> <p>Mode: Intermittent</p>	<p>Essential oil+</p> <p>Ultrasonic gel+</p> <p>Device</p>	<ol style="list-style-type: none"> 1. Left-right order: lay the arm flatwise, apply and rub oil in from lower arm to the entire arm and slide to the fingers, 3 times. 2. Push the entire arm with both palms, 3 times. 3. Push hands upwards tracing three main collateral channels: Large Intestine Meridian(LI) - Triple Energizer Meridian(TE) - Small Intestine Meridian(SI) to armpits with fingers spread and hukou pointing upwards, 3 times. 4. Caress the treatment area, 3 times. 5. Trace and rub three main collateral channels on arms with kneeling finger back and forth till them turn hot, 3 times. 6. Caress the treatment area, 3 times. 7. Lay inner arm upwards, and push hands tracing three yin channel on inner arm: Lung Meridian(LU) - Pericardium Meridian(PC) - Heart Meridian(HT) to armpit, with hukou pointing upwards, 3 times respectively. 8. Rub three channels back and forth with kneeling fingers, 3 times. 9. Caress the treatment area and slide to fingers. 10. Repeat the previous operation on another side. 11. The end of massage 	<p>Technique1</p>  <hr/> <p>Technique2,3,4,5,6</p>  <hr/> <p>Technique7</p>  <hr/> <p>Technique8,17,19</p>  <hr/> <p>Technique9</p>  <hr/> <p>Technique13,16</p>
---	--	--	--







	<p>technique.</p> <p>12. Apply gel evenly on the arm.</p> <p>13. RF device operation: lay arms flatwise and start from fat part of lower arm, tracing three channels to armpit, 3 times.</p> <p>14. Move the device in annular motion, tracing three channels to armpit, 3 times.</p> <p>15. You may double the operation on flabby arms, 3 times.</p> <p>16. Trace three channels to armpit, 3 times.</p> <p>17. Lay inner arm flatwise and push three channels on upper arm to armpit, 3 times.</p> <p>18. Trace three channels on upper arm and move the device in circles to armpit, 3 times.</p> <p>19. Trace three channels from upper arm to armpit.</p> <p>20. Repeat the previous operation with 40K operation.</p> <p>21. Do the same on the other side.</p>	 <hr/> <p>Technique14,15</p> 
--	--	---




Suggested Treatment:

10 treatments account for a full course. Once done, it will feel warm and tight, and the arms will be relaxed. A course of treatment help reduce excessive arm fat and the skin appear to be tight and plump. 2 courses of treatment began to shape the arms, strengthen the effect, dredge meridians, and enhance the physique. 3 courses help consolidate stability and prevent rebound.

Body-Shaping on Back: 60 min, once/week

<p>RF energy: 30%-80%</p> <p>Vacuum&RF energy: 30%-80%</p> <p>Suction: 0.3-1.5</p> <p>Release: 0-1.0</p>	<p>Essential oil+ Device</p>	<ol style="list-style-type: none">1. Apply oil on back and press Gb-20 and Du-16 points.2. (Starting from hairline) stroke Da Ban Jin with thumb, 3-5 times.3. Move thumb outwards from Bladder Meridian(BL) to the 8 Crevice Area and then to Gb-20 and Du-16 points.4. Starting from neck to caudal vertebra, move hands in S-shape, 3 times.5. Push Bladder Meridian(BL) with thumb in left-right order, 3 times.6. Push Bladder Meridian(BL) to the 8 Crevice Area with both thumbs, at the same time, 3 times.7. Push Bladder Meridian(BL) in three kneeling fingers, 3 times.8. Push scapula slot with both hands alternately in left-right order, 3-6 times.9. Push scapula slot with both hands horizontally, 3-6 times.10. Caress the back with both hands and push Gb-20 and Du-16 points, 3 times.11. Overlap both thumbs and push SI-11 point and slide to the arm, 3 times.12. Repeat the previous	<p>Technique1,10</p>  <hr/> <p>Technique2</p>  <hr/> <p>Technique3</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5,6,7,13,16</p>
--	--------------------------------------	--	---


		<p>action on another arm, 3 times.</p> <p>13. Rub Meridian GV and Bladder Meridian(BL) with hands until they turned hot.</p> <p>14. The end of massage technique.</p> <p>15. Apply gel evenly on back.</p> <p>16. RF Operation: Starting from Meridian GV to Bladder Meridian(BL), move the device from neck to the 8 Crevice area, 3-5 times.</p> <p>17. Move the device in circles on Du-14 point, and the 8 Crevice area , 3-5 times respectively.</p> <p>18. From Da Ban Jin to neck, move the device back and forth, 3-5 times.</p> <p>19. Lift the device on scapula slot in left-right order, 3-5 times.</p> <p>20. Move the device horizontally as if writing an "8" number to the 8 Crevice area, 3 times.</p> <p>21. Move the probe upwards in annular motion to armpit in left-right order on both sides, 3 times.</p> <p>22. Lift from sides of the waist upwards till armpit, 3-5 times.</p> <p>23. Clean with hot towel.</p> <p>24. Apply essential oil evenly on the back.</p> <p>25. Repeat the same technique with Vacuum&RF Operation.</p> <p>26. Clean with fresh towel.</p>	 <hr/> <p>Technique8,19</p>  <hr/> <p>Technique9</p>  <hr/> <p>Technique11,12</p>  <hr/> <p>Technique17</p>  <hr/> <p>Technique18</p>  <hr/> <p>Technique20</p>
--	--	---	--







			 <hr/> <p>Technique21</p>  <hr/> <p>Technique22</p> 
--	--	--	--

Suggested Treatment:

10 treatments account for a full course. Once done, back feel significantly relaxed, and it can also alleviate neck stiffness, and lift skin. A course help relax shoulder, reduce back fat, ease the lump on the nape of the back, and shape the back line. After two courses, the lump appears to be smaller, and it help dredge meridians, regulate sub-health, and improve sleep. 3 courses help with back thinning, define back line, strengthen viscera function, improve digestion, detoxification and metabolism, (a period of treatment is to reduce size, two courses is to strengthen effect, three courses help with consolidation and avoid rebound.)

Body-Shaping on Buttocks: 60 min, once/week

<p>RF Energy: 30%-80%</p> <p>Mode: PRO</p> <p>40K Energy: 30%-80%</p> <p>Mode: Consecutive</p>	<p>Essential oil+ Ultrasonic gel+ Device</p>	<ol style="list-style-type: none"> 1. Standing on the side, accumulate oil on hands and then slide to the waist and then lift up along the hips from the waist, 3 times (this is to caress the treatment area.) 2. Repeatedly push the 8 Crevice area with two thumbs. 3. Caress the treatment area 	<p>Technique1,4,10</p>  <hr/> <p>Technique2</p>
--	--	--	--



	<p>for 3 times and then point: Bl-23, 8 Crevice area, Du-1, Gb-30, Bl-36, 3 times.</p> <p>4. Caress the treatment area, 3 times.</p> <p>5. First left, then right, with both hands extending from the root of the thigh from bottom to top -- Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) - Meridian BV, 3 times each.</p> <p>6. Both hands follow from the root of thigh from bottom to top -- Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) - Meridian BV, 3 times each.</p> <p>7. Overlap both palms and push up together from thigh root to pulse (ascension) ,3 to 5 times.</p> <p>8. Push your hands from bottom to top from both sides of your hips to the top of your hips, 3 times.</p> <p>9. Repeat the previous operation.</p> <p>10. Caress the treatment area.</p> <p>11. The technique on the other side is the same as above.</p> <p>12. Apply gel evenly on buttocks.</p> <p>13. RF Operation: lift from the thigh root to Meridian BV in lines , 3 times.</p> <p>14. Lift 3 times from bottom to top on both sides of the</p>	 <hr/> <p>Technique3</p>  <hr/> <p>Technique5,13,16</p>  <hr/> <p>Technique6</p>  <hr/> <p>Technique8,14,17</p>  <hr/> <p>Technique15</p>  <hr/>
--	--	---







		<p>buttocks to the highest point of the buttocks.</p> <p>15. Circle your hips, 3 to 5 times.</p> <p>16. Lift 3 times from thigh root to vein.</p> <p>17. Lift from bottom to top on both sides of the buttocks to the highest point of the buttocks, 3 times.</p> <p>18. Do the same on the other side.</p> <p>19. 40K operation method is the same as above.</p> <p>20. Clean with hot towel.</p>	
--	--	--	--




Suggested Treatment:

10 treatments account for a full course. After one treatment, buttocks will experience lifting and heat. After a course of treatment buttocks start to show obvious curve lines, and excessive fat begin to disappear slowly, and skin appear to be tight. After 3 courses, the conditions of cold body and uterus, other gynecological problems such as irregular menstruation are improved. It also help consolidate the effect, and increase the female charm.

Body-Shaping on Legs: 70 min, once/week

<p>Vacuum&RF Energy: 30%-80%</p> <p>Suction: 0.3-1.5</p> <p>Release: 0-1.0</p> <p>Mode: M1</p> <p>40K Energy: 30%-80%</p> <p>Mode: Consecutive</p>	<p>Essential oil (massage cream)+ Ultrasonic gel+ Device</p>	<ol style="list-style-type: none"> 1. Rear leg, left-right order: apply and rub oil in from calf to thigh to heel, 3 times. 2. Alternately push the entire leg from bottom to top and wrap back into the heel, 3 times. 3. Move hands upwards to four main collateral channels: Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) with fingers spread and hukou pointing upwards. 4. Push hands alternately to popliteal fossa, 3 times. 5. Twist both hands upwards 	<p>Technique1,2,3,4,6,7,8</p>  <hr/> <p>Technique5</p>  <hr/> <p>Technique10,12</p>
--	--	--	---

		<p>alternately as if twisting a fried dough, 3 times.</p> <p>6. Caress the treatment area, 3 times.</p> <p>7. Push four main collateral channels upwards with kneeling fingers of both hands, 3 times.</p> <p>8. Caress the treatment area, 3 times.</p> <p>9. Do the same on the other side.</p> <p>10. Vacuum&RF operation: move the device from popliteal fossa to Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB), 3 times.</p> <p>11. Move the probe in small circular motion on calf, 3 times.</p> <p>12. Move the device upwards to popliteal fossa and caress the treatment area with hands, 3 times.</p> <p>13. Push from popliteal fossa tracing 4 main collateral channels to the end of thigh, 3 times.</p> <p>14. Move the device in circular motion from popliteal fossa to the end of thigh, 3 times.</p> <p>15. Push from popliteal fossa tracing main collateral channels to the end of thigh, 3 times.</p> <p>16. Clean with towel and apply gel evenly on the treatment area.</p> <p>17. Repeat the same techniques with 40K operation.</p>	 <hr/> <p>Technique11</p>  <hr/> <p>Technique13,15</p>  <hr/> <p>Technique14</p>  <hr/> <p>Technique20,21,22,23, 24</p>  <hr/> <p>Technique25,29</p> 
--	--	---	--





	<p>18. Do the same on the other side.</p> <p>19. Clean with hot towel.</p> <p>20. Foreleg: rub oil with both hands from heel to the thigh root (namely to caress the treatment area), 3 times.</p> <p>21. Alternately push the entire leg from bottom to top and wrap back into the heel, 3 times.</p> <p>22. Move hands upwards to four main collateral channels: Spleen Meridian(SP) - Stomach Meridian(ST) - Liver Meridian(LV) - Gallbladder Meridian(GB) with fingers spread and hukou pointing upwards, 3 times.</p> <p>23. Push four main collateral channels upwards with kneeling fingers of both hands, 3 times.</p> <p>24. Vacuum&RF operation: Lift the device towards knee while tracing 3 meridians, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>25. Move the device on two sides of calf in annular motion, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>26. Lift the device from knee to thigh root in lines, 3 times.</p> <p>27. Move the device in small circles on thigh, 3 times.</p> <p>28. Coupled with hand, lift from 2 sides of leg to the middle, 3 times.</p> <p>29. Lift the device from knee to thigh root in annular</p>	<p>Technique26,30</p>  <hr/> <p>Technique27</p>  <hr/> <p>Technique28</p> 
--	---	--







		<p>motion, 3 times.</p> <p>30. Lift the device from knee to thigh root in lines, 3 times.</p> <p>31. Clean with towel and apply gel evenly on the treatment area.</p> <p>32. Repeat the same techniques with 40K operation.</p> <p>33. Do the same on the other side.</p> <p>34. Clean with hot towel.</p>	
--	--	--	--







Suggested Treatment:

10 treatments account for a full course. After one treatment, leg is relaxed, and circulation is accelerated. After a course of treatment, leg becomes thin, skin collagen tissue experience contraction due to heat effect, and tightening effect can obviously be felt on the skin. After 2 courses, leg is tight and slim, and obvious results are seen. 3 courses help consolidate the effect. Super RF + super ultrasound is to stimulate the dermis of the skin, which can sustain collagen regeneration, so the treatment effect will be more and more obvious.

Body-Shaping&Slimming: 180 min, once/week

<p>RF Energy: 30%-80%</p> <p>Vacuum RF Energy: 30%-80%</p> <p>Suction time: 0.3-1.5</p> <p>Release time: 0-1.0</p> <p>Mode: M1</p> <p>Ultrasonic Energy: 30%-80%</p> <p>Mode: Intermittent</p>	<p>Essential oil (massage cream) + Gel + Device</p>	<p>Back</p> <ol style="list-style-type: none">1. Apply oil on back and press Gb-20 and Du-16 points.2. (Starting from hairline) stroke Da Ban Jin with thumb, 3-5 times.3. Move thumb outwards from Bladder Meridian(BL) to the 8 Crevice area and then to Gb-20 and Du-16 point.4. Starting from neck to caudal vertebra, move hands in S-shape, 3 times.5. Push Bladder Meridian(BL) with thumb in left-right order, 3 times.6. Push Bladder Meridian(BL) to the 8 Crevice area with both thumbs, at the same time, 3 times.7. Push Bladder Meridian(BL) in three kneeling fingers, 3 times.8. Push scapula slot with both hands alternately in left-right order, 3-6 times.9. Caress the back with both hands and push Gb-20 and Du-16 points, 3 times.10. Overlap both thumbs and push SI-11point and slide to the arm, 3 times.	<p>Technique1,9</p>  <hr/> <p>Technique2,16</p>  <hr/> <p>Technique3</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5,6,7,12,14</p>
--	---	--	--

		<p>11. Repeat the previous action on another arm, 3 times.</p> <p>12. Rub Meridian GV and Bladder Meridian(BL) with hands until they turned hot.</p> <p>13. Apply gel evenly on back.</p> <p>14. RF Operation: Starting from Meridian GV to Bladder Meridian(BL), move the device from neck to the 8 Crevice area, 3-5 times.</p> <p>15. Move the probe in circles on Du-14 point, and the 8 Crevice area, 3-5 times respectively.</p> <p>16. From Da Ban Jin to neck, move the device back and forth, 3-5 times.</p> <p>17. Lift the device on scapula slot in left-right order, 3-5 times.</p> <p>18. Move the device horizontally as if writing an "8" number to the 8 Crevice area, 3 times.</p> <p>19. Move the device upwards in annular motion to armpit in left-right order on both sides, 3 times.</p> <p>20. Lift from sides of the waist upwards till armpit, 3-5 times.</p> <p>21. Clean with hot towel.</p> <p>22. Apply essential oil evenly on the back.</p> <p>23. Repeat the same technique with Vacuum&RF Operation.</p> <p>24. Clean with fresh towel.</p> <p>25. The end of back treatment.</p> <p>Legs</p>	 <hr/> <p>Technique8,17</p>  <hr/> <p>Technique10,11</p>  <hr/> <p>Technique15</p>  <hr/> <p>Technique18</p>  <hr/> <p>Technique19</p>  <hr/> <p>Technique20</p>
--	--	--	---

	<p>26. Massage leg in left-right order: apply and rub oil in from calf to thigh to heel, 3 times.</p> <p>27. Alternately push the entire leg from bottom to top and wrap back into the heel, 3 times.</p> <p>28. Move hands upwards to four main collateral channels: Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) with fingers spread and hukou pointing upwards.</p> <p>29. Push hands alternately to popliteal fossa, 3 times.</p> <p>30. Twist both hands upwards alternately as if twisting a fried dough, 3 times.</p> <p>31. Push four main collateral channels upwards with kneeling fingers of both hands, 3 times.</p> <p>32. Caress the treatment area, 3 times.</p> <p>33. Do the same on the other side.</p> <p>34. Vacuum&RF operation: move the device from popliteal fossa to Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB), 3 times.</p> <p>35. Move the device in small circular motion on calf, 3 times.</p> <p>36. Move the device upwards to popliteal fossa and caress the treatment area with hands, 3 times.</p>	 <hr/> <p>Technique26,27,28,29,31,32</p>  <hr/> <p>Technique30</p>  <hr/> <p>Technique34,36</p>  <hr/> <p>Technique35</p>  <hr/> <p>Technique37,39</p>  <hr/> <p>Technique38</p>
--	---	---

37. Push from popliteal fossa tracing 4 main collateral channels to the end of thigh, 3 times.

38. Move the device in circular motion from popliteal fossa to the end of thigh, 3 times.

39. Push from popliteal fossa tracing main collateral channels to the end of thigh, 3 times.

40. Clean with towel and apply gel evenly on the treatment area.

41. Repeat the same techniques with 40K operation.

42. Do the same on the other side.

43. Clean with hot towel.

44. The end of leg treatment.

Arms

45. Massage arms in left-right order: lay the arm flatwise, apply and rub oil in from lower arm to the entire arm and slide to the fingers, 3 times.

46. Push the entire arm with both palms, 3 times.

47. Push hands upwards tracing three main collateral channels: Large Intestine Meridian(LI) - Triple Energizer Meridian(TE) - Small Intestine Meridian(SI) to armpits with fingers spread and hukou pointing upwards, 3 times.


48. Trace and rub three main collateral channels on arms with kneeling finger back and forth till them turn hot, 3



Technique45



Technique46,47,48,49, 55,58



Technique50



Technique51,59,61



Technique52



Technique56,57

times.

49. Caress the treatment area, 3 times.

50. Lay inner arm upwards, and push hands tracing three yin channel on inner arm: Lung Meridian(LU) - Pericardium Meridian(PC) - Heart Meridian(HT) to armpit, with hukou pointing upwards, 3 times respectively.

51. Rub three channels back and forth with kneeling fingers, 3 times.

52. Caress the treatment area and slide to fingers.

53. Repeat the previous operation on another side

54. Apply gel evenly on the arm.

55. RF device operation: lay arms flatwise and start from fat part of lower arm, tracing three channels to armpit, 3 times.

56. Move the device in annular motion, tracing three channels to armpit, 3 times.

57. You may double the operation on flabby arms, 3 times.


58. Trace three channels to armpit, 3 times.

59. Lay inner arm flatwise and push three channels on upper arm to armpit, 3 times.


60. Trace three channels on upper arm and move the device in circles to armpit, 3 times.

61. Trace three channels from upper arm to armpit.


62. Repeat the previous operation with 40K operation.




Technique65,71,74




Technique66



Technique67



Technique68,82



Technique69



Technique70

63. Do the same on the other side.

64. The end of arm treatment.

Abdomen

65. Apply and rub oil on abdomen with hands, 3 times.

66. Rub stomach back and forth with both hands, 3-5 times.

67. Knead abdomen with both hands using chiropractic techniques, 3 times.

68. Lift Meridian BV on both sides of waist with both hands alternately, 16 times.

69. Move hands as if writing an "8" number on waist, 3 times.

70. Overlap hands and message the intestinal canal clockwise, 3 times.

71. Move hands in circular motion and caress the treatment area, 3 times.


72. Push on the following acupoints: Ren-13, Ren-12, Ren-10, Ren-8, Ren-6, Ren-4, Ren-3, St-25, SP-15, Ren-2, 2 times.

73. Starting from Ren-3 point, push both thumbs to belly navel, slide them to waist and the lift upwards to groin, 3 times.

74. Caress the treatment area with hands and slide to groin.

75. Apply gel evenly on abdomen, 1 min.

76. 40K operation: starting from one side of the waist, lift the device to abdomen and




Technique72




Technique73



Technique76





Technique78



Technique79



Technique80

		<p>groin, 3 times.</p> <p>77. Repeat the previous operation on the other side, 3 times.</p> <p>78. Move the device on abdomen in small circular motions, 3 times.</p> <p>79. Move the device on abdomen in big circular motions, 3 times.</p> <p>80. Move the probe to shape the waist as if writing “8”, 3-5 times.</p> <p>81. Move the probe up and down to shape the waist, 3-5 times.</p> <p>82. Combined with hand, lift Meridian BV on both sides of the waist, 16 times.</p> <p>83. Wipe clean abdomen with hot towel, 2 min.</p> <p>84. Apply essential oil evenly on the abdomen.</p> <p>85. Repeat the same technique with vacuum&RF operation</p> <p>86. Clean with fresh towel.</p> <p>87. The end of abdomen treatment.</p>	 <hr/> <p>Technique81</p>  <hr/>
--	--	--	--







Suggested Treatment:







10 times account for a full course of treatment. After one treatment, your body will feel relaxed, blood circulation boosted, and skin tightened. After a course you will see less fat on different parts of the body, collagen regenerated, and firm skin. After 2 courses fat reduces, and curve lines begin to show. 3 courses will help with enhance the effects and shape up, and S curve will appear.







Overall facial Care: 100 minutes, once/week

<p>RF Face Energy: 30%-80%</p> <p>Mode: PRO</p>	<p>Makeup remover+ Cleanser+ Hydrating facial mask+</p>	<ol style="list-style-type: none"> 1. Remove makeup and cleanse face and neck, 5 min. 2. Apply toner, 2 min. 3. Apply cleaning base 	<p>Technique6</p>
---	---	--	-------------------

<p>RF Eye Energy: 30%-80%</p> <p>Mode: PRO BIO Mode: Face, Neck: Consecutive Eye: Intermittent</p> <p>Cold&Hot hammer energy: 30%-80%</p> <p>Mode: Cold Hammer</p>	<p>Hot&Cold Steam+ Massage cream+ Device+ Essence+ Facial mask</p>	<p>mask+hot steam, 10 min.</p> <p>4. Wash face clean, 2 min.</p> <p>5. Apply massage cream evenly on face and caress face, 3 times.</p> <p>Push the following acupoints: Ren-24, Du-26, St-4, St-6, LI-20, BI-1, BI-2, EM3, TE-23, Gb-1, the Temple, St-1, St-2, 3 times.</p> <p>6. Caress the whole face, 2 times.</p> <p>7. Alternate both hands to lift in 3 lines from chin to ear lobe, corner of mouth to ear gate, nose wing to temple, lower eyelid to temple, 3 times.</p> <p>8. Do the same on the other side.</p> <p>9. Lift towards hairlines on forehead, 3 times.</p> <p>10. Move hands in circular motion on eyes and push temples, 3 times.</p> <p>11. Move middle and ring fingers horizontally below lower eyelids and slide to back of the ear, 3 times.</p> <p>12. Lift middle and ring fingers from inner eye corner to upper eyelids and slide to back of the ear, 3 times.</p> <p>13. Lift eye corners with "scissor" hand gesture, 3-5 times.</p> <p>14. Repeat the previous operation on the other side, 3 times.</p> <p>15. Move hands in circular motion on eyes and push temples, 3-5 times.</p> <p>16. Use both hands to lift</p>	<div data-bbox="1171 226 1310 450" data-label="Image"> </div> <p>Technique8</p> <div data-bbox="1165 551 1289 748" data-label="Image"> </div> <p>Technique10</p> <div data-bbox="1163 840 1289 1039" data-label="Image"> </div> <p>Technique11,16</p> <div data-bbox="1163 1128 1289 1328" data-label="Image"> </div> <p>Technique12</p> <div data-bbox="1163 1420 1289 1619" data-label="Image"> </div> <p>Technique13</p> <div data-bbox="1163 1711 1289 1910" data-label="Image"> </div> <p>Technique14</p>
--	--	--	--

		<p>face as if plucking the string, 1 min.</p> <p>17. Move in letter Z on forehead with ring and middle fingers and push on temples, 3 times.</p> <p>18. Neck massage: apply and rub oil on the chest with both hands to the back of the neck and point to Gb-20 and Du-16 acupoints, 3 times.</p> <p>19. Alternately lift the lower jaw and double chin to pass behind the ear and under the armpit, 3 times.</p> <p>20. Push downwards with hand spread and hukou pointing down, tracing 3 meridian channels from the neck side to the armpit, 3 times.</p> <p>21. Push downwards with 4 fingers, tracing 3 meridian channels from neck side to the armpit, 3 times.</p> <p>22. Push downwards with kneeling fingers till the neck turns hot and slide to the armpit, 3 times.</p> <p>23. Push downwards with hand spread and hukou pointing down from external collarbone to armpit, 3-5 times.</p> <p>24. Alternate hands in lifting the neck, 3 times.</p> <p>25. Do the same on the other side.</p> <p>26. Face RF Operation: Move the probe from jaw to earlobe, mouth corner to ear gate, nose wing to temple in circles, 3 times.</p>	 <p>Technique18</p>  <p>Technique19</p>  <p>Technique20</p>  <p>Technique21,22,23</p>  <p>Technique24</p>  <p>Technique25</p>
--	--	---	--

	<p>27. Brace the skin with the probe and another free hand to keep the treatment surface tight and lifted; move from jaw to earlobe, mouth corner to ear gate, nose wing to temple, lower eyelid to temple to lift the treatment area, 3 times.</p> <p>28. Combined with hands to lift towards hairline on forehead, 3 times.</p> <p>29. Do the same on the other side.</p> <p>30. Eye RF operation: move the probe beneath lower eyelid in small circles and move all the way till the temple, 3-5 times.</p> <p>31. Combined with hand, move the probe from lower eyelid to temple to lift the treatment area, 3-5 times.</p> <p>32. Lift from the beginning of the eyebrows to the end, 3-5 times.</p> <p>33. Combined with hand, move the probe to lift the corner of the eyes towards the hairline, 3-5 times.</p> <p>34. Do the same on the other side.</p> <p>35. Neck RF operation: coupled with hand massage, lift double chin to pass behind the ear and under the armpit, 3 times.</p> <p>36. Move the device in circles on one side of neck to pass under armpit, 3 times.</p> <p>37. Slide the device from internal and external clavicle to the armpit, 3-5 times.</p>	 <p>Technique27</p>  <p>Technique28,44,45</p>  <p>Technique29,46</p>  <p>Technique31,49</p>  <p>Technique32,34,50</p>  <p>Technique33,51</p>
--	--	---

	<p>38. Move the device in circles on the whole neck area, 3 times.</p> <p>39. Coupled with hand, lift neck in lines, 3-5 times.</p> <p>40. Do the same on the other side.</p> <p>41. Clean face and neck, 3 min.</p> <p>42. Apply (anti-aging)essence evenly on the face, 1 min.</p> <p>43. Face BIO Operation: Move the probe from jaw to earlobe, mouth corner to ear gate, nose wing to temple in circles, 3-5 times.</p> <p>44. Combined with hand and repeat the above operation, 3-5 times.</p> <p>45. Move the probe horizontally in circles on forehead, 3 times.</p> <p>46. Do the same on the other side.</p> <p>47. Apply (anti-aging)essence evenly around the eyes, 1 min.</p> <p>48. Eye BIO operation: move the probe beneath lower eyelid in small circles and move all the way till the temple, 3-5 times.</p> <p>49. Combined with hand, move the probe from lower eyelid to temple to lift the treatment area, 3-5 times.</p> <p>50. Lift from the beginning of the eyebrows to the end, 3-5 times.</p> <p>51. Do the same on the other side.</p> <p>52. Apply (anti-aging)essence evenly</p>	 <p>Technique36,54</p>  <p>Technique37,55</p>  <p>Technique39,57</p>  <p>Technique40,58</p>  <p>Technique61</p>  <p>Technique63</p>
--	--	---

	<p>on the neck, 1 min.</p> <p>53. Neck BIO operation: coupled with hand massage, lift double chin to pass behind the ear and under the armpit, 3 times.</p> <p>54. Move the device in circles on one side of neck to pass under armpit, 3 times.</p> <p>55. Slide the device from internal and external clavicle to the armpit, 3-5 times.</p> <p>56. Move the device in circles on the whole neck area, 3 times.</p> <p>57. Coupled with hand, lift neck in lines, 3-5 times.</p> <p>58. Do the same on the other side.</p> <p>59. Apply (hydrating) essence on face, 1 min.</p> <p>60. Ultrasonic Operation: Move the probe from jaw to earlobe, mouth corner to ear gate, nose wing to temple in circles, 3 times.</p> <p>61. Do the same on the other side.</p> <p>62. Move the probe horizontally in circles on forehead, 3 times.</p> <p>63. Apply hydrating eye mask, face mask, neck mask (at the same time)+cold hammer to cover the whole face in circles, 3 times.</p> <p>64. Remove the mask and wash face and neck clean, 5 min.</p> <p>65. Apply toner, essence, cream, eye cream, neck cream and sunscreen.</p>	<div data-bbox="1161 219 1289 421" data-label="Image"> </div> <hr/> <p data-bbox="1082 448 1248 481">Technique64</p> <div data-bbox="1161 510 1289 712" data-label="Image"> </div>
--	--	--

Suggested Treatment:

Ten treatments account for one full course. After one treatment, the skin will be firm, tender, delicate and smooth, and the eyes will be elevated to accelerate the blood circulation and lymph circulation of the eyes. After one course, it lightens fine lines, dark circles under the eyes and brightens the skin. After two courses, the overall skin is smooth, firm and elastic, and the double chin is improved. Three courses help prevent and delay skin aging, relaxation and sagging, and restore skin to its youthful state.