

User Manual

NM-SN6SP

Preface

Dear Users:

We're pleased to present to you N6SP beauty machine. It's main focus is tightening and lifting of skin through massage. With the help of RF, DMA and LED technologies, it targets on facial fine lines, wrinkles around eyes and facial dehydration and brings great results of anti-aging, tightening and anti-puffiness of eyes. At the same time, it delivers a safe, comfortable and and effective skincare experience. RF machines are anti-aging machines for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore we advise any personnel read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

Table of Contents

Preface
Part I
Brief Introduction
Advantages
Usage and Efficacy
Facial Absorption
Skin Rejuvenation
Facial Soothing
Anti-puffiness around Eyes
Part II
1. Detailed Operations
2. Technical Parameters
3. Safety Precautions
4. Dos and Don'ts
5. Troubleshooting & Solutions
6.FAQs
7.Packing List
8.Operational Diagrams

Part I

Brief Introduction

N6SP can be used to massage and lift skin around eyes as well as better absorb skincare products. It complies with the theory of human engineering in design so that it can be as close to the skin and saves effort to use. When it comes to the effects, it does not only improve issues like eye bags, fine lines and dark circles, but also deals with aging, sagging and loose skin. Radio frequency and ultrasonic beauty machines are currently the most popular beauty machines that work on tightening, body-shaping and anti-aging. It's an effective replacement of all the facial and body care programs in beauty salons. Moreover, they are convenient and easy to operate. It solves both skin and body management issues for beauty-lovers. This equipment requires no injection, no medication, no operation, thus having no side-effects. It is used externally during the whole process and has an immediate effect. RF works by heating up collagen fibers with high frequency waves and DMA technology to activates muscles with high frequency massage so that they can reach deep into the hypodermis and cause muscular movements around eyes and of the whole face. Deeply heated eye masks can also help with blood circulation and metabolism so that puffiness and edema of eyes can be treated as well.

Advantages

1. RF beauty treatment has a wide range, effective results, and has been given a lot of attention by the beauty industry.
2. It's portable and easy to operate.
3. FR technology targets at issues like sagging and loose facial skin, indistinct jawline, tear trough and fine lines around eyes.
4. It will not result in unevenness, bleeding, swelling or blood stasis.
5. It is multifunctional and has numerous approaches to use.
6. The red and blue rays emitted from the device can be used together with other functions. They can sterilize and reduce inflammation, and help skin to better absorb nutrients.
7. Deeply heated eye masks can effectively improve blood circulation and puffiness of eyes.

Usage and Efficacy

1. Facial Absorption

RF

Radio frequency can alternate electrical polarities of biological tissues on the treated area millions of times per second and change directions in the same frequency concurrently. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy, which stimulates the dermis to produce more collagen and fill up the gap of the lost collagens, therefore plumping up the skin and revive elasticity.

Vibration

Through high frequency vibration, it massages the skin, thus achieving effects like relaxing muscles, speeding up metabolism, and introducing nutrients like essence deep into hypodermis. Other results are achieved as well, such as increased density, tightness, and elasticity of the skin. It also prevents sagging and loose skin conditions.

Effects

1. Deeply induce nutrients.
2. Accelerate nutrients absorption, and deeply nourish the skin.
3. Open up pores and accelerate skin metabolism.
4. Help with skin dehydration and fine lines.
5. Tighten the skin, and increase skin elasticity.

Applicable Range

1. Those with facial dullness.
2. Those with rough skin and large pores.
3. Those with dryness, dehydration and peeling of skin.
4. Those with dry and fine lines; those with a slow metabolism.

Inapplicable Range

1. Those in menstruation, pregnancy, lactation, surgical recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors; those who have just implanted metal materials.
3. Those with skin trauma, wound and inflammation; those who have just undergone

plastic surgery.

4. Those who are having an allergic reaction, or people with severe sensitive skin.
5. Those with skin diseases and infectious diseases.

Dos and Don'ts

1. Use sunscreen and facial masks to maintain hydrated skin.
2. Avoid smoking, drinking and staying up late to prevent dryness.
3. Avoid washing face with overheated water. Warm water is advised.
4. Avoid irritant skin products, such as alcohol or AHAs.
5. Remove make-up and clean face thoroughly before bed.
6. Avoid changing skin products too often.

2. Skin Rejuvenation

RF

Radio frequency is a high frequency wave, whose highest frequency can reach millions of times per second. It directly goes into the bottom layers of skin and continuously produce new collagen, and revive lost collagen so as to slow down the process of collagen loss. The effect it has achieved is more overt than applying and consuming collagen products.

It can alternates electrical polarities of biological tissues on the treated area millions of times per second and change directions in the same frequency concurrently. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy, which stimulates the dermis to produce more collagen and fill up the gap of the lost collagens, therefore plumping up the skin and revive elasticity.

Vibration

Through high frequency vibration, it massages the skin,thus achieving effects like relaxing muscles, speeding up metabolism, and introducing nutrients like essence deep into hypodermis. Other results are achieved as well, such as increased density, tightness, and elasticity of the skin. It also prevents sagging and loose skin conditions.

Red LED

Red Ray: With the wavelength of 635nm, it heats up the bottom layers of skin, therefore reactivating collagens and speeding up the micro-circulation of skin so that it appears to be more elastic, smooth, glossy and has less wrinkles. The theory behind is to heat collagen fibers in dermis to induce collagen through RF technology and improve cell bio-activity through red ray, thus achieving effects like smooth skin with less wrinkles.

Effects

1. Accelerate metabolism and skin circulation.
2. Better absorption of skincare products; faster rate of absorption.
3. Stimulate collagen in bottom layers of skin; improve skin elasticity and tightness.
4. Prevent hyperpigmentation, fade melanin and whiten skin.
5. Prevent sagging, rough skin with wrinkles; anti-aging.

Applicable Range

1. Those with dry, rough and dull skin.
2. Those with poor absorption and dehydration of skin.
3. Those with dullness and unevenness.
4. Those with sagging skin.
5. Those with aging and wrinkled skin.

Inapplicable Range

1. Those who have just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.
4. Those with skin trauma or wound.
5. Those who are overaging.
6. Those in menstruation and surgical recovery.
7. Those with skin diseases and infectious diseases.

Dos and Don'ts

1. Use sun protection products; avoid exposing under the sun after operation.
2. Avoid scratching and toughing with hands.
3. Apply more facial masks to hydrate the skin; prevent dryness; accelerate metabolism.
4. Avoid changing skin products too often.
5. Remove make-up and clean face thoroughly before bed.
6. Avoid over-heated environment.

3. Facial Soothing

RF

Radio frequency can alternate electrical polarities of biological tissues on the treated area millions of times per second and change directions in the same frequency concurrently. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy, which stimulates the dermis to produce more collagen and fill up the gap of the lost collagens, therefore plumping up the skin and revive elasticity.

Vibration

Through high frequency vibration, it massages the skin, thus achieving effects like relaxing muscles, speeding up metabolism, and introducing nutrients like essence deep into hypodermis. Other results are achieved as well, such as increased density, tightness, and elasticity of the skin. It also prevents sagging and loose skin conditions.

Blue LED

Blue LED: With a certain wavelength, blue ray helps to shrink pores, reduce oil, calm skin and treat acne. It also helps to lock moisture and essence in skin, making it more elastic. Blue ray helps with inflammation of skin and treatment of acne by fighting bacterial infection.

Effects

1. Boost skin care product absorption and permeability.
2. Calm skin.
3. Help with acne and sensitive skin.
4. Shrink pores and smooth fine lines.

5. Repair and strengthen skin barrier; after sun repair and sunburn relief.
6. Reduce oil production and regulate oil gland secretion.
7. Brighten skin; fight bacterial infection and treat acnes.

Applicable Range

1. Those with dry skin.
2. Those with large and clogged pores; those with dull, rough and dark skin.
3. Those with sensitive skin: weak skin barrier, thin stratum corneum, and couperose.
4. Those with acne and sensitive skin.
5. Those with oily, dry and combination skin.

Inapplicable Range

1. Those in menstruation, pregnancy, lactation, surgical recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors; those who have just implanted metal materials.
3. Those with skin trauma, wound and inflammation; those who have just undergone plastic surgery.
4. Those who are having an allergic reaction, or people with severe sensitive skin.
5. Those with skin diseases and infectious diseases.

Dos and Don'ts

1. Use sunscreen and facial masks to maintain hydrated skin.
2. Avoid smoking, drinking and staying up late to prevent dryness.
3. Avoid washing face with overheated water. Warm water is advised.
4. Avoid irritant skin products, such as alcohol or AHAs.
5. Remove make-up and clean face thoroughly before bed.
6. Avoid changing skin products too often.

4. Anti-puffiness around Eyes

Eye Patch

Radio frequency works through the circuits formed by electrodes in pair. Radio frequency waves penetrate skin directly and produce thermal energy created by strong vibration from tissue resistance (1 MHz). During the treatment, RF device stimulates new collagen production by heating up hypodermis around eyes, which

helps to tighten the skin, remove dark circles and bags, treat couperose and fine lines.

Effects

1. Accelerate blood circulation and metabolism
2. Stimulate collagens and slow down the process of skin aging.
3. Relieve eye fatigue, and swollen eyelids; Soothe sore eyes.

Applicable Range

1. Those with wrinkles, fine lines, eye bags, dark circles.
2. Those with eye fatigue with symptoms like dryness.
3. Those who face computer and phone for a long time.
4. Those who stay up late often.
5. Those who just ended a relationship and cries a lot.

Inapplicable Range

1. Those who have just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.
4. Those with skin trauma or wound.
5. Those who are overaging.
6. Those in pregnancy and surgical recovery.
7. Those with skin diseases and infectious diseases.

Dos and Don'ts

1. Use sun protection products and maintain dehydrated skin.
2. Apply eye masks, preferably 3 times a week.
3. Avoid spicy and greasy food.
4. Avoid smoking, drinking and staying up late; eat more vegetables, fruits, and less greasy food.
5. Avoid drinking water excessively half an hour before bed.

Part II



1. Detailed Operations




Interface Display



Energy Level Adjustment: press once to add one level, press twice to add two levels; do it at one sitting.



Start/Off/Mode (Press for 5 sec to turn it on; press once to select mode.)





Working Time Display



RF Energy Display (The maximum energy setting is 5.)




Heat Mode (Press  once to turn on the heat mode.)



Vibration Mode (Press  twice to turn on the vibration mode.)



Photon Mode (Press  once to select red LED and press twice to select blue LED.)

2. Technical Parameters

Model: NM-SN6SP

Weight: 0.18kg

Size: 17cm*5cm*7cm

Working Environment: 5-40 °C

Power supply

Input: AC100-240V 50/60HZ 0.6A

Output: DC9V 2A

Warranty: 1 Year

Radio Frequency

Frequency: 1MHz

Power Consumption: 25W

3. Safety Precautions

Please consult your doctor or professionals before using the device if you have the following symptoms:

1. Those who are pregnant or during lactating.
2. Those with heart disease or equipped with cardiac pacemaker.
3. Those with unhealed wounds or who are during convalescent period.
4. Those with epilepsy, severe diabetes and hyperthyroidism.
5. Those with malignant tumors, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.

4. Dos and Don'ts

1. Please remove the cotton pad and clean the device with water and keep it properly stored after each operation.
2. Use a grounding pin and make sure the power socket of the device is properly connected before each operation.
3. Ensure that the voltage of the device is appropriate. If the voltage of the local power supply is unstable, we recommend that users increase to the matching power voltage.
4. In order to ensure the therapeutic effect and normal service life expectancy of the device, please use the specified accessories provided or recommended by the original manufacturer.
5. The device should not be placed in wet environment or near water, nor should they be directly exposed to sunlight.
6. Do not place the device near a strong heat source, because this may affect the life and normal use of the device.
7. Before treatment, please remove all metal objects from the body first, so as to avoid unexpected conditions and influence the curative effect.
8. Do not use in eyes, thyroid gland, parathyroid gland, testis, abdomen of pregnant women and atrial pulse generator.
9. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
10. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use, so as to ensure the safety of the electrical products.
11. Use device or train device operators in strict accordance with instructions in the user manual.
12. When operating, the device should fully touch the skin to avoid uneven heat.
13. Start at the lowest level of energy at the beginning of the operation, then slowly add up.
14. When using this device, the operating parts must be kept moist and dry skin

should be avoided.

15. Clean the device with normal saline after operation to ensure its cleanliness and hygiene and prolong its service life.

16. Each pair of eye films is advised to use 3-5 times; avoid long term use; avoid sharing with other users so as to prevent bacterial infection.

5. Troubleshooting & Solutions

1. Power does not turn on and the key lamp on the back of the device doesn't work?

- A. Ensure that the power cord is connected to a working power outlet.
- B. Ensure the fuse tube on the back of the instrument is not loose or being burnt.

2. No RF output?

- A. Please check whether the plug connecting the device handle and the fuselage is tightly connected.
- B. Please make sure that the treated area is clean. Oil or oily essential oil products may cause poor contact between the handle of the device and the human body, resulting in no RF output.

3. Reduced RF output?

- A. Please check whether the plug connecting the device handle and the fuselage is tightly connected.
- B. Please check whether the handle is contaminated with non-conductive grease and other substances, which will cause poor contact and reduce the output.
- C. Please check whether the product used is an device-specified adapter.

4. The device can be activated, but the wrong information appears on the screen?

- A. Please unplug the power plug behind the device and wait about 1 minute before restarting the device.
- B. If the above method did not work, please contact the device distributor for assistance.

6. FAQs

1. Q: What is collagen?

A: Collagen is a kind of biological macromolecule substance, a kind of white, opaque and non-branched fibrous protein. It can supplement the nutrition needed by all layers of the skin, enhance the collagen activity in the skin, lock moisture, nourish the skin, delay aging, beautify, relieve facial relaxation, provide nutrients for the hair and so on. Collagen is a nutrient that human body must replenish to delay aging. With the growth of age, collagen will gradually be lost. At the age of 20, women have begun to

age, collagen gradually decrease. At the age of 25, collagen loss reaches its peak. At the age of 40, the content is less than half of which one has at age 18. It is the loss of collagen and water that causes the wrinkles on the face of the elderly. The breakage of collagen fibers and elastic mesh will cause skin tissue to be oxidized, atrophied and collapsed. The skin will be dry, wrinkled, flabby and inelastic etc,. Therefore, in order to delay aging, collagen must be supplemented.

2.Q: When can I see the effects of RF?

A: Usually on the same day and within a week of treatment, the collagen tissue of the skin is heated to produce a tightening effect. The tightening sensation of the skin can be clearly felt. As radio frequency stimulates the subcutaneous tissue and sustains the regeneration of collagen so more often the treatment is, the more obvious the effects are.

3. Q: Is RF harmful to skin?

A:RF treatment is a non-surgical project. It stimulates the regeneration and metabolism of collagen at the bottom layer of the skin, doing no damage to the skin. It is only a normal phenomenon of local skin fever and redness after operation, which belongs to the acceleration of blood circulation. It will disappear after a while so there's no need to worry.


4. Q: How long is the whole treatment?


A: Including facial mask, each treatment takes about half an hour. At the beginning of the operation, do it once a day for 3 consecutive days. After a month, do it 3 times a week. After 2 months, do it once a week to consolidate the effects.

7.Packing Lists

- 1 . Main Body × 1
- 2 . Cotton Pad × 1
- 3 . Power Cord × 1
- 4 . Eye Patch × 1

8.Operational Diagrams


Parameter Adjustment	Product	Techniques	Diagram
Facial Absorption:40 minutes, 3-4 times/week			
RF Energy: 5 Vibration Mode: On	Make-up removing oil+Facial Wash+Toner+Hydrating essence	<ol style="list-style-type: none"> 1. Remove make-up and cleanse face. 5 minutes. 2. Apply toner. 1 minute. 3. Charge the device, and fixate the cotton pad on the probe with the help of the magnetic cotton holder; apply essence on the cotton and choose RF and vibration mode after starting the machine. 4. Move the device in circles. (from chin to earlobe; from mouth corners to ear gate; from the wing of the nose to the temple; from the risorius to the temple; from lower eyelid to the temple.) Move the device horizontally in circles on forehead. Do all the above operation for 3 times. 5. Repeat step 4 on other side of the face. 6. Move the device upwards. (from chin to earlobe; from mouth corners to ear gate; from 	



		<p>the wing of the nose to the temple; from the risorius to the temple; from lower eyelid to the temple.) Move the device upwards to the hairline with hand massages and the device.</p> <p>Do all the above operation for 3 times.</p> <p>7. Repeat the previous operation on other side of the face.</p> <p>8. Apply moisturizer and make-up base.</p>	
<p>Recommended treatment: Once a day the first 3 consecutive days; 3-4 times a week in the first month; once a week after that.</p>			
<p>Skin Rejuvenation: 40 minutes, 3-4 times/week</p>			
<p>RF Energy: 5</p> <p>Vibration Mode: On</p> <p>Red LED: On</p>	<p>Make-up removing oil+Facial Wash+Toner+Firming essence</p>	<ol style="list-style-type: none"> 1. Remove make-up and cleanse face. 5 minutes. 2. Apply toner. 1 minute. 3. Charge the device, and fixate the cotton pad on the probe with the help of the magnetic cotton holder; apply essence on the cotton and choose RF, LED red mode , vibration mode after starting the machine. 4. Move the device in circles. (from chin to earlobe; from mouth corners to ear gate; from the wing of the nose to the temple; from the risorius to the temple; 	

		<p>from lower eyelid to the temple.) Move the device horizontally in circles on forehead.</p> <p>5. Repeat the previous operation on other side of the face.</p> <p>6. Move the device upwards. (from chin to earlobe; from mouth corners to ear gate; from the wing of the nose to the temple; from the risorius to the temple; from lower eyelid to the temple.) Move the device upwards to the hairline with hand massages and the device.</p> <p>Do all the above operation for 3 times.</p> <p>7. Repeat step 6 on other side of the face. (Facial mask can be applied for better effect after the device operation.)</p> <p>8. Apply moisturizer and make-up base.</p>	
--	--	--	--

Recommended treatment:
Once a day the first 3 consecutive days; 3-4 times a week in the first month; once a week after that.

Facial Soothing: 40 minutes, 3-4 time/week

<p>RF Energy: 5</p> <p>Vibration Mode: On</p> <p>Blue LED: On</p>	<p>Make-up removing oil+Facial Wash+Toner+Repair essence</p>	<ol style="list-style-type: none"> 1. Remove make-up and cleanse face. 5 minutes. 2. Apply toner. 1 minute. 3. Charge the device, and fixate the cotton pad on the probe with 	
---	--	--	---

		<p>the help of the magnetic cotton holder; apply essence on the cotton and choose RF, LED blue mode , vibration mode after starting the machine.</p> <p>4. Move the device in circles. (from chin to earlobe; from mouth corners to ear gate; from the wing of the nose to the temple; from the risorius to the temple; from lower eyelid to the temple.) Move the device horizontally in circles on forehead.</p> <p>5. Repeat the previous operation on other side of the face.</p> <p>6. Move the device upwards. (from chin to earlobe; from mouth corners to ear gate; from the wing of the nose to the temple; from the risorius to the temple; from lower eyelid to the temple.) Move the device upwards to the hairline with hand massages and the device.</p> <p>Do all the above operation for 3 times.</p> <p>7. Repeat the previous operation on other side of the face. (Facial mask can be applied for better effect after the device operation.)</p> <p>8. Apply moisturizer and make-up base.</p>	 <hr/> 
--	--	---	---

Recommended treatment:
Once a day the first 3 consecutive days; 3-4 times a week in the first month; once a week after that.

Anti-puffiness around Eyes: 20 minutes, 3-5 times/week

Energy level of eye films: 2-5	Make-up removing oil+Facial Wash+Toner+Eye films	<ol style="list-style-type: none"> 1. Remove make-up and cleanse face. 5 minutes. 2. Apply toner. 1 minute. 3. Connect eye patches with the power cord and apply it evenly beneath the lower eyelids. 4. Hang the cord up the ears. 5. Choose the energy level subjecting to your own comfort level. Starting from level 1 and adds up slowly. 6. Wait for 15 minutes. 7. Turn off the device and take off the eye patches. 8. Wipe clean the remnant essence on eyes and face. 	
--------------------------------	--	---	--


Recommended treatment:
Once a day the first 3 consecutive days; 3-4 times a week in the first month; once a week thereafter.

Aided with Facial Mask: 15 minutes, 2-3 times/week

RF Energy:3- 5 Vibration Mode: On LED Mode: Red	Make-up removing oil+Facial Wash+Toner+Facial mask	<ol style="list-style-type: none"> 1. Remove make-up and cleanse face. 5 minutes. 2. Apply toner. 1 minute. 3. Put on a facial mask and turn on the vibration mode, RF and 	No diagram
---	--	---	------------

& Blue		<p>LED mode(red light mode for normal skin type and blue light mode for sensitive skin type.)</p> <p>Wait for 15 minutes.</p> <p>4. Take off the mask and wash your face clean. 2 minutes.</p> <p>5.Apply toner, essence, moisturizer and sunscreen.</p>	
<p>Recommended treatment: 2-3 times a week; no set period of treatment; it can be used each time when applying facial masks.</p>			

Whole Body Treatment: 45 minutes, 2-3 times/week

<p>RF Energy: 4-5</p> <p>Vibration Mode: On</p> <p>LED Mode: Red & Blue</p> <p>Eye Films RF: 4-5</p>	<p>Make-up removing oil+Facial Wash+Toner+Facial Mask</p>	<ol style="list-style-type: none"> 1. Remove make-up and cleanse face. 5 minutes. 2. Apply toner. 1 minute. 3. Charge the device, and fixate the cotton pad on the probe with the help of the magnetic cotton holder; apply essence on the cotton and choose RF, LED red mode , vibration mode after starting the machine. 4. Move the device in circles. (from chin to earlobe; from mouth corners to ear gate; from the wing of the nose to the temple; from the risorius to the temple; from lower eyelid to the temple.) Move the device horizontally in 	
--	---	--	---

		<p>circles on</p> <p>5.Repeat the previous operation on other side of the face.</p> <p>6.Move the device upwards. (from chin to earlobe; from mouth corners to ear gate; from the wing of the nose to the temple; from the risorius to the temple; from lower eyelid to the temple.) Move the device upwards to the hairline with hand massages and the device.</p> <p>Do all the above operation for 3 times.</p> <p>7.Repeat the previous operation on other side of the face.</p> <p>8.Repeat the above steps after choosing RF, vibration and LED blue mode.</p> <p>9.Put on a facial mask while at the same time turn on the vibration, RF, LED blue or red mode to better absorb the essence from the mask. Wait for 15 minutes.</p> <p>10.Turn off the device and take off the mask.</p> <p>11.Wipe face clean.</p> <p>12.Apply toner, essence, moisturizer and sunscreen.</p>	
<p>Recommended treatment: Once a day in the first 3 consecutive days; 3-4 times a week in the first month; once a week after that.</p>			

