

User Manual

NM-S1



Preface

Dear Users,

We're pleased to present to you our Eye RF Massage, which is a multifunctional eye care device to improve eyestrain, dark circles, edema and fine lines. By using advanced technology, we aim at a safe and, at the same time, effective result of protecting eyes, nourishing skin, removing wrinkles and reducing eyestrain. Multifunctional Eye RF Massage is a eye care device for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore we advise any personnel to read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

Table of Contents

Preface
Part I
Introduction
Advantages
Working Principles
Effects
Applicable Range
Inapplicable Range
Notes after Operation
Part II
1. Detailed Operations
2. Technical Parameters
3. Safety Precautions
4. Dos and Don'ts
5. Troubleshooting & Solutions
6. FAQs
7. Packing List
8. Operational Diagrams

Part I

Brief Introduction

Eye RF Massage is currently the most popular eye protection and care device that can effectively solve eye skin issues and eyestrain for beauty-lovers. Moreover, it is convenient and easy to operate, that is comfortable, effort-saving and have clear effects. Eye RF Massage is believed to replace the eye care project and the traditional eye care technique. This equipment is safe and requires no injection, no medication, no operation, thus having no side-effects. The Eye RF Massage combines RF and EMS technologies. EMS can stimulate muscles to exercise and relax muscles. It has bionic eye exercises for relieving eyestrain, dry eyes and so on.

Advantages

1. It is small, easy to carry and operate.
2. Eye RF Massage combines RF and EMS technologies, meeting the needs of caring eyes and eye health.
3. It is easy to learn and operate with immediate effect.
4. Three modes can be switched freely according to your needs or feeling.
5. It is suitable for many occasions: at home, walking, taking a bus or subway or plane.
6. Automatically turning off in 10 min timing, safe and comfortable.
7. It has disposable eye patch, which is clean and sanitary, so you can share it with your friends and family.
8. 42°C RF massage help relieve fatigue and anti-aging.

Principles

EMS

EMS is short for Electronic Muscle Stimulation, referred to as muscle electrical stimulation technology, also known as motor technology. The simple explanation is that the normal muscle after receiving signal sent from the brain start activities actively, and EMS technology through current stimulation, directly convey signals to the muscles, pushing muscles for passive movement. When the movement is at eyes, it can relieve eyestrain, eye puffy, sore, dry and so on. EMS technology through the electric current to stimulate the skin layer and the capillary and nervous inside the epidermal layer cells to produce hydration, fuel the folds of skin oxygen and nutrients, restore elasticity to drooping skin.

RF

RF works through the circuits formed by electrodes in pairs. RF waves directly penetrate the skin and make use of the impedance generated by the skin to generate thermal energy by generating a strong resonant rotation of the cell molecules (1million/s). In the process of eye anti-aging treatment, RF rapidly heats up the bottom layer of the skin under the eye and promotes the growth of collagen of the skin, so as to tighten the skin around eyes, remove under eye bags, edema, improve dark circles, improve couperose skin and fade fine lines.

Main Effects

1. Massage eyes, relieve sore, bloat, ache around eyes.
2. Relax eyes, relieve eyestrain, lacrimation and blurred vision.
3. Relieve dry eyes, eye irritation and give you a healthy look.
4. Improve eye bags and eye edema.
5. Improve eye fine lines, crow's feet, eye grains.
6. Stimulate collagen regeneration of eyes, to prevent flabby skin around eyes
7. Restrain melanin growth, fade dark circles.

Applicable Range

1. Those who have overuse of eyes, eyestrain, blurred vision or lacrimation.
2. Those who often face computer or lower head to check phones.
3. Those with fine lines or crow's feet around eyes.
4. Those with flabby, drooping eye skin.
5. Those with dark circles, eye bags or eye edema.
6. Those with eye sore, dry eye or eye irritation.

Inapplicable Range

1. Those with serious "3 Hs" (hypertension, hyperglycemia, hyperlipidemia) and who have heart disease, malignant tumor, hemophilia or severe disease.
2. Those with epilepsy, serious diabetes or hyperthyroidism.
3. Those with heart disease or cardiac pacemaker.
4. Those in pregnancy, menstruation or lactation.
5. Those with skin disease or infectious disease.
6. Those who are sensitive to the electricity.
7. Those who has metal implants in the body.

Notes after Operation

1. Do not stay up and have a rest.
2. Have a rest with eyes closed for 5 min after operation.
3. Moisturize your eyes area and apply eyes mask and eye cream before sleep every night.

Part II

1. Detailed Operation

The installation is as follows:



1. First, connect the cord and the main device.
2. Then buckle the two pole buckles of the other end of the cord to the eye patches.



1 x Main Device



2 x Connecting Cord



2 x Eye Patch




Power Outlet (when operating, insert two connecting cord separately into these two outlets)



Power Cord Inlet

The detailed operation is as follows:





This concave is the keypad of Power On/Power Off and Level Adjustment Pad. Long press here for 3 seconds is Power On/Power Off, tap here is to switch levels when it is Power On. The level range is 1~3. The higher the level is, the more obvious the feeling is. Level 1 is advised for using for the first time.



The display of the first level (when the power is on, the default level is the first level). The first level is suitable for those who use it for the first time and those who want to relieve eyestrain and sore eyes.



The display of the second level (tap once after the power is on, then the level is changed to the second level). The second level is suitable for those who are adaptive and need to relieve eyestrain, sore caused by over use of eyes, improve dark circles and eye edema.



The display of the third level (tap twice after the power is on, then the level is changed to the third level). The third level is suitable for those who have sore eyes, dry eyes or eyes irritation caused by over use and need eye care, including tightening, anti-aging, eye bags removing and so on.

10 min auto-off | 1-2 times/day



1. Connect the eye patches with the main device.



2. Rest the cables behind the ears. Paste patches beneath lower eyelids. (Store away the plastic films)



3. Touch-on for 3 seconds to turn on the device. Adjust energy level. Enjoy the massage for 10 min.



4. Place the patches back on the plastic films and store away.

2. Technical Parameters

Frequency of Radio frequency: 1MHz

Temperature Range: 0-45 °C

N. W: 60g

G.W: 80g

Dimension: 80 (L)x 40 (W) x 20 (H) mm

Package Size: 150(L) x 110(W) x 30 (H) mm

3. Safety Precautions

Patients with the following symptoms should be cautious about using this instrument. Please consult a doctor or a professional before using the instrument. The details are as follows:

1. Those in pregnancy, menstruation or lactation. The pregnant woman must consult a doctor before using this device.
2. Those who have heart pacemaker, defibrillator or other metal implants in the body.
3. Those with serious “3 Hs” (hypertension, hyperglycemia, hyperlipidemia) and who have heart disease, malignant tumor, hemophilia or severe disease.
4. Those with skin disease or infectious disease.
5. Those who are sensitive to the electricity.
6. Those who are under 16 years old.

4. Dos and Dont's

1. In order to ensure the therapeutic effect and normal service life expectancy of the device, please use the specified accessories provided or recommended by the original manufacturer.
2. The device should not be placed in wet environment or near water, nor should they be directly exposed to sunlight.
3. Do not place the device near a strong heat source, because this may affect the life and normal use of the device.
4. Please restrain from using the device on eyes, thyroid, parathyroid, testicles, abdomen of pregnant women, and pacemakers.
5. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
6. Use device or train device operators in strict accordance with instructions in the user manual.
7. For the first time of use, the minimum level is advised, adjust the level when adapted.
8. When the power is low, please remember to charge it in time to avoid power failure, which will affect the effect.
9. Do not use the instrument while charging to avoid electric shock or affecting the service life of the instrument.
10. Please use this instrument with caution if you are sensitive to current.
11. Clean your eyes area with wet towel before using this instrument, to avoid the cross infection from the bacteria around eyes. Each pair of eye patches are advised to use 5~6 times at most.
12. Do not use this instrument at heart, front of the neck and temples of both sides directly.
13. Do not use it while driving.
14. Do not use this instrument in bathtub or shower.
15. Do not immerse the instrument in any liquid.
16. Please make sure that both eye patches are tightly attached to the skin before starting the device. The device will not work normally if the eye patches are not tightly attached to the skin or only single eye patch is in use.
17. Tightly attach both eye patches to the skin first, then start the device.

5. Troubleshooting & Solutions

1. The device cannot be started?

- A. Check out the battery is enough or not.
- B. Check out whether you are long pressing the keypad for 3 seconds, tap cannot start the device.

2. There is no any feeling after the device starts?

- A. Check out the connecting cord buckles the eye patches.
- B. Check out the eye patches attach tightly to the skin.

3. The device is turned off automatically after it starts?

- A. The battery may be low and you need to charge it.
- B. Make sure you are using a pair of eye patches rather than one single eye patch.
- C. Make sure both eye patches attach tightly to the skin.
- D. If the above methods cannot be handled, please contact the instrument dealer for assistance.

6. FAQs

1. Q: How often do I use this device? How long does it take for one use?

A: This device can be used every day. It is usually used for 1-2 times a day, 10 min/time. It is easy to be carried that can be used at home or outside.

2. Q: Can I use this device if I have fine lines around my eyes?

A: Yes. This device combines RF and EMS technologies. RF high-frequency waves can heat up the dermis, promote the skin's metabolism and collagen regeneration, thus lifting skin, removing wrinkles, anti-aging and making skin tightened and tender. The EMS micro-current stimulates the microvessels and nerve triggering cells in the skin layer and epidermis layer to produce hydration, replenish the nutrients and oxygen that the wrinkled skin lacks, and restore the elasticity of sagging skin.

3. Q: Dose this device work to eyestrain, dry eyes and eye irritation?

A: Yes. The EMS function of this device can stimulate muscles exercising, relax muscles and has bionic eye exercises to reach the effects of relieving eyestrain, dry eyes and eye irritation.

7. Packing List

- 2*eye patches
- 1*main device

8. Operational Diagrams

Parameters	Product	Techniques	Diagrams
Eyes Health Treatment: 15 min, 1-2 times/day			
Level: 1~3 (adjust according to one's feeling)	Disposable Wipes+ Eye Patches	<ol style="list-style-type: none"> 1. First install the connecting cord to the main device, then buckle the connecting cord to the eye patches. 2. Clean the eyes area with wipes, 1 min. 3. Attach the eye patch to the area from lower eyelid to eye corner. 4. Do the same on the other eye. 5. Long press for 3 seconds to turn it on. 6. Select proper level to start operation. 7. Wait for around 10 min. 8. Take off the eye patches. 9. Close eyes and have a rest for 5 min. 	No diagrams
<p>Suggested Treatment: One month of treatments account for one full course. One-time of use can relieve eyestrain and ache eyes. One-month of continuous use can improve ache, sore, bloat, dry and irritation of eyes. Two-month of continuous use can cure ache, sore, bloat of eyes, eyestrain and glazed eyes caused by over use of eyes. Long term of use can care your eyes, protect your eyes, which can prevent dry eyes, eye irritation, blurred vision and other eye diseases.</p>			
Eyes Care: 25 min, 3-5 times/week			
Level: 2~3 (adjust according to one's feeling)	Makeup Remover+ Face Wash+ Essence+	<ol style="list-style-type: none"> 1. First install the connecting cord to the main device, then buckle the connecting cord to the eye patches. 	No diagrams

	Eye Patches+ Eye Mask	<ol style="list-style-type: none"> 2. Remove makeups and clean the face. 3. Apply toner. 4. Apply eye essence evenly to the eyes area until the essence is absorbed. 5. Attach the eye patch to the area from lower eyelid to eye corner. 6. Do the same on the other eye. 7. Long press for 3 seconds to turn it on. 8. Select proper level to start operation. 9. Wait for around 10 min. 10. Take off the eye patches. 11. Apply the eye mask, 10 min. 12. Take off the eye mask and clean the area around eyes. 13. Apply toner, essence, eye cream, face cream and sunscreen. 	
--	-----------------------	--	--

Suggested Treatment:

One month of treatments account for one full course. In the early stage, treatments are advised for every day. After 15 days of consecutive care, you can have treatments for 3~5 times a week. One-time of care can relieve eyestrain, eyes edema. One month of care can promote blood circulation, improve eyes edema and fade dark circles. Two months of care can improve eye bags, fine lines around eyes and increase the elasticity of skin. Long term of care can strengthen the effect, prevent anti-aging, flabby skin and wrinkles around eyes.