

# R79-3S

VER:1.0

## USER MANUAL



---

## Warning

Dear user:

Thanks for using our company's products, in order to fully make use of the products, we sincerely suggests follow:

- ◆◆ Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆ Please follow the instruction to fit and operate the machine correctly.
- ◆◆ Please do not remove and change any accessories of the machine.
- ◆◆ It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆ It forbids others to change and maintain the machine except the professional people.
- ◆◆ Please contact us if there is any trouble during operation.

**This product does not belong to medical,  
please don't use for disease treatment!**

## Parts Introduction



- 1.Touch screen
- 2.Rotary encoder adjustment knob/confirming button
- 3.Output indicator
- 4.Accessories frame
- 5.Upgrade socket
- 6.Finger jack
- 7.General power switch
- 8.Power supply lines
- 9.Distilled water inlet
- 10.Cryo head handle jack
- 11.Distilled water level window
- 12.Distilled water outlet
- 13.Overflow




**The host does not have a cart by default, if you need a cart, please purchase separately.**

**Any part of instrument made of PMMA material (for example: accessories shelf, panel, transparent handles,etc.). Do not use alcohol or corrosive solvent to wipe clean.**

## Accessories



cryo head handle



Set box




Finger(XL L M S)

## Transport Preparation

1. The distilled water should be cleared before transportation.

2. Method:


Insert the plug of silicone tube with seal ring into socket [12] distilled water outlet. Discharge the distilled water and shake the apparatus to help the liquid clean out completely. And then use thumb to press the spring, pull out the silicone tube.



**Notice: While discharging the distilled water, the handle should be placed at the handle support.**

**3. One the plug of silicone tube insert into inlet and the other plug of it insert into overflow during transportation. Remove the silicone tube when in normal use in case that effect the vacuum of cryolipolysis handle.**

**4. For transport ,please pull out ice sculpture head handle, to plug with a hose Ice sculpture head handle jack.**



---

---

## Add Distilled Water Steps

1. Use distilled water or deionized water;
2. Make sure the distilled water outlet[12] is closed;
3. Use special funnel to connect to Distilled water inlet [9];
4. Add distilled water or deionized water from the funnel;
5. After filling it with water, take off the funnel.

### **NOTE:**

1. Add about 500~700ml distilled water or deionized, change the water after one week.
2. When could not add into the distilled water, please open the distilled water outlet[12], Add the distilled water and observe the valve. when the liquid is smooth outflow, blocked it and continue to fill...
3. If do not use the apparatus 15 days or above, we are recommended to put the distilled water clean and clear.

## Basal Operation

1. Put the accessories away and connect them to the mainframe.

(1) Insert the cryo head handle into [10] jack;

(2) Insert the finger into [6] jack;

**Notice: In normal use, Do not insert frequently ice sculpture head handle.**


2. Insert the plug of Power supply line into [8] socket, turn on Power switch [7]. input default password "000000" to (Fig 1) enter interface as shown in (Fig 1).

3. In the interface as (Fig 1):

(1) Click on "Basic setting" to select Language, set Date and Time.



**Attention: The "Advanced setting" is only for manufacturer, Users do not use.**

(2) In the interface as (Fig 1) Click on "Items" to enter operating interface.


4. This machine can click the keys on the touchscreen to set the working parameters, it also can use [2] Function selecting knob & Confirming button to set. The detail operation method as follows: Turn the knob when you choose the function, the icon "  " move to the right place, then press the knob [2] once to confirm.

5. This use's manual use touchscreen for example.

6. **Use the finger sleeve attention: do not cross use.**



Correct




Wrong




## RV FINGER Treatment

### Operation steps

1. Click "RV finger" in the interface as shown in (fig 1) to enter the nursing interface as shown in (fig 2).
2. Clean and care parts.
3. The default channel of the system is "R" finger, Click on L (left hand) and R (right hand) finger are selected according to the nursing site. Both can be done at the same time.



(fig 2)

4. Set work time. The system default is 30 minutes. Click on “-” and “+” adjust. The range is 1 to 120 minutes. Each part about 10 to 20 minutes.
5. Set intensity. You can click on the part of the person icon to adjust the strength. Click on “-” and “+” adjust. The range is 1(weak) to 10(strong).
  - a. Face selected default strength “1” ,The range is 1(weak) to 10(strong).
  - b. Abdominal selection default strength “4” ,The range is 1(weak) to 15(strong).
  - c. The back default strength “2” ,The range is 1(weak) to 10(strong).
6. Set Q-IR mode. The system default is ON. Click on “ Q-IR ” icon. The icon is marked “ON”, The icon is darkened to “Off”.
7. Set the vibration. The system default is ON. Click on “ VIBRATION ” icon. The icon is marked “ON”, The icon is darkened to “Off”.
8. Set depth. The system default is Shallow. Click on the “ Shallow ” “Deep” conversion.
9. Apply gels on operating parts. Click on “  ” , If you need to pause, Click on “  ” .
10. The beautician holds the handle and moves slowly in circle. Do not stop at the same parts without operating.
11. The system will pause automatically if setting time is over.
12. Clean the leftovers on skin with warm water.
13. Please turn off the general power switch and unplug the power supply line. Then clean the finger covers with a clean towel.
14. Click on “  ” to return to (Fig 1) .

### Cautions:

1. Do not wear any aglet, both for the beautician and the users. (Contraindications for other items)
2. The electrode probe must be fully in contact with skin.
3. Do not make the electrode stay on any part without moving. Keep it moving on skin slowly.
4. Skin around eyes and forehead is very thin, the output intensity should be relatively weakened.
5. Do not let the Q-IR light straightly shine upon the eyes in case of eyes hurt.


## EMS FINGER Treatment

### Operation steps

1. Click "EMS finger" in the interface as shown in (fig 1) to enter the nursing interface as shown in (fig 3).
2. Clean and care parts.
3. The default channel of the system is "R" finger, Click on L (left hand) and R (right hand) finger are selected according to the nursing site. Both can be done at the same time.
4. Set work time. The system default is 30 minutes. Click on “-” and “+” adjust. The range is 1 to 120 minutes. Each part about 10 to 20 minutes.
5. Set intensity. You can click on the part of the person icon to adjust the strength. Click on “-” and “+” adjust. The range is 1(weak) to 10(strong).
  - a. Face selected default strength “1” ,The range is 1(weak) to 10(strong).
  - b. Abdominal selection default strength “4” ,The range is 1(weak) to 15(strong).
  - c. The back default strength “2” ,The range is 1(weak) to 10(strong).
6. Set Q-IR mode. The system default is ON. Click on “ Q-IR ” icon. The icon is marked “ON”, The icon is darkened to “Off”.
7. Set the vibration. The system default is ON. Click on “ VIBRATION ” icon. The icon is marked “ON”, The icon is darkened to “Off”.
8. Set mode. The system default is P1, Click on “-” and “+” adjust. The range is P1(slow) to P3(fast).
9. Apply gels on operating parts. Click on “⏸” , If you need to pause, Click on “▶” .
10. The beautician holds the handle and moves slowly in circle. Do not stop at the same parts without operating.
11. The system will pause automatically if setting time is over.
12. Clean the leftovers on skin with warm water.
13. Please turn off the general power switch and unplug the power supply line. Then clean the finger covers with a clean towel.
14. Click on “🏠” to return to (Fig 1).

### Cautions: ⚠

1. Do not wear any aglet, both for the beautician and the users.(Contraindications for other items)
2. The electrode probe must be fully in contact with skin.
3. Do not make the electrode stay on any part without moving. Keep it moving on skin slowly.
4. Skin around eyes and forehead is very thin, the output intensity should be relatively weakened.
5. Do not let the Q-IR light straightly shine upon the eyes in case of eyes hurt.




(fig 3)



## CRYO & Heat Treatment

### Cryo and Heat Treatment operation steps

1. Click "CRYO" in the interface as shown in figure 1 to enter the nursing interface as shown in (fig 4).
2. Clean and care parts.
3. The system default is "CRYO" therapeutics. Click "HEAT" to select heat therapy.
4. Set work time. The system default is 30 minutes. Click on "-" and "+" adjust. The range is 1 to 120 minutes. Each part about 10 to 20 minutes.



(fig 4)

5. Set cool temperature. You can click on the person icon to adjust the temperature, click on "-" and "+" adjust, The ICE range is  $-5^{\circ}\text{C}$  to  $20^{\circ}\text{C}$ . The heat temperature range is  $37^{\circ}\text{C}$  to  $55^{\circ}\text{C}$ .
  - a. Face selected default refrigeration intensity  $-1^{\circ}\text{C}$ , Heat temperature  $42^{\circ}\text{C}$ .
  - b. Abdominal selection default refrigeration intensity  $-5^{\circ}\text{C}$ , Heat temperature  $45^{\circ}\text{C}$ .
  - c. The back default refrigeration intensity  $-0^{\circ}\text{C}$ , Heat temperature  $40^{\circ}\text{C}$ .
6. Do cool treatment according to zone, and adjust appropriate required temperature.
7. Apply gels on operating parts. , If you need to pause, Click on "▶" key.
8. The beautician holds the handle and moves slowly in circle. Do not stop at the same parts without operating.
9. The system will pause automatically if setting time is over.
10. Clean the leftovers on skin with warm water.
11. Please turn off the general power switch and unplug the power supply line. Then clean the finger covers with a clean towel.
12. Click on "🏠" to return to (Fig 1).

### Cautions: ⚠️

1. Do not wear any aglet, both for the beautician and the users. (Contraindications for other items)
2. The electrode probe must be fully in contact with skin.
3. Do not make the electrode stay on any part without moving. Keep it moving on skin slowly.
4. Skin around eyes and forehead is very thin, the output intensity should be relatively weakened.

---

---

## Cautions During Treatment

- 1.The gamepad function can only be started when the gamepad is correctly connected to the host computer.
- 2.Don't operate the machine on the electriferous floor.
- 3.Other people DO NOT touch the electrode flat when the machine is operating.
- 4.The operator and customer DO NOT be with the metal decoration during in treatment time.
- 5.The handle should be in full contact with the skin during nursing operation, as improper contact may cause irritation and burn to the skin.
- 6.The electrode probe is moved on the skin by the circulation way, but don't stop on one part .
- 7.Body care can choose greater energy, but not too high to make people feel uncomfortable.

## Forbidden Group

- 1.People with heart disease or high blood pressure, or configured cardiac pacemaker.
- 2.Patients with acute inflammation, asthma, deep vein thrombosis, thyroncus, cancer.
- 3.People with hemorrhagic disease, trauma or who is bleeding.
- 4.Pregnant women and child.
- 5.Medical Plastic parts, or parts with artificial in filling inside.
- 6.People with metal inside the body(include metallic tooth).
- 7.Patient with skin inflammation or with edema.
- 8.People with an abnormal immune system.
- 9.Numb or insensitive to heat.

---

---

## Technical Parameters


Input Voltage: AC220~230V 50/60Hz □ AC110V 50/60Hz □

Output Power: 135VA






Fuse: FL5X20\3.15A with buffer (with spare fuse)

## Packing List

1.Host	1SET	<input type="checkbox"/>
2.Large finger	2PCS	<input type="checkbox"/>
3.Middle finger	2PCS	<input type="checkbox"/>
4.Small finger	2PCS	<input type="checkbox"/>
5.Enlarge finger	2PCS	<input type="checkbox"/>
6.User manual	1COPY	<input type="checkbox"/>
7.Accessory frame	1PCS	<input type="checkbox"/>
8.Hexagonal screw(M4X8)	2PCS	<input type="checkbox"/>
9.Hexagonal screw wrench(M3)	1PCS	<input type="checkbox"/>
10.Set box	2PCS	<input type="checkbox"/>
11.Two pieces of liquid addition accessory	1PCS	<input type="checkbox"/>
12.Silicone tube	1PCS	<input type="checkbox"/>
13.cryo head handle	1PCS	<input type="checkbox"/>
14.Slider gel	1PCS	<input type="checkbox"/>
15.Sealing ring	10PCS	<input type="checkbox"/>
16. Power supply line	1PCS	<input type="checkbox"/>
17.Set box armllet	2PCS	<input type="checkbox"/>



## Symbols and Meanings

Warning Symbol	Meaning
 Danger	If used incorrectly will bring the risk of death or serious injury!
 Warning	If used incorrectly will cause possibility ill of death or serious injury!
 Attention	If used incorrectly will cause personnel injury or damage to the goods!
 Danger	Device exists high pressure risk of high voltage output terminals!
 Warning	If used incorrectly will cause personnel injury!

\*Things damage refers to housing, property, livestock and pet damage.

## Contraindications

1. The following symptoms were forbidden to use

- (1) Patients who used by Heart rate regulator, cardiac pacemaker, ECS;
- (2) The use of artificial bone; Patients with silicosis;
- (3) Drunk or high fever patients;
- (4) Patients with higher sensitivity mental disease and epilepsy patients
- (5) Patients with no independent expression consciousness and can not take care of themselves.
- (6) Patients with metal materials in body (including metal teeth), pregnant women, cachexia, malignant tumor.
- (7) Severe hypertension and severe patients with cardiovascular and cerebrovascular disease.
- (8) Patients with renal failure.
- (9) Patients with severe infection



2. The following please use under the guidance of professionals

- (1) Women during menstruation;
- (2) Children;
- (3) Hemorrhagic disease or trauma, inflammation, skin disease, skin infection;
- (4) Patients with feeling numb, no response to cold and heat;
- (5) The immune system has defects, scar formation history, physical disease;
- (6) The people who has taken face-lifting or her body has artificial thing.



## Maintenance

1. Appliances and medical equipment close to the machine may be resulting in operational failure due to interference.
2. Prohibit use of the equipment in the following environment
  - Hot and humid environment
  - Dusty environment
  - Environmental with much smoke or steam
3. If there comes a failure, please shut down the power immediately, and send it to the distributor or manufacturer for maintenance. Do not repair by yourself.

#### 4. Exclusion of liability.

**The Company is not responsible for the following circumstances:**

- (1) If fire, earthquake, war occurred;
  - (2) The fault of the user, operator error;
  - (3) Damage caused by using other method (not included in the use manual).
  - (4) Operate without the manual's guide.
5. Do not make the machine close to the wall, keeping a 30CM space around for cooling.

## Safety and Environmental Considerations

### Toxic and hazardous substances reference list Electronic information product pollution control and management measures Table of contents of toxic and hazardous substances or elements

Part Name	Toxic and Hazardous Substance or Element						Remark
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr(VI))	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)	
Plastic Board	○	○	○	○	○	○	
Printed Circuit Board	×	○	○	○	○	○	
Power Supply	×	○	○	○	○	○	
Accessory	×	○	○	○	○	○	

#### Declaration:

- The content of the restricted substances in all homogeneous materials of this part is below the limits set by the SJ/T11363-2006 standard.
  - ×
- The content of a homogeneous material in the part of the material is beyond the SJ/T11363-2006 standard requirements. The user's normal use of this product within the time limit of environmental protection, these substances or elements will not occur leakage or sudden change, will not cause the personal, property damage. For such substances or elements, the user shall not dispose of it yourself, please according to government regulations to the relevant departments of the government designated recycling.