

User Manual

SR-AE8782



Preface

Dear users:

We're pleased to present to you *5 IN 1 Multifunction Electric Face Facial Cleansing Brush Spa Skin Care Massage*. It cleans and massages face. By adopting dual frequency + rotation technology, it realizes effects of cleansing and massaging. *5 IN 1 Multifunction Electric Face Facial Cleansing Brush Spa Skin Care Massage* is a specialized skin care device for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore we advise any personnel read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

Table of Contents

Preface
Part I
Brief Introduction
Advantages
Working Principles
Main Effects
Applicable Range
Inapplicable Range
Notes after Treatment
Part II
1. Detailed Operations
2. Safety Precautions
3. Dos and Don'ts
4. Troubleshooting & Solutions
5. FAQs
6. Packing List
7. Operational Diagrams

Part I

Brief Introduction

5 IN 1 Multifunction Electric Face Facial Cleansing Brush Spa Skin Care Massage can cleanse skin, unclog pores and at the same time massage face and eyes, boost skin absorption. Ergonomic design brings you a comfortable, effort-saving experience. With its powerful effect, it can replace traditional cleansing, removing blackheads, exfoliating and massage devices. Moreover, it's convenient and easy to operate. This device requires no injection, no medication, no operation, thus having no side-effects. It is used externally throughout the entire process and comes with immediate result. Its super soft bristles can reach deep into pores and clean debris and gunks, helping skin metabolize so as to reduce skin issues of breakouts, blackheads, enlarged pores and so on.

Advantages

1. Different operation heads can be changed according to the desired effect.
2. Small, convenient, portable and easy to operate.
3. Simple operation, easy to learn, remarkable effect and no side effects.
4. High-quality materials, soft and not harmful to the skin, safe and environmental-friendly.
5. Streamlined, ergonomic design, more convenient to hold, more comfortable.
6. The speed of the two gears can be adjusted and controlled freely, which is very convenient.
7. Adopt cutting-edge technology -- to solve four facial problems: dull sallow skin, brightening, acne removal, exfoliating, and blackhead removal.

Working Principle

The high frequency vibration principle of ultrasonic wave is mainly used to make the brush bristles vibrate, so as to achieve efficient cleaning and massage effect without over stimulation. It can help cleanser form high density even foam, thoroughly clean the face, and achieve the beauty result of skin care.

Nylon brush head:

Function: Thoroughly clean the dirt that sticks to the face

Usage: Clean your face with facial cleanser

Sponge brush head:

Function: Remove cosmetic residues, absorb oil and make it easy to remove makeup, clean skin, and achieve the effect of massage and cleaning at the same time

Usage: Clean your face with facial cleanser

Latex sponge head:

Function: Apply to around the eyes and cheeks, massage and clean, effectively relieving eye fatigue, reducing dark circles, under eye bags

Usage: Clean your face with facial cleanser

Exfoliating massage head:

Function: For hands and feet, exfoliate and exfoliate dead skin to make skin smooth and delicate

Usage: For aged cutin of hands and feet

Rolling massage head:

Function: Promote blood circulation, making skin more elastic. Long-term use has certain effect face toning

Usage: It can be used on the face, hands and legs. Increase level leads to obvious effect on face slimming. It can be used with anti-wrinkle cream, whitening cream and essential oil to boost skin absorption.

Main Effects

1. Deeply clean pores and promote skin absorption.
2. Improve dark yellow and dull skin.
3. Exfoliate and smooth the skin.
4. Reduce acne and blackheads.
5. Minimize large pores and improve coarse skin.
6. Massage your face to promote blood circulation and metabolism.
7. Adjust skin oil secretion and improve skin quality.

Applicable Range

1. Those with exuberant oil secretion, breakout and blackheads.
2. Those who don't clean for a long time and often wear makeup.
3. Those whose needs exfoliating.
4. Those with poor skin absorption, metabolism and lack of nutrition.
5. Those with large pores, coarse skin and dull complexion.
6. Those with eye strain, sore eyes and dark circles.

Inapplicable Range

1. Use with caution for those have serious heart disease.

2. Those with severe skin allergies, pustules and ulcers.
3. Those with skin diseases and infectious diseases, such as eczema and psoriasis.
4. Those with skin lesions, inflammation and other local signs of disease.
5. Those who have just had microneedling, intradermal injection, plastic surgery or implants.

Notes after Treatment

1. Always hydrate and use sun protection.
2. It is advised to wear a mask at least 3 times a week.


Part II

1. Detailed Operation




Pull and switch different brush heads for different skin needs.



Install battery;  AA battery



Take off the battery cover and install; have the battery - end 

pushes against the spring 



ON/OFF & Frequency +/-

There's low frequency and high frequency.

Low frequency is suitable for first-time users and those who have thin horny skin layers.

Low frequency can be used for facial massage.

High frequency is suitable for those who adapt quickly and have thick horny skin layers.

High frequency can be used for body massage.



Nylon brush head:

Function: Thoroughly clean the dirt that sticks to the face

Usage: Clean your face with facial cleanser

**Sponge brush head:**

Function: Remove cosmetic residues, absorb oil and make it easy to remove makeup, clean skin, and achieve the effect of massage and cleaning at the same time

Usage: Clean your face with facial cleanser

**Latex sponge head:**

Function: Apply to around the eyes and cheeks, massage and clean, effectively relieving eye fatigue, reducing dark circles, under eye bags

Usage: Clean your face with facial cleanser

**Exfoliating massage head:**

Function: For hands and feet, exfoliate and exfoliate dead skin to make skin smooth and delicate

Usage: For aged cutin of hands and feet

**Rolling massage head:**

Function: Promote blood circulation, making skin more elastic. Long-term use has certain effect face toning

Usage: It can be used on the face, hands and legs. Increase level leads to obvious effect on face slimming. It can be used with anti-wrinkle cream, whitening cream and essential oil to boost skin absorption.

How To Use

1. Put in 2 AA batteries into the device



2. Moisten face with warm water; squeeze some cleanser on palm, and add a bit warm water; lather up and evenly apply on facial area



3. Turn on the device and start cleaning; cleanse downwards, tracing T zone; cleanse from cheeks on two sides inward; foamy cleanser brings better cleansing effect



4. Cleanse and massage for 1-3 min; wash face clean with warm water, which is more gentle



5. Gently wipe off the remaining water on the face. Gently pat. Do not rub with strong force.



6. Apply skin care products. Gently pat to help absorb. Switch to the rolling massage head. Turn on the device and massage with low frequency to help absorb.

7. Turn off the device and detach the massage head. Rinse it off with clean water. Be gently so as to not damage the massage head.

2. Safety Precautions

Patients with the following symptoms should be cautious about using this instrument. Please consult a doctor or a professional before using the instrument. The details are as follows:

1. Use with caution for those have serious heart disease.
2. Use with caution for those who are pregnant or sick.
3. Those with severe skin allergies, pustules and ulcers.
4. Those with skin diseases and infectious diseases, such as eczema and psoriasis.
5. Those with skin lesions, inflammation and other local signs of disease.

3. Dos and Dont's

1. The device should not be placed in wet environment or near water, nor should they be directly exposed to sunlight.
2. Do not place the device near a strong heat source, because this may affect the life and normal use of the device.
3. Take the batteries out when not using the device for a long time.
4. When damaged, contact professionals for maintenance before opening the device up yourself.
5. Pregnant ladies or patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
6. Use device or train device operators in strict accordance with instructions in the user manual.
7. Do not soak the device in water.
8. Keep the treatment area moist when operating. Avoid cleaning or massaging on dry skin.
9. Start with low frequency for first-time users and change to high frequency after adaptation. Adjust based on one's own comfort level.
10. Do not use exfoliating massage head on face. It will does damage to facial skin when it does not need exfoliating.
11. Do not press too hard or stay at one spot for long to avoid hurting skin when using exfoliating massage head.
12. For those who have just undergone microneedling, intradermal injection, plastic surgeries, etc., wait for skin to repair before using the device.
13. Please cease immediately if experiencing discomfort.
14. Frequent use of the device will not make you reliant on it.

4. Troubleshooting & Solutions

1. I cannot start the device?

- A. Check if you've put in batteries.
- B. Check if you've place the batteries in wrong directions.
- C. Check if your batteries have run out of power.

2. How come the device stopped working in the middle of the treatment

- A. Check if your batteries have run out of power.
- B. Restart the device.

5. FAQs

1. Q: Can I use the device everyday?

A: Yes. It can be used everyday. The operation head is made of premium materials,

which is soft and does not hurt the skin. When using the brush head, you can first wet the brush head and gently brush on the skin surface. You must use the strength well and do not use too much force.

2.Q: Can I use the device to remove makeup?

A: It is not recommended to use this device to remove makeup. In general, makeup contains a lot of oil. If you use the device to remove makeup, it is easy to accelerate the aging of the brush head. It's best to remove makeup before cleansing skin with sponge brush head.

3. Q: How often should I switch brush head?

A: Under normal circumstances, the brush head will wear off after being used for a period of time, and the ability to clean the skin's oil and dirt will also be reduced. In order to ensure the cleaning ability of the device and the hygiene of the brush, it is recommended to replace the brush head for three months, so as to ensure the cleanliness and the strength of cleaning.

4. Q: Does it applicable to all skin types?

A: No, this device is a mild cleansing and massage device, but if the skin is very sensitive, often grow acne, and has cuts, then it is not recommended to use this device to wash your face.

6. Packing List

1 x Main machine

1 x Nylon brush head

1 x Sponge brush head

1 x Latex sponge head

1 x Exfoliating massage head

1 x Rolling massage head

7.Operational Diagram

Parameter Adjustment	Product	Techniques	Diagrams
Facial Cleaning: 20 min, 5-7 sessions/week			
<p>Work Head: Nylon brush head</p> <p>Frequency: Cheeks: LOW T Zone: HIGH (including chin, nose and forehead)</p>	<p>Makeup remover+ Cleanser+ Toner+ Sheet mask</p>	<ol style="list-style-type: none"> 1. Install nylon brush head. 2. Wet the surface with water. 3. Remove makeup and dust off face. 4. Squeeze pea sized cleanser on the brush head. 5. Squeeze certain amount of cleanser on palm and add warm water. Lather up and apply foam evenly on forehead, nose, chin and two cheeks. Gently massage face for 3-5 seconds and make sure foam is evenly spread all over the facial area. 6. Adjust frequency and start the treatment. 7. Move slowly up and down the nose, 3-5 times. 8. Move slowly from mouth corner to around the mouth, 3-5 times. 9. Move slowly up and down the cheek, 2-3 times. 10. Do the same on the other side. 11. Move in circular motions on forehead, 3-5 times. 12. Wash face clean with warm water. 	<p>No Diagrams</p>

		13. Apply toner. 14. Apply sheet mask, 15 min. 15. Wash face clean. 16. Apply toner, lotion, cream and sunscreen.	
Suggested Treatment: Basic facial cleaning requires long-term use. After the first session, no visible change can be seen but skin feels smoother. After a month of persistent use, skin is less congested. 2 months of persistent use can reduce blackheads, shrink pores and reduce excessive oil secretion. 3 months of persistent use can help skin absorb better, and metabolize better. Long-term persistent use can help improve skin texture, [revent blackhead, breakouts and pigmentation.			
Note: Sponge brush head shares the same techniques with nylon brush head.			
Facial Care: 30 min, 5-7 sessions/week			
Work Head: Cleaning: Nylon brush head Massage: Rolling massage head Frequency: Cleaning: Cheeks: LOW T Zone: HIGH (including chin, nose and forehead) Massage: LOW Adjust based on one's own comfort level	Makeup remover+ Cleanser+ Toner+ Massage cream+ Sheet mask	1. Install latex sponge head. 2. Wet the surface with water. 3. Remove makeup and dust off face. 4. Squeeze pea sized cleanser on the brush head. 5. Squeeze certain amount of cleanser on palm and add warm water. Lather up and apply foam evenly on forehead, nose, chin and two cheeks. Gently massage face for 3-5 seconds and make sure foam is evenly spread all over the facial area. 6. Adjust frequency and start the treatment. 7. Move slowly from between eyebrows towards the end of the eyebrow, 3-5 times. 8. Move slowly up and down the nose, 3-5 times. 9. Move slowly from mouth	No Diagrams

		<p>corner to around the mouth, 3-5 times.</p> <p>10. Move in circular motions around cheek, 203 times.</p> <p>11. Move slowly from lower eyelid towards the temple, 2-3 times.</p> <p>12. Move slowly from the beginning of the eyebrow towards the end, 2-3 times.</p> <p>13. Do the same on the other side.</p> <p>14. Move in circular motions on forehead, 3-5 times.</p> <p>15. Turn off the device and replace the brush head with rolling massage head.</p> <p>16. Wash face clean with warm water.</p> <p>17. Apply toner.</p> <p>18. Apply anti-wrinkle cream evenly all over the face.</p> <p>19. Adjust frequency and start the treatment.</p> <p>20. Move in parallel lines from chin to earlobe, mouth corner to ear gate, nose wing to temple to lift the treatment site, 3 times.</p> <p>21. Do the same on the other side.</p> <p>22. Lift towards hairline on forehead, 3 times.</p> <p>23. Apply sheet mask, 15 min.</p> <p>24. Wash face clean.</p> <p>25. Apply toner, lotion, cream, eye cream and sunscreen.</p>	
--	--	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--

Suggested Treatment:

Basic facial care requires long-term use. After the first session, no visible change can be

seen but skin feels smoother. After a month of persistent use, skin is less congested, skin absorbs better and facial fine lines faded. 2 months of persistent use can reduce blackheads, shrink pores, fade fine lines and brighten skin. 3 months of persistent use can help smooth out skin, and fade fine lines. Long-term persistent use can help improve skin texture, increase skin elasticity, prevent blackheads, breakouts and pigmentation.

Hand Exfoliating: 30 min, once/week

<p>Work Head: Exfoliating massage head</p> <p>Frequency: LOW</p>	<p>Exfoliating cream+ Hand mask+ Hand cream+ Towel</p>	<ol style="list-style-type: none"> 1. Install exfoliating massage head. 2. Soak hands in warm water for 2-3 min. 3. Apply certain amount of exfoliating cream on the back of the hand and wrist. 4. Adjust frequency and start the treatment. 5. Move slowly on wrist, 3 times. 6. Move in circular motions on the back of the hand, 3 times. 7. Move slowly back and forth on fingers. Start from thumb and move to the pinkie in sequence, 2-3 times. 8. Do the same on the other side. 9. Wash hands clean with warm water. 10. Apply hand mask, 10-15 min. 11. Wash hands clean. 12. Dry the treated area with towel. 13. Apply moisturizing hand cream evenly on both hands. 	<p>No Diagrams</p>
------------------------------------------------------------------------------	------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------

Suggested Treatment:
10 treatments account for one full course. Skin metabolizes fast in summer so exfoliate once a week. Skin metabolizes slower in winter so exfoliate once every two weeks.

After one session, skin is visible less coarse. After a full course of persistent use, skin is smooth and brightened. Two months of persistent use help brighten and smooth out skin. Three courses help enhance effect and avoid rebound. Long-term use help get ideal hand skin.

Foot Exfoliating: 30 min, once/week

<p>Work Head: Exfoliating massage head</p> <p>Frequency: HIGH</p>	<p>Exfoliating cream+ Lotion+ Towel</p>	<ol style="list-style-type: none"> 1. Install exfoliating massage head. 2. Soak feet in warm water for 5-8 min. 3. Apply certain amount of exfoliating cream on the feet. 4. Adjust frequency and start the treatment. 5. Move in circular motions around ankles, 3 times. 6. Move slowly back and forth on instep, 3 times. 7. Move in circular motions on toes. Start from thumb toe and move to the pinkie toe in sequence, 2-3 times. 8. Pause the device. 9. Apply certain amount of exfoliating cream on the sole. 10. Adjust frequency and start the treatment. 11. Move in circular motions on the sole, 2-3 times. 12. Move slowly back and forth along the sole, 2-3 times. 13. Do the same on the other foot. 14. Wash feet clean with warm water. 15. Dry the treated area with towel. 16. Apply moisturizing foot 	<p>No Diagrams</p>
-------------------------------------------------------------------------------	-----------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------

		cream evenly on both feet.	
--	--	----------------------------	--

Suggested Treatment:

10 treatments account for one full course. Skin metabolizes fast in summer so exfoliate once a week. Skin metabolizes slower in winter so exfoliate once every two weeks. After one session, skin is visible less coarse. After a full course of persistent use, skin is much tender smooth and refined. Two months of persistent use help brighten and smooth out skin. Three courses help enhance effect and avoid rebound. Long-term use help prevent accumulation of dead skin.