

# User Manual

SR-JK673



# Preface

Dear users:

We're pleased to present to you our Fractional RF, a highly efficient anti-aging beauty machine. It targets at both face and body, mainly focusing on lifting and tightening of face, removing wrinkles, lifting skin around eyes and neck, and face line shaping and body shaping. We aim at a safe and, at the same time, effective beauty result by using high-end technologies. Fractional RF beauty machines are anti-aging, machines for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse consequences. Therefore we advise any personnel read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

# Table of Contents

Preface
Part I
Brief Introduction
Advantages
Principles & Treatments
Anti-Aging: Face
Anti-Aging: Eyes
Anti-Aging: Neck
Repair Post-Pregnancy Stretch Marks
Part II
1. Detailed Operations
2. Technical Parameters
3. Safety Precautions
4. Dos and Don'ts
5. Troubleshooting & Solutions
6.FAQs
7.Packing List
8.Operational Diagrams

# Part I

## Brief Introduction

Fractional RF beauty machines are currently the most popular beauty machines that work on tightening, body-shaping and anti-aging. It's an effective replacement of all the facial and body care programs in beauty salons. Moreover, they are convenient and easy to operate. It solves both skin and body management issues for beauty-lovers. This equipment requires no injection, no medication, no operation, thus having no side-effects. It is used externally during the whole process and has an immediate effect. The unique treatment head transmits high frequency energy waves to the skin layer, causing collagen to contract, resulting in immediate skin firming effect and stimulating collagen regeneration in the skin to achieve long-lasting skin lifting and firming effect.

## Advantages

1. Multiple functions, wide range of use and significant effects.
2. High energy and high frequency, maximum effect, accurate positioning, balanced and comprehensive usage.
3. There is no recovery period after treatment, so it can be done as you go, which will not affect your work and life.
4. Unique patented technology -- the key to cluster heat source, bonding with surrounding collagen to form a stable support.
5. Effective thermal damage of heat in dermal fiber layer, thus truly stimulating collagen regeneration.
6. No surgery, no injection, high safety, no side effects.
7. Once treated, deep skin collagen can last for up to six months to a year with long-lasting effects.
8. The cool hammer equipped with the instrument has sedative repair effect, which will prevent skin from redness or swelling.

## Anti-aging: Face

### Principles

## Fractional RF : Face (Probe B)

Through transferring numerous cluster heat, this unique probe can send the matrix molecular energy waves of up to 6 million times per second, repair skin in deep layers, from the focal facade of RF in hypodermis, make a strong impact on dermal tissues, which produces an aggregation effect of electric field in the dermal tissues, make the skin temperature up to 60 °C or so, produce immediately shrinkage, using the electromagnetic wave in the skin resistance to generate heat, and heat energy into deep skin, direct affect dermal tissue, stimulate collagen in the skin, as as to achieve long-term skin lifting and firming effect.

## Cool Hammer

Use the principle of air conditioner. Special wafer ceramics are used to absorb energy and refrigerate at - 5 degrees in an instant. The method of physical heat accumulation produced by wafer energy absorption makes one side produce heat up to 40 degrees. The fan discharges excess heat energy and achieves the effect of constant temperature and freezing alternately. It also promotes metabolism, alleviates allergic reaction, shrinks pore, calms skin, whitens and tenders skin.

Biological effects: During the cooling process, tissue temperature cools down, blood vessels and veins contract, metabolic rate reduces, cell activity inhibits, peripheral nerve sensitivity reduces, thus having a sedative and soothing effect, counteracting or reducing the symptoms of allergic reaction, such as redness, swelling, and pain, therefore achieving certain effects, such as calming the skin, reducing redness and swelling, shrinking pores and so on.

## Effects

1. Tighten loose skin.
2. Improve large pores and skin tones.
3. Elevate nasolabial fold.
4. Remove wrinkles and fine lines; improve acne scars.
5. Stimulate collagen regeneration.

## Applicable Range

1. Those with facial wrinkles, forehead wrinkles.
2. Those with fine lines and wrinkles around lips.

3. Those with loose, sagging, inelastic facial skin.
4. Those with large pores and facial dullness.
5. Those with undefined contours.

## **Inapplicable Range**

1. Those who have been implanted with a pacemaker or other similar electronic device.
2. Those who are pregnant.
3. Those who have been taking isotretinoin in the past 9 months (Accutane/Roaccutane/Oratane).
4. It is not recommended to do it for those who have filler inside one's face, such as gold wire implantation.
5. Those with severe heart disease, diabetes, and hyperthyroidism.
6. Those whose surgical wound is not healed.
7. Those with late-stage tumor.
8. Those who have unrealistic expectations of results.

## **Notes after Operation**

1. Do not wash your face with overheated water within three days (lukewarm or cold water can be used).
2. Strengthen moisturizing and sunscreen because skin might be dry after treatment.
3. Do not soak in hot springs, steam saunas or exercise strenuously, etc. within three days.
4. Use facial masks every day for a straight week.
5. Spicy, greasy food should be avoided. Staying up late, smoking and drinking are prohibited. Eat more vegetables, fruits.
6. Avoid food that will cause 'three highs' (high blood pressure, high cholesterol and high blood sugar), less greasy food is advised.

## **Anti-aging: Eyes**

### **Principles**

### **Fractional RF: Eyes (Probe A)**

Through transferring numerous cluster heat, this unique probe can send the matrix

molecular energy waves of up to 6 million times per second, repair skin in deep layers, from the focal facade of RF in hypodermis, make a strong impact on dermal tissues, which produces an aggregation effect of electric field in the dermal tissues, make the skin temperature up to 60 °C or so, produce immediately shrinkage, using the electromagnetic wave in the skin resistance to generate heat, and heat energy into deep skin, direct affect dermal tissue, stimulate collagen in the skin, as as to achieve long-term skin lifting and firming effect.

## Cool Hammer

Use the principle of air conditioner. Special wafer ceramics are used to absorb energy and refrigerate at - 5 degrees in an instant. The method of physical heat accumulation produced by wafer energy absorption makes one side produce heat up to 40 degrees. The fan discharges excess heat energy and achieves the effect of constant temperature and freezing alternately. It also promotes metabolism, alleviates allergic reaction, shrinks pore, calms skin, whitens and tenders skin.

Biological effects: During the cooling process, tissue temperature cools down, blood vessels and veins contract, metabolic rate reduces, cell activity inhibits, peripheral nerve sensitivity reduces, thus having a sedative and soothing effect, counteracting or reducing the symptoms of allergic reaction, such as redness, swelling, and pain, therefore achieving certain effects, such as calming the skin, reducing redness and swelling, shrinking pores and so on.

## Effects

1. Remove wrinkles around the eyes, crow's feet and eye bags.
2. Tighten loose skin around the eyes.
3. Speed up metabolism and prevent hyperpigmentation.
4. Firm and delicate skin and lift the canthus.
5. Improves eye blood circulation.

## Applicable Range

1. Those with wrinkles, fine lines, eye bags , and dark circles.
2. Those with dry lines and dynamic wrinkles caused by skin dehydration.
3. Those with drooping canthus and eyelids.
4. Those with loose skin around eyes.

5. Those with dynamic fine lines around eyes.

## **Inapplicable Range**

1. Those who have been implanted with a pacemaker or other similar electronic device.
2. Those who are pregnant.
3. Those who have been taking isotretinoin in the past 9 months (Accutane/Roaccutane/Oratane).
4. It is not recommended to do it for those who have filler inside one's face, such as gold wire implantation.
5. Those with severe heart disease, diabetes, and hyperthyroidism.
6. Those whose surgical wound is not healed.
7. Those with late-stage tumor.
8. Those who have unrealistic expectations of results.

## **Notes after Operation**

1. Do not wash your face with overheated water within three days (lukewarm or cold water can be used).
2. Strengthen moisturizing and sunscreen because skin might be dry after treatment.
3. Do not soak in hot springs, steam saunas or exercise strenuously, etc. within three days.
4. Use facial masks every day for a straight week.
5. Spicy, greasy food should be avoided. Staying up late, smoking and drinking are prohibited. Eat more vegetables, fruits.
6. Avoid food that will cause 'three highs' (high blood pressure, high cholesterol and high blood sugar), less greasy food is advised.

## **Anti-aging: Neck**

### **Principles**

### **Fractional RF: Neck (Probe C/D)**

Through transferring numerous cluster heat, this unique probe can send the matrix molecular energy waves of up to 6 million times per second, repair skin in deep layers,



from the focal facade of RF in hypodermis, make a strong impact on dermal tissues, which produces an aggregation effect of electric field in the dermal tissues, make the skin temperature up to 60 °C or so, produce immediately shrinkage, using the electromagnetic wave in the skin resistance to generate heat, and heat energy into deep skin, direct affect dermal tissue, stimulate collagen in the skin, as as to achieve long-term skin lifting and firming effect.

## Cool Hammer

Use the principle of air conditioner. Special wafer ceramics are used to absorb energy and refrigerate at - 5 degrees in an instant. The method of physical heat accumulation produced by wafer energy absorption makes one side produce heat up to 40 degrees. The fan discharges excess heat energy and achieves the effect of constant temperature and freezing alternately. It also promotes metabolism, alleviates allergic reaction, shrinks pore, calms skin, whitens and tenders skin.

Biological effects: During the cooling process, tissue temperature cools down, blood vessels and veins contract, metabolic rate reduces, cell activity inhibits, peripheral nerve sensitivity reduces, thus having a sedative and soothing effect, counteracting or reducing the symptoms of allergic reaction, such as redness, swelling, and pain, therefore achieving certain effects, such as calming the skin, reducing redness and swelling, shrinking pores and so on.

## Effects

1. Improve fine lines and wrinkles around neck.
2. Tighten loose, rough and dull neck skin.
3. Improve your double chin.
4. Accelerate lymphatic detoxification and improve facial skin quality.

## Applicable Range

1. Those with fine lines or wrinkles on the neck.
2. Those with loose and inelastic neck skin.
3. Those with skin dullness.
4. Those who bow their heads frequently.

## Inapplicable Range

1. Those who has been implanted with a pacemaker or other similar electronic device.
2. Those who are pregnant.
3. Those who have been taking isotretinoin in the past 9 months (Accutane/Roaccutane/Oratane).
4. It is not recommended to do it for those who has filler inside one's face, such as gold wire implantation.
5. Those with severe heart disease, diabetes, and hyperthyroidism.
6. Those whose surgical wound is not healed.
7. Those with late-stage tumor.
8. Those who have unrealistic expectations of results.

## Notes after Operation

1. Do not wash your face with overheated water within three days (lukewarm or cold water can be used).
2. Strengthen moisturizing and sunscreen because skin might be dry after treatment.
3. Do not soak in hot springs, steam saunas or exercise strenuously, etc. within three days.
4. Use facial masks every day for a straight week.
5. Spicy, greasy food should be avoided. Staying up late, smoking and drinking are prohibited. Eat more vegetables, fruits.
6. Avoid food that will cause 'three highs' (high blood pressure, high cholesterol and high blood sugar), less greasy food is advised.

# Repair Post-Pregnancy Stretch Marks

## Principles

### Fractional RF: Body (Probe D/E)

Through transferring numerous cluster heat, this unique probe can send the matrix molecular energy waves of up to 6 million times per second, repair skin in deep layers, from the focal facade of RF in hypodermis, make a strong impact on dermal tissues, which produces an aggregation effect of electric field in the dermal tissues, make the skin temperature up to 60 °C or so, produce immediately shrinkage, using the electromagnetic wave in the skin resistance to generate heat, and heat energy into deep skin, direct affect dermal tissue, stimulate collagen in the skin, as as to achieve long-term skin lifting and firming effect.

## Cool Hammer

Use the principle of air conditioner. Special wafer ceramics are used to absorb energy and refrigerate at - 5 degrees in an instant. The method of physical heat accumulation produced by wafer energy absorption makes one side produce heat up to 40 degrees. The fan discharges excess heat energy and achieves the effect of constant temperature and freezing alternately. It also promotes metabolism, alleviates allergic reaction, shrinks pore, calms skin, whitens and tenders skin.

Biological effects: During the cooling process, tissue temperature cools down, blood vessels and veins contract, metabolic rate reduces, cell activity inhibits, peripheral nerve sensitivity reduces, thus having a sedative and soothing effect, counteracting or reducing the symptoms of allergic reaction, such as redness, swelling, and pain, therefore achieving certain effects, such as calming the skin, reducing redness and swelling, shrinking pores and so on.

## Effects

1. Repair mild post-pregnancy stretch marks.
2. Stimulate collagen regeneration in hypodermis.
3. Tighten loose, saggy skin.
4. Repair light stretch marks.

## Applicable Range

1. Those with loose abdominal skin after delivery.
2. Those with loose waist and abdomen; Those who wish to look better in clothes.
3. Those with stretch marks on the waist and abdomen after birth.
4. Those with stretch marks.

## Inapplicable Range

1. Those who have been implanted with a pacemaker or other similar electronic device.
2. Those who are pregnant.
3. Those who have been taking isotretinoin in the past 9 months (Accutane/Roaccutane/Oratane).

4. It is not recommended to do it for those who has filler inside one's face, such as gold wire implantation.
5. Those with severe heart disease, diabetes, and hyperthyroidism.
6. Those whose surgical wound is not healed.
7. Those with late-stage tumor.
8. Those who have unrealistic expectations of results.


## **Notes after Operation**

1. Do not wash your face with overheated water within three days (lukewarm or cold water can be used).
2. Strengthen moisturizing and sunscreen because skin might be dry after treatment.
3. Do not soak in hot springs, steam saunas or exercise strenuously, etc. within three days.
4. Use facial masks every day for a straight week.
5. Spicy, greasy food should be avoided. Staying up late, smoking and drinking are prohibited. Eat more vegetables, fruits.
6. Avoid food that will cause 'three highs'(high blood pressure, high cholesterol and high blood sugar), less greasy food is advised.


## **Part II**

### **1. Detailed Operation**

After checking that the instrument is connected properly, the following starting interface will appear when the power supply is turned on.



### 1.1 Function Selection



Time Selection



Time: Decrease



Time: Increase



Pause




Start

## 1.2 Detailed Operation: Fractional RF



Click [Thermage](#) to go to the interface below



Intensity Adjustment



Pulse Time Adjustment



Pulse Time: Increase




Pulse Time: Decrease




Probe A Display(It lights up when using probe A. )




Probe B Display(It lights up when using probe B. )



Probe C Display(It lights up when using probe C. )



Probe D Display(It lights up when using probe D. )




Probe E Display(It lights up when using probe E. )

### 1.3 Detailed Operation: Cool Hammer



Click [Cool hammer](#) to go to the interface below



Cool Hammer

Energy Intensity



Cool Hammer Energy Intensity: Decrease



Cool Hammer Energy Intensity: Increase



Intensity Level Display



## 2. Technical Parameters

Fractional RF Frequency: 10MHZ

Fractional RF Energy Level: 1-20J/CM

Fractional RF Head : 5 Types

Voltage: 100V-240V

RF Tips Pins: 16 pins, 49pins, 100pins, 169pins, 576pins

Fractional RF Size:  $\Phi$ 8mm,  $\Phi$ 20mm, 4.0cm<sup>2</sup>

Control System: 8 inch color touch screen

Thermal RF Types: Interval pulse, continuous pulse

Output Power: 350W

G.W: 15 KG

MEAS: 53\*48\*53 CM

## 3. Safety Precautions

**Please consult your doctor or professionals before using the device if you have the following symptoms:**

1. Those who has been implanted with a pacemaker or other similar electronic device.
2. Those who are pregnant.
3. Those who have been taking isotretinoin in the past 9 months (Accutane/Roaccutane/Oratane).
4. It is not recommended to do it for those who has filler inside one's face, such as gold wire implantation.
5. Those with severe heart disease, diabetes, and hyperthyroidism.
6. Those whose surgical wound is not healed.
7. Those with late-stage tumor.
8. Those who have unrealistic expectations of results.

## 4. Dos and Dont's

1. Use a grounding pin and make sure the power socket of the device is properly

connected before each operation.

2. Ensure that the voltage of the device is appropriate. If the voltage of the local power supply is unstable, we recommend that users increase to the matching power voltage.
3. In order to ensure the therapeutic effect and normal service life expectancy of the device, please use the specified accessories provided or recommended by the original manufacturer.
4. The device should not be placed in wet environment or near water, nor should they be directly exposed to sunlight.
5. Do not place the device near a strong heat source, because this may affect the life and normal use of the device.
6. Before treatment, please remove all metal objects from the body first, so as to avoid unexpected conditions and influence the curative effect.
7. Please restrain from using the device on eyes, thyroid, parathyroid, testicles, abdomen of pregnant women, and pacemakers.
8. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
9. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use, so as to ensure the safety of the electrical products.
10. Use device or train device operators in strict accordance with instructions in the user manual.
11. In the process of operation, the adjustment of energy and output time should not be in direct proportion. When the adjustment of energy is high, the output time should be low.
12. Before treatment, special marking paper should be used to precisely locate the marking, so as to ensure the balanced distribution of treatment energy.
13. Paint the temporary marking grid with 95% alcohol.
14. Make sure you have enough gel on your skin.
15. Start with the lowest energy at the beginning of operation, and gradually add up after adaptation.
16. When using this instrument, the operating area must be kept moist and avoid operation on dry skin.
17. After the operation, scrub the instrument with normal saline to ensure its cleanliness and extend its service life.
18. During operation, the contact surface of the operating head must fully touch the skin, and then press the button to avoid discomfort or scalding.

## 5. Troubleshooting & Solutions

### 1. The instrument can't start and the key lamp on the back of the instrument doesn't work?

- A. Ensure that the power cord is connected to a working power outlet.

B. Ensure the fuse tube on the back of the instrument is not loose or being burnt.

## **2. No energy output of the instrument?**

A. Please check whether the plug connecting the instrument handle and the body is tightly connected.

B. Please check whether the treatment area has been cleaned. Grease or oil essential oil products may cause poor contact between the treatment handle of the instrument and human body, resulting in no RF output.

## **3.The cool hammer is not cold?**

A. Please check whether the plug connecting the instrument and the handle is tightly connected.

B. Please check whether the continuous use time exceeds 20 minutes.

## **4.The device can be activated, but the wrong information appears on the screen?**

A. Please unplug the power plug behind the device and wait about 1 minute before restarting the device.

B. If the above method did not work, please contact the device distributor for assistance.

## **6. FAQs**

### **1. Q: Why does it take 6-10 times to tighten skin and remove wrinkles with regular means while only once is needed with Fractional RF?**

A: Fractional RF works through high-frequency energy wave conduction technology, stimulate hypodermal collagen proliferation. It's an high-tech non-surgical anti-aging instrument. Only one treatment can make loose skin tighten, bags, hairline become lighter, jaw line is clear, the overall skin is smooth and compact, facial contour is clear, the whole person more youthful and beautiful.

### **2. Q: How long can a treatment last?**

A: After one treatment, the skin will be firm and lifted.What's more, continuous improvement can be achieved in the following three months, and the effect can be maintained for about half a year. Results differ from person to person and the repair factor after the operation exists. Skin that absorbs well can maintain relatively longer.

### **3. Q: How do you feel during the treatment?**

A: step 1: when the instrument is turned on and the probe touches the skin, the skin will feel comfortable and cool. Step 2: then the heat goes deep into the skin. When it begins to release heat, the experiencer can feel a brief probe emitting heat in the deeper tissues of the skin. At this point, RF heat is already effectively stimulating

the collagen, causing it to instantly contract and gradually tighten the skin. Step 3: cool down again. When the RF heat energy decreases to the lowest point, the experimenter can feel the cool feeling of the skin again.

**4. Q: What's the difference between Fractional RF and traditional facelift?**

A: Compared with the traditional electric wave skin tightening surgery, on the probe of Fractional RF, the heat energy goes deep into the skin through the way of cluster heat transfer of countless dot arrays, thus the fascia layer forms a dot matrix of collagen adhesion, making the tightening skin firm and reposition on the collagen scaffold, so as to achieve the purpose of lifting the face.


**7. Packing List**

- 1. Main Body x 1
- 2. Bracket x 1
- 3. Fuse x 2
- 4. Handle x 1
- 5. Cool Hammer x 1
- 6. Probe x 5
- 7. Power Cord x 1

**8. Operational Diagrams**

Parameter Adjustment	Product	Techniques	Diagrams
<b>Anti-Aging on Face: 40 min, once a week</b>			
Fractional RF Energy Intensity: 5-15  Fractional RF	Makeup remover oil+ cleanser+ facial marking	1. Remove makeup and clean face, 5 min. 2. Place marking paper on face and apply alcohol on paper, 2-3 sec. 3. Remove marking paper.	No Diagrams

Pulse Time: 0.5-2.0  Cool Hammer Energy Intensity: 30%-80%	paper+ gel +facial mask+95% alcohol	4. Apply gel evenly on one side (thickness:coin), 1 min. 5. Device operation: install face probe on handle and operate on marking cubes accordingly, 3-5 times. (Make sure the device is close to the skin.) 6. Apply gel evenly on another side (thickness:coin), 1 min. 7. Device operation: install face probe on handle and operate on marking cubes accordingly, 3-5 times. (Make sure the device is close to the skin.) 8. Do the same on forehead, 8 minutes. 9. Scrape off the rest gel with mask spatula, 2 min. 10. Wash face clean, 2 min. 11. Apply restorative face mask+cool hammer operation; cool hammer operation: move the probe in circles and lift the face in lines from chin to forehead; cover the whole face, 15 min. 12. Wash face clean, 2 min. 13. Apply skincare products(toner, essence, cream, sunscreen)	
<p><b>Suggested Treatment :</b>  5 treatments account for 1 full course. After one course of treatment, you can feel the heat in the bottom of the skin and feel the facial improvement. Also,the skin is significantly firmer and lifted, and the skin color gradually turns ruddy. Two courses of treatment delay the aging of the skin and restore the skin to its youthful state.</p>			
<p><b>Anti-Aging on Eyes: 40 min, once a week</b></p>			
Fractional RF	Makeup	1. Remove makeup and	

<p>Energy Intensity: 5-15</p> <p>Fractional RF Pulse Time: 0.5-2.0</p> <p>Cool Hammer Energy Intensity: 30%-80%</p>	<p>remover oil+ cleanser+ facial marking paper+ gel +eye mask+95% alcohol</p>	<p>clean face, 5 min.</p> <ol style="list-style-type: none"> <li>2. Place marking paper on lower eye area and apply alcohol on paper, 2-3 sec.</li> <li>3. Remove marking paper.</li> <li>4. Apply gel evenly on one side (thickness:coin), 1 min.</li> <li>5. Device operation: install eye probe on handle and operate on marking cubes accordingly, 3-5 times. (Make sure the device is close to the skin.)</li> <li>6. Do the same on another side, 5 min.</li> <li>7. Scrape off the rest gel with mask spatula, 2 min.</li> <li>8. Wash eyes clean, 2 min.</li> <li>9. Cool Hammer Operation: move slowly from inner canthus to brow ridge; start with one side and finish with another, 10 min.</li> <li>10. Apply eye mask, 15 min.</li> <li>11. Wash eyes clean, 2 min.</li> <li>12. Apply skincare products(toner, essence, eye cream, sunscreen)</li> </ol>	<p>Technique 9</p> 
<p><b>Suggested Treatment :</b></p> <p>5 treatments account for 1 full course. After one treatment, it tighten the skin of eyes and increase blood circulation to the eyes. After one course of treatment, fine lines become lighter and wrinkles are reduced. After two courses of treatment, eye skin becomes firmer, tender and shiny, improving and preventing eye aging.</p>			
<p><b>Anti-Aging on Neck: 40 min, once a week</b></p>			
<p>Fractional RF Energy Intensity: 5-15</p> <p>Fractional RF Pulse Time:</p>	<p>Makeup remover oil+ cleanser+ facial marking paper+ gel</p>	<ol style="list-style-type: none"> <li>1. Remove makeup and clean neck, 5 min.</li> <li>2. Place marking paper on neck and apply alcohol on paper, 2-3 sec.</li> <li>3. Remove marking paper.</li> <li>4. Apply gel evenly on one</li> </ol>	<p>Technique 7</p>




<p>0.5-2.0</p> <p>Cool Hammer Energy Intensity: 30%-80%</p>	<p>+neck mask+95% alcohol</p>	<p>side (thickness:coin), 1 min.</p> <p>5. Device operation: install neck probe on handle and operate on marking cubes accordingly, 3-5 times. (Make sure the device is close to the skin and avoid throat.)</p> <p>6. Scrape off the rest gel with mask spatula, 2 min.</p> <p>7. Wash neck clean, 2 min.</p> <p>8. Cool Hammer Operation: move from two sides of throat in circles; cover the whole neck, 15 min.</p> <p>9. Apply neck mask, 15 min.</p> <p>10. Wash neck clean, 2 min.</p> <p>11. Apply skincare products(toner, essence, neck cream, sunscreen).</p>	
<p><b>Suggested Treatment :</b></p> <p>5 treatments account for 1 full course. After one treatment, you can feel the heat at the bottom of the skin, and the skin color changes. After one treatment, the neck skin gradually becomes tight and the neck wrinkles become light. After two treatment, the skin becomes tight and plump, delicate and shiny.</p>			
<p><b>Repair Post-Pregnancy Stretch Marks: 40 min, once a week</b></p>			
<p>Fractional RF Energy Intensity: 5-15</p> <p>Fractional RF Pulse Time: 0.5-2.0</p> <p>Cool Hammer Energy Intensity: 30%-80%</p>	<p>Cleanser+ facial marking paper+ gel +restorative mask+95% alcohol</p>	<p>1. Clean waist and abdomen with cleanser, 5 min.</p> <p>2. Place marking paper on abdomen and apply alcohol on paper, 2-3 sec.</p> <p>3. Remove marking paper.</p> <p>4. Apply gel evenly on one side (thickness:coin), 1 min.</p> <p>5. Device operation: install abdomen probe on handle and operate on marking cubes accordingly, 3-5 times. (Make sure the device is close to the skin and avoid</p>	<p>No Diagrams</p>

		<p>navel.) Add gel during the process if needed.</p> <p>6. Scrape off the rest gel with mask spatula, 2 min.</p> <p>7. Wash the treatment area clean, 2 min.</p> <p>8. Cool Hammer Operation: move from left to right in circles; cover the whole abdomen, 15 min.</p> <p>9. Apply restorative mask, 15 min.</p> <p>10. Remove mask and wash abdomen clean, 2 min.</p>	
--	--	--	--

**Suggested Treatment :**

5 treatments account for 1 full course. After a course of treatment you will feel heat on and tightened skin, which will be smooth and delicate with the help of restorative facial mask. It also help prevents loose and sagging skin. After two courses of treatment, stretch marks will be diluted, fine lines and laxity will be improved. After three courses of treatment, stretch marks will be repaired to tighten and smooth skin.

**Anti-Aging: The Whole Body: 70 min, once a week**

<p>Fractional RF Energy Intensity: 5-15</p> <p>Fractional RF Pulse Time: 0.5-2.0</p> <p>Cool Hammer Energy Intensity: 30%-80%</p>	<p>Cleanser+ facial marking paper+ gel +restorative mask+95% alcohol</p>	<p>1. Remove makeup and wash face, eyes and neck.</p> <p>2. Place marking paper on face and eyes and apply alcohol on paper, 2-3 sec.</p> <p>3. Remove face and eye marking paper.</p>	<p>Technique 14.15</p> 
		<p>4. Place marking paper on neck (avoid throat and aorta) and apply alcohol on paper, 2-3 sec.</p> <p>5. Remove marking paper.</p>	<p>Technique 16</p> 
		<p>6. Apply gel evenly on one side of face, eyes and neck(thickness:coin), 2 min.</p> <p>7. Device operation: install face probe on handle and operate on marking cubes of neck accordingly, 3-5 times. (Make sure the device is</p>	<p>Technique 27</p> 



	<p>close to the skin.)</p> <p>8. Operation on Face: operate on marking cubes accordingly, 3-5 times.</p> <p>9. Do the same on forehead, 8 min.</p> <p>10. Repeat the operation on another side.</p> <p>11. Switch to eye probe and operate on eyes.</p> <p>12. Operate on marking cubes accordingly, 3-5 times.</p> <p>13. Repeat the operation on another side.</p> <p>14. Cool Hammer Operation: move in circles on neck, 3 times.</p> <p>15. Move from neck to face in circular and linear lines, 3 times.</p> <p>16. Move horizontally in circles on forehead, 3 times.</p> <p>17. Scrape off the rest gel with mask spatula, 2 min.</p> <p>18. Wash face, eyes and neck clean, 2 min.</p> <p>19. Apply face and neck mask, 15 min.</p> <p>20. Wash face and neck clean, 2 min.</p> <p>21. Apply skincare products(toner, essence, eye cream, sunscreen)</p> <p>22. Waist and Abdomen: Wash it with cleanser, 5 min.</p> <p>23. Place marking paper on abdomen and apply alcohol on paper, 2-3 sec.</p> <p>24. Remove marking paper.</p> <p>25. Apply gel evenly on one side (thickness:coin), 2 min.</p> <p>26. Device operation: install abdomen probe on handle and operate on marking</p>	 <p>The diagram shows a female torso from the waist up. On the abdomen, there are several red dots and small circles indicating marking points. On the forehead, there are also red dots and small circles. The diagram is positioned above a horizontal line.</p>
--	---	---

		<p>cubes accordingly, 3-5 times. (Make sure the device is close to the skin and avoid navel.) Add gel during the process if needed.</p> <p>27. Cool Hammer Operation: move from left to right in circles; cover the whole abdomen, 15 min.</p> <p>28. Scrape off the rest gel with mask spatula, 2 min.</p> <p>29. Wash the treatment area clean, 2 min.</p> <p>30. Apply restorative mask, 15 min.</p> <p>31. Remove mask and wash abdomen clean, 2 min.</p>	
<p><b>Suggested Treatment :</b></p> <p>5 treatments account for 1 full course. After one course, the treatment area will experience heat and tightening feeling. The effect of Skin firming and lifting become gradually obvious. After two courses, the collagen in the skin layer increases, the abdominal stretch marks fade, and the skin is firm and elastic.</p>			