

# User Manual

SR-MJ1809P



## Preface

To our valued users,

Thanks for choosing our latest integrated beauty equipment.

This is the latest skin cleansing and rejuvenating equipment that combines together skin scrubber, high frequency, electric face washing brush, cold & hot hammer, roller ion, facial ultrasound, eye ultrasound, diamond dermabrasion and oxygen injector. SR-MJ1809P is a professional beauty equipment, which needs to be operated by people with professional training. Any improper use of this equipment may bring adverse consequences to human body, thus we advise all the people to read this manual thoroughly and follow the instructions specified hereby strictly before operation.

We believe our quality product will give you a good return, and our perfect after-sales service will make you worry-free.

Thank you.

# Table of Contents

Preface
Part I
Brief Introduction
Advantages
Working Principle
Classification of Items
Facial Deep Cleansing
Facial Deep Infusion
Facial Skin Cleansing & Acne Removing
Eye Tightening & Anti-aging
Neck Refining & Anti-aging
Contraindications
Matters Needing Attention After Treatment
Part II
Detailed Operations
Technical Parameters
Contraindications
Precautions for Use
Troubleshooting & Solutions
FAQs
Packing List
Operational Diagrams

# Part I

## Brief Introduction

SR-MJ1809P is an integrated beauty equipment that applies high frequency, cold & hot effect, positive and negative ion, ultrasound and vacuum. All these technologies are currently the most popular and widely used ones for beauty and anti-aging. It can effectively replace most of the facial care programs in beauty salon, and can be operated easily, conveniently and swiftly, and solves skin problems for people pursuing beauty. Moreover, this equipment works safely and effectively, and requires no skin breaking, no injection, no medication and no operation. It has no side effects, and is operated externally throughout, and has an instant effect.

## Advantages

1. 9-in-1 multifunctional beauty equipment that rapidly fixes facial care problems.
2. It has multiple cleansing probes that effectively fix facial clogged pores problem.
3. It has multiple operating handles, which can be changed according to different demands.
4. Ozone with high frequency brings good news to people with acne.
5. It suits all kinds of skin, and can improve skin quality and flatten fine lines.
6. It's comfortable, painless, noninvasive and injury-free during treatment, and there is no convalescence, people can leave immediately, and without affecting people's normal work and life after treatment.
7. No consumption, low cost, but with quick return.
8. It has a wider treatment range, and has an instant and evident effect.
9. Unevenness, bleeding and swelling will not occur.

## Working Principle

### Ultrasonic Skin Scrubber

It applies the principle of ultrasound with high frequency ultrasound. Energy produced by mechanical vibration of 28000 times per second has an obvious permeation effect on skin. When ultrasonic cavitation effect encounters water, it atomizes water into small bubbles. The small bubbles, through constant vibration, go deep into skin's pores breaking acne, blackheads and dark acne into granulum which will be vibrated onto skin's surface through high frequency and ultrasonic vibration. Low frequency emulsifies sebum and excreta that clogged pores, and meanwhile it can also do sterilization, atomizes water molecules, removes redundant thick cutin, purifies pores and decomposes grease through mechanical vibration after entering into pores, so as to make the pores unobstructed, contract the pores, and make the skin bright, clear, soft and smooth.

**Physical Effect:** Sound wave of high frequency vibration, with the matched use of softened products or exfoliating cream, can soften the impaired and rough stratum corneum, open pores, and deeply remove the wastes from pores. Ultrasonic Skin Scrubber breaks nutrient into molecular particles with high frequency ultrasonic vibration of 3MHz, which can effectively enhance skin absorption, infuse nutrition deeply into dermis, and make skin shiny and elastic.

## Transparent Tube Head with High Frequency

Electro-therapeutic apparatus with high frequency contacts skin through high frequency micro-current generated by glass electrode tube and after having the power on. Micro-current in glass electrode tube ionizes oxygen in the air producing ozone, which will improve blood circulation and exert physical therapy on the skin.

## Electric Face Washing Brush

The face washing brush adopts skin penetration technology with ultrasonic pulse. It acts directly on tiny brush by using ultrasonic wave, making it vibrate and exerting the effect on tiny pores on face so as to impose high frequency vibration on pores and achieve the effects of pores cleansing and facial muscles massaging.

## Cold & Hot Hammer

It uses the principle of expansion caused by heat and contraction caused by cold. When the temperature of our skin stays at 5 degrees, metabolic rate reduces, cell activity inhibits and slows down, and sensitivity of nerve endings lowers, which plays a role in inhibiting pain, swelling, and ache caused by skin allergy. Due to the same principle, a temperature of 5 degrees contracts blood vessels, and causes pores shrinkage which lessens large pores and firms skin, thus flattening wrinkles and crow's-feet on face.

Physiological functions and effects of cold treatment with a temperature of 5 degrees:

1. It shrinks blood vessels, contracts capillary pores, and tremendously minify large capillary pores.
2. It strengthens collagenous fiber, makes skin tightened and elastic, and flattens wrinkles and crow's-feet.
3. It contracts blood vessels, effectively alleviates varicose veins and spider veins bulging. Moreover, harmful wastes that backflow to veins at any time can be easily purified and detoxed. It also has an obvious effect on spots fading and removing.
4. It reduces metabolic rate, inhibits and slows down cell activity, lowers the sensitivity of nerve endings with a sedative and soothing effect, offsets skin pain, swelling and ache, and it, especially, has a good curative effect for allergic skin.

Physiological functions and effects of hot treatment with a temperature of 40 degrees:

1. Tepidity makes people feel warm and comfortable, and relax muscles, balances autonomic nerve, and dispel pain and soreness.
2. It accelerates blood circulation, increases oxygen content in blood, opens pores,
3. It boosts skin's metabolic rate, increases heat energy in tissue, discharges wastes and toxin naturally, and it has effect on eye bags and black eyes removal.

## EMS Roller

The full name for EMS is Electrical Muscle Stimulation, which is called Electric Wave Resurfacing by beauty salon. The so called EMS micro-current refers to a low level electric current, which is used to transmit a mild and gentle electric wave. And this electric wave penetrates epidermis reaching directly to facial dermis and stimulates and massages facial muscle. It's prove that micro-current can stimulate the production of ATP, which is the key material for generating collagen and elastin. And this collagen and elastin make facial contour clearer, thus achieving the effect of lifting, tightening and face-lifting. Meanwhile, EMS micro-current can also fully exercise local fat layer, improve body's basic metabolism, and achieve the effect of fat and proud fresh removing.

## Facial & Eye Ultrasound

By utilizing the principle of high frequency vibration, ultrasonic beauty equipment, with a high frequency vibration wave of 28,000 Hz per second, can powerfully reach 4 to 6 millimeter of subcutaneous tissue. Ultrasound is of powerful energy and strong permeability and can cause cells to vibrate with it, which can not only remove grease from pores but also change cell's volume so as to improve blood circulation and achieve the purposes of stimulating cells' regeneration capacity and making skin glossy and elastic. In addition, tepidity effect generated by high frequency vibration also has the function of rapidly importing beauty nutrient solution into skin, which can make skin absorb nutrient content from skin care products better.

## Diamond Dermabrasion

Embedded with tiny diamond particles of different sizes, the specialized diamond toned wand, matched with vacuum suction, grinds off a shallow layer of epidermis by rubbing back and forth. It's a noninvasive and physical dermabrasion. It removes dead skin cells of epidermis to regenerate them and to improve skin quality. After treatment, skin's fine lines lessen, and skin becomes smooth and flat, and pores turn more delicate. Meanwhile, matched with hands techniques, it can achieve the effect of lifting.

**Physical Effects:** Diamond dermabrasion adopts vacuum suction technology. Through decreasing internal air pressure, it effectively removes old and useless stratum

corneum and sucks out wastes in pores like blackheads and excessive grease, thus achieving the effect of free breath and circulation for skin, dredging pores, improving large pores, and brightening skin color.

Dermabrasion adopts diamond crystals that have been processed specially to ill skin, and use vacuum suction to remove necrotic cells from skin as well as remove the deep waste and dirty of pore, help skin double cleansing and promote growth of new skin.



## Spray Bottle

Extracting oxygen with a concentration of only 21% of the air by using the principle of low pressure jetting. The oxygen, passing through a molecular filter made with special technology, will be compressed into oxygen active factor with a purity over 95% which then will be released under a pressure of 0.2 Pa. It turns aqueous nutrients into nebulous oxygen which will be transmitted to skin stratum basale, which quickly opens the channel for skin absorption, kills anaerobion in the skin, activates cells, metabolizes old and useless cutin, deeply and densely replenishes moisture, shrinks pores, and clean the deep skin.

## Classification of Items

### Facial Deep Cleansing

#### Matched Probe

Electric Brush, Ultrasonic Skin Scrubber, Diamond Dermabrasion, Spray Bottle, EMS Roller and Cold & Hot Hammer

### Effects

1. Remove old and useless stratum corneum.
2. Cleanse deep skin.
3. Dredge pores.

4. Balance the secretion of grease.
5. Lighten skin color and refine pores.
6. Improve large pores and fade fine lines.
7. Strongly remove strawberry-like nose, acne and blackheads.

## Indications

1. People with clogged pores, or relatively thick stratum corneum.
2. People with dull or lackluster skin.
3. People with coarse skin or large pores.
4. People with strong grease secretion.
5. People who are prone to have pimples, closed comedones, acne, or with inflammation on underlying skin.
6. People with thin or impaired stratum corneum.
7. People with poor skin absorption, fine lines, wrinkles or color spots.

## Facial Deep Infusion

### Matched Probe

Electric Brush, EMS Roller, Facial Ultrasound, Spray Bottle and Cold & Hot Hammer

### Effects

1. Tighten face and infuse nutrition into deep skin.
2. Increase skin's elasticity and make skin delicate.
3. Moisten skin and increase skin's absorptivity.
4. Accelerate blood circulation and metabolism.
5. Stimulate collagen hyperplasia and delay aging.

## Indications

1. People with dark or lacklustre face.
2. People with slack or saggy skin.
3. People with fine lines, nasolabial folds or periorbital lines.
4. People with vague facial contour.
5. People with coarse skin or large pores.
6. People who are under long-term exposure of UV radiation in workplace.
7. People with slack skin, edema or puffiness after child delivery.

## Facial Skin Cleansing & Acne Removing

### Matched Probe



Electric Brush, High Frequency Handle, EMS Roller, Spray Bottle and Cold & Hot Hammer. (To avoid cross-infection, it's not advised to use skin scrubber and diamond dermabrasion on acne skin.)

## Effects

1. Diminish inflammation and do disinfection for deep skin.
2. Remove acne, fade marks, rejuvenate and brighten skin.
3. Improve coarse skin and large pores.
4. Remove acne and repair acne skin.
5. Do disinfection with ozone, remove acne and mites.
6. Accelerate blood circulation and increase skin metabolism.
7. Lighten skin color and refine pores.

## Indications

1. People with strong grease secretion, or large pores.
2. People with acne skin, or who always have acne.
3. People with extensive acne on the skin.
4. People who have dark acne after having pimples.
5. People with clogged pores, or accumulated mites.

## Eye Tightening & Anti-aging

### Matched Probe

Electric Brush, Eye Ultrasound and Cold & Hot Hammer

## Effects

1. Relieve eye fatigue, black eyes, eye bags and edema.
2. Improve periorbital lines and crow's feet.
3. Fade pigment and accelerate blood circulation.
4. Accelerate metabolism and prevent pigment from sediment.
5. Replenish nutrition for eyes and moisten skin.
6. Tighten skin and make it tender, and lift the corner of eyes.
7. Accelerate blood circulation of the eyes, and help with the efficient absorption.

## Indications

1. People with wrinkles, fine lines, eye bags or black eyes.
2. People with dry skin, dry lines or false wrinkles.
3. People whose eyes are prone to fatigues or dryness. People
4. People who always face computer or cellphone.

5. People who always stay up late or with black eyes.
6. People who always exposure to dry or high temperature environment.

## Neck Refining & Anti-aging

### Matched Probe

Electric Brush and Facial Ultrasound(For neck, it's not advised to use cold & hot hammer, EMS Roller and Spray bottle.)

### Effects

1. Improve neck with slack, coarse or lacklustre skin.
2. Tighten skin and increase skin's elasticity.
3. Relieve double chin.
4. Accelerate lymphatic detox and improve facial skin.
5. Prevent neck and lymphatic diseases.

### Indications

1. People whose neck have coarse skin or dim skin color.
2. People whose neck have clogged lymph.
3. People whose neck have slack or inelastic skin.
4. People who always lower their heads.
5. People who are not satisfied with their necks' skin color.

## Contraindications

1. People who just had plastic surgery, or had prostheses or metal objects implanted inside the treated parts.
2. People who recently had injections, such as hyaluronic acid, mesotherapy, wrinkles removal, or plastic surgery.
3. People with hypertension, heart disease, diabetes, severe thyroid, malignant tumors, etc.
4. People in allergic period, with severely sensitive skin, or who are allergic to metals.
5. People with skin trauma or wound.
6. People who are overly aging.
7. People in menstruation, pregnancy, lactation or surgical recovery.
8. People with skin diseases or infectious diseases.
9. People who have unrealistic illusion about the effect.

## Matters Needing Attention After Treatment

1. Do not wash the face with overheated water within three days(but can wash with

warm and cold water).

2. Keep hydrating, and protect yourself from the sun.
3. Do not enjoy hot springs, saunas, or do strenuous exercise, etc., within three days.
4. It's advised to apply facial mask at least three times in a week.
5. Avoid eating spicy and greasy food, staying up late, smoking and drinking alcohol. Instead, eat more vegetables, fruits and light food.
6. Do not eat food causing Three Highs(hypertension, hyperlipidemia and hyperglycemia), and mostly eat light food.
7. Don't use irritant skin care products that contain scrub, exfoliant, AHA(Alpha Hydroxyl Acid), etc.

# Part II

## 1. Detailed Operations

Time, energy level and mode can be set according to individual skin conditions. There is no fixed values, and it varies from region to region, and individual tolerance is also different. Normally, energy level is between 3 and 7, while time and mode can be adjusted in line with individual habits. Information provided in this user manual is for reference only. Set the parameters in accordance with actual conditions.

Check the equipment and make sure it's properly connected. The following initial page will appear after having the power switch on.



### 1.1 Functions Selection



Time Setting & Display/ Time Range: 0 to 60 minutes



Time Setting(decrease)



Time Setting(increase)



Energy Level Setting



Energy Level Setting(decrease)



Energy Level Setting(increase)



Energy Level Display/ Energy Level Range: 1 to 9



Mode Switching



Mode A



Mode B



Mode C



Start/ Pause

## 1.2 Detailed Operations for Cold & Hot Hammer



Select Cold & Hot Hammer icon(When the icon is color red, then it's selected.).



Time Setting



Mode Switching



There are two modes available. A is cold hammer and B is hot hammer.



Energy Level Setting. When in cold hammer, the higher the energy level,

the colder the hammer is. When in hot hammer, the higher the energy level, the hotter the hammer is.



To operate after clicking Start.

### 1.3 Detailed Operations for EMS Roller



Select EMS icon.



Time Setting



Mode Switching



There are three modes available. A is micro-current fluctuation with relatively slower frequency. B is micro-current fluctuation with relatively faster frequency. C is micro-current fluctuation with fast frequency.



Energy Level Setting



To operate after clicking Start.

### 1.4 Detailed Operation for Ultrasound(This icon is for starting both face and eye ultrasound.)



Select Ultrasound icon.



Time Setting



Mode Switching



There are two modes available.

Mode A is continuous working, namely, ultrasound works all the time after having the equipment on, suitable for people with relatively dry skin, poor skin absorption, or who can bear with the strong energy output.

Mode B is discontinuous working, namely, ultrasound works discontinuously after having the equipment on, suitable for people with normal skin condition, poor skin absorption, or first-time user.



Energy Level Setting



To operate after clicking Start.

## 1.5 Detailed Operations for Diamond Dermabrasion & Spray Bottle

### 1) Diamond Dermabrasion



Select Diamond Dermabrasion/ Spray Bottle icon.



Time Setting



Rotating Knob for vacuum suction setting. Rotate to the left to decrease suction and to the right to increase.(advises for suction: 20% to 80%, which can be adjusted according to actual situation. Don't set the suction too high, otherwise it may cause redness.)



This is air inlet(air inlet), insert one end of rubber hose into air inlet and the other end into handle.



Black Rubber Ring (This rubber ring must be installed when operating, otherwise there is no suction or suction will be weakened.)



Filter Cotton(It must be installed when operating, otherwise dirt may enter into the equipment causing damages to it. Meanwhile, the cleaning effect can be seen.)



Unscrew the lid from the handle tail, encase it into white rubber hose, then insert the handle into the rubber hose, and last screw the lid up.



Choose the proper operating head according to the size of the skin.



Big diamond heads are suitable for treating part with relatively large size, like the periphery of eyes, the two sides of nose, cheek, facial wrinkles removal and body parts.



Round heads are suitable for cleansing forehead, chin, nose wing and cheek, and facial wrinkles removal.



Sharp ones are suitable for nasal fossa, obstinate blackheads, mites, clogged pores and facial wrinkles removal.



To operate after clicking Start.

## 2) Spray Bottle



Select Diamond Dermabrasion/ Spray Bottle icon.



This is air outlet(across from the equipment and on the left side is air outlet), insert one end of rubber hose into air outlet and the other end into the spray





bottle as shown by red arrow .

Open the black lid, and pour essence, toner, pure water and normal saline into the bottle.


Essence/ Toner: can deeply replenish moisture and brighten skin color.

Pure Water: can deeply replenish moisture and do deep cleansing.

Normal Saline: suits people with acne skin and sensitive skin.



To operate after clicking Start.



When operating, the part shown by red arrow must be sealed with index finger to prevent water from coming out, otherwise no water come out after having the equipment on.

## 1.6 Electric Face Washing Brush



Switch is on the handle. Press the switch on the handle to operate.

## 1.7 Skin Scrubber



Switch is on the handle. Press the switch as shown by the red arrow to operate.

## 1.8 Transparent Tube Head with High Frequency



Switch is on the handle. Twist the part of connecting wire as shown by the red arrow to operate.

## 2. Technical Parameters

Power Supply Input: 100V-240V  
Power: 20W 190W

Ultrasonic Handles for Face  
Frequency: 1MHz 3MHZ  
Power: 20W 50W

Cold Therapy Hammer  
Temperature: 0-5 Degree Centigrade  
Power: 20W 10W

Ultrasonic Shape Head  
Frequency: 1MHz 3MHZ  
Power: 20W 50W

Double Deep Cleanse-dermabrasion  
Total power: 14W  
Voltage: 12V  
Diameter of the Micro-carving Head: 22mm, 20mm, 16mm, 16mm, 15mm, 15mm,  
15mm, 08mm, 06mm  
Length of Micro-carved Rod: 110mm  
Suction Range: 0 to 80Kpa

Ion Probe:  
Power: 20W 5W

High Frequency  
Frequency: 50HZ

## 3. Contraindications

People with the following situations should use the equipment with caution. Please consult a doctor or professional before using this equipment. The details are as follow:

- (1) Women in pregnancy, or lactation.
- (2) People with heart disease, or heart pacemaker;

- (3) People with unhealed surgical wound, or in surgical recovery;
- (4) People with epilepsy, severe diabetes or hyperthyroidism.
- (5) People with malignant tumors, hemophiliac or severe bleeding;
- (6) People with skin diseases or infectious diseases should use with caution.
- (7) People with metal implants inside body such as stent, pacemaker, etc.
- (8) People who are sensitive to metals and electric current should use with caution.

#### 4. Precautions for Use

- (1) Plug with ground pin must be used, and power socket that already factually grounded must be ensured before using the equipment.
- (2) To make sure equipment's voltage is adaptive. If voltage of local power supply is unstable, we suggest users add voltage regulator with matching power between the mains supply and the equipment.
- (3) To guarantee curative effect and normal service life of the equipment, please uniformly use specified parts provided or suggested by original manufacturer.
- (4) The equipment can't be placed in damp places or near water, and also can't be exposed to direct sunlight.
- (5) Don't place the equipment near a strong heat source since this may affect its service life and normal use.
- (6) Please remove all the metal objects from body before treatment so as to avoid unexpected situations which may affect curative effect.
- (7) Please restrain from using the equipment aiming at eyes, thyroid, parathyroid, testicles, pregnant woman's abdomen, pacemaker, etc.
- (8) People who are suffering from illness should use with caution unless he gets permission from a doctor.
- (9) Please turn off power switch of the equipment if no one uses it, and ensure main power is off after someone used it and before he left, thus safety of the equipment can be guaranteed.
- (10) Start with the lowest energy level when just started, and increase the level gradually after client got used to it.
- (11) It's prohibited to use it on the rim of eyes, ears, the periphery of parts that have skin diseases, and parts that have plastic surgery.
- (12) Please set vacuum suction in line with client's comfort level when operating diamond dermabrasion.
- (13) Rubber ring must be mounted when operating diamond dermabrasion, otherwise there is no suction or suction will be weakened.
- (14) Filter cotton must be set when operating diamond dermabrasion, otherwise there is dirt after having skin filtered.
- (15) Please maintain the wetness of skin, and avoid dry skin operation, which make people feel more comfortable.
- (16) Clean the equipment with normal saline after operation so as to ensure its cleanliness and hygiene, as well as extending its service life..
- (17) Using the equipment and training the operators in strictly accordance with the

instructions specified in the manual.

## 5. Troubleshooting & Solutions

(1) The equipment can't be started?

- A. Make sure the power cord is plugged into a good power socket.
- B. Check whether the power cord in the back is fully inserted.

(2) The equipment can be started, but there is error message on the monitor?

- A. Take out the plug in the back, and wait for about 1 minute and re-plug and restart it.
- B. If the above-mentioned solutions are not working, please contact the equipment dealer for assistance.

(3) Some handle is not working?

- A. Check whether the handle's switch or the equipment's button is started.
- B. If the connection is good but it's still not working, probably the handle get damaged, then you need to contact after-services for repair or re-order.

(4) High frequency handle or EMS roller doesn't have output after starting?

- A. Check whether the handle is properly installed. If no, install it properly and then start to operate.
- B. If there is back-up operating head, replace a new one and set the energy level to the highest to have a test.

(5) The diamond dermabrasion doesn't have suction or with small suction?

- A. Please check whether the rubber ring of the diamond head is twisted too tight that causes air leakage.
- B. Please check whether the rubber ring is installed.
- C. Please check whether the rubber ring is properly connected to main machine, or whether there is damage to the rubber ring.

## 6. FAQs

(1) Q: Why the cutin must be removed?

A: Stratum corneum is the top layer of epidermis, on which cells are completely dead. Depending on body parts, the stratum corneum is consists of 15 to 20 layers of deal cells. Human's skin has a metabolic cycle of 28 days, which is exactly the cycle for people in adolescence. For adults, the cycle is the actual age plus 10. However, due to various reasons, the aging and dead cutin cells may not shed timely, which causes acanthokeratoderma, causes no transparency to skin, affects skin color and skin's metabolism and refresh. Therefore, if someone can't carry out the metabolism himself under normal conditions, then an artificial method must be used to achieve the purpose.

(2) Q: Can I wear makeup after using this equipment?

A: Absolutely. After the treatment, the skin turns moisturized and plump. It's more appropriate to wear makeup at this time. This treatment requires no skin breaking and is injury-free, safe, healthy, rapid and effective.

(3) Q: How long does it take to see the effect after doing facial infusion?

A: The skin becomes clean and bright obviously after one treatment, which will become more obvious if you stick to it. This treatment can fix problems like hypopenia in deep skin, dull skin, fine lines and large pores.

(4) Q: How long does it take to see the effect after doing acne removal?

A: The reasons for the formation of acne vary from person to person. There are both internal and external causes. The effect is related to the severity of acne, skin's self-repairing capabilities, regular management and daily care habit. Acne removing a long process and should be done in a gradual way.

(5) Q: Does the acne will come out again after acne removal?

A: Acne are caused by various reasons, and it related to lifestyle, dietary habits and skin care habits. If you don't want to have acne, you should do from the following aspects: remove excessive grease from hair follicle, keep hair follicle unobstructed, and regulate the balance between water and oil.

(6) Q: What is collagen?

A: Collagen is a biological macromolecule substance, and a fibrous protein that is white, opaque and non-branched. It can supplement the nutrition needed by all the skin layers, and enhance collagen activity in the skin, and with the effects of locking moisture, nourishing skin, delaying aging, beautifying looks, relieving facial relaxation, and caring hair. Collagen is a nutrient that the human body must be supplemented for delaying aging. With increasing age, the collagen will drain gradually. Women have started to aging, and their collagen have drained, and the content of collagen have decreased little by little after the age of 20. At the age of 25, collagen loss reaches its peak. At the age of 40, the content of collagen is less than half of that of the age of 18. Collagen and water loss leads to the breakage of collagen fibers and elastic mesh that support the skin, which is the exact cause for the wrinkles on the face of the elderly. It causes the oxidation, atrophy and collapse of skin tissues, which makes the skin dry, wrinkled, slack and inelastic and etc. Therefore, collagen must be supplemented for aging delaying.

(7) Q: Does skin will become thin and sensitive after diamond dermabrasion treatment?

A: It will not happen. Diamond dermabrasion only gets rid of aged cuticle during treatment. Our skin refreshes every 28 days, thus diamond dermabrasion only makes skin become more and more smooth, tender, transparent and causes no sensitive

symptoms.

(8) Q: Is there any dependency for this treatment?






A: There is no dependency. It's target at skin problems, help it improve, and cleanse deep skin, and maintains and prevents skin from having problems caused by other factors, and helps skin absorb nutrition, solve yellowish skin and lack of moisture, as well as rough skin. If quit in later stage, the skin will only back to previous conditions of secretion and senility, and it will not cause dependency.


## 7. Packing List

Main Machine x1  
Power Cord x1  
Skin Scrubber x1  
Electric Face Washing Brush x1  
High Frequency Handle x1  
Transparent Tube Probe Head with High Frequency x1  
Cold & Hot Hammer x1  
EMS Roller Handle x1  
Eye Ultrasound Handle x1  
Facial Ultrasound Handle x1  
Transparent Glass Tube x3  
Diamond Head x9  
Spray Bottle x2  
Black Rubber Ring x several  
Filter Cotton x1 pack  
Fuse x2

## 8. Operational Diagrams

Parameter	Product	Technique	Diagram
Setting			
Facial Deep Cleansing: 90 Minutes/ Once a Week			

<p><b>Electric Face Washing Brush:</b> 2 minutes</p> <p><b>Skin Scrubber:</b> 2 to 5 minutes</p> <p><b>Diamond Dermabrasion:</b> 5 to 8 minutes</p> <p><b>Spray Bottle:</b> 3 to 5 minutes</p> <p><b>EMS Roller:</b> Advised Time: 10 to 15 minutes</p> <p>Energy Level: 3 to 7</p> <p>Mode: A</p> <p><b>Cold &amp; Hot Hammer:</b> Advised Time: 5 to 8 minutes</p> <p>Energy Level: 3 to 7</p>	<p>Makeup Remover + Cleansing Mousse + Hydromask + Hot Steam Machine + Toner, Lotion, Facial Cream + Hyaluronic Acid + Essence + Equipment</p>	<ol style="list-style-type: none"> <li>1. Remove makeup, 1 minute.</li> <li>2. Do cleansing with <b>electric face washing brush</b>. Press cleansing mousse onto face, turn the brush on, and treat the face with the brush moving circlewise, 2 minutes.</li> <li>3. Clean the face, 1 minute.</li> <li>4. Apply toner, 1 minute.</li> <li>5. Apply hydromask and use hot steam, and wait for 10 minutes.</li> <li>6. Cleanse pores of skin's shallow layer with <b>skin scrubber</b>. Treat from the bottom up and one line after another on the face, 2 to 5 minutes</li> <li>7. Deep Cleansing with <b>diamond dermabrasion</b>: Treat the whole forehead and nose with diamond dermabrasion, appress the device onto forehead and leave it there for 1 second, then gently revolve the device with fingers and bring it up(chin→ nose→ forehead), 5 minutes.</li> <li>8. Treat nasal fossa with <b>diamond dermabrasion</b>. Slide back and forth on nasal fossa, 1 minute.</li> <li>9. Treat the cheek. Start with chin, aligning three points and treating it again after having every two to three points be treated, and move one line after another till Temple, 3</li> </ol>	<p>Technique 2, 15, 17</p>  <hr/> <p>Technique 16</p>  <hr/> <p>Technique 6, 9, 10, 18, 19, 28</p>  <hr/> <p>Technique 11, 20, 29</p>  <hr/> <p>Technique 25</p>  <hr/> <p>Technique 31</p>
--	--	---	---

	<p>minutes.</p> <p>10. Do the same on the other side.</p> <p>11. Have more treatment to forehead, and slide towards hairline.</p> <p>12. Clean the face, 2 minutes.</p> <p>13. Do hydrating with <b>spray bottle</b>. Add diluted hyaluronic acid stoste, replenish water for skin, keep away from eyes and nostril. Spray over the whole face for 3 to 5 minutes.</p> <p>14. Rub and pat the hyaluronic acid that on the face until it's absorbed, 1 minute.</p> <p>15. Apply massage cream evenly to the face, and caress it for 3 times.</p> <p>16. Press Chengjiang(REN-24), Renzhong(DU-26), Dicang(ST-4), Jiache(ST-6), Yingxiang(LI-20), Jingming(BI-1), Cuanzhu(BI-2), Yuyao(EM3), Sizhukong(TE-23), Temple, Tongziliao(GB-1), Chengqi(ST-1), and Sibai(ST-2). Repeat 3 times.</p> <p>17. Caress the whole face for 3 times.</p> <p>18. Lift in three lines respectively with two hands doing it alternately, from chin to earlobe, from the corner of mouth to Ermen(SJ21), from nose</p>	
--	--	---





		<p>wing to Temple, from lower eyelid to Temple, and lift the corner of eyes. Repeat 3 times.</p> <p>19. Do the same on the other side.</p> <p>20. Lift from forehead towards hairline. Repeat 3 times.</p> <p>21. Flip and lift the face with rotating fingers. Repeat 3 to 5 times.</p> <p>22. Flip from forehead towards hairline with one hand. Repeat 3 to 5 times.</p> <p>23. One side, flip and lift the face with two hands doing it alternately. Do it for 2 to 3 minutes.</p> <p>24. Do the same on the other side.</p> <p>25. Move zigzag on the whole forehead with middle and ring finger, slide to the front of ear and lift to and fro for 3 times, and pass over the back of ear.</p> <p>26. Treatment done.</p> <p>27. Clean the face, and apply essence, 3 minute.</p> <p>28. Techniques for <b>EMS Roller</b>: Lift upwards from chin to the back of ear, from the corner of mouth to Ermen(SJ21), from nose operating to Temple, from lower eyelid to Temple.</p> <p>29. Move from forehead towards hairline.</p> <p>30. Apply facial mask, and wait for 15 minutes(can use <b>cold &amp; hot hammer</b> while applying facial mask,</p>	
--	--	---	--






		<p>or use facial mask and essence that not yet fully absorbed after removing the mask.)</p> <p>31. If operate on facial mask, please use fixed-point stamp method. If operate on essence, please treat one line after another and circlewise.</p> <p>32. Remove the mask and clean the face.</p> <p>33. Apply toner, lotion and facial cream for daily skin care.</p> <p>34. Treatment done.</p>	
--	--	--	--

**Course of Treatment Recommended:**

One Treatment a week, and ten treatments is a course. After one treatment, skin becomes clean, skin and flat and unrough. After a course, skin gets smooth, delicate and glossy. After two courses, skin turns tightened, bright and elastic. After three courses, an pores, shiny and ruddy skin can be maintained.

**Facial Deep Infusion: 90 Minutes/ Once a Week**

<p><b>Electric Face Washing Brush:</b> 2 minutes</p> <p><b>EMS Roller:</b> Advised Time: 10 to 15 minutes</p> <p>Energy Level: 3 to 7</p> <p>Mode: A</p> <p><b>Facial Ultrasound:</b> Advised Time: 10 to 15 minutes</p>	<p>Makeup Remover + Cleansing Mousse + Hydromask + Hot Steam Machine + Toner, Lotion, Facial Cream + Hyaluronic Acid + Essence + Equipment</p>	<ol style="list-style-type: none"> <li>1. Remove makeup, 1 minute.</li> <li>2. Do cleansing with <b>electric face washing brush</b>. Press facial cleansing mousse onto face, turn the brush on, and treat the face with the brush moving circlewise, 2 minutes.</li> <li>3. Clean the face, 1 minute.</li> <li>4. Apply toner, 1 minute.</li> <li>5. Apply hydromask and hot steam, and wait for 10 minutes.</li> <li>6. Apply massage cream evenly to the face, and caress it for 3 times.</li> </ol>	<p>Technique 2, 6, 8</p>  <hr/> <p>Technique 7</p>  <hr/> <p>Technique 9, 19</p>
--	--	---	--

<p>Mode: A</p> <p><b>Spray Bottle:</b> 3 to 5 minutes</p> <p><b>Cold &amp; Hot Hammer:</b> Advised Time: 5 to 8 minutes</p> <p>Energy Level: 3 to 7</p>		<p>7. Press Chengjiang(REN-24), Renzhong(DU-26), Dicang(ST-4), Jiache(ST-6), Yingxiang(LI-20), Jingming(BI-1), Cuanzhu(BI-2), Yuyao(EM3), Sizhukong(TE-23), Temple, Tongziliao(GB-1), Chengqi(ST-1), and Sibai(ST-2). Repeat 3 times.</p> <p>8. Caress the whole face for 3 times.</p> <p>9. Lift in three lines respectively with two hands doing it alternately, from chin to earlobe, from the corner of mouth to Ermen(SJ21), from nose wing to Temple, from lower eyelid to Temple, and lift the corner of eyes. Repeat 3 times.</p> <p>10. Do the same on the other side.</p> <p>11. Lift from forehead towards hairline. Repeat 3 times.</p> <p>12. Flip and lift the face with rotating fingers. Repeat 3 to 5 times.</p> <p>13. Flip from forehead towards hairline with one hand. Repeat 3 to 5 times.</p> <p>14. One side, flip and lift the face with two hands doing it</p>	 <p>Technique 11, 20</p>  <p>Technique 16</p>  <p>Technique 21, 26</p>  <p>Technique 22, 26</p> 
---	--	--	--






		<p>alternately. Do it for 2 to 3 minutes.</p> <p>15. Do the same on the other side.</p> <p>16. Move zigzag on the whole forehead with middle and ring finger, slide to the front of ear and lift to and fro for 3 times, and pass over the back of ear.</p> <p>17. Treatment done.</p> <p>18. Clean the face, and apply essence, 3 minute.</p> <p>19. Techniques for <b>EMS Roller</b>: Lift upwards from chin to the back of ear, from the corner of mouth to the middle of ear, from philtrum to the middle of ear, from nose wing to Temple, from lower eyelid to Temple.</p> <p>20. Move from forehead towards hairline.</p> <p>21. Techniques for <b>Facial Ultrasound</b>: Lift from chin to the back of ear, from the corner of mouth to the middle of ear, from nose wing to Temple, from lower eyelid to Temple.</p> <p>22. Move from forehead towards hairline.</p> <p>23. Do hydrating with <b>spray bottle</b>. Add diluted hyaluronic acid stoste, replenish</p>	
--	--	---	--

		<p>water for skin, keep away from eyes and nostril. Spray over the whole face for 3 to 5 minutes.</p> <p>24. Rub and pat the hyaluronic acid that on the face until it's absorbed, 1 minute.</p> <p>25. Apply facial mask, and wait for 15 minutes(can use cold &amp; hot hammer while applying facial mask, or use facial mask and essence that not yet fully absorbed after removing the mask.)</p> <p>26. If operate on facial mask, please use fixed-point stamp method. If operate on essence, please treat one line after another and circlewise.</p> <p>27. Remove the mask and clean the face.</p> <p>28. Apply toner, lotion and facial cream for daily skin care.</p> <p>29. Treatment done.</p>	
--	--	--	--

**Course of Treatment Recommended:**

A course consists of ten treatments. After one treatment, skin becomes tightened and rejuvenated, and be supplemented with the needed nutrition. After a course, skin improves and turns glossy. After two courses, skin becomes smooth and plump, and its outline gets clearer gradually. After 3 courses, skin metabolism and detox enhance. And it helps to fade pigment and dark skin, and contributes to prevent and delay skin aging, loosening and sagging, and helps with collagen hyperplasia in underlying skin, which makes skin tender and smooth, firmed and glossy.

Facial Skin Cleansing & Acne Removing: 45 Minutes/ Once  
a Week






<p><b>Electric Face Washing Brush:</b> 2 minutes</p> <p><b>High Frequency Handle:</b> 10 to 15 minutes</p> <p><b>Spray Bottle:</b> 3 to 5 minutes</p> <p><b>EMS Roller:</b> Advised Time: 10 to 15 minutes</p> <p>Energy Level: 3 to 7</p> <p>Mode: A</p> <p><b>Cold &amp; Hot Hammer:</b> Advised Time: 5 to 8 minutes</p> <p>Energy Level: 3 to 7</p>	<p>Makeup Remover + Cleansing Mousse + Hydromask + Hot Steam Machine + Toner, Lotion, Facial Cream + Repairing Solution + Essence + Equipment</p>	<ol style="list-style-type: none"> <li>1. Remove makeup, 1 minute.</li> <li>2. Do cleansing with <b>electric face washing brush</b>. Press facial cleansing mousse onto face, turn the brush on, and treat the face with the brush moving circlewise, 2 minutes.</li> <li>3. Clean the face, 1 minute.</li> <li>4. Apply toner, 1 minute.</li> <li>5. Apply hydromask and hot steam, and wait for 10 minutes.</li> <li>6. Remove acne with acne pin, 5 to 10 minutes.</li> <li>7. Technique for <b>transparent tube handle with high frequency</b>: Diminish inflammation and do disinfection to acne part with fixed-point stamp method. Diminish inflammation to the whole face with moving circlewise. 10 to 15 minutes</li> <li>8. Clean the whole face.</li> <li>9. Anti-allergy with <b>spray bottle</b>. Add diluted repairing solution or freeze-dried power, inject to the whole face with low pressure, avoid hands touching, and do anti-allergy.</li> <li>10. Apply essence to the whole face. 1 minute</li> <li>11. Techniques for <b>EMS Roller</b>: Lift upwards from chin to the back</li> </ol>	<p>Technique 2</p>  <hr/> <p>Technique 11</p>  <hr/> <p>Technique 12</p>  <hr/> <p>Technique 7, 15</p>  <hr/> <p>Technique 7, 15</p>  <hr/>
---	---	--	---

		<p>of ear, from the corner of mouth to the middle of ear, from philtrum to Temple, from nose wing to Temple, from lower eyelid to Temple.</p> <p>12. Treat from forehead towards hairline.</p> <p>13. Clean it, 2 minutes.</p> <p>14. Apply facial mask, and wait for 15 minutes(can use cold &amp; hot hammer while applying facial mask, or use facial mask and essence that not yet fully absorbed after removing the mask.)</p> <p>15. If operate on facial mask, please use fixed-point stamp method. If operate on essence, please treat one line after another and circlewise.</p> <p>16. Remove the mask and clean the face.</p> <p>17. Apply toner, lotion and facial cream for daily skin care.</p> <p>18. Treatment done.</p>	
--	--	--	--

**Course of Treatment Recommended:**

Ten treatments is a course. After one treatment, pores become visibly clean, acne decrease. After one course, skin color turns even and bright, pores improve, acne lessen, and acne marks fade. After three courses, acne will come out again and acne marks fade until disappear if you control your dietary, and develop a good lifestyle and skin care habits.

**Eye Tightening & Anti-aging: 45 Minutes/ Once a Week**

<p><b>Electric Face Washing Brush:</b> 2 minutes</p> <p><b>Eye Ultrasound:</b> Advised Time: 10 to 15 minutes</p> <p>Mode: A</p> <p><b>Cold &amp; Hot Hammer:</b> Advised Time: 5 to 8 minutes</p> <p><b>Energy Level:</b> 3 to 7</p>	<p>Makeup Remover + Cleansing Mousse + Essence/ Eye Cream + Instrument + Eye Mask</p>	<ol style="list-style-type: none"> <li>1. Remove makeup, 1 minute.</li> <li>2. Do cleansing with <b>electric face washing brush</b>. Press facial cleansing mousse onto face, turn the brush on, and treat the face with the brush moving circlewise, 2 minutes.</li> <li>3. Apply toner, 1 minute.</li> <li>4. Apply essence(eye cream) evenly to the eyes with hands moving in circles. Repeat 3 times.</li> <li>5. <b>Eye Ultrasound</b> Technique, Set energy parameters, mode and time(about 10 to 15 minutes)</li> <li>6. Lift the device from lower eyelid to the corner of eye. Repeat 3 to 6 times.</li> <li>7. Lift the device from lower eyelid to Temple. Repeat 3 to 6 times.</li> <li>8. Lift the device from lower eyelid moving in small circles to Temple. Repeat 3 to 6 times.</li> <li>9. Lift the device from lower eyelid to Temple. Repeat 3 to 6 times.</li> <li>10. Lift the device from brow ridge to hairline. Repeat 3 to 6 times.</li> <li>11. Do the same on the other side.</li> <li>12. Clean the face, 1 minutes.</li> <li>13. Apply facial mask, and wait for 15 minutes.</li> <li>14. Remove the mask and clean the face, 2</li> </ol>	<p>Technique 4</p>  <hr/> <p>Technique 6</p>  <hr/> <p>Technique 7</p>  <hr/> <p>Technique 8</p>  <hr/> <p>Technique 9, 10</p>  <hr/>
---	---	--	---







		minutes. 15. Apply eye essence. Treatment done.	
--	--	---	--

**Course of Treatment Recommended:**

Ten treatments is a course. After one treatment, metabolism speeds up, eyes tighten, fine lines fade. After one course, fine lines fade obviously, and the effect maintains. After two courses, the effect becomes obvious gradually, skin tightening improve, collagen in underlying skin increase, eyes tail lift evidently, and eyes become more confident and sharp. After three courses, the effect consolidates, eyes' fine line disappear, and wrinkles will not rebound if you do daily maintenance.

**Neck Refining & Anti-aging: 45 Minutes/ Once a Week**

<p><b>Electric Face Washing Brush:</b> 2 minutes</p> <p><b>Facial Ultrasound:</b> Advised Time: 10 to 15 minutes</p> <p>Mode: A</p>	<p>Makeup Remover + Cleansing Mousse + Massage Cream + Essence + Neck Mask</p>	<ol style="list-style-type: none"> <li>1. Remove makeup.</li> <li>2. Do cleansing with <b>electric face washing brush</b>. Press facial cleansing mousse onto face, turn the brush on, and treat the necks with the brush moving circlewise, 2 minutes.</li> <li>3. Apply toner, 1 minute.</li> <li>4. Massage neck, and lift in the front of chest and apply oil to it with both hands, then move to the back of neck, and press Fengchi(GB20) and Fengfu(DU16). Repeat 3 times.</li> <li>5. Lift jaw and double chin with both hands doing it alternately, and lift to armpit via the back of ear, and pass under armpit. Repeat 3 times.</li> <li>6. Stroke the three channels and collaterals on the side of neck with the V-shaped webbed area</li> </ol>	<p>Technique 4</p>  <hr/> <p>Technique 5</p>  <hr/> <p>Technique 6, 7, 8, 13</p>  <hr/> <p>Technique 9, 10, 12, 14</p> 
---	--	--	--

	<p>between thumb and index finger doing it circlewise, and till armpit and pass under there. Repeat 3 times.</p> <p>7. Stroke the three channels and collaterals on the side of neck with four fingers, and till armpit and pass under there. Repeat 3 times.</p> <p>8. Rub the side of neck with kneeling fingers till it turns hot, and rub till armpit and pass under there. Repeat 3 times.</p> <p>9. Stroke external collarbone with the V-shaped webbed area between thumb and index finger, and till armpit. Repeat 3 to 5 times.</p> <p>10. Do the same on the other side.</p> <p>11. Treatment done.</p> <p><b>12. Facial Ultrasound</b> Technique: Coupled with hands, lifting from double chin to armpit via the back of ear, and pass under the armpit. Repeat 3 times</p> <p>13. Move the equipment circlewise on the side of neck till armpit in three lines respectively, and pass under there. Repeat 3 times.</p> <p>14. Slide the equipment from internal and external collarbones to armpit, and pass under there. Repeat 3 to 5 times.</p> <p>15. Move the equipment circlewise around the neck.</p>	<p>Technique 15</p> 
--	--	---

		Repeat 3 times. 16.Do the same on the other side. 17.Clean the neck with hot towel, 2 minutes. 18.Apply neck mask, and wait for 15 minutes. 19.Clean the neck with hot towel, and use hot compress for 5 minutes. 20.Apply toner, neck essence, neck cream and sunscreen. 21.Treatment done.	
--	--	--	--

**Course of Treatment Recommended:**

A course consists of ten treatments. After one treatment, skin becomes rejuvenated and delicate. After a course, neck skin turns smooth, and circulation of neck lymph accelerates. After two courses, skin becomes tightened, cervical stripe fade, and double chin improves. After 3 courses, skin turns delicate, glossy, firmed and plump, and lymphatic detox speeds up, dark face and acnes improves. It also delays skin aging and rejuvenate it.