

# User Manual

WL-9391



## Preface

Dear Users,

We're pleased to present to you our latest 5-in-1 comprehensive beauty machine WL-9391, which includes multipolar RF, vacuum&RF and 40K technologies for body slimming and skin firming. It is a professional beauty machine, which needs to be operated by professionals with proper training. Any improper use will result in adverse outcomes. Therefore, we advise any personnel to read this manual thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

# Table of Contents

Preface	
Part I	
Introduction	
1. Advantages	
2. Working Principles	
3. Face Anti-ageing	
4. Eye Anti-wrinkle	
5. Neck Care	
6. Body Management	
(1). Waist&Abdomen Shaping	
(2). Arm Shaping	
(3). Back Shaping	
(4). Hip Shaping	
(5). Leg Shaping	
Part II	
1. Detailed Operations	
2. Technical Parameters	
3. User Contraindications	
4. Dos and Don'ts	
5. Troubleshooting & Solutions	
6. FAQs	
7. Packing List	
8. Operational Diagrams	

# Part I

## Brief Introduction

Multi-polar RF, Vacuum&RF and 40K ultrasound fat-burning handpieces are currently among the most popular tools for skin care, skin firming, body shaping and anti-ageing. So, this machine is an effective replacement of all the facial and body care programs in beauty salons. Moreover, it is convenient and easy to operate. It solves both skin and body management problems for beauty-seekers. This equipment requires no injection, no medication, and no operation, for which it has no side effects. It is used externally during the whole process and has an immediate effect. RF heats deeper skin layers, stimulates collagen to regenerate and reorganize in dermis, and promotes blood circulation at skin bottom, which can achieve the results of skin firming and anti-aging. Ultrasounds with specific frequencies target those parts in which there are stubborn fat to bring customers amazing results.


## Advantages

1. This 5-in-1 beauty machine is used for facial and body care.
2. The ultrasonic wave can powerfully burn and remove fat, which surpasses liposuction to solve the problems for obese people.
3. There are several hand pieces for treatment. You can use corresponding one based on different needs and different treatment areas.
4. The perfect combination of ultrasound and radio frequency solves major problems in losing weight. It helps with weight loss and meanwhile firms skin, increases skin elasticity, promotes metabolism, and ultimately strengthen people's constitutions.
5. It adopts the advanced technology of ultrasonic wave fat-burning.
6. It is painless and non-invasive treatment process. Once the treatment is over, customer can go back home and there is no downtime, which means it will not affect customer's life and work.
7. No consumption, low cost and quick returns.
8. Wider treatment range for body and face.
9. Unevenness, bleeding, and swelling will not appear after treatment.

## Face Firming and Anti-ageing

### 1. Working Principles

**RF:** high-frequency waves are used to heat up dermis of face during the treatment to promote skin's metabolism and collagen regeneration. In this way, it can lift skin, reduce wrinkles, slow down aging process, firm and rejuvenate skin.




## 40K

The principle of ultrasonic wave is used to make human adipose cells take a strong impact and produce friction movement among them after strong sound waves enter into the human body, which can effectively consume heat and water of cells and shrink adipose cells. In addition, when the sound wave vibrates, it can produce a strong impact on cells. And then the adipose cells will burst instantly, and they will be reduced so as to achieve the effect of fat removal.

**Ultrasonic Cavitation:** Tens of thousands of tiny bubbles, namely cavitation bubbles, are produced by vibration of liquid. These bubbles grow in the negative pressure region formed by the longitudinal transmission of ultrasound, which will close rapidly in the positive pressure region to be compressed and stretched under alternating positive and negative pressures. When the bubbles burst, there will be huge instantaneous pressure of up to tens of MPa to hundreds of MPa, and strong vibration and noise will be produced.

**Advantages:** It only aims for low-density adipose tissue for its specific frequency, and protects high-density tissue such as vascular and nerve tissue.

Unioisetion cavitation 2.0 vibration: Ultrasound waves penetrate into the fat layer of body to produce cavitation effect, effectively making fat cells burst and emulsify into glyceride



## Vacuum&RF

Massaging skin and muscles with a special vacuum suction handpiece can effectively improve the humoral fluidity of human cells, and increase the movement of cells to improve skin elasticity. Meanwhile, it can accelerate the blood circulation of microvessels, expel toxins from the body through the normal circulation of the lymphatic system, and reduce the probability of occurrence of adverse conditions such as spots, pigmentation and congestion in various parts of the body. The kneading sensation generated by vacuum suction can increase the tissue activity of skin and muscle, and therefore help to reduce tough tissue and increase the elasticity of skin tissue, so that both body slimming and shaping can be realized at the same time. In addition, the movement of vacuum cup can also stimulate sympathetic nervous system to lower skin sensitivity. The suction and release process of air pressure can not only improve capillary system but also promote the flow between deeper veins and lymphoid gland network to strengthen blood vessels and alleviate the condition of varicose veins.

**Advantages:** it is superior to the previous vacuum working mode which is single and simple. Using different modes for different parts can be more efficient and more effective for slimming and shaping body, and it will never harm human body. Uniquely designed RF tips integrate RF and vacuum as a whole, in which RF and vacuum can exist independently and also are complementary to each other. Compared to the ordinary one-polar RF, it can ideally shape body in a more effective, faster and evener way

## Vacuum Physical Effects

### 1. Skin layer

The suction and release of air pressure:

(1) It can improve the fluidity of cells, thus increasing the movement of cells, and can treat

diseases related to blood stasis and blood stagnation.

(2) It can drain excess water of lymph glands and veins.

Effects: redundant water is discharged from tissue fibers.

## 2. Vascular layer

The suction and release of air pressure:

(1) Improve the blood circulation of microvessels for the capillary system, and also promote the flow between the deeper veins and lymphatic networks.

(2) Remove harmful substance in the body.

Effects: Strengthen blood vessels, eliminate toxins and relieve the condition of varicose veins.

## 3. Fibrous layer

The suction and release of air pressure:

It can deliver two kinds of effects in different tissues (including skin, muscle and so on).

(1) Repair cells and increase their activity.

(2) Restore and improve skin elasticity.

(3) Stimulate the production of bone collagen and improve skin plumpness.

(4) Improve the oxygen absorption of skin and increase the consumption of carbon dioxide.

Effects: Break the hard fibers of cellulite, and make it more elastic so as to achieve the effect of body-shaping.

## 4. Nerve layer


The suction and release of air pressure:

It stimulates the surface and deep layer of sympathetic nervous system.

(1) Lower the sensitivity of skin.

(2) Restore skin elasticity and help resist tissue fibrosis.

Effects: Repair cells and lower skin sensitivity.



## 3. Face Anti-ageing

### Effects

1. Firm face and lead nutrients into deep skin.
2. Increase skin elasticity and exquisiteness.
3. Moisturize skin and increase the absorption rate of the skin.
4. Accelerate blood circulation and metabolism.
5. Relieve the condition of double chin and tighten skin.
6. Stimulate collagen regeneration and slow down aging.

### Applicable Range

1. Those whose faces are dark and dull.

2. Those whose skin is flabby and sagging.
3. Those who have fine lines, nasolabial folds, wrinkles at eye corners.
4. Those whose face contours are not clear.
5. Those who have rough skin and large pores.
6. Those who work for a long time in the office and are exposed to ultraviolet rays.
7. Those who have sagging skin, edema or puffiness after giving birth

### **Inapplicable Range**

1. Those who have just undergone plastic surgery with prostheses and metal materials, etc., implanted at the operation area.
2. Those who have used injection products recently, such as hyaluronic acid, skin booster or wrinkle remover.
3. Those who have hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
4. Those who are having allergic reactions, who have severely sensitive skin or metal allergies.
5. Those who are of advanced years.
6. Those who are in menstruation, pregnancy, lactation, or operation recovery period.
7. Those who have skin diseases and infectious diseases.
8. Those who have unrealistic expectations about the effects.

### **Do's and Don'ts after Operation**

1. Do not wash your face with overheated water for three days (you can wash your face with warm or cold water).
2. Strengthen hydration, moisturizing and avoid sunburn.
3. Do not go for hot springs, saunas, violent exercise, etc., for three days.
4. It is recommended to apply the mask at least three times a week.
5. Avoid eating spicy and greasy food, avoid staying up late, smoking and drinking, and eat more fruits and vegetables and food with less grease.
6. Avoid those food which leads to hypertension, hyperglycemia, and hyperlipemia).

## **4. Eye Anti-wrinkle**

### **Effects**

1. Relieve eye fatigue, dark circles, eye bags and edema.
2. Reduce lines around eyes and crow's feet.
3. Reduce dark spots and accelerate blood circulation.
4. Accelerate metabolism and prevent pigmentation.
5. Supplement eyes with nutrition and moisturize skin.
6. Firm and soften skin, lift the corners of the eyes.
7. Accelerate blood circulation around eyes and help skin absorb nutrients efficiently.

### **Applicable Range**

1. Those with wrinkles, fine lines, bags and dark circles around the eyes.
2. Those with dry skin, dry lines and dynamic wrinkles.
3. Those with eyes prone to fatigue and dryness.
4. Those who frequently face computer screens and mobile phones.
5. Those who often stay up late and have dark circles.
6. Those who often stay in a dry or hot environment.

## **Inapplicable Range**

1. Those who had just undergone plastic surgery.
2. Those who have hypertension, heart diseases, diabetes, severe thyroid diseases and malignant tumors.
3. Those who are having an allergic reaction and have severely sensitive skin or are allergic to metals.
4. Those who have skin traumas or wounds.
5. Those who are of advanced years.
6. Those who are pregnant or who are on the road to recovery from operations.
7. Those who have skin diseases and contagious diseases.

## **Do's and Don'ts after Operation**

1. Avoid excessive exposure to the sun. Suncream is encouraged to use.
2. After operation, the treated area is relatively dry and so being hydrated and protection against sunburn are encouraged.
3. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
4. Avoid facial washing with overheated water, taking a sauna, going for hot springs or other violent exercise within 7 days after operation.
5. Moisturizers and eye masks are encouraged to use. It will be preferable to use eye masks 3 times or more a week.

## **5. Neck Care**

### **Effects**

1. Reduce the flabbiness, roughness and dullness of the neck skin.
2. Firm skin and increase skin elasticity.
3. Remove double chin.
4. Accelerate lymph detox and improve facial skin.
5. Prevent neck and lymph diseases.

### **Applicable Range**

1. Those with rough and dull skin on the neck.
2. Those with lymphatic blocking
3. Those with sagging and inelastic neck skin.
4. Those who often bow their heads
5. Those who are not satisfied with their neck tone.

### **Inapplicable Range**

1. Those who have just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin.
4. Those with skin trauma or wound.
5. Those who are of advanced years.
6. Those who are pregnant or convalescents.
7. Those with skin disease and other contagious diseases.



## **Do's and Don'ts after Operation**

1. Avoid sunburn. Keep the neck warm.
2. After operation, the treated area is relatively dry, moisturizing and sun protection is encouraged.
3. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
4. Avoid facial washing with overheated water, taking a sauna, go for hot springs or violent exercise within 7 days after operation.
5. Moisturizing and mask for the neck is encouraged, Use masks at least 3 times a week. Essence or neck cream is encouraged.
6. Avoid bowing head for a long time.

## **6. Body Management**

### **(1). Waist&Abdomen Shaping**

#### **Effects**

1. Relieve cold hands and cold feet, cold womb or cold body of women.
2. Tighten the skin on the waist and abdomen.
3. Reduce lumbar and abdominal fat.
4. Tighten skin, relieve stretch marks, obesity lines and increase skin elasticity.
5. Accelerate metabolism, relieve constipation and improve intestinal peristalsis.

#### **Applicable Range**

1. Those with cold feet, cold hands, cold womb and cold body.
2. Those with lumbar and abdominal fat or who have sagging skin after birth.
3. Those who often sit for a long time, or have unsatisfactory waistlines.
4. Those with stretch marks.
5. Those with constipation or obstruction of abdominal meridians.

#### **Inapplicable Range**

1. Those with "3 Hs"(hypertension, hyperglycemia, hyperlipidemia) and who have heart disease.
2. Those in pregnancy, menstruation or lactation.
3. Those whose incisions are healing or who are on the road to recovery.
4. Those with epilepsy and severe diabetes and hyperthyroidism.
5. Those with malignant tumors, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those with severe gynecological diseases.
8. Those whose gynecological diseases are being treated.

## **Do's and don'ts after Operation**

1. Avoid crop top or anything exposing waist and abdomen as well as cold and windy environment.
2. Avoid binge eating or drinking. No alcohol, no spicy and greasy food. Staying up late should be avoided. Drink warm water.
3. Keep abdomen warm. Take a bath after 4-6 hours.
4. Steam sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.

5. Before going to bed at night, you can rub your abdomen clockwise with your hands, so as to help with metabolism and weight loss.

## **(2). Arm Shaping**

### **Effects**

1. Stimulate collagen growth in the bottom of the skin to shape and firm it.
2. Improve sagging skin.
3. Improve flabby arms and thick arms.
4. Soothe and firm sagging skin.
5. Tighten skin.
6. Accelerate blood circulation and dredge meridians and collaterals.

### **Applicable Range**

1. Those with thick arms and want to look better in clothes.
2. Those with flabby arms.
3. Those with sagging arm skin.
4. Those whose arms are prone to pain and numbness.
5. Those who often hold baby or children

### **Inapplicable Range**

1. Those who have just had plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
3. Those who are during allergic period and has severely sensitive skin.
4. Those with skin trauma or cut.
5. Those who are of advanced years.
6. Those who are pregnant and who are recovering from surgery.
7. Those with skin diseases and infectious diseases.

## **Do's and Don'ts after Operation**

1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours.
3. Drink plenty of warm water to stay hydrated.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise within one week after operation.

## **(3). Back Shaping**

### **Effects**

1. Relieve shoulder and back soreness, and improve the lump on back of the neck.
2. Dredge channels and collaterals, and improve channels and collaterals blocking.
3. Increase blood circulation and metabolism.
4. Improve blood supply to the head and sleep.
5. Regulate the functions of viscera and strengthen the body.
6. Firm skin and prevent skin from sagging.
7. Reduce excess fat on the back and shape back.

### **Applicable Range**

1. Those with sore shoulders and backs and stiff necks.

2. Those who suffer from insomnia, dream a lot or have a worse memory.
3. Those who are prone to fatigue, drowsiness and poor circulation of qi and blood.
4. Those with thick back and want to look better in clothes.
5. Those who have a lump on the back of the neck.

### **Inapplicable Range**

1. Those with metal implants in the body, such as stents and pacemakers, are allergic to metal.
2. Those are during pregnancy, menstruation and lactation.
1. Those whose surgical wound is healing or convalescence
2. Those who has “three Hs” , heart disease, epilepsy and severe diabetes and hyperthyroidism.
3. Those with malignant tumor and hemophilia or severe bleeding.
4. Those with skin diseases and infectious diseases.
5. Those who are too weak.
6. Those who are drunk, thirsty, overworked, on a full or empty stomach.

### **Do's and Don'ts after Operation**

1. Keep warm, avoid cold air and drink plenty of hot water.
2. Take shower in 4-6 hours after operation.
3. Avoid staying up late, drinking and overeating.
4. Avoid raw, cold and spicy foods and get enough sleep.
5. Avoid shoulder-baring and backless clothes

## **(4). Hip Shaping**

Hip is located in the middle of human body, and are the key hub of qi and blood circulation of meridians. It is the main switch of six meridians and also the bridge connecting the qi and blood circulation of upper and lower parts of human body, which plays an important role in building S-shaped curve and promote female attractiveness.

### **Effects**

1. Improve blood circulation and speed up metabolism.
2. Improve menstrual pain, irregular menstruation, abnormal leucorrhea and other female diseases.
3. Improve sleep quality, improve female sexual function and tighten the vagina.
4. Activate the warm nest function, stimulate the secretion of glands, and increase the feelings between husband and wife.
5. Make complexion ruddy, fade color spots and return to young state.
6. Shape hips, improve sagging and outward expansion of hips, tighten skin and increase elasticity.

### **Applicable Range**

1. Those with saggy hips and fat accumulation.
2. Those with stretch marks induced by obesity and obesity.
3. Those whose butt shape is not good-looking, flat and soft outward expansion.
4. Those with cold and cool hips with low hip temperature.
5. Those with cold womb, dysmenorrhea, irregular menstruation, gynecological inflammation and other problems.
6. Those with decreased estrogen levels and poor sex lives.

### **Inapplicable Range**

1. Those who are in menstruation, pregnancy, lactation and operation recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc
3. Those with skin diseases, patients with infectious diseases and skin sensitive period.
4. Those with wounds during recovery from surgery.
5. Those with allergic and severely sensitive skin.
6. Those who have just had liposuction.
7. Those who are of advanced years.

## **Do's and Don'ts after Operation**

1. Keep hips warm and avoid wearing miniskirts and shorts.
2. Take Shower in 4-6 hours after operation.
3. Drink plenty of hot water to avoid cold air.
4. Avoid staying up late, drinking and overeating.
5. Avoid eating raw, cold and spicy foods and get enough sleep.
6. Sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.

## **(5). Leg Shaping**

### **Effects**

1. Tighten skin and prevent sagging.
2. Stimulate collagen regeneration and improve fat lines.
3. Increase leg circulation and detoxification, metabolism.
4. Promote blood circulation and remove blood stasis, dredge channels and collaterals, and prevent varicose veins.
5. Tighten excessive protruding fat on legs and get rid of fat on thick thighs.

### **Applicable Range**

1. Those with poor circulation of lower limbs, edema and obese people.
2. Those who have low immunity and feel discomfort and pain all over the body are prone to colds.
3. Those with constipation and has coarse, flabby skin.
4. Those whose channels and collaterals of the legs are blocked, and the proportion of the legs is poor.

### **Inapplicable Range**

1. Those who are in menstruation, pregnancy or lactation.
2. Those with hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc
3. Those with skin diseases or infectious diseases.
4. Those who have incisions on their bodies.
5. Patients with severe varicose veins and tumors.
6. Those who are having allergic reactions or who have severely sensitive skin.
7. People who have just had liposuction.
8. Those who are of advanced years.
9. Pregnant women or who are on the road to recovery.


## **Do's and Don'ts after Operation**

1. Keep warm after operation. Do not eat cold food.
2. Take shower in 4-6 hours later.

3. Drink more warm water to keep hydrated and speed up metabolism.
4. Do not eat and drink too much or stay up late.
5. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
6. Wear pants after operation. Avoid wearing miniskirts and shorts.

## Part II

### 1. Detailed Operation



Connect the handpieces to the machine as the picture above shows and press power switch, and then you will see the following interface



TIME 30:00 TIME time display and adjustment(range: 00: 00-60: 00)

TIME button for decreasing time

TIME button for increasing time

ENERGY ENERGY energy level display and adjustment(range: 1-9)

ENERGY button for lowering energy level

ENERGY button for raising energy level



button for switching mode



start/stop button

## 1. 40K Function Interface



choose to enter the following interface



time adjustment



buttons for adjusting energy level(users may experience tinnitus during the treatment. And the higher the energy level is, the stronger the sensation will be. The energy level should be above zero and no more than 7. Adjust it based on individual endurance)



button for switching mode

Mode A: constant working modem(suitable for those use this to reduce fat in the early stage)

Mode B: discontinuous working mode(suitable for those who want to consolidate treatment effects in the later stage)

Apply gel evenly to the area, set time, energy level and choose working mode(default mode is

mode A) and then press  to start the treatment.

## 2. Vacuum&RF Function Interface



choose to enter following interface





time adjustment



energy level adjustment(1-7 is recommended)



suction time adjustment(right button is for increasing time, left button is to decreasing time)



release time adjustment(right button is for increasing time, left button is to decreasing time)



button for switching mode

Mode A: the energy of the RF handpiece directly reaches preset level after start button is pressed.

Mode B: the energy of the RF handpiece slowly reaches preset level after start button is pressed and the handpiece touches skin.



knob for adjusting suction level(Turning clockwise is to increase suction and turning anticlockwise is to decrease suction.)

When release time is not zero and suction time is longer than release time, the handpiece is at Suction&Release mode. When release time is 0.0, the handpiece is at Constant Suction mode.

Apply essential oil evenly to the area, set time, energy level, suction and release time and



suction level and choose working mode(default mode is mode A) and then press to start the treatment.

### 3. Sextupole RF Function Interface



choose to enter the following interface



time adjustment



energy level adjustment(1-7 is recommended)




button for switching mode


Mode A: the energy of the RF handpiece directly reaches preset level after start button is pressed.

Mode B: the energy of the RF handpiece slowly reaches preset level after start button is pressed

and the handpiece touches skin.


Apply RF cream evenly to the area, set time, energy level and choose working mode(default mode is mode A) and then press  to start the treatment.


#### 4. Quadrupole RF Function Interface

choose  to enter the following interface




 time adjustment

 energy level adjustment(1-7 is recommended)


 button for switching mode

Mode A: the energy of the RF handpiece directly reaches preset level after start button is pressed.

Mode B: the energy of the RF handpiece slowly reaches preset level after start button is pressed and the handpiece touches skin.

Apply RF cream evenly to the area, set time, energy level and choose working mode(default mode is mode A) and then press  to start the treatment.

#### 5. Tripolar RF Function Interface

choose  to enter the following interface



 time adjustment





energy level adjustment(1-7 is recommended)




button for switching mode

Mode A: the energy of the RF handpiece directly reaches preset level after start button is pressed.

Mode B: the energy of the RF handpiece slowly reaches preset level after start button is pressed and the handpiece touches skin.

Apply massage cream evenly to the area, set time, energy level and choose working

mode(default mode is mode A) and then press  to start the treatment.

## 2. Technical Specifications

Multipolar RF Head for Body

Frequency:1MHz

Power:50W

40KHz Cavitation Head

Frequency:40KHz

Power:20W

Vacuum Bipolar RF Head

Quadrupole RF with Vacuum

Frequency:1MHz

Power:50W

Vacuum: <-80kPa

Pressure: >250kPa

Air flow: >10L/minute

Noise level: <70dB (30cm away)

Tripolar RF Head for Face and eyes' around

Frequency:1MHz

Power:10W

Quadrupole RF Head for Face and Body

Frequency:1MHz

Power:30W

Dimensions: 400x370x290(mm)

Net Weight: 3.8kg

Gross Weight : 6.55kg

Input Power: AC 110V/220V 50/60Hz

Power Consumption:MAX 150W

Display: Touch Screen

## 3. User Contraindications

People described as below should be cautious about using this instrument. Please consult a doctor or a professional before using the instrument. The details are as follows:

1.Those who are in pregnancy or lactation.

2. Those with heart disease or who are equipped with pacemaker.
3. Those whose wounds have not healed and who are recovering from the operation.
4. Those with epilepsy, severe diabetes and hyperthyroidism.
5. Those with malignant tumor, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those who have undergone plastic surgery or used injections for beauty purpose in recent three months.

#### **4. Dos and Dont's**

1. The plug with grounding pin must be used before the instrument is used and the power socket of the instrument is properly grounded.
2. Ensure that the voltage of the instrument is suitable. If the voltage of the local power supply is not stable, we suggest that the user add a voltage regulator with matching power between the mains and the instrument.
3. In order to ensure the treatment effect and normal service life of the instrument, please use the specified parts provided or suggested by the original manufacturer.
4. The instrument should not be placed in damp places or near water sources, and should not be directly exposed to sunlight.
5. Do not place the instrument near a strong heat source, as this may affect the life and normal use of the instrument.
6. Before treatment, please remove all metal objects from the body first, so as to avoid unexpected conditions and influence the curative effect.
7. Please restrain from using the device on eyes, thyroid, parathyroid, testicles, abdomen of pregnant women, and pacemakers.
8. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
9. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use, so as to ensure the safety of the electrical products.
10. When the treatment receiver is taking other weight-loss drugs, it is recommended to stop taking them for 1 to 2 months before treatment. Extension of the course of treatment is recommended if the treatment receiver want to use this machine immediately.
11. Don't use the instrument on an empty stomach. After you are full, you should wait for at least 1 hour before taking the treatment.
12. During operation, the instrument should be in full contact with the skin to avoid uneven heating.
13. Set the energy at lowest level at the beginning and slowly raise it after the customer adapts to it.
14. The treatment area should be kept moisturized during the operation and do not perform the treatment on dry skin.
15. The suction level should not be too high when the vacuum&RF handpiece is used on face to avoid bad skin reactions.
16. Clean the instrument with normal saline after operation to ensure its cleanliness to prolong its service life.
17. Use machine or train machine operators in strict accordance with instructions in the user manual.

## 5. Troubleshooting & Solutions

**1. The instrument cannot be started, and the button light on the back of the instrument does not work?**

- A. Make sure the power cord is connected to the socket with normal power supply.
- B. Whether the fuse on the back of the instrument is loose or blown.

**2. No RF output from the instrument?**

- A. Please check whether the cord connects the handpiece to the machine properly
- B. Please check whether the treatment area has been cleaned. Grease or essential oil products may cause poor contact between the handpiece and human body to cause this problem.

**3. Weakened RF output?**

- A. Please check whether the cord connects the handpiece to the machine properly
- B. Please check whether non-conductive grease and other substances are on the treatment handpiece, which can cause poor contact between the handpiece and skin to weaken the output.
- C. Please check whether the skin product used is specified by the instrument.

**4. There is no suction or weakened suction?**

- A: Please check whether the filter is inserted tight. If not, air may get in and there will be no suction or weakened suction
- B: Please check whether the hose is connected to the machine and the Vacuum&RF handpiece properly.
- C: If the above methods cannot solve the problem, please contact the instrument dealer for assistance.

**5. I can start the instrument, but the screen shows erroneous message?**

- A. Please unplug the power cord behind the instrument and wait for about 1 minute before restarting the instrument.
- B. If the above method cannot solve the problem, please contact the instrument dealer for assistance.

## 6. FAQs

**1. Q: How long can I start seeing results of RF treatment?**

A: Normally, you can see the treatment results that day when you take the operation or within one week. Skin collagen tissue is heated and then contracts, and you can obviously feel the sensation of skin being tightened. Because RF works to stimulate subcutaneous tissue to promote collagen regeneration, so the more you take the treatment, the more noticeable the effects will be.

**2. Q: Is RF harmful to skin?**

A: RF treatment is non-surgical operation. It stimulates collagen regeneration at the skin bottom and boosts metabolism, and therefore has no damage to the skin. Local redness and hotness are the normal conditions of accelerated blood circulation, which will be relieved after a moment. So there's no need to worry about that.

**3. Q: How long does the body-firming treatment take?**

A: One treatment session is 60 minutes. We will combine professional techniques and instruments to achieve remarkable results.

**4. Q: What functions does this instrument have?**

A: It has the functions of dissolving fat and shaping body, smoothing&rejuvenating skin, firming &lifting skin and anti-aging. It can be operated on every part of the body. In terms of facial care, it can form face contour, alleviate wrinkles and saggy skin. In terms of body care, it can reduce regional fatness and build S-shaped curve, and promote metabolism and detoxification of the whole body. Meanwhile, it can also enhance the functions of viscera and solve sub-health problems.

**5. Q: Which one is better in terms of weight loss, liposuction or this machine?**

A: The principle of liposuction is to suck out the excess fat in one part of the body through negative pressure suction so as to achieve the purpose of rapid weight loss. It's a surgical treatment requiring anesthesia and so there is recovery period, which has risks despite its immediate results. However, our instrument is risk-free without side effects. It can not only help you reduce fat but also stimulate collagen regeneration to lift and tighten skin. It can also build perfect curves by lifting skin, and accentuate women's sexual charm.

**6. Q: Will I regain weight after operation?**

A: Ultrasound works to dissolve fat to reduce weight, and users are not prone to regain their weight. Because what's dissolved through ultrasonic treatment is fat, not water. Fat needs a longer period of time to accumulate, so normally people will not regain weight.

**7. Q: Do I need to be on a diet?**

A: Pay attention to what you eat. Because radio frequency and ultrasonic treatment work to dissolve fat and accelerate metabolism. It's best to avoid eating spicy, oily, or fried food in case they affect regular metabolism. Proper workout for sweating helps with fat metabolism process. In such way, the results will be more noticeable.

**8. Q: Does it have side effects on the body?**

A: RF beauty treatment is non-invasive and is one of the most effective and the safest approach to remove wrinkles and shape body. Generally speaking, side effects won't appear. A small number of people may experience transient redness or swelling that will disappear after a few hours. Those with dry skin may experience atrophy after treatment, which will disappear after three days because collagen regeneration will make skin plump up though RF heat can cause moisture loss of skin. Overall, there are no side effects on people's health.

**9. Q: Can ultrasound treatment be performed on every part of the body?**

A: Ultrasound treatment works to rupture fat cells, by which fat is released from cells and then metabolized through lymphatic system. So its power is very high. Heart is sensitive to high frequency sound waves since it is a hollow organ. After being vibrated by sound waves, the heart will reflect the waves due to the different transmission of sound waves by the muscle tissue and the blood in the heart. This force can pull valves away from the heart muscle. If you use the instrument on eyes, it can lead to retinal detachment, so the ultrasound treatment should not be performed on the areas of eyes and heart. (not on the areas of waist, back and chest)

**10. Q: Does ultrasound treatment have side effects on human body?**

A: Ultrasonic beauty treatment is non-surgical and non-invasive program without anaesthesia. It relies on ultrasound (mechanical sound wave) - dissolve fat through cavitation - ultrasound focusing. And the results are fat smash - fat dissolving - skin firming and body sculpting respectively. It only targets low-density fat cells but leaves high-density tissue like blood vessels intact, and therefore has no side effects on human body. Slight tinnitus will probably occur

during the treatment, which is within normal range. So there's no need to worry about that.

**11. Q: Why does tinnitus occur?**

A: It's because ultrasound comes with a high frequency vibration, which could go up to 20kHz. It penetrates into 20mm-deep skin fat layers through ultrasound focusing. Highly intensified ultrasound can cause high-speed friction among fat cells, which can heat and dissolve them. And during the process, sound may be produced so that customer will have the sensation of tinnitus.

**12. Q: What is collagen?**

A: Collagen is a kind of biological macro molecule substance, a kind of white, opaque and non-branched fibrous protein. It can supplement the nutrients needed by all layers of the skin, enhance the collagen activity in the skin, retain moisture, nourish skin, slow down ageing, relieve facial relaxation, provide nutrients for the hair and so on. Collagen is a nutrient that human body must replenish to delay ageing. With the growth of age, collagen will gradually be lost. For women, collagen will begin to decrease since the age of 20. At the age of 25, collagen loss will reach its peak. At the age of 40, the collagen content is less than half of that one has at age of 18. It is the loss of collagen and water that causes the wrinkles on the face of the elderly. The breakage of collagen fibers and elastic fibers will cause skin tissue to be oxidized, atrophied and collapsed. The skin will be dry, wrinkled, flabby and inelastic. Therefore, in order to delay ageing, there must be plenty of collagen.

**13. Q: Why do I need hip treatment?**

A: Because it can help with detox of lymphatic system of body, relieve gynecological diseases, and protect women from harm brought by buttock blockage. Hip blockage will bring about gynecological conditions for the fact that hips are the areas where many important organs of human body exist. If hips suffer from the squeeze of upper and lower energizer meridians, they will be vulnerable to cold, damp and blood stasis. When hips are affected by cold, hip meridian contraction will occur, and then there will be some symptoms such as dysmenorrhea, irregular menstruation, dark menses, blood clot, blood and meridian blockage. It also affects women's fertilization.






**14. Q: How does vacuum detoxify?**

A: Vacuum can cause capillary congestion, stimulate cell to increase vitality. The suction and release of pressure during operation can make local pores continue to open and close, which can promote skin breathing, increase the amount of skin oxygen, and speed up waste removal. So it achieves the effects of promoting qi and blood circulation, dispelling cold and damp, dredging meridians and collaterals, clearing toxins, dissipating heat and so on.

## **7. Packing List**

- 1 x Machine Body
- 1 x Tripolar RF Handpiece
- 1 x Quadropole RF Handpiece
- 1 x Sextupole RF Handpiece
- 1 x Vacuum&RF Handpiece
- 1 x 40K Handpiece
- 1 x Power Cord
- 1 x Filter
- 2 x Fuse

## 8. Operational Diagrams


Parameter Adjustment	Product	Techniques	Diagrams
<h3>Face Firming and Anti-ageing: 70 min once a week</h3>			
<p><b>Quadruple RF</b>            Recommended energy level: 1~7            Time: 15 min            Mode: A(default mode)</p>	<p>Makeup Remover+Face Cleanser+Cold&amp;Hot Steam+Massage Cream+Hydro Mask+Essence+Face Mask</p>	<ol style="list-style-type: none"> <li>1. Remove makeup and clean face, 5 minutes.</li> <li>2. Apply water mask and hot steam, 10 minutes</li> <li>3. Clean the face, 2 minutes.</li> <li>4. Apply toner , 1 minute</li> <li>5. Apply massage cream evenly to face and soothe face, 3 times.</li> <li>6. Apply digital pressing technique on the following acupoints (RN-24, DU-26, ST-4, ST-6, LI-20, BL-1, BL-2, EM3, TE-23, the temple, GB-1, ST-1, ST-2), 3 times.</li> <li>7. Soothe face, 3 times</li> <li>8. Alternately apply lifting technique from chin to ear lobe, mouth corner to Ermen(SJ21), and nose wing to temple, lift eye corner from lower eyelid to temple with both hands, 3 times.</li> <li>9. Repeat the technique on the other side.</li> <li>10. Apply lifting technique towards hairline on forehead, 3 times</li> <li>11. Apply plucking technique with fingers alternately to lift face.</li> <li>12. Apply plucking technique with one hand towards hairline on forehead.</li> <li>13. One side, apply plucking technique with both hands alternately to lift face, 2-3 minutes</li> <li>14. Repeat the technique on the other side.</li> <li>15. Move in a zigzag motion</li> </ol>	<p>Technique5, 7</p>  <hr/> <p>Technique6</p>  <hr/> <p>Technique8, 17</p>  <hr/> <p>Technique10, 12, 19</p>  <hr/> <p>Technique15</p> 





		<p>on forehead with middle and ring fingers, slide to the front side of ears and lift back and forth for 3 times and then slide to the back of ears.</p> <p>16. The end.</p> <p>17. <b>Quadruple RF</b> operation: lift from chin to earlobe, mouth corner to ear center and nose wing to temple line by line, 3 times.</p> <p>18. Repeat the techniques above on the other side.</p> <p>19. Move the handpiece towards hairline on forehead, 3 times.</p> <p>20. Wash face clean, 2 minutes.</p> <p>21. Apply face mask, 15 minutes.</p> <p>22. Clean face, 2 minutes.</p> <p>23. Apply toner, essence, face cream and sunscreen.</p> <p>24. The end.</p>	
--	--	--	--

**Treatment Suggestions**

One course of treatment consists of ten treatment sessions. After one treatment session, skin will be firmed and rejuvenated as the nutrients it needs are absorbed. After one course of treatment, skin will lifted with glow. Two courses will make skin smooth and plump, leaving face contour clearer. Three courses will promote skin metabolism and detox processes, help fade pigmentation and skin darkness, slow down skin ageing and prevent skin from sagging, and stimulate collagen regeneration to leave skin smooth and firm with glow.

**Eye Anti-wrinkle: 45 min 2-3 times a week**



<p><b>Tripolar RF</b>          Recommended energy level: 1~7          Time: 15 min          Mode:          A(default mode)</p>	<p>Make-up remover+face cleanser+essence /eye cream+eye mask+WL-9391</p>	<ol style="list-style-type: none"> <li>1. Remove makeup and clean face, 5 minutes</li> <li>2. Apply toner, 1 minute</li> <li>3. Apply essence(eye cream) evenly to eye area in circular motion, 3 times.</li> <li>4. <b>Tripolar RF</b> operation: adjust energy level, mode and time: 10 minutes.</li> <li>5. Move the instrument to lift skin from lower eyelid to eye corner, 3-6 times.</li> <li>6. Move the instrument to lift skin from lower eyelid to</li> </ol>	<p>Technique3</p>  <p>Technique5</p>
--	--	--	---

		<p>temple, 3-6 times.</p> <p>7. Move the instrument to lift skin from lower eyelid to temple in small circles, 3-6 times.</p> <p>8. Move the instrument to lift skin from lower eyelid to temple, 3-6 times.</p> <p>9. Lift from brow bone to hairline.</p> <p>10. Repeat the operation on the other side.</p> <p>11. Wash eye area clean, 1 minute</p> <p>12. Apply eye mask, 15 minutes</p> <p>13. Remove mask and wash the area clean, 2 minutes</p> <p>14. Apply eye essence and cream.</p> <p>15. The end</p>	 <p>Technique6</p>  <p>Technique7</p>  <p>Technique8, 9</p> 
--	--	--	--





**Treatment Suggestions**

It is recommended to take this treatment 2-3 times per week. After one treatment, eye area is lifted and firmed and blood circulation is accelerated. After one month, fine wrinkles and black circles are reduced and skin colour is brightened. After three months, skin in eye area is totally firmed and rejuvenated with glossiness. If the customers keep using this instrument, eye aging can be relieved and prevented.

**Neck Care: 60 min once a week**

<p><b>Quadruple RF</b>          Recommended energy level: 1~7          Time: 15 min          Mode: A(default mode)</p>	<p>Makeup Remover+Face Cleanser+Massage Cream+Essence+Neck Mask</p>	<ol style="list-style-type: none"> <li>1. Remove makeup and clean neck, 5 min.</li> <li>2. Apply toner, 1 min.</li> <li>3. Neck massage, rub oil into chest with both hands towards the back of the neck and press Gb-20 and Du-16 acupoint, 3 times.</li> <li>4. Alternately lift the lower jaw and double chin area to pass the back of ear to the armpit, 3 times.</li> <li>5. Push downwards with thumb and index finger</li> </ol>	<p>Technique3</p>  <p>Technique4</p>  <p>Technique5, 6, 7</p>
--	---	---	---















		<p>splaying, tracing 3 meridians from the neck side to the armpit, 3 times.</p> <p>6. Push downwards with 4 fingers, tracing 3 meridians from neck side to the armpit, 3 times.</p> <p>7. Push downwards with bending fingers till the neck turns hot and slide to the armpit, 3 times.</p> <p>8. Push downwards with thumb and index finger splaying from external collarbone to armpit, 3-5 times.</p> <p>9. Do the same on the other side.</p> <p>10. The end</p> <p>11. <b>Quadruple RF</b> operation, use one hand and the handpiece to lift double chin area and then pass the back of ear to armpit, 3 times.</p> <p>12. Move the handpiece in circles along three lines on one side of neck to armpit, 3 times.</p> <p>13. Slide the handpiece from inner and outer sides of clavicle to the armpit, 3-5 times.</p> <p>14. Move the handpiece in circles on the whole neck area, 3 times.</p> <p>15. Do the same on the other side.</p> <p>16. Wipe neck clean with warm towel, 2 minutes</p> <p>17. Apply neck mask, 15 minutes</p> <p>18. Wipe neck clean with warm towel, and put the warm towel on neck for 5 minutes.</p> <p>19. Apply toner, neck essence, neck cream and suncream</p> <p>20. The end.</p>	 <p>Technique11</p>  <p>Technique12</p>  <p>Technique14</p> 
--	--	--	---

### Treatment Suggestions

A course of treatment consists of ten treatment sessions. After one treatment, skin becomes rejuvenated and smooth. After a full course, neck skin becomes smooth, and the neck lymphatic circulation is accelerated. After 2 courses of treatment, skin will be tightened, neck lines will be reduced, and the condition of double chin will be relieved. After 3 courses of treatment, skin will be firmed with glow, lymphatic detox will be accelerated and some problems will be resolved such as acne and skin dullness. Moreover, skin aging will be prevented and skin will remain in a young state.

## Waist&Abdomen Shaping: 45 min once a week

<p><b>40K</b>          Recommend ed energy level: 1~7          Time: 15 min          Mode: A(default mode)</p> <p><b>Vacuum&amp;RF</b>          Time: 15 min          Suction time: above zero          Release time: zero          Recommend ed energy level: 1~7          Mode: A(default mode)</p>	<p>Massage Cream(Essential Oil)+Gel+WL-9391</p>	<ol style="list-style-type: none"> <li>1. Massage techniques</li> <li>2. Rub oil into abdomen with hands moving in circular motion, 3 times.</li> <li>3. Rub abdomen back and forth with both hands, 3-5 times.</li> <li>4. Knead abdomen with both hands using chiropractic technique, 3 times.</li> <li>5. Lift Meridian BV(Belt Vessel) on both sides of waist with both hands alternately, 16 times.</li> <li>6. Move hands in 8-shaped motion with right hand and left hand moving clockwise and anti-clockwise respectively on waist, 3 times.</li> <li>7. Overlap hands and message the intestinal canal clockwise, 3 times.</li> <li>8. Move hands in circular motion to soothe the treatment area, 3 times.</li> <li>9. Apply digital pressing technique on the following acupoints: RN13, RN12, RN10, RN8, RN6, RN4, RN3, ST25, SP15. 2 times.</li> <li>10. Starting from RN3 point, push to belly navel, slide to waist and then lift upwards to groin with both thumbs, 3 times.</li> <li>11. Soothe the treatment area with hands and then</li> </ol>	<p>Technique2, 8, 11</p>  <hr/> <p>Technique3, 19</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5</p>  <hr/> <p>Technique6</p>  <hr/> <p>Technique7</p>  <hr/> <p>Technique9</p>
---	---	--	---

		<p>slide to groin.  12. The end  13. <b>40K</b> handpiece operation: starting from one side of the waist, lift the handpiece to abdomen and groin, 3 times.  14. Repeat the previous operation on the other side, 3 times.  15. Move the handpiece on abdomen in small circular motion, 3 times.  16. Move the handpiece on abdomen in big circular motion, 3 times.  <b>17. Vacuum&amp;RF</b> handpiece operation, Constant Suction mode, one side, lift from waist side to abdomen line by line, 3 times  18. Do the same on the other side  19. Lift Meridian BV(Belt Vessel) horizontally back and forth, 3 times  20. Move the cup up and down to shape waist, 3 times  21. Move cup in big circles on abdomen, 3 times  22. One side, lift from waist side to abdomen line by line, 3 times  23. The end.</p>	 <p>Technique10</p>  <p>Technique 13, 17, 22</p>  <p>Technique15</p>  <p>Technique16, 21</p>  <p>Technique20</p> 
--	--	---	---

**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment session, metabolism will be accelerated. After one course, fat will be slowly metabolized, body will start to slim down. Constipation and other issues will be relieved at the same time. After two courses, the effects will gradually become obvious, skin will be tightened and lifted, and skin collagen will be increased. Stretch marks and flabby skin are also alleviated. Three courses will help consolidate the effects. Unwanted fat will gradually disappear, and waist and abdomen curves will start to take shape.

## Arm Shaping: 60 min once a week

### 40K

Time: 20 min(10 min on each side)  
 Recommended energy level: 1~7  
 Mode: A(default mode)

### Sextupole RF


Time: 20 min(10 min on each side)  
 Recommended energy level: 1~7  
 Mode: A(default mode)

Massage Cream(Essential Oil)+Gel+WL-9391

1. Massage techniques
2. Left-right order: lay the arm flat, rub oil into skin from lower arm to the entire arm and slide to the fingers(soothing technique), 3 times.
3. Push the entire arm with both palms, 3 times.
4. Push hands upwards tracing three meridians: Large Intestine Meridian(LI) - Triple Energizer Meridian(TE) - Small Intestine Meridian(SI) to armpits with thumbs and index fingers splaying, 3 times.
5. Soothe the area, 3 times
6. Rub three meridians on arms with bending fingers back and forth to warm them up, 3 times.
7. Soothe the area, 3 times
8. Lay inner arm upwards, and push hands tracing three yin meridians on inner arm: Lung Meridian(LU) - Pericardium Meridian(PC) - Heart Meridian(HT) to armpit with thumbs and index fingers splaying, 3 times respectively.
9. Rub three meridians back and forth with bending fingers to warm them up, 3 times.
10. Soothe the area and slide to fingers
11. Repeat above techniques on the other arm.
12. Wipe arms clean with towel, and then apply gel



Technique 2, 3, 4, 5, 6, 7




Technique 8, 9



Technique 10



Technique 13, 16, 21, 24




Technique 14, 15, 22, 23



Technique 17, 25

Technique 19, 27





	<p>evenly to the entire arm.</p> <p>13. <b>40K</b> operation: lay customer's arm flat, push three meridians from fat part of lower arm to armpit, 3 times.</p> <p>14. Move the handpiece in annular motion along three meridians to armpit, 3 times.</p> <p>15. Redouble the operation on flabby arms, 3 times.</p> <p>16. Move the handpiece along three meridians to armpit, 3 times.</p> <p>17. Lay inner arm upwards and push three meridians on upper arm to armpit, 3 times.</p> <p>18. Move the handpiece in circles along three meridians on upper arm to armpit, 3 times.</p> <p>19. Move the handpiece along three meridians from upper arm to armpit.</p> <p>20. Repeat those techniques on the other arm.</p> <p>21. <b>Sextupole RF</b> operation: lay customer's arm flat, push three meridians from fat part of lower arm to armpit, 3 times.</p> <p>22. Move the handpiece in annular motion along three meridians to armpit, 3 times.</p> <p>23. Redouble the operation on flabby arms, 3 times.</p> <p>24. Move the handpiece along three meridians to armpit, 3 times.</p> <p>25. Lay inner arm upwards and push three meridians on upper arm to armpit, 3 times.</p> <p>26. Move the handpiece in</p>	 <hr data-bbox="1050 383 1305 387"/>
--	---	---









		<p>circles along three meridians on upper arm to armpit, 3 times.</p> <p>27. Move the handpiece along three meridians from upper arm to armpit.</p> <p>28. Repeat those techniques on the other arm.</p> <p>29. The end.</p>	
--	--	--	--

**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment session, it will feel warm and tight, and the arms will be relaxed. After one course of treatment, excessive fat on arms will be reduced and the skin will appear to be tight and plump. Two courses of treatment will help shape arms, dredge meridians, and enhance constitution to strengthen the effects. 3 courses will help consolidate effects and prevent arms from regaining fat.

**Back Shaping: 45 min once a week**

<p><b>Body RF:</b> Time: 15 min Recommended energy level: 1~7 Mode: A(default mode)</p> <p><b>Vacuum&amp;RF</b> Time: 15 min Suction time: above zero Release time: zero Recommended energy level: 1~7 Mode: A(default mode)</p>	<p>Massage Cream(Essential Oil)+Gel+WL-9391</p>	<ol style="list-style-type: none"> <li>1. Massage techniques</li> <li>2. Apply oil to back and press Gb-20 and Du-16 acupoints.</li> <li>3. Pluck Dabanjing(the area connecting neck and shoulder) (from hairline), 3-5 times</li> <li>4. Move thumb outwards from Bladder Meridian(BL) to the sacral region and then to Gb-20 and Du-16 points.</li> <li>5. Starting from neck to caudal vertebra, move hands in S-shape, 3 times.</li> <li>6. Push Bladder Meridian(BL) with thumb in left-right order, 3 times.</li> <li>7. Push Bladder Meridian(BL) to the sacral region with both thumbs at the same time, 3 times.</li> <li>8. Push Bladder Meridian(BL) in three bending fingers, 3 times.</li> <li>9. Push along the medial border of scapula with both hands alternately in left-right order, 3-6 times.</li> <li>10. Push medial border of scapula with both hands</li> </ol>	<p>Technique2, 11</p>  <hr/> <p>Technique 3, 18, 25</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique5</p>  <hr/> <p>Technique6, 7, 8, 14, 16, 23</p>
--	---	--	--

		<p>horizontally, 3-6 times.</p> <p>11. Soothe the back with both hands and push Gb-20 and Du-16 points, 3 times.</p> <p>12. Overlap both thumbs and push SI-11 point and slide to the arm, 3 times.</p> <p>13. Do the same on the other side, 3 times</p> <p>14. Rub Meridian GV and Bladder Meridian (BL) with hands to warm them up.</p> <p>15. The end</p> <p>16. <b>Body RF</b> handpiece operation, move the handpiece along Meridian GV and Bladder Meridian (BL) respectively from neck to the sacral region, 3-5 times.</p> <p>17. Move the handpiece in circles on Du-14 point and sacral region, 3-5 times respectively.</p> <p>18. From neck, move the handpiece on Dabanjing (the area connecting neck and shoulder) back and forth, 3-5 times.</p> <p>19. Move the handpiece to lift the medial border of scapula back and forth in left-right order, 3-5 times.</p> <p>20. Move the handpiece in 8-shaped motion downwards to sacral region, 3 times</p> <p>21. Starting from the left side, move the handpiece in annular motion upwards to armpit, 3 times.</p> <p>22. Apply lifting technique by moving the handpiece from one side of the waist upwards to armpit, 3-5 times.</p> <p>23. <b>Vacuum &amp; RF</b> handpiece operation, move the handpiece along Meridian GV and Bladder Meridian (BL) respectively from neck to the sacral region, 3-5 times.</p> <p>24. Move the handpiece in</p>	 <hr/> <p>Technique 9, 19, 26</p>  <hr/> <p>Technique 10</p>  <hr/> <p>Technique 12, 13</p>  <hr/> <p>Technique 17, 24</p>  <hr/> <p>Technique 20, 27</p>  <hr/> <p>Technique 21, 28</p>  <hr/> <p>Technique 22, 29</p>  <hr/>
--	--	--	---







		<p>circles on Du-14 point and sacral region, 3-5 times respectively.</p> <p>25. From neck, move the handpiece on Dabanjing(the area connecting neck and shoulder) back and forth, 3-5 times.</p> <p>26. Move the handpiece to lift the medial border of scapula back and forth in left-right order, 3-5 times.</p> <p>27. Move the handpiece in 8-shaped motion downwards to sacral region, 3 times</p> <p>28. Starting from the left side, move the handpiece in annular motion upwards to armpit, 3 times.</p> <p>29. Apply lifting technique by moving the handpiece from one side of the waist upwards to armpit, 3-5 times.</p> <p>30. The end</p>	
--	--	---	--

**Treatment suggestions**

A course of treatment consists of 10 treatment sessions. After one treatment session, back will be significantly relaxed, neck stiffness will also be alleviated, and skin will be lifted. A course will help relax shoulder, reduce back fat, relieve the lump on the nape of the back, and shape the back line. After two courses, the lump will appear to be smaller, meridians will be dredged, sub-health problems will be solved, and sleep quality will be improved. 3 courses will help with back thinning, back line refinement, strengthening visceral function, improving digestion, detoxification and metabolism(one course of treatment is to reduce fat, two courses are to strengthen effects, and three courses help consolidate effect and prevent the reappearance of above condition.)

**Hip Shaping: 45 min once a week**





<p><b>Vacuum&amp;RF</b>  Time: 30 min  Suction time: above zero  Release time: zero  Recommended energy level: 1~7  Mode: A(default mode)</p>	<p>Massage Cream(Essential Oil)+WL-9391</p>	<ol style="list-style-type: none"> <li>1. Massage techniques.</li> <li>2. Standing on the side, rub oil into skin by moving both hands clockwise and anticlockwise respectively from sacral region on hips , 3 times (this is soothing technique.)</li> <li>3. Apply pushing technique on the sacral region with two thumbs, 3 times.</li> <li>4. Soothe the treatment area for 3 times and then apply digital pressing technique on the following points: Bl-23, Baliao points(BL31-34), Du-1, Gb-30, BL-36, 3 times.</li> <li>5. Soothe the treatment area, 3 times.</li> <li>6. Left-right order, apply pushing technique from the top of thigh along the four meridians -- Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) to Meridian BV respectively with both hands, 3 times each.</li> <li>7. Alternately apply pushing technique with both hands along Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) from the top of thigh to Meridian BV respectively, 3 times each.</li> <li>8. Overlap both palms and then apply pushing technique from the top of thigh to Meridian BV (lifting), 3 to 5 times.</li> <li>9. Apply pushing technique from both sides of hips to the top of hips(shaping) back and forth, 3 times.</li> <li>10. Repeat technique 7</li> <li>11. Soothe the treatment area.</li> <li>12. Repeat the techniques above on the other side.</li> </ol>	<p>Technique2,5,11</p>  <hr/> <p>Technique3</p>  <hr/> <p>Technique4</p>  <hr/> <p>Technique6, 14, 17</p>  <hr/> <p>Technique7, 10</p>  <hr/> <p>Technique9, 15, 18</p> 
---	---	--	---








		<p>13. The end.</p> <p>14. <b>Vacuum&amp;RF</b> operation: lift from the top of thigh to Meridian BV line by line, 3 times.</p> <p>15. Lift upwards from both sides of the buttocks to the highest point of the buttocks 3 times.</p> <p>16. Move the handpiece in small spiral curve on hip, 3-5 times</p> <p>17. Lift from the top of thigh to Meridian BV line by line, 3 times.</p> <p>18. Lift upwards from both sides of the buttocks to the highest point of the buttocks 3 times.</p> <p>19. Do the same on the other side.</p> <p>20. The end.</p>	
--	--	---	--

**Treatment Suggestions**

One course of treatment consists of 10 treatment sessions. After one treatment, buttocks will experience the sensations of being heated and lifted. After one course of treatment, buttocks start to show obvious curve lines, excessive fat begin to disappear slowly, and skin appear to be tight. After 3 courses, the conditions of cold body and uterus, other gynecological problems such as irregular menstruation are relieved, which can help consolidate the effect to make females more attractive.

**Leg Shaping: 100 min once a week**

<p><b>40K</b> Time: 40 min(10 min on each side) Recommended energy level: 1~7 Mode: A(default mode)</p> <p><b>Vacuum&amp;RF</b> Time: 40 min(10 min on each side) Suction time: above zero Release time:</p>	<p>Massage Cream(Essential Oil)+Gel+WL-9391</p>	<ol style="list-style-type: none"> <li>1. Massage techniques</li> <li>2. Customer lies on his front, back side of the leg, left-right order: rub oil into the skin from calf to thigh and then move back to heel, 3 times</li> <li>3. Alternate the heels of both palms to push the entire leg upwards and move back to heel, 3 times</li> <li>4. Move hands upwards to push four main meridians: Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) with thumbs and index fingers splaying.</li> <li>5. Push popliteal fossa with</li> </ol>	<p>Technique2, 3, 4, 5, 7, 8, 9</p>  <p>Technique6</p>  <p>Technique 12, 14, 21</p>
--	---	---	---

<p>zero  Recommended energy level: 1~7  Mode:  A(default mode)</p>		<p>both hands alternately, 3 times.</p> <p>6. Move both hands upwards alternately in an undulating motion, 3 times.</p> <p>7. Soothe the area, 3 times</p> <p>8. Push four main meridians upwards with bending fingers of both hands, 3 times.</p> <p>9. Soothe the area, 3 times</p> <p>10. Repeat the operation on the other side.</p> <p>11. The end.</p> <p>12. <b>4OK</b> operation: push Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) upwards to popliteal fossa respectively, 3 times.</p> <p>13. Move in small circles on calf to dissolve fat, 3 times</p> <p>14. Move the handpiece upwards to popliteal fossa and soothe the treatment area with hand, 3 times.</p> <p>15. Push 4 meridians from popliteal fossa to the top of thigh respectively, 3 times.</p> <p>16. Move in small circles from popliteal fossa to the top of thigh to dissolve fat, 3 times.</p> <p>17. Push meridians from popliteal fossa to the top of thigh, 3 times</p> <p>18. Move the handpiece in small circles on thigh to dissolve fat, 3-5 times</p> <p>19. Repeat the techniques above on the other leg.</p> <p>20. <b>Vacuum&amp;RF</b> operation: lift from heel to popliteal fossa line by line with one hand and the cup, 3 times</p> <p>21. Move cup in small circles on calf, 3 times</p> <p>22. Move the cup on calf from top to bottom back and forth,3-5 times.</p>	 <hr/> <p>Technique 13, 22, 23</p>  <hr/> <p>Technique 15, 17, 24</p>  <hr/> <p>Technique 16, 18, 25</p>  <hr/> <p>Technique26</p>  <hr/> <p>Technique27</p>  <hr/> <p>Technique 29, 30, 31, 37, 40, 43</p>  <hr/> <p>Technique39, 46</p>
--	--	--	--


23. Lift from popliteal fossa to the top of thigh, 3 times.  
 24. Move cup in small circles on thigh, 3 times.  
 25. Lift from two sides of leg to the middle, 3 times.  
 26. Lift upwards from two sides of thigh in annular motion  
 27. Do the same on the other side.  
 28. Massage the front side of legs: rub oil into skin from feet to thigh root(soothe the treatment area), 3 times  
 29. Push four meridians: Bladder Meridian(BL) - Kidney Meridian(KI) - Liver Meridian(LV) - Gallbladder Meridian(GB) to the top of thigh with thumb and index finger splaying.  
 30. Push four meridians with bending fingers of both hands,3 times  
 31. The end  
 32. Do the same on the other side.  
 33. **40K** handpiece operation: lift the handpiece towards knee from calf while tracing 4 meridians, 3 times(start from thigh if there's not enough unwanted fat on calf.)  
 34. Move the handpiece on two sides of calf in annular motion, 3 times(start from thigh if there's not enough unwanted fat on calf.)  
 35. Lift the handpiece from knee to the top of thigh line by line, 3 times.  
 36. Move the handpiece in small circles on thigh, 3 times.  
 37. Lift the handpiece from knee to the top of thigh in annular motion, 3 times.  
 38. Lift the handpiece from knee to the top of thigh line



Technique38, 44



Technique45



		<p>by line, 3 times.</p> <p>39. <b>Vacuum RF</b> operation: lift the cup towards knee from calf while tracing 4 meridians, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>40. Lift from knee to the top of thigh,3 times</p> <p>41. Move in small circles on thigh, 3 times.</p> <p>42. Lift from two sides of thighs to the middle, 3-5 times.</p> <p>43. Lift upwards on two sides of thigh in annular motion.</p> <p>44. Do the same on the other thigh.</p> <p>45. The end.</p>	
--	--	---	--

**Treatment Suggestions**

A course of treatment consists of ten treatment sessions. After one treatment, leg will be relaxed, and circulation will be accelerated. After one course of treatment, leg will become thin, there will be contraction of skin collagen due to heat effect, and tightening effect can obviously be felt on the skin. After 2 courses of treatment, there will be obvious results with the slimness and firmness of legs promoted. 3 courses will help consolidate the effects.