

User Manual

WL-VS18



Preface

Dear users:

We're pleased to present to you our WL-VS18 infrared slimming blanket featured with weight loss, figure slimming, and body coldness dispelling and so on.

As the name suggested, it is a slimming blanket, targeted at body fat, mainly focusing on dissolving fat, detoxification, dredging meridians and collaterals and body shaping. We aim at a safe and, at the same time, effective beauty result by using high-end technologies. Infrared slimming blankets are for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore we advise any personnel read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

Table of Contents

Preface
Part I
Brief Introduction
Advantages
Working Principles
Effects
Applicable Range
Inapplicable Range
Notes after Operation
Part II
1.Detailed Operations
2.Technical Parameters
3. Safety Precautions
4. Dos and Don'ts
5. Troubleshooting & Solutions
6.FAQs
7.Packing List
8.Operational Diagrams

Part I

Brief Introduction

Infrared slimming blankets are currently the most popular beauty machines that work on weight-losing, body-shaping and detoxifying. It's an effective replacement of all the traditional weight loss programs in beauty salons such as sauna, steam and detox. Moreover, they are convenient and easy to operate. It solves body management issues for beauty-lovers. This blanket requires no injection, no medication, no operation, thus having no side-effects. It is used externally during the whole process and has an immediate effect. Infrared slimming blanket is a blanket-shape weight loss beauty instrument, which works through heat energy effect that make human body sweat. The coverage area is large, which can cover all the parts below the customer's neck. 3 steps of warm up, perspire, drop in temperature help with detoxification and pass out redundant moisture with sweat, thereby achieving the result of reducing weight and slimming body.

Advantages

1. The thermal energy effect of infrared slimming blanket accelerates metabolism, thus promoting the rapid decomposition of fat cells; Complete body shaping trilogy quickly: loosen adipose, fast combustion of adipose, release adipose, thus achieving the result of reducing weight and slimming body.
2. The operation is simple, convenient, time-saving, labor-saving, safe and effective.
3. No medicine, no injection, no diet, no exercise, no side effects.
4. Far infrared ray is a kind of reproductive light, which is invisible. It can restore the diameter of blood vessels, improve the circulation of qi and blood, and strengthen metabolism.
5. In winter, it can unclog the body qi and blood, warm hands and feet, improve sleep quality and help you fall asleep more effectively.

Working Principles

Using infrared dry steam, it accelerates the blood circulation, promotes the metabolism of the body, achieve health care physical therapy effect to the body. Equipped with adjustable digital control, far-infrared can dissolve fat through heat, according to the situation of the customer design temperature, promote blood circulation, promote metabolism, reduce fat accumulation, to promote sweat and subcutaneous fat, so as to achieve weight loss, the effect of beautifying build. At the same time can physiotherapy rheumatism and other symptoms. With the use of essential oil, slimming cream, body skin care products, etc., it can achieve the effects of aromatherapy, slimming, detoxification and lubricating the skin.

Effects

1. Activated the activity of biological macromolecules.
2. Promote and improve blood circulation.
3. Accelerate metabolism and increase muscle nutrients.
4. Reduce inflammation and swelling.
5. Accelerate metabolism and promote the rapid decomposition of fat cells.
6. Loosen fat and burn it fast.
7. Improve insomnia, constipation, endocrine disorders, arthritis and other diseases, improve human immunity.
8. Promote the absorption of the product, so that the nutrients in the product can penetrate into the skin more easily.
9. Improve human microcirculation and enhance cell vitality.
10. Replenish the vitality of the human body, eliminate fatigue and increase physical strength.

Applicable Range

1. Those with systemic obesity and partial obesity.
2. Those who want to lose weight but don't want to exercise.
3. Those who have much body fat and protruding fat accumulation.
4. Those with poor sleep, insomnia, dreaminess and endocrine disorders.
5. Those with poor metabolism.
6. Those who often feel tired owing to frequent staying up late.
7. Those office workers who face computer for long hours.

Inapplicable Range

1. Those with hypertension, heart disease, epilepsy, diabetes and severe hyperthyroidism.
2. Those with malignant tumor, hemophilia or severe bleeding.
3. Those whose wounds have not healed and who are recovering from the operation.
4. Those with skin diseases and infectious diseases.
5. Those who are having an allergic reaction or those with super sensitive skin.
6. Those with skin trauma or wound.
7. Those who are drunk, full, empty, thirsty and overworked.

Notes after Operation


1. Shower 4-6 hours after operation.
2. Dress in time after operation, keep warm.


3. Do not eat immediately after the operation to avoid indigestion and upset stomach.
4. Drink more hot water to replenish water and accelerate metabolism.

Part II

1. Detailed Operation

After checking that the instrument is connected properly, the following starting interface will appear when the power supply is turned on.





2. Technical Parameters

Material: PVC

Power: 650w

Frequency: 50HZ-60Hz

Voltage: 110V or 220-240V

3. Safety Precautions

Patients with the following symptoms should be cautious about using this instrument. Please consult a doctor or a professional before using the instrument. The details are as follows:

1. Those with hypertension, heart disease, epilepsy, diabetes and severe hyperthyroidism.
2. Those with malignant tumor, hemophilia or severe bleeding.
3. Those whose wounds have not healed and who are recovering from the operation.
4. Those with skin diseases and infectious diseases.
5. Those who are having an allergic reaction or those with super sensitive skin.
6. Those with skin trauma or wound.
7. Those who are drunk, full, empty, thirsty and overworked.

4. Dos and Dont's

1. Use a grounding pin and make sure the power socket of the device is properly connected before each operation.
2. Ensure that the voltage of the device is appropriate. If the voltage of the local power supply is unstable, we recommend that users increase to the matching power voltage.
3. In order to ensure the therapeutic effect and normal service life expectancy of the device, please use the specified accessories provided or recommended by the original manufacturer.
4. The device should not be placed in wet environment or near water, nor should they be directly exposed to sunlight.
5. Do not place the device near a strong heat source, because this may affect the life and normal use of the device.
6. Patients who are suffering from illness should consult with a doctor and gain permission from the doctor before using.
7. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use, so as to ensure the safety of the electrical products.
8. Use device or train device operators in strict accordance with instructions in the user manual.
9. When taking other weight-loss drugs, it is recommended to stop taking them for 1 to 2 months before treatment. Extension of the course of treatment is recommended if you wish to lose weight immediately.
10. Don't use the instrument with an empty stomach. After a full meal, wait at least 1 hour before a course of weight loss treatment.
11. Clean the instrument with normal saline after operation to ensure its cleanliness and hygiene and prolong its service life.
12. Stop using if the leather damage occurs on the blanket surface.
13. Do not turn on the power when folded; do not glue the upper and lower blankets together and use in high temperature.
14. During operation, the user shall remove metal objects on the body and enter naked (or wear disposable underwear), and the head shall be exposed.

5. Troubleshooting & Solutions

1. **The instrument cannot be started, but the key light on the back of the instrument does not work?**
 - A. Make sure the power cord is connected to a valid power socket.
 - B. Whether the fuse tube on the back of the instrument is loose or burnt out.

2. No thermal sensation after turning on the instrument?

- A. Please check whether the slimming blanket is connected with the connecting wire.
- B. Please check whether there is any damage to the wire.
- C. If the above methods cannot resolve the issue, please contact the local dealer of the instrument for assistance.

6.FAQs

1. Q: Is infrared slimming blanket harmful to the body?

A: The main effect of infrared fat dissolving blanket is dredging blood veins and meridians and collaterals, dissolve body fat, and detoxify. Far infrared ray medical care principle scientists discover any object can emit infrared ray with only the difference of wavelength. The human body is also an infrared emitter, which emits 5. The 6-15 micron far infrared ray, which accounts for 50% of the total energy of the whole human body, and the materials used in the far infrared ray fiber products effectively emit 5.6-15 microns of far infrared, accounting for the overall wavelength of 90%. As is known to all, when two waves interact with each other, resonance phenomenon will be generated. The human body is an organism. Over 70% to 80% of the human body is composed of water molecules. You can rest assured that infrared fat dissolving blanket is harmless to the body.

2. Q:How long does the it take and how often should I do it?

A: It takes about 60 minutes for a treatment to operate normally. There is no complicated procedure before and after operation. The best time for operation is 45~60 minutes. Under normal circumstance you can do it 3-5 times per week.

7. Packing List

1 set mainframe

1 pc power wire

1 set blanket

2 pcs fuse

8.Operational Diagrams

Parameter Adjustment	Product	Techniques	Diagrams
Body Slimming: 40 min, 2-3 times/week			
Operation Time: 40-60 min Temperature: 40-70°C	Fresh towel+ Infrared slimming blanket	<ol style="list-style-type: none"> 1. Spread out the infrared slimming blanket, attach the two connecting wires with the main machine outlet ports and make sure they stay connected. 2. Have the customer lay inside the blanket and zip up. 3. Set time, temperature and turn on the machine, 45 min. 4. Turn off the machine. 5. Disconnect wires with the main machine. 6. Zip down the blanket. 7. Wipe sweat away with fresh towel. 	No Diagram
<p>Suggested Treatment: Ten treatments account for a full course. At earlier stages you can do it 2-3 times a week. After 2 courses you can switch to 1-2 times per week. After one treatment, body will feel relaxed and blood circulation is improved. After a full course, body fat is gradually reducing and at the same time skin feels more smooth and fine. After 2 course, curve line starts to show, and clothes look better on you. After 3 courses redundant adipose gradually disappear, S curve and slim waist start to take shape.</p>			