

User Manual

WL-WD676A



Preface

Dear users:

We're pleased to present to you our 6-in-1 beauty machine featured with weight loss, figure slimming, skin lifting, anti-aging and so on. It is a multifunctional beauty machine, targeted at both face and body, mainly focusing on lifting and tightening of face, removing wrinkles, and body shaping. We aim at a safe and, at the same time, effective beauty result by using high-end technologies. RF, ultrasonic and lipo laser beauty machines are anti-aging machines for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore we advise any personnel read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

Table of Contents

Preface
Part I
Introduction
Advantages
Working Principles & Treatments
Facial Anti-aging
Anti-aging around Eyes
Anti-aging around Neck
Body-shaping: Waist & Abdomen
Warm Womb & Detoxification
Body-shaping: Breasts
Body-shaping
Arms
Back
Buttocks
Legs
Body Sculpting & Slimming
Overall Facial Care
Part II
1.Detailed Operations
2.Technical Parameters
3. Safety Precautions
4. Dos and Don'ts
5. Troubleshooting & Solutions
6.FAQs
7.Packing List
8.Operational Diagrams

Part I

Brief Introduction

RF and ultrasonic beauty machines are currently the most popular beauty machines that work on tightening, body-shaping and anti-aging. Adding on the effective fat dissolving effect of laser paddles, it's an effective replacement of all the facial and body care programs in beauty salons. Moreover, they are convenient and easy to operate. It solves both skin and body management issues for beauty-lovers. This equipment requires no injection, no medication, no operation, thus having no side-effects. It is used externally during the whole process and has an immediate effect. RF machines heats deeper skin layers, stimulates collagen cells to reform in dermis, increases blood circulation in hypodermis, thus achieving skin tightening and anti-aging. The ultrasound of special frequency is aimed at the parts that are difficult to lose weight in daily life, bringing amazing results.

Advantages

1. 6-in-1 beauty machine that works on both of your facial and body care.
2. Ultrasonic fat burning replaces and transcends liposuction, bringing good news for treating obesity.
3. Various work heads can be replaced in different parts according to different requirements.
4. The perfect combination of ultrasound and radio frequency solves the dilemma of losing weight. It helps with weight loss while at the same time maintains firm skin, increases skin elasticity, metabolism, and ultimately better health.
5. Adopting the cutting-edge technology of ultrasonic fat burning.
6. Painless, non-invasive during the process. Skipping any recovery periods without affecting any work and life plans.
7. No consumption, low cost and quick returns.
8. Wider treatment range and faster and more visible effects.
9. Unevenness, bleed, swelling and stasis will not appear after treatment.
10. Ultrasound, vacuum and 40K all have the energy release of red light respectively. They can dissolve fat and massage while at the same time accelerate metabolism, and sterilize skin. The effect is remarkable with a higher comfort level.
11. One handle is equipped with multiple functions hence saves time and energy and comes with more powerful results.
12. HIFU is with a higher energy level and evenness.
13. Red light RF is with a more effective result.
14. Lipo laser: Safe, painless, even dissolving of fat cells, leaving skin smooth and scar-free.

15. Ultrasound, vacuum and RF all have the energy release of red light respectively. They can dissolve fat and massage while at the same time accelerate metabolism, and sterilize skin. The effect is remarkable with a higher comfort level.

Facial Anti-aging

Principles

RF Face

Radio frequency works through the circuits formed by electrodes in pair. The optimum frequency of radio frequency biothermal effect is 3 MHz, and the maximum speed of extreme change is 8 million charged particles in tissues per second. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy. During the process, RF promotes the metabolism and proliferation of skin collagen through high electromagnetic wave change, therefore heating up dermis layer of the skin, thus achieving effects like tightening, lifting, wrinkle removal, anti-aging and smooth skin.

Biological Effect

Radio frequency can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, making the skin more compact and elastic.

Every time energy comes into the skin, it feels a temporary sense of heat. This means that collagen has been stimulated to produce fever, resulting in tightening of the skin. While operating sliding, the radio frequency wave emitted by the therapeutic head of the super-frequency wave system, under the cold protection of the epidermis, has acted on the deep dermis through the epidermis, rapidly heating water molecules to produce safe Bio-thermal energy. When the natural friction of collagen tissue is heated up to 45°C~60°C, it will shrink immediately and stimulate the collagen to grow continuously. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, and release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

BIO

BIO Cell Activation Instrument, which imitates human biological current, enters muscle cells through skin, stimulates ATP contained in cells, and restores normal operation and function of cells. It collects the latest technology to shape facial contour, thin face, reduce double chin and wrinkles and shrink pores for customers in the shortest time possible. BIO makes facial muscles tight and delicate, and rosy cheeks

will appear naturally . It is more effective on loose skin, so as to achieve the goal of shaping a beautiful face.

Effects

1. Tighten skin and channel nutrients deep into skin layers.
2. Whiten skin and reduce pigmentation and dark spots.
3. Nourish skin and improve skin absorption.
4. Accelerate blood circulation and metabolism.
5. Help relieve double chin and sagging skin.
6. Delay aging and stimulate collagen cells.

Applicable Range

1. Those with facial dullness.
2. Those with skin relaxation and sagging.
3. Those with fine lines, nasolabial folds, crow's feet.
4. Those with indistinct facial contour.
5. Those with rough skin and large pores.
6. Those who are under long-term exposure of UV radiation at workplace.
7. Those with lax skin, edema or obesity after birth.

Inapplicable Range

1. Those who have just undergone plastic surgery which have implanted prostheses, metal materials and so on.
2. Those who have recently received certain injections, such as hyaluronic acid, intradermal injection, wrinkle removal, or plastic surgery.
3. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
4. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.
5. Those with skin trauma or wound.
6. Those who are overaging.
7. Those in menstruation, pregnancy, lactation, surgical recovery.
8. Those with skin diseases and infectious diseases.
9. Those who are unrealistic about the effects.

Notes after Operation

1. Do not wash your face with overheated water within three days (lukewarm or cold water is recommended).
2. Enhance moisturizing and sunscreen.

3. Do not soak in hot springs, steam saunas or exercise strenuously, etc. within three days.
4. Use facial masks at least three times a week.
5. Spicy, greasy food should be avoided. Staying up late, smoking and drinking are prohibited. Eat more vegetables, fruits, and less greasy food.
6. Avoid food that will cause '3 Hs'(high blood pressure, high cholesterol and high blood sugar), less greasy food is advised.

Anti-aging around Eyes

Principles

BIO

BIO Cell Activation Instrument, which imitates human biological current, enters muscle cells through skin, stimulates ATP contained in cells, and restores normal operation and function of cells. It collects the latest technology to shape facial contour, thin face, reduce double chin and wrinkles and shrink pores for customers in the shortest time possible. BIO makes facial muscles tight and delicate, and rosy cheeks will appear naturally . It is more effective on loose skin, so as to achieve the goal of shaping a beautiful face.

Effects

1. Relieve eye fatigue, dark circles, bags and edema.
2. Reduce lines around eyes and crow's feet.
3. Reduce dark spots and accelerate blood circulation.
4. Accelerate metabolism and prevent hyperpigmentation.
5. Supplement nutrition to the eyes and moisturize skin.
6. Tighten and soften skin, lift the corners of the eyes.
7. Accelerate blood circulation around eyes and help skin absorb nutrients efficiently.

Applicable Range

1. Those with wrinkles, fine lines, bags and dark circles around the eyes.
2. Those with dry skin, dry lines and dynamic wrinkles.
3. Those with eyes prone to fatigue and dryness.
4. Those frequently face computer screens and mobile phones.
5. Those who often stay up late and have dark circles.
6. Those often stay in a dry or hot environment.

Inapplicable Range

1. Those who had just undergone plastic surgery.

2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin or allergic to metals.
4. Those with skin trauma or wound.
5. Those who are overaging.
6. Those who are pregnant or convalescents.
7. Those with skin disease and contagious diseases.

Notes after Operation

1. Avoid excessive exposure to the sun. Sun protection is encouraged.
2. After operation, the treated area is relatively dry, moisturizing and sun protection is encouraged.
3. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
4. Avoid facial washing with overheated water, sauna steaming, hot springs or other strenuous exercise within 7 days after operation.
5. Moisturizers and eye masks are encouraged to use, preferably 3 times or more on eye masks a week.

Anti-aging around Neck

Principles

RF

RF instrument can achieve the effect of tightening skin, lifting skin and shaping face. The optimum frequency of radio frequency biothermal effect is 3 MHz, and the maximum speed of extreme change is 8 million charged particles in tissues per second. In the rapid change of electrodes, the dermal tissue naturally resists the radio frequency current and generates thermal energy.

RF wave directly penetrates the skin and generates energy by the resistance action of the skin, which makes the temperature of the skin bottom rise, rapidly and continuously heats the tissue and promotes the growth of collagen in the skin. The principle is tightening collagen and stimulating the proliferation of collagen in the dermis and then is utilized to enhance the tightness and durability in time. After 2-6 months, collagen will gradually proliferate and recombine, making sagging or relaxed skin tight and elevated.

Biological Effect

Radio frequency can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, making the skin more compact and elastic.

Every time energy comes into the skin, it feels a temporary sense of heat. This means that collagen has been stimulated to produce fever, resulting in tightening of the skin. While operating sliding, the radio frequency wave emitted by the therapeutic head of the super-frequency wave system, under the cold protection of the epidermis, has acted on the deep dermis through the epidermis, rapidly heating water molecules to produce safe Bio-thermal energy. When the natural friction of collagen tissue is heated up to 45°C~60°C, it will shrink immediately and stimulate the collagen to continue to grow. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

Effects

1. Reduce fine lines and wrinkles on the neck.
2. Improve the flabbiness, roughness and dullness of the neck skin.
3. Tighten skin and increase skin elasticity.
4. Remove double chin.
5. Accelerate lymphatic detoxification and improve facial skin.
6. Prevent cervical and lymphatic diseases.

Applicable Range

1. Those with fine lines and wrinkles on the neck.
2. Those with loose and inelastic neck skin.
3. Those with dark skin.
4. Those who often bow their heads.
5. Those with clogged lymph system.

Inapplicable Range

1. Those who have just undergone plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe thyroid and malignant tumors.
3. Those who are having an allergic reaction, or people with severe sensitive skin.
4. Those with skin trauma or wound.
5. Those who are overaging.
6. Those who are pregnant or convalescents.
7. Those with skin disease and other contagious diseases.

Notes after Operation

1. Sunscreen is encouraged . Keep the neck warm.

2. After operation, the treated area is relatively dry, moisturizing and sun protection is encouraged.
3. Within 1-3 days after operation, it is better not to use alcohol, AHA or other exfoliating products.
4. Avoid facial washing with overheated water, sauna steaming, hot springs or strenuous exercise within 7 days after operation.
5. Moisturizing and mask for the neck is encouraged, at least 3 times masks a week. Essence or neck cream is encouraged.
6. Avoid bowing head for a long time.

Body-shaping: Waist & Abdomen

Principle

40K&RF

40K: When strong sound wave enters human body, human adipocytes can produce a strong impact and friction movement between adipocytes, which can effectively consume heat, water and shrink adipocytes, so as to achieve the effect of fat removal.

RF: RF wave directly penetrates through the skin and utilize the resistance skin formed to generate energy, heats up dermis layer and promote collagen regeneration, thus achieving effects like tightening, lifting, and refining of skin.

Red LED: With the wavelength of **635nm**, it heats up the bottom layers of skin, therefore reactivating collagens and speeding up the micro-circulation of skin so that it appears to be more elastic, smooth, glossy and has less wrinkles. The theory behind is to heat collagen fibers in dermis to induce collagen through RF technology and improve cell bio-activity through red ray, thus achieving effects like smooth skin with less wrinkles.

Vacuum&RF

Massage skin and muscles with the special negative pressure suction head, which can effectively improve body fluid flow and increase cell activity to achieve effects like improving skin elasticity, accelerating blood circulation of capillary blood vessels, clearing out toxins through lymphatic drainage, and improving and lowering the occurrence of dark spots, hyperpigmentation and blood stasis. The kneading action generated by vacuum suction can increase the tissue activity of skin and muscle, thus helping to reduce the cellular tissue and increase the elasticity of skin tissue, so that both slimming and shaping can be completed at the same time. In addition, vacuum negative pressure movement can also stimulate the surface and deep sympathetic nervous system, improve skin sensitivity. The suction and release process of air pressure can not only improve the capillary system but also improve the flow between

deeper veins and lymph network, strengthen blood vessels, and improve varicose veins.

Advantages: breaking through the previous single and simple working mode of negative air pressure, using different modes for different parts can be more efficient and more effective, slimming, shaping, and will never harm human body. Unique RF function head design, RF and negative pressure are relatively independent but complementary metabolic system as a whole, RF has negative pressure, negative pressure has RF, each other as a whole. Compared to the ordinary single-stage RF can be more effective, faster and more uniform perfect remodeling body.

Vacuum Physical Effects

1. Skin layer

The suction and release of air pressure:

(1) It can improve the fluidity of cells, thus increasing the movement of cells, and can treat diseases related to blood stasis and blood stagnation.

(2) Draining lymph glands and veins.

Effects: redundant water is discharged from tissue fibers.

2. Vascular layer

The suction and release of air pressure:

(1) Improve the blood circulation of microvessels, not only for the capillary system, but also the flow between the deeper veins and lymphatic networks.

(2) Remove excess toxins in the body.

Effects: Strengthen blood vessels, eliminate toxins and relieve varicose veins.

3. Fibrous layer

The suction and release of air pressure:

(1) Repair cell tissue, increase its activity.

(2) Lift and revive skin elasticity.

(3) Stimulate the production of bone collagen and improve skin plumpness.

(4) Improve the oxygen support capacity of skin and increase the consumption of carbon dioxide.

Effects: Break the hard fibers of cellulite, make it more elastic, so as to achieve the effect of body-shaping.

4. Nerve layer

The suction and release of air pressure:

(1) Stimulates the surface and deep sympathetic nervous system.

(2) Improves skin sensitivity.

(3) Repair skin elasticity and resist tissue fibrosis.


Effects: Repair and improve skin sensitivity.

Laser Fat Removal

Laser fat-dissolving device is a red laser with a wavelength of 635nm-650nm, which is a visible spectrum. This wavelength of light has a strong penetrating power, which can effectively activate and repair fat cells, penetrate the fat layer, and dissolve subcutaneous fat by heat. Its treatment method is safe, without side effects or pain. Low levels of laser energy chemical signals are sent to fat cells to break down triglycerides stored in the body's fat layer into free fatty acids and glycerol and release them through membrane channels. Fatty acids and glycerol are transported to tissues that produce metabolic energy. The release of fatty acids is a natural response when the body needs to use stored energy reserves, so there is no unnatural reaction in the body, nor does it affect or damage surrounding structures such as skin, blood vessels and peripheral nerves. After a period of exercise therapy, it will ensure that the free fatty acids of the body are eliminated through complete metabolism.

Biological Effect

It utilizes the body's metabolism to absorb and discharge naturally. In the process of operation, the laser seal small blood vessels, significantly reduce the injuries, introduce the heat treatment area, stimulates the production of collagen, to make it smooth and tight and also achieve the effect of body-sculpting.



Effects

1. Relieve women's cold hands and cold feet, cold womb or cold body.
2. Tighten the skin on the waist and abdomen.
3. Reduce lumbar and abdominal fat, like beer belly and love handles.
4. Tighten skin, relieve stretch marks, obesity lines and increase skin elasticity.
5. Accelerate metabolism, relieve constipation and encourage intestinal peristalsis.
6. Keep fit and prevent obesity.

Applicable Range

1. Those with cold hands and cold feet and cold womb.
2. Those with lumbar and abdominal fat or who have sagging skin after birth.
3. Those who sit for a long time, or with unideal waistlines.
4. Those with striae due to obesity and pregnancy.
5. Those with constipation or obstruction of abdominal meridians.

Inapplicable Range

1. Those with "3 Hs"(hypertension, hyperglycemia, hyperlipidemia) and who have heart disease.
2. Those in pregnancy, menstruation or lactation.
3. Those who have just undergone surgical wounds or convalescents.
4. Those with epilepsy and severe diabetes and hyperthyroidism.
5. Those with malignant tumors, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those with severe gynecological diseases.
8. Those whose gynecological diseases are being treated.

Notes after Operation

1. Avoid crop top or anything exposing waist and abdomen as well as cold and windy environment.
2. Avoid binge eating or drinking. No alcohol, no spicy and greasy food. Staying up late should be avoided. Drink warm water.
3. Keep abdomen warm. Take a bath after 4-6 hours.
4. Steam sauna, hot spring or strenuous exercise should be avoided within 7 days after operation.
5. Before going to bed at night, you can rub your abdomen clockwise with your hands, so as to better the effect of weight loss and metabolism.

Warm womb & detoxification

Ovary is an important organ of women, and can stimulate the hormone secretion of women again, increase the level of hormones in the body, also known as the source of life of women. Good ovarian maintenance can make facial skin fine, smooth and fresh, always maintain toughness and elasticity. It can also promote reproduction and body health, regulate and secrete female hormones, and improve the quality of life of both sexes. Ovarian health can make a woman's bosom plump, compact, and mellow. Bad ovarian function will affect estrogen secretion, sexual function, quality and color of skin, and women BWH(bust,waist,hip) measurements, which results in yellowish face, bloated body, vaginal dryness, and advance entering to anile condition. So ovarian maintenance is very important for women. Women now have a lot of bad habits, such as dressing for grace and does not care about temperature, in addition to that they often eat some cold food. If continuous for a long time they are prone to have cold womb, which can cause symptoms such as dysmenorrhea, menstrual edema and even amenorrhoea. In a word cold womb does great harm to the health of women, so be sure to take means to warm the womb,which can also achieve the effect of body-shaping.

Principle

RF

RF energy works directly on the dermis through the skin epidermis and heats up the body, with its temperature reaching 40 °C-60 °C. Through biothermal effect it enhances the blood circulation and lymph circulation of the cervical region, fast decompose adipose tissue metabolism, thereby to regulate the uterus, increase the period flow, which is good for inhibiting cold womb.

Negative Pressure RF

Through special negative pressure suction head on the skin, and muscle massage, can effectively improve the cells in the body fluids liquidity, increase cell movement can activate cells to improve the effect of skin elasticity, at the same time, accelerate the blood circulation of imperceptible blood-vessel, excess toxins through the lymphatic system the normal cycle of body outside, to reduce and improve the dark spots and pigment, as well as blood stasis.

Advantages: breaking through the previous single and simple working mode of negative air pressure, using different modes for different parts can be more efficient and more effective, slimming, shaping, and will never harm human body. Unique RF functional head design, RF and negative pressure are relatively independent but complementary to each other as one of the metabolic system, deep heating while helping massage and metabolism, so that the effect is twice the result with half the effort.

Effects

1. Improve women's gynecological problems.
2. Improve women's cold hands and feet, cold womb and body.
3. Regulate the function of female reproductive system and enhance the ovarian function.
4. Promote follicular development and regulate the secretion of sex hormones (estrogen, progesterone, androgen).
5. Maintain youthful appearance, delay the arrival of menopause and delay aging.
6. Tighten skin, improve stretch marks and increase skin elasticity.

Applicable Range

1. Those with cold hands and feet, cold womb and body.
2. Those with gynecological inflammatory menstruation and dysmenorrhea.
3. Those whose faces begin to sag.
4. Those who are prone to sore and soft waist.
5. Those with cold abdomen and protruding lower abdomen.

Inapplicable Range

1. Those with 'three Hs' (hypertension, hyperlipidemia, hyperglycemia) and heart disease.
2. Those who are during pregnancy, menstruation and lactation.
3. Those whose surgical wound is healing or convalescence.
4. Those with epilepsy, severe diabetes and hyperthyroidism.
5. Those with malignant tumor and hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Those with severe gynecological diseases.
8. Those whose gynecological diseases are being treated.

Notes after Operation

1. Avoid wearing crop tops, and exposure to wind and cold.
2. Avoid overeating, staying up late, drinking, eating raw, cold, spicy and greasy food, and drinking more hot water.
3. Keep your abdomen warm and bathe 4-6 hours after the treatment.
4. Avoid sauna, hot spring or strenuous exercise within 7 days after operation.

Body-shaping: Breasts

Principle

RF Body

RF heat effect produces specific radio frequency waves in a specific depth of the subcutaneous skin, which directly acts on the dermis through the skin epidermis, thus heating the tissue, promoting the catabolism of subcutaneous fat, and stimulating the proliferation and reorganization of collagen and elastic fibers. With the continuous regeneration and rearrangement of collagen in the dermis, it can be achieved. Let the flabby skin feel the effect of lifting upward and compacting immediately after the treatment. At the same time, massage can promote breast blood circulation, accelerate blood circulation, dredge congested blood clots and silt, and free fat to the chest under the condition of vibration. At the same time, more nutrients are absorbed into the breasts, which not only prevents atypical lobular hyperplasia and breast cancer, but also enables you to have natural, healthy, beautiful and strong breasts, and makes beautiful female friends more confident.

Effects

1. Adjust chest type and reduce accessory breast.

2. Relieve thoracic nodules, slight hyperplasia, alleviate menstrual breast pain.
3. Reduce chest expansion.
4. Alleviate irregular menstruation and reduce facial spots or skin without elasticity.
5. Alleviate postpartum breast atrophy, breast relaxation, breast duct blockage.

Applicable Range

1. Those with bad breast shape and accessory breasts.
2. Those with nodules in the chest, slight hyperplasia and painful breast during menstruation.
3. Those with extended chest fat extended.
4. Those with low immunity.
5. Those with irregular menstruation and facial spots or whose skin is inelastic.
6. Those with improperly developed mammary glands, postpartum breast atrophy, breast relaxation, breast duct blockage.

Inapplicable Range

1. Those with cardiovascular and cerebrovascular diseases, diabetes mellitus, abnormal coagulation function and impaired vital organ function.
2. Those who are long-term or ongoing users of anticoagulants, vasodilators, corticosteroids and other drugs.
3. Those with infectious lesions on the chest skin.
4. Patients with severe thoracic proliferative or fibroma or cyst.
5. Those in pregnancy and lactation.

Notes after Operation

1. Keep warm.
2. Wear styled and comfortable underwear. Don't press your chest hard.
3. Avoid getting cold. Take a bath after 4-6 hours.

Body-shaping

Principle

RF Body

Through RF thermal heating, it produces specific RF waves in specific subcutaneous depth. And it travels through skin epidermis and works directly on the dermis, thus heating up skin cells, promoting the catabolism of subcutaneous fat, at the same time stimulating collagen and elastic fibers hyperplasia restructuring. With the

regeneration and rearrangement of dermis collagen, sagging skin can feel lifted, tightened after treatment.

Biological Effect

Radio frequency can shake millions of high frequency radio waves per second, penetrate the epidermis, directly act on the dermis, deep heat collagen cells, stimulate the growth of collagen fibers, so that it has a large number of new collagen support, making the skin more compact and elastic. RV waves travel through epidermis and directly work on dermis in depth, and quickly heat up water molecules and generate safe thermal energy. When the natural friction of collagen tissue is heated up to 45°C~60°C, it will shrink immediately and stimulate the collagen to continue to grow. At the same time, bio-heat can effectively accelerate blood flow in adipocytes, release free fatty acids (ffa), thereby increasing the dissolution of superficial fat, and achieve the purpose of tightening and lifting the sagging part of the body.

40K&RF

40K: When strong sound wave enter human body, human adipocytes can produce a strong impact and friction movement between adipocytes, which can effectively consume heat, water and shrink adipocytes, so as to achieve the effect of fat removal.

RF: RF radio frequency wave directly penetrates through the skin and utilize the resistance skin formed to generate energy, heats up dermis layer and promote collagen regeneration, thus achieving effects like tightening, lifting, and refining of skin.

Red LED: With the wavelength of **635nm**, it heats up the bottom layers of skin, therefore reactivating collagens and speeding up the micro-circulation of skin so that it appears to be more elastic, smooth, glossy and has less wrinkles. The theory behind is to heat collagen fibers in dermis to induce collagen through RF technology and improve cell bio-activity through red ray, thus achieving effects like smooth skin with less wrinkles.

Vacuum&RF

Massage skin and muscles with the special negative pressure suction head, which can effectively improve body fluid flow and increase cell activity to achieve effects like improving skin elasticity, accelerating blood circulation of capillary blood vessels, clearing out toxins through lymphatic drainage, and improving and lowering the occurrence of dark spots, hyperpigmentation and blood stasis. The kneading action generated by vacuum suction can increase the tissue activity of skin and muscle, thus helping to reduce the cellular tissue and increase the elasticity of skin tissue, so that both slimming and shaping can be completed at the same time. In addition, vacuum negative pressure movement can also stimulate the surface and deep sympathetic

nervous system, improve skin sensitivity. The suction and release process of air pressure can not only improve the capillary system but also improve the flow between deeper veins and lymph network, strengthen blood vessels, and improve varicose veins.

Advantages: breaking through the previous single and simple working mode of negative air pressure, using different modes for different parts can be more efficient and more effective, slimming, shaping, and will never harm human body. Unique RF function head design, RF and negative pressure are relatively independent but complementary metabolic system as a whole, RF has negative pressure, negative pressure has RF, each other as a whole. Compared to the ordinary single-stage RF can be more effective, faster and more uniform perfect remodeling body.

Vacuum Physical Effects

1. Skin layer

The suction and release of air pressure:

(1) It can improve the fluidity of cells, thus increasing the movement of cells, and can treat diseases related to blood stasis and blood stagnation.

(2) Draining lymph glands and veins.

Effects: redundant water is discharged from tissue fibers.

2. Vascular layer

The suction and release of air pressure:

(1) Improve the blood circulation of microvessels, not only for the capillary system, but also the flow between the deeper veins and lymphatic networks.

(2) Remove excess toxins in the body.

Effects: Strengthen blood vessels, eliminate toxins and relieve varicose veins.

3. Fibrous layer

The suction and release of air pressure:

(1) Repair cell tissue, increase its activity.

(2) Lift and revive skin elasticity.

(3) Stimulate the production of bone collagen and improve skin plumpness.

(4) Improve the oxygen support capacity of skin and increase the consumption of carbon dioxide.

Effects: Break the hard fibers of cellulite, make it more elastic, so as to achieve the effect of body-shaping.

4. Nerve layer

The suction and release of air pressure:

(1) Stimulates the surface and deep sympathetic nervous system.

(2) Improves skin sensitivity.

(3) Repair skin elasticity and resist tissue fibrosis.

Effects: Repair and improve skin sensitivity.

Laser Fat Removal

Low levels of laser energy chemical signals are sent to fat cells to break down triglycerides stored in the body's fat layer into free fatty acids and glycerol and release them through membrane channels. Fatty acids and glycerol are transported to tissues that produce metabolic energy. The release of fatty acids is a natural response when the body needs to use stored energy reserves, so there is no unnatural reaction in the body, nor does it affect or damage surrounding structures such as skin, blood vessels and peripheral nerves. After a period of exercise therapy, it will ensure that the free fatty acids of the body are eliminated through complete metabolism.

Biological Effect

Laser mainly acts as a low-energy laser (biological stimulation) to reach fat in depth. It stimulates biological cells and induces or strengthens a series of physiological reactions by giving appropriate energy, including promoting local blood circulation, regulating cell function, enhancing immune function and promoting cell metabolism. It can at the same time tighten skin, wrinkles, reduce fat, and shape body. Laser penetrating power is very strong, and can activate or repair fat cells more effectively. It can also penetrate the fat layer, heat up and dissolve subcutaneous fat. The treatment method is safe, and has no side effects, no pain.

Body-shaping: Arms

Effects

1. Stimulate collagen hyperplasia in the bottom of the skin to shape and firm it.
2. Improve loose skin.
3. Improve flabby arms and thick arms.
4. Soothe and firm loose skin.
5. Dissolve the fat, and say goodbye to flabby arms.
6. Accelerate blood circulation and dredge channels and collaterals.

Applicable Range

1. Those with thick arms and want to look better in clothes.
2. Those with jiggly and flabby arms.
3. Those with loose arm skin.
4. Those whose arms are prone to pain and numbness.

Inapplicable Range

1. Those who have just had plastic surgery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc.

3. Those who are during allergic period and has severely sensitive skin.
4. Those with skin trauma or cut.
5. Those who are overly old.
6. Those who are pregnant and who are recovering from surgery.
7. Those with skin diseases and infectious diseases.

Notes after Operation

1. Keep warm after operation. Do not eat cold food.
2. 4-6 hours later in the shower.
3. Drink plenty of warm water to stay hydrated.
4. Refuse to overeat or stay up late.
5. Avoid sauna, hot spring or strenuous exercise within one week after operation.

Body-shaping: Back

Effects

1. Relieve shoulder and back soreness, and improve the lump on back of the neck.
2. Dredge channels and collaterals, and improve channels and collaterals blocking.
3. Increase blood circulation and metabolism.
4. Improve blood supply to the head and sleep.
5. Regulate the functions of viscera and strengthen the body.
6. Firm skin and prevent sagging.
7. Improve excess fat on the back and shape the back.

Applicable Range

1. Those with sore shoulders and backs and stiff necks.
2. Those with insomnia, dreaminess and fading memory.
3. Those who are prone to fatigue, drowsiness and poor circulation of qi and blood.
4. Those with thick back and want to look better in clothes.
5. Those who has a lump on the back of the neck.

Inapplicable Range

1. Those with metal implants in the body, such as stents and pacemakers, are allergic to metal.
2. Those are during pregnancy, menstruation and lactation.
3. Those whose surgical wound is healing or convalescence
4. Those who has "three Hs" , heart disease, epilepsy and severe diabetes and hyperthyroidism.
5. Those with malignant tumor and hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.

7. Those who are too weak.
8. Those who are drunk, full, empty, thirsty and overworked.

Notes after Operation

1. Keep warm, avoid cold air and drink plenty of hot water.
2. Shower 4-6 hours after operation.
3. Avoid staying up late, drinking and overeating.
4. Avoid raw, cold and spicy foods and get enough sleep.
5. Avoid shoulder-baring and backless clothing.

Body-shaping: Buttocks

Buttocks are located in the middle of human body, and are the key hub of meridian qi and blood operation. They are the main switch of six meridian and also the bridge connecting the upper jiao qi and lower jiao qi blood operation of human body. They are the important factor of the body management whose S curve increases the feminine charm.

Effects

1. Improve blood circulation and speed up metabolism.
2. Improve menstrual pain, irregular menstruation, abnormal leucorrhea and other female diseases.
3. Improve sleep quality, improve female sexual function and tighten the vagina.
4. Activate the warm nest function, stimulate the secretion of glands, and increase the feelings between husband and wife.
5. Make complexion ruddy, fade color spots and return to young state.
6. Shape hips, improve sagging and outward expansion of hips, tighten skin and increase elasticity.

Applicable Range

1. Those with saggy hips and fat accumulation.
2. Those whose butt shape is not good-looking, flat and soft outward expansion.
3. Those with cold and cool hips with low hip temperature.
4. Those with dysmenorrhea, dysmenorrhea, irregular menstruation, gynecological inflammation and other problems.
5. Those with decreased estrogen levels and poor sex lives.

Inapplicable Range

1. Those who are in menstruation, pregnancy, lactation and operation recovery.

2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc
3. Those with skin diseases, patients with infectious diseases and skin sensitive period.
4. Those with wounds during recovery from surgery.
5. Those with allergic and severely sensitive skin.
6. Those who have just had liposuction.
7. Those who are overly old.

Notes after Operation

1. Keep hips warm and avoid wearing miniskirts and shorts.
2. Shower 4-6 hours after operation.
3. Drink plenty of hot water to avoid cold air.
4. Avoid staying up late, drinking and overeating.
5. Avoid eating raw, cold and spicy foods and get enough sleep.
6. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.

Body-shaping: Legs

Effects

1. Tighten skin and prevent sagging.
2. Stimulate collagen regeneration and improve fat lines.
3. Increase leg circulation and detoxification, metabolism.
4. Promote blood circulation and remove blood stasis, dredge channels and collaterals, and prevent varicose veins.
5. Tighten excessive protruding fat on legs and rid thick thighs.

Applicable Range

1. Those with poor circulation of lower limbs, edema and obese people.
2. Those who have low immunity and feel discomfort and pain all over the body are prone to colds.
3. Those with constipation and has coarse, flabby skin.
4. Those whose channels and collaterals of the legs are blocked, and the proportion of the legs is poor.

Inapplicable Range

1. Those who are in menstruation, pregnancy, lactation and operation recovery.
2. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc
3. Those with skin diseases, patients with infectious diseases and skin sensitive period.
4. Those with wounds during recovery from surgery.

5. Patients with severe varicose veins and tumors.
6. Those with allergic and severely sensitive skin.
7. People who have just had liposuction.
8. Those who are overly old.
9. Those who are pregnant or recovering from surgery.

Notes after Operation

1. Keep warm after operation. Do not eat cold food.
2. 4-6 hours later in the shower.
3. Drink more warm water to replenish water and speed up metabolism.
4. Refuse to overeat or stay up late.
5. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
6. Wear pants when you're done. Avoid miniskirts and shorts.

Body Sculpting&Slimming

Effects

1. Relieve shoulder and back soreness, and improve the lump on back of the neck.
2. Dredge channels and collaterals, and improve channels and collaterals blocking.
3. Accelerate blood circulation and dredge channels and collaterals.
4. Promote blood circulation and remove blood stasis, dredge channels and collaterals, and prevent varicose veins.
5. Regulate the functions of viscera and strengthen the body.
7. Stimulate collagen regeneration, firming and shaping.
8. Improve loose and soft waist and abdomen skin.
9. Improve the jiggly and thick arm.
10. Improve unwanted fat on the back, waist and abdomen.
11. Tighten skin and prevent sagging.
12. Improve stretch marks, fat marks, and increase skin elasticity.
13. Accelerate metabolism, improve constipation and intestinal peristalsis.

Applicable Range

1. Those with cold hands, feet, womb and body.
2. Those who are sedentary and don't look good at the waist.
3. Those with unwanted fat on waist and abdomen, prominent belly, loose skin after delivery.
4. Those with fat marks, stretch marks.
5. Those with constipation, abdominal meridian blocking.
6. Those with big arms, thick back, poor leg shape ratio and want to look better in clothes.
7. Those with jiggly and flabby arms.

8. Those with loose skin on arms and legs.
9. Those whose hands and legs are easily sore and numb.
10. Those with sore shoulders and backs and stiff necks.
11. Those with insomnia and dream, memory decline of the population.
12. Those who are prone to fatigue, drowsiness and poor circulation of qi and blood.
13. Those with thick back and unattractive clothes.
14. Those with poor circulation of lower limbs, edema and obese people.
15. Those who have low immunity and feel discomfort and pain all over the body are susceptible to colds.

Inapplicable Range

1. Those who are during pregnancy, menstruation and lactation.
2. Those with hypertension, heart disease, diabetes, thyroid, epilepsy and severe hyperthyroidism.
3. Those with malignant tumor and hemophilia or severe bleeding.
4. Those who have just had liposuction.
5. Those with severe varicose veins and tumors.
6. Those whose surgical wound is healing or recovering.
7. Those with skin diseases and infectious diseases.
8. Those with severe gynecological diseases.
9. Those with gynecological diseases being treated.
10. Those with allergic and severely sensitive skin.
11. Those with skin trauma or cut.
12. Those who are overly old.
13. Those with metal implants in the body, such as stents and pacemakers, are allergic to metal.
14. Those who are too weak.
15. Those who over-drink, full, empty, thirsty and overworked.

Notes after Operation

1. Avoid wearing clothes that show navel, shoulder or back, miniskirts or shorts, and avoid getting cold from the wind.
2. Avoid overeating, staying up late, drinking and eating raw, cold, spicy and greasy food.
3. Keep warm after operation, and do not bathe until 4-6 hours apart.
4. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
5. Before going to bed at night, you can knead your abdomen clockwise with your hands to make the metabolic effect of weight loss better.
6. Drink plenty of hot water to hydrate and speed up metabolism.

Overall Facial Care

Effects

1. Improve facial wrinkles, and crow's feet.
2. Deeply introduce nutrition, fade pigment, spots and whiten skin.
3. Increase skin elasticity, firm and delicate skin, and lift canthus.
4. Relieve eye fatigue, dark circles, bags and edema.
5. Increase blood circulation and metabolism.
6. Improve sagging and relieve sagging.
7. Stimulate collagen hyperplasia and delay aging.
8. Accelerate lymphatic detoxification and improve facial skin quality.

Applicable Range

1. Those with dark complexion, rough skin, large pores and dry skin.
2. Those whose skin is flabby, saggy and inelastic.
3. Those with fine lines, decree lines, canthus lines.
4. Those with wrinkles, fine lines, bags under the eyes, and dark circles under the eyes.
5. Those who face computer screens and cell phones for long hours.
6. Those who often stay up late with dark circles under their eyes.
7. Those who always be in a dry or hot environment.
- 8 Those who after giving birth to a baby has loose skin, edema or fat.
9. Those who work long hours in the office , facing ultraviolet.

Inapplicable Range

1. Those who have just undergone plastic surgery, the treatment site is implanted with prostheses, metal materials, etc.
2. Those who has recent injection products, such as hyaluronic acid, hydro-gloss, wrinkle removal, or plastic surgery.
3. Those with hypertension, heart disease, diabetes, severe thyroid gland, malignant tumor, etc.
4. Those who are in the allergic period, severely sensitive skin and allergic to metal.
5. Those who are in menstrual period, pregnancy period, lactation period, operation recovery period.
6. Those with skin trauma or cut.
7. Those who are overly old.
8. Those with skin diseases and infectious diseases.
9. Those who have unrealistic expectations about results.


Notes after Operation

1. Do not wash your face with overheated water (either warm or cold) for seven days.
2. Strengthen hydrating and moisturizing, avoid insolation, and pay attention to sunscreen.
4. It is recommended to apply facial mask and eye mask at least 3 times a week.
5. Within 1-3 days after operation, it is best not to use products such as alcohol, fruit acid and scrub.
6. Avoid sauna, hot spring or strenuous exercise for 7 days after operation.
7. Avoid spicy and greasy food, stay up late, smoke and drink, eat more fruits and vegetables, light food.
8. Avoid three high food, mainly less greasy food.
9. Apply facial essence and eye essence.

Part II

1. Detailed Operation

After checking that the instrument is connected properly, the following starting interface will appear when the power supply is turned on.





1.1 Function Selection




Start/Pause



When it's on



Energy Adjustment



Time Adjustment



Decrease of Time



Increase of Time



Restore



Home




Mode Selection

1.2 Detailed Operation: 40K & RF



Select to go to the interface below



40K Mode Selection



M1 is consecutive working mode




M2 is intermittent working mode

1.3 Detailed Operation: BODY VACUUM & RF



Select **BODY VACUUM & RF** to go to the interface below



Vacuum suction time adjustment



Vacuum release time adjustment



Decrease of vacuum time



Increase of vacuum time



Time display of suction and release

1.4 Detailed Operation: BODY RF



Select **BODY RF** to go to the interface below



Energy adjustment of body RF

1.5 Detailed Operation: FACE RF



Select to go to the interface below



Energy adjustment of face RF

1.6 Detailed Operation: BIO



Select  to go to the interface below



M1 is intermittent mode, suitable for beginners and first-time user.




M2 is accelerated intermittent mode, suitable for those have slight dark and dull skin.




M3 is consecutive mode, suitable for those have acne skin, coarse and dull skin.

1.7 Detailed Operation: LIPO LASER



Select  to go to the interface below



Energy Display of Lipo

Laser



Decrease of Energy Level



Increase of Energy Level



M1 is the consecutive mode (light stays on the entire time)



M2 is the intermittent mode with relatively slower rate of light change



M3 is the intermittent mode with relatively faster rate of light change

2. Technical Parameters

Power supply Input:110V/220V 60HZ/50HZ

Power Output: Max 110W

40K RF Cavitation Head

Frequency: 40KHz

Power: 45W

Vacuum+RF head:

Power supply output: 66V

Frequency: 3MHz

Power:90W

12-polar RF head:

Power supply output:66V

Frequency:5MHz

Power:50W

9- Polar RF for Body

Frequency:1MHz

Power:40W

RF Head For Eye and Face:

Power supply output:66V

Frequency:1MHz

Power:30W

LED Laser

Wavelength –635nm

Energy output – 52 x5mW

Totaled output:52 x5mw=260mw

Power output consumption– 100VA

Red wavelength:650 - 730nm

3. Safety Precautions

Patients with the following symptoms should be cautious about using this instrument. Please consult a doctor or a professional before using the instrument. The details are as follows:

1. Those who are pregnant or lactating.
2. Those with heart disease or pacemaker equipped.

3. Those whose wounds have not healed and who are recovering from the operation.
4. Those with epilepsy, severe diabetes and hyperthyroidism.
5. Those with malignant tumor, hemophilia or severe bleeding.
6. Those with skin diseases and infectious diseases.
7. Avoid wearing metal jewelry during operation.
8. It is not allowed to be use around metal, plastic, silicon and other parts embedded in the body.

4.Dos and Dont's

1. After each use of the instrument, please clean the treatment head with water and keep it properly.
2. The plug with grounding pin must be used before the instrument is used and the power socket of the instrument is properly grounded.
3. Ensure that the voltage of the instrument is suitable. If the voltage of the local power supply is not stable, we suggest that the user add a voltage regulator with matching power between the mains and the instrument.
4. In order to ensure the treatment effect and normal service life of the instrument, please use the specified parts provided or suggested by the original manufacturer.
5. The instrument should not be placed in damp places or near water sources, and should not be directly exposed to sunlight.
6. Do not place the instrument near a strong heat source, as this may affect the life and normal use of the instrument.
7. Please remove all metal objects from your body before treatment, so as not to affect the curative effect due to unexpected conditions.
8. Do not use in the eyes, thyroid gland, parathyroid gland, testicles, pregnant woman's abdomen and atrial pulse generator.
9. Patients who are currently ill should also use this instrument with caution unless approved by a physician.
10. Please turn off the power switch of the instrument when no guest is using it, and make sure to turn off the main power of the instrument before leaving after daily use to ensure the safety of electrical products.
11. Use the instrument or train the operator strictly according to the instruction of the manual.
12. If you are taking other weight-loss drugs, it is recommended to stop taking the drugs for 1 to 2 months before you lose weight. If you want to lose weight immediately, you should extend the course of treatment
13. Do not use the instrument on an empty stomach. Only after at least 1 hour after a full meal can you carry out the weight-loss treatment.
14. During operation, the instrument should be in full contact with the skin to avoid uneven heating.
15. Start with the lowest energy at the beginning of operation, and gradually add energy after adaptation.

16. When using this instrument, the operating area must be kept moist and avoid dry skin operation.
17. After the operation, scrub the instrument with normal saline to ensure its cleanliness and extend its service life.
18. During operation, the laser paddle must be fixed on the customer first, and then the instrument should be started to adjust the energy. The energy should be adjusted slowly from low to high to make the customer feel warm and comfortable.
19. Must apply products with high oil content, and keep it moist.
20. During operation, the whole surface of the laser paddle should be in full contact with the skin, not only with a small area.
21. Laser paddle can be used alone, directly fixed where the customer wants to lose fat, and adjust energy level. Operation time of one treatment area should be about 20 minutes. Combined with massage, RF and ultrasonic instruments it achieves better results.

5. Troubleshooting & Solutions

1. The instrument cannot be started, and the key light on the back of the instrument does not work.

- A. Make sure the power cord is connected to a valid power socket.
- B. Whether the fuse tube on the back of the instrument is loose or burnt out.

2. No RF output of the instrument?

- A. Please check whether the plug connecting the instrument handle and the body is tightly connected.
- B. Please check whether the treatment area has been cleaned. Grease or oil essential oil products may cause poor contact between the treatment handle of the instrument and human body, resulting in no RF output.

3. Weakened RF output?

- A. Please check whether the plug connecting the instrument handle and the body is tightly connected.
- B. Please check whether non-conductive grease and other substances are stained on the treatment handle, which will cause poor contact and weaken the output.
- C. Please check whether the product used is the adaptive product specified by the instrument.

4. No suction or very little suction?

- A. Please close the instrument and check the filter element specially used for the instrument, which may need to be replaced.
- B. Please check the rubber ring of the handle plug, because the rubber ring may be worn and cause air leakage.

C. Please check whether the oil filter cup outside the filter is tightened and whether the rubber ring is worn. It is possible that air leakage at this position may lead to poor air pressure.

D. If the above methods cannot be handled, please contact the instrument dealer for assistance.

5. I can start the instrument, but the screen show's error message?

A. Please unplug the power plug behind the instrument and wait for about 1 minute before restarting the instrument.

B. If the above methods cannot be handled, please contact the instrument dealer for assistance.

6.FAQs

1. Q: How long can I start seeing results of RF treatment?

A: Under normal circumstances within that day or a week. Skin collagen tissue is heated and produces contractions, and you can obviously feel skin tightening. Because RF works by stimulating subcutaneous tissue and bring about sustainable collagen regeneration, so the more you do it, the more visible effects you shall receive.

2. Q: Is RF harmful to skin?

A: RF treatment falls under the category of non-surgery operation. It stimulates collagen regeneration at the skin bottom and boosts metabolism, therefore has no damage to the skin. Partial redness and fever is the normal phenomenon of accelerated blood circulation, and will be resolved by oneself after a moment, so there's no need to worry.

3.Q: How long does the body firming treatment take?

A: One operation time is 60 minutes. We will combine professional techniques and instruments to achieve remarkable results.

4.Q: What are all the functions of this instrument?

A: Dissolve fat, body shaping, smooth skin, firming, lifting, and anti-aging. It can be operated all over the body. In terms of facial care, it can contour facial outline, alleviate wrinkles and saggy skin. In terms of body care, it can reduce regional fatness and build S curve, promote the metabolism and detoxification of the whole body. Meanwhile, it can also enhance the functions of viscera and regulate the sub-health of the body.

5.Q:Which one is better in terms of weight loss, liposuction or this machine?

A: The principle of liposuction is to suck out the excess fat in one part of the body through the method of negative pressure suction, so as to achieve the purpose of rapid weight loss on parts of the body. It achieves fast results, requires anesthesia, recovery period. It falls under the category of surgery operation and comes with risks.

However, our machine is risk-free and has no side effects. While losing weight, it can stimulate collagen regeneration to lift and tighten skin. By lifting, it can also shape perfect curves and accentuate women's natural charm.

6.Q: Will I experience rebound after operation?

A: RF works by dissolving fat to reduce weight, and are not prone to rebound. Because what's dissolved through RF treatment is fat, not water. The formation of fat needs a longer period of time to accumulate, so it is not prone to rebound.

7.Q: Do I need to be on a diet?

A: Pay attention to what you eat. Because radio frequency and ultrasonic treatment end with blasting adipose and accelerated metabolism. It's best to avoid eating spicy, oily, or fried food, to avoid affecting regular metabolism. Also proper workout and sweating help with the process, in a way weight loss will be more visible as well.

8.Q: Does it have any side effects on the body?

A: RF beauty treatment is non-invasive and so far the most safe and effective method to remove wrinkles and shape body. Generally speaking side effects won't appear. A small number of people may experience transient redness or swelling that will disappear after a few hours. Those with dry skin may experience atrophy after the initial treatment. The skin will lose moisture due to RF heat. However, the skin will become plump at the beginning of collagen regeneration and these symptoms will disappear after three days. There are no side effects on health.

9.Q: Can RF instrument be used for chest shaping?

A: As the growth of age and female physiology period, sub-health state, gravity, stimulation during lactation, a large number of nutrition loss... All the above will cause collagen loss, decreased blood flow, damaged elastic tissue, flabby ligament gland, and prolapsed breasts. Through RF energy conduction, collagen protease can be promoted to repair elastic tissue and ligament glands, thus achieving the effect of chest lifting and shaping.

10. Q: Can ultrasound be operated all over the body?

A: Ultrasonic cavitation ruptures fat cells and causes mechanical disruption of fat cell membrane. Liquefied fat cells then cleared away by metabolic pathways. Our heart is a hollow organ, and sensitive to high frequency sound waves. Back and forth reflex occur due to the fact that muscle tissue and blood in the heart do not conduct sound waves in the same way. The reflex force can detach cardiac valves from cardiac muscles. If directed at the eye, it can cause retinal detachment. In a word, avoid eye and heart area when using ultrasonic beauty device. (It is also advised not to use it on waist, back and chest.)

11.Q: Does ultrasound have side effects on human body?

A: Ultrasonic beauty machine requires no surgery, no opening up and no anesthesia, therefore is indeed non-invasive. It relies on ultrasound (mechanical sound wave) - blast fat through cavitation - ultrasound focusing. And the results are fat smash - fat dissolving - skin firming and body sculpting respectively. It only targets at low density fat cells while leaving high density tissues like blood vessels intact, therefore having no side effects on human body. Slight tinnitus might appear to happen during the treatment, which is within normal range so there's no need to worry.

12.Q: Why does tinnitus occur?

A: It's because ultrasound comes with a high frequency vibration, which could go up to 20kHz. It also goes deep into skin fat layers, which could go down to 20mm below epidermis. Highly intensifies ultrasound cause high speed friction among fat cells, resulting in them bursting and heating up, and eventually emulsified. And during the process, ultrasound can cause microvibration to the tissue that make us feel like as "tinnitus".

13.Q: What is collagen?

A: Collagen is a kind of biological macromolecule substance, a kind of white, opaque and non-branched fibrous protein. It can supplement the nutrition needed by all layers of the skin, enhance the collagen activity in the skin, lock moisture, nourish the skin, delay aging, beautify, relieve facial relaxation, provide nutrients for the hair and so on. Collagen is a nutrient that human body must replenish to delay aging. With the growth of age, collagen will gradually be lost. At the age of 20, women have begun to age, collagen gradually decrease. At the age of 25, collagen loss reaches its peak. At the age of 40, the content is less than half of which one has at age 18. It is the loss of collagen and water that causes the wrinkles on the face of the elderly. The breakage of collagen fibers and elastic mesh will cause skin tissue to be oxidized, atrophied and collapsed. The skin will be dry, wrinkled, flabby and inelastic etc,. Therefore, in order to delay aging, collagen must be supplemented.

14. Q: Why do you need breast maintenance?

A: The lymph nodes of bosom rank the top among other parts of the body, which also lead to toxin accumulation very easily, and urbanite life rhythm is very fast now, working pressure is very big, which can bring about the hyperplasia of different level. The saying goes that bosom is the cradle of children, the garden of the woman, the symbol of female identity! Breast maintenance can make your figure more beautiful, add scores for your whole temperament, thus turn head rate increases 10 times certainly! Good figure and temperament will also ensure husband's loyalty!

The evolution of breast cancer: breast milk residual, secretions, toxins -- nodules and blockages -- lobular hyperplasia, ductal hyperplasia, breast hyperplasia -- duct blockage and adhesion -- fibroma, breast cyst -- breast cancer.

15.Q: Why do you need hip maintenance?

A: Because it can help the body lymphatic system discharge poison, improve the gynaecological disease, and avoid harm brought by buttock blockage. Hip impassability surely bring about gynaecological issues. What's inside buttocks is pelvic cavity, bowel. It connects daimai, lumbar vertebra, sciatic nerve from above, uterus ovary, uterine adenexa in front, anus, vagina, inguinal lymphatic from below.

If suffering from squeeze of shangjiao and xiaojiao meridian channels, hip would be most vulnerable to cold, wet and blood stasis. When hip is affected by cold, it can cause meridian contraction of hips, dysmenorrhea, irregular menstruation, dark menstruation color, clot, deterred blood flow and meridian channels. It also affects fertility.






16.Q: How does negative pressure detoxify?




A: Negative pressure can congest capillary, stimulate cell to increase vitality. The suction and release of pressure during operation can make local pores continue to open and close, therefore promote skin breathing, increase the amount of skin oxygen, speed up the elimination of waste. So it achieves the effects of promoting qi and blood circulation, dispelling coldness and dampness, dredging channels and collaterals, clearing toxins, heat and so on.

7. Packing List

1 x 40KHz RF Cavitation Head
1 x 9-Polar Bipolar RF Head for Body
1 x 12-Polar RF Head for Body
1 x Vacuum RF Head
1 x 5-Polar RF Head for Face and Eye Around
6 x Big paddle with 8 Light//Each 5mw
2 x Small paddle with 2 light//Each 5mw
2 x holders for light paddles
1 x Power Supply Cord (Both 110v and 220v voltage (frequency: 50/60Hz) version are available. Also we stock the USA, UK, AU, EU plug to suit for the country all over the world.)

8. Operational Diagrams

Parameter Adjustment	Product	Techniques	Diagrams
Facial Anti-aging: 60 min, once/week			
RF Energy: 3-8 BIO Energy: 2-8	Makeup remover+ Cleanser+ Toner+ Massage cream+ Essence+ Device+ Facial mask	<ol style="list-style-type: none"> 1. Remove makeups and cleanse face, 5 min. 2. Apply toner, 2 min. 3. Apply message cream evenly on face, 3 times. 4. Push the following acupoints (chengjiang point, dicang point, renzhong point, yingxiang point, jiache point, temple) , 3 times. 5. Lift 4 facial areas alternately with both hands in left-right order (chin-back of the ear, corner of mouth-ear gate, wing of nose-temple, corner of eye-hairline), 3 times. 6. Lift on 2 sides of face at the same time with both hands, 3 times. 7. Lift hands on one side alternately, 3 times. 8. Start from forehead, lift from middle of the eyebrows to hairline with middle and ring fingers, 3 times. 9. Overlap middle and ring fingers and move them in circular motions on forehead, 3 times. 10. Lift and caress treatment area alternately with both hands, 3 times. 	<p>Technique 4</p>  <hr/> <p>Technique 5</p>  <hr/> <p>Technique 8</p>  <hr/> <p>Technique 9</p>  <hr/> <p>Technique 10</p> 




	<p>11. Wash face clean, 2 min.</p> <p>12. Apply firming essence evenly on face, 2 min.</p> <p>13. BIO operation: start from chin, move the device in circular motions while middle and ring fingers move upwards in circles simultaneously on another side, 3 times.</p> <p>14. Exchange sides and repeat the previous operation, 3 times.</p> <p>15. Move the device in circles while move fingers in circles simultaneously on another side, 3 times.</p> <p>16. Exchange sides and repeat the previous operation.</p> <p>17. Lift face alternately with device on one side and hand on another side of the face, 3 times.</p> <p>18. Exchange side and repeat the previous operation.</p> <p>19. Lift towards hairline with both hands alternately on forehead, 3 times.</p> <p>20. RF operation: start from jaw and move in circles and make sure one line is next to another, 3 times.</p> <p>21. From jaw to temple, lift in lines and make sure one line is next to another, 3-5 times.</p> <p>22. Lift vertically towards hairline on forehead, 3-5 times.</p> <p>23. Use both hands to caress in four lines.</p>	<p>Technique 13,14, 15, 16, 20</p>  <hr/> <p>Technique 17, 18, 21</p>  <hr/> <p>Technique 19, 22</p>  <hr/>
--	---	---







		<p>24. Wash face clean, 2 min.</p> <p>25. Apply tightening essence evenly, 2 min.</p> <p>26. Apply hydrating&tightening mask, 15 min.</p> <p>27. Take off the mask and wash face clean, 2 min.</p> <p>28. Apply toner, essence, cream and sunscreen.</p>	
--	--	--	--

Suggested Treatment:

10 treatments account for a full course. Once done, skin feels tightened. After one course, skin is brightened, and skin barrier is repaired. 2 courses help increase elasticity, boost skin metabolism and detoxification, fade dark circles, brighten dull skin, and increase skin resistance and lower the sensitivity for sensitive skin. 3 courses help maintain skin condition, prevent skin sagging and aging, and define skin contour.

Anti-aging around Eyes: 60 min, once/week

<p>BIO Energy: 2-8</p>	<p>Makeup remover oil+ Face wash+ Toner+ Eye essence+ Device+ Eye mask</p>	<p>1. Remove makeup and clean face, 5 min.</p> <p>2. Apply toner, 1 min.</p> <p>3. Apply essence evenly around eyes and gently massage to rub it in, 1 min.</p> <p>4. Push the following acupoints(jingming point, cuanzhu point, yuyao point, sizhukong point, temple, tongziliao point, chengqi point) using middle and ring fingers, 3 times.</p> <p>5. Move hands in circular motion on eyes and push temples, 3 times.</p> <p>6. Move middle and ring fingers horizontally below lower eyelids and slide to back of the ear, 3 times.</p> <p>7. Lift middle and ring fingers from inner eye corner to upper eyelids and slide to back of the ear, 3</p>	<p>Technique 3, 5, 10, 11</p> 
			<p>Technique 4</p> 
			<p>Technique 6</p> 
			<p>Technique 7</p>






		<p>times.</p> <p>8. Lift eye corners with “scissor” hand gesture, 3-5 times.</p> <p>9. Repeat the previous operation on the other side, 3 times.</p> <p>10. Move hands in circular motion on eyes and push temples, 3-5 times.</p> <p>11. Move hands on eyes as if writing an “8” number, 3 times.</p> <p>12. Eye RF device operation: move the device in small circles from lower eyelids to temple, 3-5 times.</p> <p>13. Combined with hands: lift from lower eyelids to temple, 3-5 times.</p> <p>14. Lift the device from beginning till the end of eyebrow, 3-5 times.</p> <p>15. Move the device as if writing an “8” number on temple.</p> <p>16. Combined with hands: Lift from lower eyelids to temple, 3-5 times.</p> <p>17. Apply eye mask, 15 min.</p> <p>18. Take off masks and wash eye and face clean, 2 minutes.</p> <p>19. Apply toner, essence, eye essence, facial cream, and sunscreen.</p>	 <p>Technique 8, 9</p>  <p>Technique 12</p>  <p>Technique 13, 16</p>  <p>Technique 14</p>  <p>Technique 15</p> 
--	--	--	--




Suggested Treatment:

10 treatments account for one full course. Once done, the eyes are lifted and tightened, increasing blood circulation to the eyes. After one course, it lightens fine

lines, dark circles under the eyes and brightens the skin. After 3 courses, the skin around the eyes is firm, tender and shiny, improving and preventing the aging of the eyes.

Anti-aging around Neck: 60 min, once/week






<p>RF Energy: 3-8</p>	<p>Face Makeup remover oil+ Face wash+ Massage cream+ Essence+ Device+ Neck mask</p>	<ol style="list-style-type: none"> 1. Remove makeup and cleanse face, 5 min. 2. Apply toner, 1 min. 3. Neck massage: apply and rub oil on the chest with both hands to the back of the neck and point to feng chi and feng fu acupoint, 3 times. 4. Alternately lift the lower jaw and double chin to pass behind the ear and under the armpit, 3 times. 5. Push downwards with hand spread and hukou pointing down, tracing 3 meridian channels from the neck side to the armpit, 3 times. 6. Push downwards with 4 fingers, tracing 3 meridian channels from neck side to the armpit, 3 times. 7. Push downwards with kneeling fingers till the neck turns hot and slide to the armpit, 3 times. 8. Push downwards with hand spread and hukou pointing down from external collarbone to armpit, 3-5 times. 9. Do the same on the other side. 10. Wash neck clean. 2 min. 11. Apply essence evenly on neck, 1 min. 12. RF operation: coupled 	<p>Technique 3</p>  <hr/> <p>Technique 4</p>  <hr/> <p>Technique 5, 6, 7</p>  <hr/> <p>Technique 8</p>  <hr/> <p>Technique 12</p>  <hr/> <p>Technique 13</p>
-------------------------------	--	---	--







		<p>with hand massage, lift double chin to pass behind the ear and under the armpit, 3 times.</p> <p>13. Move the device in circles on one side of neck to pass under armpit, 3 times.</p> <p>14. Slide the device from external clavicle to the armpit, 3-5 times.</p> <p>15. Move the device in circles on the whole neck area, 3 times.</p> <p>16. Coupled with hand, lift neck in lines, 3-5 times.</p> <p>17. Do the same on the other side.</p> <p>18. Apply neck mask, 15 min.</p> <p>19. Wash neck clean, 2 min.</p> <p>20. Apply essence and neck cream.</p>	 <p>Technique 15</p>  <p>Technique 16</p> 
--	--	--	--




Suggested Treatment:

10 treatments account for one full course. After one treatment, skin becomes tender, smooth. After a full course, neck skin becomes smooth, speed up the neck lymphatic circulation. After 2 courses it tightens skin, fade neck lines, improve the double chin. After 3 course, skin is fine, firm and glowing, which also promote lymph drainage, improve facial dullness, some problems such as acne, and delay the aging of the skin, make skin restore young state.

Body-shaping on Waist&Abdomen: 60 min, once/week

<p>40K&RF Energy: 40K Mode: M1</p> <p>RF Energy Level: 3-8</p> <p>Vacuum&RF Energy: RF Energy Level: 3-8</p> <p>Suction Time: 0.3-1.5</p> <p>Release Time: 0-1.5</p> <p>Laser Paddle Energy Level: 30%-80%</p>	<p>Massage cream+ Ultrasonic gel+ Device</p>	<ol style="list-style-type: none"> 1. Apply and rub oil on abdomen with hands, 3 times. 2. Rub stomach back and forth with both hands, 3-5 times. 3. Knead abdomen with both hands using chiropractic techniques, 3 times. 4. Lift daimai on both sides of waist with both hands alternately, 16 times. 5. Move hands as if writing an "8" number on waist, 3 times. 6. Overlap hands and message the intestinal canal clockwise, 3 times. 7. Move hands in circular motion and caress the treatment area, 3 times. 8. Push on the following acupoints: shangwan zhongwan xiawan shenque qihai guanyuan zhongji tianshu dasheng qugu, 2 times. 9. Starting from zhongji point, push both thumbs to belly navel, slide them to waist and the lift upwards to groin, 3 times. 10. Caress the treatment area with hands and slide to groin. 11. 40K&RF device operation: starting from one side of the waist, lift the device to abdomen and 	<p>Technique 1, 7</p>  <hr/> <p>Technique 2</p>  <hr/> <p>Technique 3</p>  <hr/> <p>Technique 4</p>  <hr/> <p>Technique 5</p>  <hr/> <p>Technique 6</p>
--	--	--	---


		<p>groin, 3 times.</p> <p>12. Repeat the previous operation on another side, 3 times.</p> <p>13. Move the device on abdomen in small circular motions, 3 times.</p> <p>14. Move the device on abdomen in big circular motions, 3 times.</p> <p>15. Starting from daimai on one side of the waist, lift the device to another side of the waist back and forth, 3 times.</p> <p>16. Move the device on waist and abdomen vertically back and forth, 3 times.</p> <p>17. Vacuum&RF operation: set it to direct suction mode, lift the device on one side of the waist to abdomen, 3 times.</p> <p>18. Repeat the previous operation on another side.</p> <p>19. Lift daimai horizontally back and forth, 3 times.</p> <p>20. Lift the device on sides of the waist vertically back and forth to shape the body.</p> <p>21. Move the device on abdomen in big circular motion, 3 times.</p> <p>22. Lift the device from one side of the waist to abdomen, 3 times.</p> <p>23. Laser paddle operation: Strap on waist&abdomen, and fixate laser paddles on waist and abdomen, around 20min.</p> <p>24. Take off the laser</p>	 <p>Technique 8</p>  <p>Technique 9</p>  <p>Technique 10, 21</p>  <p>Technique 11, 12, 17, 18, 22</p>  <p>Technique 13</p>  <p>Technique 14</p>
--	--	---	--







		paddles and straps.	 <hr/> <p>Technique 15, 19</p>  <hr/> <p>Technique 16, 20</p> 
--	--	---------------------	--

Suggested Treatment:

Ten treatments account for one full course of treatment. Once done, abdomen is heating up, metabolism is accelerated. After one course, fat slowly reduce through metabolism, body starts to slim down. Constipation and other issues are addressed at the same time. After 2 courses the effect gradually becomes obvious, skin is tightened and lifted, skin collagen increased. Stretch marks and flabby skin are also improved. 3 courses help reinforce the effect. Redundant adipose gradually disappear, lumbar abdomen curve start to take in shape, and slim waist starts to show.

Warm Womb & Detoxification: 60 min, once/week

<p>RF Energy level: 3-8</p> <p>Vacuum&RF Energy: RF Energy Level: 3-8</p> <p>Suction Time: 0.3-1.5</p>	<p>Massage cream+ Device</p>	<ol style="list-style-type: none"> 1. Apply and rub oil on treatment area with hands, 3 times. 2. Rub stomach back and forth with both hands, 3-5 times. 3. Lift daimai on both sides of waist with both hands alternately, 16 times. 4. Move hands as if writing an "8" number on waist, 3 times. 	<p>Technique 1, 6</p>  <hr/> <p>Technique 2</p>
--	------------------------------	--	--

<p>Release Time: 0-1.5</p>		<p>5. Push channels and collaterals from top to bottom with hands alternately (ren mai -- open 2 inches of kidney channels -- open 4 inches of stomach channels -- open 6 inches of spleen channels -- open 8 inches of liver channels -- open gallbladder channels on both sides) for 3 times)</p> <p>6. Move hands in circular motion and caress the treatment area, 3 times.</p> <p>7. Rub hands together until they're warm and put them on lower abdomen.</p> <p>8. RF device operation: starting from one side of the waist, lift the device to abdomen and groin, 3 times.</p> <p>9. Move the device on abdomen in small circular motions, 3 times.</p> <p>10. Move the device on abdomen in big circular motions, 3 times.</p> <p>11. Push the device on abdomen through channels in bottom-up order to groin, 3 times.</p> <p>12. Combined with hands: lift daimai on two sides of the waist, 16 times.</p> <p>13. Move the device on abdomen in small circles, 3-5 times.</p> <p>14. Vacuum&RF operation is the the same as the above.</p>	<div style="text-align: center;">  </div> <hr/> <p style="text-align: center;">Technique 3, 12</p> <div style="text-align: center;">  </div> <hr/> <p style="text-align: center;">Technique 4</p> <div style="text-align: center;">  </div> <hr/> <p style="text-align: center;">Technique 5, 11</p> <div style="text-align: center;">  </div> <hr/> <p style="text-align: center;">Technique 8</p> <div style="text-align: center;">  </div> <hr/> <p style="text-align: center;">Technique 9, 13</p> <div style="text-align: center;">  </div> <hr/> <p style="text-align: center;">Technique 10</p>
--------------------------------	--	--	---



Suggested Treatment:

10 treatments account for a full course. Once done, abdominal temperature increases, meridian dredges. After a course, metabolism speeds up, constipation and dysmenorrhea improves. After 2 courses, menstrual flow returns to normal, cold feet improves. 3 courses help strengthens the function of warm womb, regulate the female reproductive system problems, reduce leucorrhea, increase vaginal secretion, delay menopause, delay aging, and skin giving out a healthy complexion.


Body-shaping on Breasts: 45 min, once/week

RF Energy:
3-8


Massage
cream+
Device

1. Stand at the head of the bed, both hands from shan zhong point began to rub oil to armpit up at the same time lift the suspension ligament (this is to caress treatment area) 3 times.
2. Both hands thumb point: shanzhong, rugen, dabao, yingchuang, zhongfu, yunmen, 3 times.
3. Caress the treatment area, 3 times.
4. Overlap both palms from shan zhong point to chest both sides back and forth as if writing an "8" number, 3 times.
5. Lift the accessory breast to the suspensory ligament alternately with both hands, first left, then right, 10 times.
6. Caress the treatment area, 3 times.
7. Sit beside the customer, first left, then right, with


Technique 1, 3, 6, 10









Technique 2





Technique 4



Technique 5

		<p>both hands to soothe and lift the chest, 3-5 times.</p> <p>8. Alternately push the breast duct around the chest (one circle of the chest), 3-5 times.</p> <p>9. Use the thumb to dredge the nodules in the chest, 3 times.</p> <p>10. Lift the whole chest with hands, 3-5 times.</p> <p>11. The technique on the other side is the same as above.</p> <p>12. RF operation: first left, then right, from the bottom to the top along the direction of the chest to the nipple, combined with hand pacification, 5-8 times.</p> <p>13. Draw half a circle along the chest to lift and pull from the bottom to the collarbone, combining with hand comfort, 5-8 times.</p> <p>14. Make small circles to clear the nodule site for 3-5 times.</p> <p>15. The customer lies on his side and stands up with her arms. The device burns fat in circles at the position of auxiliary breast, 5-8 times.</p> <p>16. The instrument is pushed from the position of underarm auxiliary breast to the chest (shaping and receiving auxiliary breast), 5-8 times</p> <p>17. The technique on the other side is the same as above.</p>	 <hr/> <p>Technique 7, 13</p>  <hr/> <p>Technique 8</p>  <hr/> <p>Technique 9</p>  <hr/> <p>Technique 12</p>  <hr/> <p>Technique 14</p>  <hr/> <p>Technique 15</p>
--	--	--	--

			 <hr/> <p>Technique 16</p>  <hr/>
--	--	--	---


Suggested Treatment:

10 times account for a course of treatment. After a full course, breasts will experience heat, blood circulation is boosted, and a certain improvement effect will be seen. After a course of treatment the effect is more obvious, nodules are improved. 3 courses help tighten the skin and consolidate shaping effect, and increase the elasticity of the breasts. Regulating endocrine makes women more attractive.







Body-Shaping on Arms: 60 min, once/week



RF Level: 3-8	Energy	Massage cream+ Ultrasonic gel+	1. Left-right order: lay the arm flatwise, apply and rub oil in from lower arm to the entire arm and slide to the	Technique 1, 2, 3, 4, 5, 6
------------------	--------	-----------------------------------	---	----------------------------

<p>40K&RF Energy: 40K Mode: M1</p> <p>RF Energy Level: 3-8</p> <p>Laser paddle energy level: 30%-80%</p>	<p>Device</p>	<p>fingers, 3 times.</p> <p>2. Push the entire arm with both palms, 3 times.</p> <p>3. Push hands upwards tracing three main collateral channels(dachangjing-sanjiaojing-xiaochangjing) to armpits with fingers spread and hukou pointing upwards, 3 times.</p> <p>4. Caress the treatment area, 3 times.</p> <p>5. Trace and rub three main collateral channels on arms with kneeling finger back and forth till them turn hot, 3 times.</p> <p>6. Caress the treatment area, 3 times.</p> <p>7. Lay inner arm upwards, and push hands tracing three yin channel on inner arm(feijing-xinbaojing-xinjing) to armpit, with hukou pointing upwards, 3 times respectively.</p> <p>8. Rub three channels back and forth with kneeling fingers, 3 times.</p> <p>9. Caress the treatment area and slide to fingers.</p> <p>10. The end of massage techniques.</p> <p>11. Repeat the previous operation on another side.</p> <p>12. RF device operation: lay arms flatwise and start from fat part of lower arm, tracing three channels to armpit, 3 times.</p> <p>13. Move the device in annular motion, tracing three channels to armpit, 3 times.</p>	<div data-bbox="1129 212 1273 421" data-label="Image"> </div> <p>Technique 7, 8</p> <div data-bbox="1109 504 1292 712" data-label="Image"> </div> <p>Technique 9</p> <div data-bbox="1109 795 1292 1003" data-label="Image"> </div> <p>Technique 12, 15</p> <div data-bbox="1129 1086 1273 1294" data-label="Image"> </div> <p>Technique 13, 14</p> <div data-bbox="1129 1377 1273 1585" data-label="Image"> </div> <p>Technique 16</p> <div data-bbox="1109 1668 1292 1877" data-label="Image"> </div> <p>Technique 18</p>
--	---------------	--	---

		<p>14. You may double the operation on flabby arms, 3 times.</p> <p>15. Trace three channels to armpit, 3 times.</p> <p>16. Lay inner arm flatwise and push three channels on upper arm to armpit, 3 times.</p> <p>17. Trace three channels on upper arm and move the device in circles to armpit, 3 times.</p> <p>18. Trace three channels from upper arm to armpit.</p> <p>19. Do the same on the other side.</p> <p>20. Repeat the previous operation with 40K&RF operation.</p> <p>21. Do the same on the other side.</p> <p>22. Laser paddle operation: Strap on arms, and fixate laser paddles on fat part of arms, around 20 min.</p> <p>23. Take off the laser paddles and straps.</p>	
<p>Suggested Treatment: 10 treatments account for a full course. Once done, it will feel warm and tight, and the arms will be relaxed. A course of treatment help reduce excessive arm fat and the skin appear to be tight and plump. 2 courses of treatment began to shape the arms, strengthen the effect, dredge meridians, and enhance the physique. 3 courses help consolidate stability and prevent rebound.</p>			
<p>Body-Shaping on Back: 60 min, once/week</p>			
<p>RF energy: 3-8</p> <p>Vacuum RF energy: RF Energy</p>	<p>Essential oil+ Device</p>	<p>1. Apply oil on back and press Fengchi and fengfu points.</p> <p>2. (Starting from hairline) stroke badanjin with thumb, 3-5 times.</p>	<p>Technique 1, 10</p>

<p>Level: 3-8</p> <p>Suction: 0.3-1.5</p> <p>Release: 0-1.5</p> <p>Laser paddle energy level: 30%-80%</p>		<p>3. Move thumb outwards from pangguangjin to baliao area and then to fengchi and fengfu point.</p> <p>4. Starting from neck to caudal vertebra, move hands in S-shape, 3 times.</p> <p>5. Push pangguangjin with thumb in left-right order, 3 times.</p> <p>6. Push pangguangjin to baliao area with both thumbs, at the same time, 3 times.</p> <p>7. Push pangguangjin in three kneeling fingers, 3 times.</p> <p>8. Push scapula slot with both hands alternately in left-right order, 3-6 times.</p> <p>9. Push scapula slot with both hands horizontally, 3-6 times.</p> <p>10. Caress the back with both hands and push fengchi and fengfu points, 3 times.</p> <p>11. Overlap both thumbs and push tianzong point and slide to the arm, 3 times.</p> <p>12. Repeat the previous action on another arm, 3 times.</p> <p>13. RF Operation: Starting from dumai to pangguangjin, move the device from neck to baliao area, 3-5 times.</p> <p>14. Move the device in circles on dazhui point, and baliao area, 3-5 times respectively.</p> <p>15. From dabanjin to neck,</p>	<div data-bbox="1129 215 1273 421" data-label="Image"> </div> <p>Technique 2, 15, 25</p> <div data-bbox="1129 506 1273 712" data-label="Image"> </div> <p>Technique 3</p> <div data-bbox="1129 797 1273 1003" data-label="Image"> </div> <p>Technique 4</p> <div data-bbox="1129 1088 1273 1294" data-label="Image"> </div> <p>Technique 5, 6, 7</p> <div data-bbox="1129 1379 1273 1585" data-label="Image"> </div> <p>Technique 8, 16, 26</p> <div data-bbox="1129 1671 1273 1877" data-label="Image"> </div> <p>Technique 9</p>
---	--	---	--



		<p>move the device back and forth, 3-5 times.</p> <p>16.Lift the device on scapula slot in left-right order, 3-5 times.</p> <p>17.Move the device horizontally as if writing an “8” number to baliao area, 3 times.</p> <p>18.Move the device upwards in annular motion to armpit in left-right order on both sides, 3 times.</p> <p>19.Starting from neck, slide the device from dumai and pangguangjin to baliao area, 3-5 times.</p> <p>20.Negative Pressure RF Operation: starting from neck, slide the device from dumai and pangguangjin to baliao area.</p> <p>21.Move the device in circular movetion on dazhui point and baliao area, 3 times.</p> <p>22.Combine device with hands, lift the device and hands from armpits to dumai and waist, 3 times.</p> <p>23.Move the device from bottom back to armpits in annular motion and left-right order, 3 times.</p> <p>24.Lift the device upwards to armpits on sides of the waist, 3-5 times.</p> <p>25.Lift dabanjin horizontally back and forth, 3-5 times.</p> <p>26.Lift scapula slots back and forth, 3-5 times.</p> <p>27.Laser paddle operation: Strap on back, and fixate laser paddles on fat part of</p>	 <hr/> <p>Technique 11, 12</p>  <hr/> <p>Technique 13, 19, 20</p>  <hr/> <p>Technique 14, 21</p>  <hr/> <p>Technique 17</p>  <hr/> <p>Technique 18, 23</p>  <hr/> <p>Technique 24</p>
--	--	---	--






		<p>the back, around 20 min. 28. Take off the laser paddles and straps.</p>	 <hr/> <p>Technique 22</p> 
--	--	--	---





Suggested Treatment:







10 times account for a full course. Once done, back feels significantly relaxed, neck stiffness alleviated, and skin feels lifted. After a full course, shoulder feels relaxed, back fat lessens, the lump on the nape of the neck improves, and the back line starts to show. 2 courses help dredge meridians, regulate sub-health, and improve sleep. 3 courses help achieve thin back, curvy back line, better viscera function and digestion, detoxification and metabolism. (one course helps reduce size, two courses help strengthen the effects, three courses help with consolidation and avoid rebound.)




Body-Shaping on Buttocks: 60 min, once a week

<p>RF Energy Level: 3-8 40K&RF Energy: 40K Mode: M1 RF Energy Level: 3-8</p>	<p>Massage cream+ Ultrasonic gel+ Device</p>	<ol style="list-style-type: none"> 1. Standing on the side, accumulate oil on hands and then slide to the waist and then lift up along the hips from the waist, 3 times (this is to caress the treatment area.) 2. Repeatedly push ba liao zone with two thumbs. 3. Caress the treatment area for 3 times and then point: "shenyu", "ba liao", "chang qiang", "huan tiao" and "cheng fu", 3 times. 4. Caress the treatment area, 3 times. 5. First left, then right, with both hands extending from the root of the thigh from 	<p>Technique 1, 4, 10</p>  <hr/> <p>Technique 2</p>  <hr/> <p>Technique 3</p>
--	--	---	---

		<p>bottom to top -- bladder tendon -- kidney -- liver -- gallbladder -- vein, 3 times each.</p> <p>6. Both hands follow from the root of thigh from bottom to top -- bladder channel, kidney channel, liver channel and gallbladder channel to the vein, 3 times each.</p> <p>7. Overlap both palms and push up together from thigh root to pulse (ascension) ,3 to 5 times.</p> <p>8. Push your hands from bottom to top from both sides of your hips to the top of your hips, 3 times.</p> <p>9. Repeat the previous operation.</p> <p>10. Caress the treatment area.</p> <p>11. The technique on the other side is the same as above.</p> <p>12. Operation of radio frequency instrument: lift from the thigh root to the vein line by line, 3 times.</p> <p>13. Lift 3 times from bottom to top on both sides of the buttocks to the highest point of the buttocks.</p> <p>14. Circle your hips, 3 to 5 times.</p> <p>15. Lift 3 times from thigh root to vein.</p> <p>16. Lift from bottom to top on both sides of the buttocks to the highest point of the buttocks, 3 times.</p> <p>17. The technique on the</p>	 <hr/> <p>Technique 5, 12, 15</p>  <hr/> <p>Technique 6</p>  <hr/> <p>Technique 8, 13, 16</p>  <hr/> <p>Technique 14</p>  <hr/>
--	--	---	---

		<p>other side is the same as above.</p> <p>18. Ultrasonic operation method is the same as above.</p>	
<p>Suggested Treatment:</p> <p>10 times account for a full course. After one treatment, buttocks feel lifted and heated. After a course, hip curve starts to show, excess fat slowly disappears, skin becomes tight. 3 courses help improve gynecological problems like cold body and cold uterus. Menstruation returns to normal. It helps consolidate the effect and increase the female sex appeal.</p>			
<p>Body-Shaping on Legs: 70 min, once/week</p>			
<p>Vacuum&RF Energy: RF Energy Level: 3-8 Suction: 0.3-8 Release: 0-1.5</p> <p>40K&RF Energy: 40K Mode: M1</p> <p>RF Energy Level: 3-8</p> <p>Laser paddle energy level: 30%-80%</p>	<p>Essential oil+ Ultrasonic gel+ Device</p>	<ol style="list-style-type: none"> 1. Rear leg, left-right order: apply and rub oil in from calf to thigh to heel, 3 times. 2. Move hands upwards to four main collateral channels(panguangjin-she njin-ganjin-danjin) with fingers spread and hukou pointing upwards, 3 times. 3. Push hands alternately to popliteal fossa, 3 times. 4. Twist both hands upwards alternately as if twisting a fried dough, 3 times. 5. Caress the treatment area, 3 times. 6. Push four main collateral channels upwards with kneeling fingers of both hands, 3 times. 7. Caress the treatment area, 3 times. 8. Do the same on the other side. 9. 40K operation: move the device from popliteal fossa to 	<p>Technique 1, 2, 3, 5, 6, 7, 8, 9</p>  <hr/> <p>Technique 4</p>  <hr/> <p>Technique 10, 17</p>  <hr/> <p>Technique 11, 16, 18</p> 

		<p>pangguangjin-shenjin-ganjin-danjin, 3 times.</p> <p>10. Move the device in circular motion on calf, 3 times.</p> <p>11. Move the device upwards to popliteal fossa and caress the treatment area with hands, 3 times.</p> <p>12. Push from popliteal fossa tracing 4 main collateral channels to the end of thigh, 3 times.</p> <p>13. Move the device in circular motion from popliteal fossa to the end of thigh, 3 times.</p> <p>14. Push from popliteal fossa tracing main collateral channels to the end of thigh, 3 times.</p> <p>15. Move the device in circular motion on thigh, 3-5 times.</p> <p>16. Vacuum&RF operation(combined with hands): lift the device from heel to popliteal fossa, 3 times.</p> <p>17. Move the device on calf in circular motions, 3 times.</p> <p>18. Move the device upwards and downwards on calf, 3-5 times.</p> <p>19. Lift the device from popliteal fossa to the end of thigh, 3 times.</p> <p>20. Lift the device from popliteal fossa to the end of thigh, 3 times.</p> <p>21. Lift the device from sides of leg to the middle, 3-5 times.</p> <p>22. Lift the device in</p>	<p>Technique 12, 14, 19</p>  <hr/> <p>Technique 13, 15, 20</p>  <hr/> <p>Technique 21</p>  <hr/> <p>Technique 22</p>  <hr/> <p>Technique 24, 25, 26, 29, 32, 34</p>  <hr/> <p>Technique 30, 35</p>  <hr/> <p>Technique 31</p>
--	--	---	--


		<p>annular motion on both sides of the leg, 3 times.</p> <p>23. Repeat the previous operation on another leg.</p> <p>24. Foreleg: rub oil with both hands from heel to the thigh root (namely to caress the treatment area), 3 times.</p> <p>25. Move hands upwards to four main collateral channels(pijing-weijing-ganjing-danjing) with fingers spread and hukou pointing upwards, 3 times.</p> <p>26. Push four main collateral channels upwards with kneeling fingers of both hands, 3 times.</p> <p>27. 40K operation: Lift the device towards knee while tracing 3 collaterals, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>28. Move the device on two sides of calf in annular motion, 3 times(start from thigh if there's not enough unwanted fat on calf.)</p> <p>29. Lift the device from knee to thigh root in lines, 3 times.</p> <p>30. Move the device in small circles on thigh, 3 times.</p> <p>31. Lift the device from knee to thigh root in annular motion, 3 times.</p> <p>32. Lift the device from knee to thigh root in lines, 3 times.</p> <p>33. Vacuum&RF operation: calf operation is the same</p>	 <hr/> <p>Technique 36</p>  <hr/> <p>Technique 37</p>  <hr/>
--	--	---	--







		<p>as the above. (start from thigh if there's not enough unwanted fat on calf.)</p> <p>34. Lift the device from knee to thigh root in lines, 3 times.</p> <p>35. Move the device in small circles on thigh, 3 times.</p> <p>36. Lift the device from knee to thigh root in lines, 3 times.</p> <p>37. Lift the device from knee to thigh root in annular motion, 3 times.</p> <p>38. Repeat the previous operation on the other side.</p> <p>39. Laser paddle operation: Strap on legs, and fixate laser paddles on fat part of the legs, around 20 min.</p> <p>40. Take off the laser paddles and straps.</p>	
--	--	--	--







Suggested Treatment:







10 treatments account for a full course. After one treatment, leg is relaxed, and circulation is accelerated. After a course of treatment, leg becomes thin, skin collagen tissue experience contraction due to heat effect, and tightening effect can obviously be felt on the skin. After 2 courses, leg is tight and slim, and obvious results are seen. 3 courses help consolidate the effect. Super RF + super ultrasound is to stimulate the dermis of the skin, which can sustain collagen regeneration, so the treatment effect will be more and more obvious.







Body-Sculpting&Slimming: 180 min, once/week

RF Energy Level: 3-8	Essential oil (massage cream) + Gel + Device	<p>1. Apply oil on back and press Fengchi and fengfu points.</p> <p>2. (Starting from hairline) stroke badanjin with thumb, 3-5 times.</p> <p>3. Move thumb outwards from panguangjin to baliao area and then to fengchi</p>	<p>Technique 1,10</p> 
Vacuum RF Energy:			Technique 2,15,25
RF Energy Level:			

<p>3-8</p> <p>Suction time: 0.3-1.5</p> <p>Release time: 0-1.5</p> <p>40K&RF Energy: 40K Mode: M1</p> <p>RF Energy Level: 3-8</p> <p>Laser paddle energy level: 30%-80%</p>		<p>and fengfu point, 3 times.</p> <p>4. Starting from neck to caudal vertebra, move hands in S-shape, 3 times.</p> <p>5. Push pangguangjin with thumb in left-right order, 3 times.</p> <p>6. Push pangguangjin to baliao area with both thumbs, at the same time, 3 times.</p> <p>7. Push pangguangjin in three kneeling fingers, 3 times.</p> <p>8. Push scapula slot with both hands alternately in left-right order, 3-6 times.</p> <p>9. Push scapula slot with both hands horizontally, 3-6 times.</p> <p>10. Caress the back with both hands and push fengchi and fengfu points, 3 times.</p> <p>11. Overlap both thumbs and push tianzong point and slide to the arm, 3 times.</p> <p>12. Repeat the previous action on another arm, 3 times.</p> <p>13. RF Operation: Starting from dumai to pangguangjin, move the device from neck to baliao area, 3-5 times.</p> <p>14. Move the device in circles on dazhui point, and baliao area, 3-5 times respectively.</p> <p>15. From dabanjin to neck, move the device back and forth, 3-5 times.</p> <p>16. Lift the device on</p>	 <p>Technique 3</p>  <p>Technique 4</p>  <p>Technique 5,6,7</p>  <p>Technique 8,16,26</p>  <p>Technique 9</p>  <p>Technique 11,12</p>
---	--	---	--

		<p>scapula slot in left-right order, 3-5 times.</p> <p>17. Move the device horizontally as if writing an "8" number to baliao area, 3 times.</p> <p>18. Move the device upwards in annular motion to armpit in left-right order on both sides, 3 times.</p> <p>19. Starting from neck, slide the device from dumai and panguangjin to baliao area, 3-5 times.</p> <p>20. Vacuum&RF Operation: starting from neck, slide the device from dumai and panguangjin to baliao area, 3-5 times.</p> <p>21. Move the device in circular movetion on dazhui point and baliao area, 3 times.</p> <p>22. Combine device with hands, lift the device and hands from armpits to dumai and waist, 3 times.</p> <p>23. Move the device from bottom back to armpits in annular motion and left-right order, 3 times.</p> <p>24. Lift the device upwards to armpits on sides of the waist, 3-5 times.</p> <p>25. Lift dabanjin horizontally back and forth, 3-5 times.</p> <p>26. Lift scapula slots back and forth, 3-5 times.</p> <p>27. Laser paddle operation: Strap on the back and fixate the laser paddle on the fat part of the back, around 20 min.</p>	 <hr/> <p>Technique 13,19,20</p>  <hr/> <p>Technique 14,21</p>  <hr/> <p>Technique 17</p>  <hr/> <p>Technique 18,23</p>  <hr/> <p>Technique 24</p>  <hr/> <p>Technique 29,31,32,33,34,35,36,39</p>
--	--	---	---

		<p>28. Take off the paddles and straps.</p> <p>29. Leg massage, left-right order: apply and rub oil in from calf to thigh to heel, 3 times.</p> <p>30. Move hands upwards to four main collateral channels(pangguangjin-she njin-ganjin-danjin) with fingers spread and hukou pointing upwards, 3 times.</p> <p>31. Push hands alternately to popliteal fossa, 3 times.</p> <p>32. Twist both hands upwards alternately as if twisting a fired dough, 3 times.</p> <p>33. Caress the treatment area, 3 times.</p> <p>34. Push four main collateral channels upwards with kneeling fingers of both hands, 3 times.</p> <p>35. Caress the treatment area, 3 times.</p> <p>36. Repeat the previous operation on another side.</p> <p>37. 40K&RF device operation: move the device from popliteal fossa to pangguangjin-shenjin-ganjin-danjin, 3 times.</p> <p>38. Move the device in circular motion on calf, 3 times.</p> <p>39. Move the device upwards to popliteal fossa and caress the treatment area with hands, 3 times.</p> <p>40. Push from popliteal fossa tracing 4 main collateral channels to the end of thigh, 3 times.</p>	 <hr/> <p>Technique 4</p>  <hr/> <p>Technique 38,45</p>  <hr/> <p>Technique 37,40,44</p>  <hr/> <p>Technique 42,44,47</p>  <hr/> <p>Technique 41,43. 48</p>  <hr/> <p>Technique 49</p>
--	--	--	---

	<p>41. Move the device in circular motion from popliteal fossa to the end of thigh, 3 times.</p> <p>42. Push from popliteal fossa tracing main collateral channels to the end of thigh, 3 times.</p> <p>43. Move the device on thigh in circular motions, 3-5 times.</p> <p>44. Vacuum&RF operation(combined with hands): lift the device from heel to popliteal fossa, 3 times.</p> <p>45. Move the device on calf in circular motions, 3 times.</p> <p>46. Move the device upwards and downwards on calf, 3-5 times.</p> <p>47. Lift the device from popliteal fossa to the end of thigh, 3 times.</p> <p>48. Move the device in circular motions on thigh, 3 times.</p> <p>49. Lift the device from sides of leg to the middle, 3-5 times.</p> <p>50. Lift the device in annular motion on both sides of the leg, 3 times.</p> <p>51. Repeat the previous operation on another leg.</p> <p>52. Laser paddle operation: Strap on legs, and fixate the paddles on the fat part of the legs, around 20 min.</p> <p>53. Take off the paddles and straps.</p> <p>54. Arms massage, left-right order: lay the arm</p>	 <hr/> <p>Technique 50</p>  <hr/> <p>Technique 54-59</p>  <hr/> <p>Technique 60,61</p>  <hr/> <p>Technique 62</p>  <hr/> <p>Technique 65,68</p>  <hr/> <p>Technique 66</p>
--	--	---

flatwise, apply and rub oil in from lower arm to the entire arm and slide to the fingers, 3 times.

55. Push the entire arm with both palms, 3 times.

56. Push hands upwards tracing three main collateral channels(dachangjing-sanjiaojing-xiaochangjing) to armpits with fingers spread and hukou pointing upwards, 3 times.

57. Caress the treatment area, 3 times.

58. Tracing and rubbing three main collateral channels on arms with kneeling finger back and forth till them turn hot, 3 times.

59. Caress the treatment area, 3 times.

60. Lay inner arm upwards, and push hands tracing three yin channel on inner arm(feijing-xinbaojing-xinjing) to armpit, with hukou pointing upwards, 3 times respectively.


61. Rub three channels back and forth with kneeling fingers, 3 times.

62. Caress the treatment area and slide to fingers.

63. The end of massage techniques.

64. Repeat the previous operation on another side.


65. RF device operation: lay arms flatwise and start from fat part of lower arm, tracing three channels to armpit, 3 times.




Technique 69,70,71




Technique 76




Technique 77









Technique 78







Technique 79



Technique 80

		<p>66. Move the device in annular motion, tracing three channels to armpit, 3 times.</p> <p>67. You may double the operation on flabby arms, 3 times.</p> <p>68. Trace three channels to armpit, 3 times.</p> <p>69. Lay inner arm flatwise and push three channels on upper arm to armpit, 3 times.</p> <p>70. Trace three channels on upper arm and move the device in circles to armpit, 3 times.</p> <p>71. Trace three channels from upper arm to armpit.</p> <p>72. Repeat the previous operation with 40K&RF operation.</p> <p>73. Do the same on the other side.</p> <p>74. Laser paddle operation: Strap on the arms and fixate the paddles on the fat part of the arms, around 20 min.</p> <p>75. Take off the paddles and straps.</p> <p>76. Waist and abdomen: apply and rub oil on abdomen with hands, 3 times.</p> <p>77. Rub stomach back and forth with both hands, 3-5 times.</p> <p>78. Knead abdomen with both hands using chiropractic techniques, 3 times.</p> <p>79. Lift daimai on both sides of waist with both hands alternately, 16 times.</p>	 <hr/> <p>Technique 81</p>  <hr/> <p>Technique 83</p>  <hr/> <p>Technique 84</p>  <hr/> <p>Technique 85</p>  <hr/> <p>Technique 86,92,97</p>  <hr/> <p>Technique 88</p>
--	--	--	--







		<p>80. Move hands as if writing an “8” number on waist, 3 times.</p> <p>81. Overlap hands and message the intestinal canal clockwise, 3 times.</p> <p>82. Move hands in circular motion and caress the treatment area, 3 times.</p> <p>83. Push on the following acupoints: shangwan zhongwan xiawan shenque qihai guanyuan zhongji tianshu dasheng qugu, 2 times.</p> <p>84. Starting from zhongji point, push both thumbs to belly navel, slide them to waist and the lift upwards to groin, 3 times.</p> <p>85. Caress the treatment area with hands and slide to groin.</p> <p>86. 40K&RF device operation: starting from one side of the waist, lift the device to abdomen and groin, 3 times.</p> <p>87. Repeat the previous operation on another side, 3 times.</p> <p>88. Move the device on abdomen in small circular motions, 3 times.</p> <p>89. Move the device on abdomen in big circular motions, 3 times.</p> <p>90. Starting from daimai on one side of the waist, lift the device to another side of the waist back and forth, 3 times.</p> <p>91. Move the device on waist and abdomen</p>	 <hr/> <p>Technique 89,96</p>  <hr/> <p>Technique 90,94</p>  <hr/> <p>Technique 91,95</p>  <hr/>
--	--	--	---







		<p>vertically back and forth.</p> <p>92. Vacuum&RF operation: set it to direct suction mode, lift the device on one side of the waist to abdomen, 3 times.</p> <p>93. Repeat the previous operation on another side.</p> <p>94. Lift daimai horizontally back and forth, 3 times.</p> <p>95. Lift the device on sides of the waist vertically back and forth to shape the body.</p> <p>96. Move the device on abdomen in big circular motion, 3 times.</p> <p>97. Lift the device from one side of the waist to abdomen, 3 times.</p> <p>98. Laser paddle operation: Strap on waist and abdomen, and fixate the paddles on waist and abdomen, around 20 min.</p> <p>99. Take off the paddles and straps.</p>	
--	--	--	--







Suggested Treatment:

10 times account for a full course of treatment. After one treatment, your body will feel relaxed, blood circulation boosted, and skin tightened. After a course you will see less fat on different parts of the body, collagen regenerated, and firm skin. After 2 courses fat reduces, and curve lines begin to show. 3 courses will help with enhance the effects and shape up, and S curve will appear.

Overall facial Care: 100 minutes, once/week

<p>Face Energy: 3-8</p>	<p>RF Makeup remover oil+</p>	<p>1. Remove makeup and clean, 5 min. 2. Apply toner, 2 min. 3. Apply hydrating facial mask, 10 min.</p>	<p>Technique 6</p> 
<p>Eye Energy: 3-8</p>	<p>RF Hydrating facial mask+ Device+ Essence+ Facial mask</p>	<p>4. Wash face clean, 2 min. 5. Apply message cream evenly on face, 3 times. 6. Push the following acupoints (chengjiang point, dicang point, renzhong point, yingxiang point, jiache point, temple) , 3 times.</p>	<p>Technique 7,12,36,37</p> 
		<p>7. Lift 4 facial areas alternately with both hands in left-right order (chin-back of the ear, corner of mouth-ear gate, wing of nose-temple, corner of eye-hairline), 3 times.</p>	<p>Technique 10,38</p> 
		<p>8. Lift on 2 sides of face at the same time with both hands, 3 times.</p>	<p>Technique 11</p> 
		<p>9. Lift hands on one side alternately, 3 times.</p>	<p>Technique 14</p> 
		<p>10. Start from forehead, lift from middle of the eyebrows to hairline with middle and ring fingers, 3 times.</p>	<p>Technique 15,20,21</p> 
		<p>11. Overlap middle and ring fingers and move them in circular motions on forehead, 3 times.</p>	
		<p>12. Lift and caress treatment area alternately with both hands, 3 times. 13. Caress the whole face. 14. Push the following acupoints(jingming point, cuanzhu point, yuyao point,</p>	

		<p>sizhukong point, temple, tongziliao point, chengqi point) using middle and ring fingers, 3 times.</p> <p>15. Move hands in circular motion on eyes and push temples, 3 times.</p> <p>16. Move middle and ring fingers horizontally below lower eyelids and slide to back of the ear, 3 times.</p> <p>17. Lift middle and ring fingers from inner eye corner to upper eyelids and slide to back of the ear, 3 times.</p> <p>18. Lift eye corners with “scissor” hand gesture, 3-5 times.</p> <p>19. Do the same on the other side, 3 times.</p> <p>20. Move hands in circular motion on eyes and push temples, 3-5 times.</p> <p>21. Move hands on eyes as if writing an “8” number, 3 times.</p> <p>22. Caress the whole and slide to the back of the ear, 3 times.</p> <p>23. Neck massage: apply and rub oil on the chest with both hands to the back of the neck and point to feng chi and feng fu acupoint, 3 times.</p> <p>24. Alternately lift the lower jaw and double chin to pass behind the ear and under the armpit, 3 times.</p> <p>25. Push downwards with hand spread and hukou pointing down, tracing 3 meridian channels from the neck side to the armpit, 3 times.</p>	 <p>Technique 16</p>  <p>Technique 17</p>  <p>Technique 18</p>  <p>Technique 23</p>  <p>Technique 24,46</p>  <p>Technique 25,26,27,47</p>
--	--	--	--

		<p>26. Push downwards with 4 fingers, tracing 3 meridian channels from neck side to the armpit, 3 times.</p> <p>27. Push downwards with kneeling fingers till the neck turns hot and slide to the armpit, 3 times.</p> <p>28. Push downwards with hand spread and hukou pointing down from external collarbone to armpit, 3-5 times.</p> <p>29. Do the same on the other side.</p> <p>30. Wash face and neck clean, 5 min.</p> <p>31. Apply (anti-aging)essence evenly on face and eyes, 1 min.</p> <p>32. Face RF operation: move the probe upwards from jaw to forehead in circles; move in 4 lines and at the same time circle face with middle and ring fingers, 3 times.</p> <p>33. Exchange hands and do the same on the other side, 3 times.</p> <p>34. Circle in 4 lines with the probe and use the other free hand to fixate the facial area, 3 times.</p> <p>35. Do the same on the other side.</p> <p>36. Alternate hands in lifting face in four lines; use one hand to lift face while the other hold the machine to lift the treated area, 3 times.</p> <p>37. Do the same on the other side.</p> <p>38. Alternate hands in lifting vertically towards hairline to</p>	 <hr/> <p>Technique 32,33,34,35</p>  <hr/> <p>Technique 39</p>  <hr/> <p>Technique 40</p>  <hr/> <p>Technique 41</p>  <hr/> <p>Technique 43</p>  <hr/>
--	--	--	---

		<p>lift the treatment area, 3 times.</p> <p>39. Eye RF operation: move the device in small circles from lower eyelids to temple, 3-5 times.</p> <p>40. Coupled with hand, lift canthus towards hairline, 3-5 times.</p> <p>41. Lift the device from beginning till the end of eyebrow, 3-5 times.</p> <p>42. Move the probe on temple as if drawing number "8".</p> <p>43. Coupled with hand, lift canthus towards hairline, 3-5 times.</p> <p>44. Do the same on the other side.</p> <p>45. Apply (anti-aging)essence evenly on neck, 1 min.</p> <p>46. Face RF operation: Combined with hand, move the device to lift double chin, pass by the back of the ear, then slide to armpit, 3 times.</p> <p>47. Move the probe in 3 circular lines on the side of the neck and slide to armpit, 3 times.</p> <p>48. Slide the probe on inner and outer collarbone, then to armpit, 3-5 times.</p> <p>49. Move the probe on the neck in circular motion, 3 times.</p> <p>50. Combined with hand, move the probe in lines with one line next to another to lift the whole neck, 3-5 times.</p> <p>51. Do the same on the other side.</p> <p>52. Apply (whitening,</p>	
--	--	---	--

		hydrating) essence on neck evenly, 1 min. 53. Apply eye, face, neck masks (at the same time), 15 min. 54. Wash clean face and neck, 5 min. 55. Apply toner, essence, cream, eye cream, neck cream, and sunscreen.	
--	--	--	--

Suggested Treatment:

Ten treatments account for one full course. After one treatment, the skin will be firm, tender, delicate and smooth, and the eyes will be elevated to accelerate the blood circulation and lymph circulation of the eyes. After one course, it lightens fine lines, dark circles under the eyes and brightens the skin. After two courses, the overall skin is smooth, firm and elastic, and the double chin is improved. Three courses help prevent and delay skin aging, relaxation and sagging, and restore skin to its youthful state.