

40K Ultrasonic Cavitation:


Using the principle of ultrasound, through the gathering of strong sound waves into the human body to cause a strong impact on the body's fat cells and the friction movement between fat cells, can effectively consume heat, consume cell water, so that fat cells shrink, in addition when the sound waves vibration can cause a strong impact between the pressure of the cells, the cells instantly burst, fat cell reduction, so as to achieve the effect of fat removal.

The cavitation principle of ultrasound: the liquid vibrates and produces tens of thousands of tiny bubbles, namely cavitation bubbles. These bubbles in the ultrasonic longitudinal propagation of the formation of negative pressure zone growth, and in the positive pressure zone quickly closed, so that the alternating positive and negative pressure under compression and stretching. The bubble is compressed until the moment of collapse, making it burst, which will produce a huge instantaneous pressure, generally up to tens of megapascals to hundreds of megapascals, and produce strong vibration and noise.

Advantages: selective tissue destruction, specific frequency, only low-density fatty tissue destruction, and protection of high-density tissue such as blood vessels and nerves.

RF:

Radiofrequency waves directly penetrate the skin, using the skin's electrical resistance to generate energy, making the skin's underlying temperature rise, rapid and sustained heating tissue and promoting the skin's collagen growth, thereby achieving the effect of skin tightening, skin lifting, and skin refinement.



BIO


By outputting biological currents that mimic the human body, it enters muscle cells through the skin, stimulating the energy ATP contained in the cells, allowing the cells to resume normal operation and function, providing energy to the muscles and

body; improving cell permeability by 30%–40% higher than usual, helping the body to better absorb nutrients and water, blood oxygenation, etc., while discharging waste, toxins and stagnant fluid as soon as possible; promoting muscle fiber synthesis and stimulating collagen regeneration, which can reduce wrinkles and make the skin softer, more elastic and lustrous.

Skin facial lifting

Microcurrent therapy -acne treatment & skin firming

Microcurrent-safe and painless facial helps erase fine lines and wrinkles, while firming your skin and defining your features. It improves muscle tone, reduces puffiness, increases cellular activity, and tightens pores.the overall result is a healthier and younger looking skin, no matter your age.



Laser lipo

The laser with a wavelength of 635nm–650nm is a visible spectrum. This light wavelength penetrates extremely well, effectively activating and repairing fat cells, penetrating the fat layer, and causing the subcutaneous fat to dissolve by heat, and its treatment is safe, without side effects or pain.

The laser uses the latest non-surgical, non-invasive laser lipolysis technology, releasing low levels of laser energy. The

low-energy laser (biostimulation) reaches the deepest part of the fat, stimulating biological cells and inducing or enhancing a series of physiological responses by giving appropriate energy, including promoting local blood circulation, regulating cellular function, enhancing immune function, and promoting cellular metabolism. It can tighten skin, remove wrinkles, reduce fat, and shape at the same time. The laser penetrates extremely well and can activate or repair fat cells more effectively. It can penetrate the fat layer and dissolve the subcutaneous fat by heat, and its treatment is safe, side-effect-free, and painless.



Hot and cold hammer

The use of the principle of thermal expansion and contraction, the use of special chip ceramic, absorbing energy, can reach $-5\text{ }^{\circ}\text{C}$ cooling effect in an instant, promote metabolism,

alleviate allergic reactions, shrink pores, calm skin, skin whitening, and rejuvenation; chip absorbing energy and generate heat, with the physical aggregation of heat generated to reach 40 °C heat effect, such a thermal effect can not only quickly into the skin function of beauty nutrients, so that the skin better absorb nutrients in skin care products, but also stimulate skin blood circulation, increase the metabolic rate of the skin.

Physiological effects and effects of the cold hammer.

1. Constrict blood vessels, produce pore contraction, and have great room for improvement on enlarged pores.
2. Strengthen muscle fibers, make the skin become firm and elastic, and smooth wrinkles, and crow's feet.
3. Shrink blood vessels, effectively reduce varicose veins, and red tendons upward, while the harmful waste flowing back into the veins is more easily purified and detoxified, and has a significant effect on the lightening and elimination of spots.
4. Reduce the metabolic rate, inhibit and slow down cellular activities, make the peripheral nerves less sensitive to play a calming and soothing role, and reduce the swelling, swelling, and pain of the skin, especially for allergic skin has a very good effect.

Physiological responses and effects of the hot hammer.

1. Warmth, making people feel warm and comfortable, relaxing muscles, balancing autonomic nerves and relieving pain and soreness. Promote the acceleration of blood circulation, increase the oxygen content in the blood, and make the pores open and blood vessels expand so that the nutrients on the skin are quickly absorbed, which is a gentle, comfortable, and effective method of introduction.

2. Increase the rate of skin metabolism, so that the tissue contains increased heat energy, natural elimination of waste toxins, the elimination of bags under the eyes, dark circles have a certain effect

Alternating hot and cold, the effect of adding and multiplying on beauty, can naturally revitalize cells, and improve and purify skin quality, is a health and beauty regimen.

