

Female Breast Pump Enhancement

Use manual



Thank you for you purchasing BEMT-002 massager,Be sure to read this instruction manual carefully before use so that you can make full and correct use of your massager.Please keep this manual available for future reference after you have read it through.

1. Introduction



A Nonsurgical, Scientifically-Proven, Breast Enlargement System Professional, effective, breast enlarging vacuum pump. Our breast pump is for firming and enlarging of underdevelopment breasts with immediate results, quick, and effective vacuum suction for females and males who want to enlarge their breasts. Or for those who want their breast to be firmer.

The device is intended for women. This system is a nonsurgical, scientifically-proven, breast enlargement system. It uses negative

pressure to apply a gentle three-dimensional pulling force that causes the breasts to slowly and naturally grow an average of about one cup size. The method lifts, shapes, and fill the breasts. The procedure creates no scar.tissue, no stretch marks, no asymmetry of breasts, and the results are permanent. (Measurements 22 months after use was stopped show no size loss.

2. How does it work

This device is based on a medical principle known as tension-induced growth: The vacuum stretches the breast tissue gently over a sustained period of time. When cells are stretched, they respond both by swelling temporarily and by replicating until, eventually, the cells divide to accommodate stretch force. The new tissue is identical to the woman's original breast tissue, and feels the same both to her and from the "outside."





absorb blackhead

usage: use cosmetic soften
blackhead , then use it without bra ,
stick to the skin of face.

absorb fat of lumbar vertebra




slimming ham



slimming crus



55kpa Suction + Two cups with different size.



3. Correct Use

Preparing & setting the electrodes

When used for the first time, please fully complete discharge then recharging it.

- (1) Firstly contact the charge head to the product first, then insert the charger to socket.
- (2) When charging, please make sure the power switch is in a state of open, and check the charger showed light is on or not.
- (3) Charging time not more than eight hours at a time
- (4) Each fully charging can be continuous use of more than half an hour .
This product can be used when charging, the slow charge while using .

- (5) More than a week the massager is not used, please make sure the charger is in good condition to charge again .

The best Treatments

- (1) Before use, please first bath, let the breast skin fully wash, promote metabolism, the best effect .
- (2) Listen to the music or the time to rest, let the body and mind and relax muscles fully
- (3) Switch on the suction, and then adjust the intensity switch, transferred to the suction acceptable range, then to small one back (A suck a put ,



the best effect)

- (4) After using this product, please fully clean water or alcohol, air-dried

Properly kept

4. FAQ

Q. Is it Safe?

A. There have been a number of studies and tests, and no damage or risks have been found. Tension-induced growth has been used for over 30 years in other areas of the body, so the science on this is very solid. There are no side effects or health risks.

Q. Is it painful when using the Breast Enlargement Pump?

A. No, If you follow the complete directions, there is no pain or discomfort when using the Breast Pump. Q. Are there any side effects?

A. There are no known side effects, besides breast growth.

Q. How often do I use the Breast Enlargement Pump, and for how long at a time?

A. Using it longer and more often- does produce faster results. Pump to a comfortable pressure and hold that for 10 to 15 minutes on each breast twice a day. Do not over pump, if there is pain or discomfort then release the pressure.

Q. Does your Breast Pump work on males and females?

A. Yes, it works on any breasts, the pump stretches the breast tissue over and over again. It is called soft tissue expansion. It is a lot like working out your muscles repeatedly

5. Note

- (1) When the charger is not used, please put the charger unplug. The fork, please be sure to seize the charger, please don't pull wire, otherwise may damage the charger and cause a fire or accident.
- (2) If the wire damaged or overheating, please stop using immediately.
- (3) Please Don't use it nearby Benzene, thinner, or spray nearby
otherwise may cause explosions or fire.
- (4) Wet hands or not dry, please do not use the charger, otherwise may cause electric shock or short circuit .
- (5) Please do not place the product and charger or use it nearby Damp places, such as the bathroom, the bath, wash basin or other high damp place, otherwise it may cause electric shock or fire.
- (6) If skin rash or skin irritation, do not use this product, In case of causing bacterial infection or infection
- (7) Please put this product out the child's reach, do not let the children to use or play .

6. Problem check

- (1) Do not repair itself or is removed or modified in this product, you may be injured by improper operation or cause a fire .
- (2) If the wire or plug is damaged, please immediately stop using this product

7. Feature

- (1) Suction is 55kpa
- (2) Two cups with different size.
- (3) Suction can be tuned.
- (4) Charger can work for 110v-250v

8. Package included

Main Device X 1

Adapter X 1

Cup X 2

Manual X 1