


Led Light PDT Skin Rejuvenation Beauty Lamp

Model:DWLED122



SAFETY INSTRUCTION

Thank you for purchasing this photon led laser skin care machine

Please read instruction carefully prior to use, and use correctly according to the instruction

Please do not use when you have the following skin problem:

- 1 The inflammation caused by acne and erythema
- 2 Any metal/plastic/colloidal silica part used in your face skin
- 3 Wound
- 4 Cutaneous,pigmentation caused by cutaneous pruritus /burn/stimulation
- 5 The inflammation caused by cosmetics
- 6 Oral interior and other mucoustissue

SPECIFICATION


Input voltage: AC100-240V(50/60Hz)

Output Power :25W

Working temperature: -5 centigrade-+40 centigrade

Storage temperature :-5 centigrade-+40 centigrade

PARTS INSTRUCTION



QUICK USER GUID

- 1 Thoroughly remove make-up and impurities first,The photon led laser treatment must be performed on the cleansed skin
- 2 Connect the machine with the electric adaptor ,Turn power on and then the

indicator light is on

3 Set Program 1 or 2 ,choose red or blue or yellow color according to your requirement

Program 1

Spring-Summer Treatment Mode

The treatment mode are special used in spring and summer,It will be overall nursing in short time,it is called program 1 as a short time nursing mode

Nursing mode for Acne/Red spot/hemangiectasis

Improve the skin problem of skin sensibility caused by imbalance diet and weary

Mode:Blue&Red Program 1(Safety goggles must be used)

Nursing mode for pigmentation/acarocecidium

Improve the pigmentation/wrinkle caused by ultraviolet radiation and age

Mode:Blue&Yellow,Program 1(Safety goggles must be used)

Program 2

Autumn-Winter Treatment Mode

The treatment mode are special used in autumn and winter ,It will be overall nursing in short time,it is called program 1 as a short time nursing mode

Nursing mode for hypersteatosis

Improve the troublesome skin problem

Mode:Yellow/Red Program 2(Safety goggles must be used)

Nursing mode for pachylosis/wrinkles/sag/pore

Improve the skin problem caused by under sleep and ruleless life

Mode:Yellow&Red Program 2(Safety goggles must be used)

Nursing mode for xerosis cutis/dry skin

Improve the skin problem caused by winter and airconditioner

Mode:Yellow. Program 2 (Safety goggles must be used)

4 Wear safety goggles, keep the distance between LED light head and skin above 15cm

5 Time: The skin cell which makes you more beautiful is produced from your body,it should take time,If you select the nursing mode of speeding skin cell reproduction,we suggestion to use LED light treatment by a little time everyday,Irradiation nurse on partial should be keep for over 15minutes

6 Care after treatment

Please donot make you in the sunlight directly after treatment

FAQ

What is LED? How is the beauty function of different color?

LED is the abbreviation of Light Emitting Diode., it is specialized in special wavelength(color), it also be diode used for traffic signal and Domestic electric appliances . As a beauty system, we always pay good attention to the selection of wavelength and power of each LED light

Is there any influence to the eye?

The LED light has not affected to the direct vision. If one instantaneously looked that the strong LED light, he/she will feel temporarily change dark to surroundings, please notice we have to wear the safety goggles when we select the treatment mode which must be used with wearing goggles. Also we suggest to close the eye when we select the treatment mode which need not wear the goggles. Please terminate use the LED system if you are not suitable.

Could the whole body may carry on the LED light treatment?

Yes, it could. When we select the goggles wearing treatment mode, Please note to avoid looking straight ahead the treatment area. The LED light system is the single plane large-size beauty system, so it could be nursing on Pate, back, buttocks and any other body part.

How about the treatment time?

The skin cell which makes you more beautiful is produced from your body, it should take time. If you select the nursing mode of speeding skin cell reproduction, we suggest to use LED light treatment by a little time everyday.

I am troubled with mode unsuitable?

If the acne has the inflammation, we only suggest to use Acne treatment mode.

Which treatment mode is better for everyday's treatment? Is there any problem to skin if I change 2 or 3 modes in one day's treatment?

We suggest to use program 1 for Spring-Summer treatment, program 2 for Autumn-Winter treatment,it is the fast way to automatic multiple treatment mode in one program. Certainly, each light is out of Ultraviolet ray, if you select the different mode for treatment, please take it easy to use. only we suggest to treatment no more than 30 minutes for same area.

Is there any preparation piror to treatment?

The patients should clean suitably to the treatment area. For example, if the face is nursing, it is needed to clean the face at first.

Is there any slight pain feeling when treatment?

The LED light system is Non-contact auxiliary beauty system without any electrolyte or micro electric current, etc. The LED light system is safe enough to treatment.

Since it will cause the skin to be somewhat dry, please note the full moisturizer