

Gold Laser Spot Removal Pen/freckle Pen/mole Removal Device



Operation Theory

Operating theory


Control by chip low temperature high frequency plasma micro discharge technology, instantly generate over 2000 degrees Celsius in epidermis spot, then carbonized and evaporate spot, 2000 degrees Celsius can kill all viruses and bacteria. Blood vessels will be frozen, meanwhile, oil layer cover wound automatically, so it has disinfection, sterilization, stopping bleeding functions.



Photothermal laser act on skin dermis.



Penetrate deep skin, absorbed by melanin granule.



Blasting instantly, decompose to small granular.



Discharge through metabolism to smooth skin.

Easy operation

GuangZhou BPP
Beauty Equipment Co.,Ltd
www.chinabpp.com/en



GuangZhou BPP
Beauty Equipment Co.,Ltd
www.chinabpp.com/en

Mobile/WhatsApp:0086 18102525729
Facebook:bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com



Mobile/WhatsApp:0086 18102525729
Facebook:bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com



Open switch , start to operate.

Beauty Equipment Co.,Ltd
www.chinabpp.com/en



Mobile/WhatsApp:0086 18102525729
Facebook:bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com



GuangZhou BPP
Beauty Equipment Co.,Ltd
www.chinabpp.com/en

GuangZhou BPP
Beauty Equipment Co.,Ltd
www.chinabpp.com/en

Mobile/WhatsApp:0086 18102525729
Facebook:bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com



LED light design

Anytime anywhere can
use,not affected by light.

Mobile/WhatsApp:0086 18102525729
Facebook:bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com

Operation



Tungsten steel needle

Needle insert into the top of the pen, needle tip outward during using, when not in use, it is inward.

GuangZhou BPP

Beauty Equipment Co., Ltd.
www.chinabpp.com/en

Mobile/WhatsApp: 0086 18102525729

Facebook: bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com


Gear adjustment

There are low flame, close, high flame three gear. Needle order is low flame, close, high flame.



Mobile/WhatsApp: 0086

Facebook: bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com



GuangZhou BPP
Beauty Equipment Co., Ltd.
www.chinabpp.com/en

GuangZhou BPP
Beauty Equipment Co., Ltd.
www.chinabpp.com/en

Power switch

When not in use, please turn off it.

Mobile/WhatsApp: 0086 18102525729
Facebook: bppbeauty010@gmail.com
Email/skype: sale010@chinabpp.com

Charging connection hole

When the red light turn into green light on charger, it means 3mm diameter fleck is about over 100, when light dims, means need to charge.



5,Function & Effect

Function:

- 1.Nevus removal

2,Pigment removal

3.Verrucous nevus removal

6,Features:

1).It can use in hospital to have an injection.

2).It can use in shallow skin and deep skin ,make it smooth and remove the wrinkle ,and so on.

3).The effect of this machine is very fast,no need to use it very often

4).We are the sole supplier in china,we can offer very competitive prices for all the customers

Technical Advantages

1). Protect the integrity of the skin; deliver the nutrition to the optimal position.

2). Easy to operate, no side effects, but significant effect.

3). Stimulate the skin for self-healing, promote skin metabolism .

4). Eliminate the deep skin's toxins and waste activate cells and repair the damaged tissue.

Packing list:

Laser pen: 1pc

Power charger: 1pc

Needle: 12pcs/set