

2in1 Far Infrared Heating Negative Ion Light Eye Massage

SR-269

Use Manual

Stunning Small for Gift



- Far-infrared heating
- Negative ion
- 360° Rotate roller ball
- LED photon therapy

PRODUCT SUPERIORITY

- 3in1 multifuncton
- Advanced technology
- 360 rotational ball
- Easy carry & use



NO.2

BEAUTY BEHIND TECHNOLOGY

combine with three technologys together



Far-infrared heating

accelerate blood circulation, improve human body physical magnetic field,Relieve fatigue and help to relax muscle and nerve

Negative ion

Promote absororption, accelerate improvement and strengthen effects.

LED Light Therapy

Stimulates collagen production which smoothes eyes' wrinkles and saggas

NO.3

28 DAYS TO WITNESS THE MIRACLE

To improve the remarkable result by using this device,we have choose people with different ages and different career,following is the investigation indicates:

1 Day

Enjoy spa comfortable at home

15-20mins/time,1-2times/day
Promote absororption,accelerate improvement and strengthen effects

7 Days

Skin become more smoother

skin turn out to be more smoother,more elasteric

28 Days

Miracle comes true

wrinkle removal 33%
skin smooth 28%
skin elasticity 15%



CATALOG

Preface	3
I Principle	4
II Functions	4
III Operation Procedure	5
IV Specification	6
V Package Listing	6
VI Warming	7
VII Attention	7

Preface

Dear User:

Thank you for purchasing !Be sure to read this instruction manual carefully before use so that you can make full and correct use of your massager.Please keep this manual available for future reference after you have read it through.

This product does not belong to medical,please don't use for disease treatment

I Principle

This product adopts far infrared heat effect principle. Thermal radiation can penetrate into deep skin, and absorbed by skin organization.

It can make atom and molecule of human body cell vibrate together to form thermal effect. It can accelerate blood circulation and metabolism and organism regeneration capacity.



II Functions

Far-infrared heating:

With 360 rotate roller ball heating up to accelerate blood circulation, improve human body physical magnetic field, relieve fatigue and help to relax muscle and nerve.

Negative ion:

Promote absorption, accelerate improvement and strengthen effects.

III Operation Procedure

Note:press any key second boot

a.Heating function

(Setting 20 minutes on high and low temperature respectively.Heating will automatically turn off after five “Di” sound.)

- (1)First press the temperature key for the 1st gear of heating up,green lamp for 38℃;
- (2)Second Press for 2ed gear of heating up,red lamp for 45℃;
- (3)Third press for turning off heating.

b.Light therapy function

- (1)First press the light therapy key for turning on light therapy,lighting on three red lights on the massage head;
- (2)Second press for three red light flashing;
- (3)Third press for the three red light s flashing in circle;
- (4)Fourth press for the light therapy function.

c.Negative ion function

(Holding conductive slice in the back of the pen to feel negative ion from the roll massage head)

- (1)First press the key for the turning on the negative ion function,lighting on the orange indication light;
- (2) Second press for turning off the negative ion function.

Attention:

d.Low voltage alarm function

When voltage is low,the alarm will be on with “Di Di Di”automatically.In this condition,please charge battery immedicately for operating.

e.Charge method (The pen could keep operating when under charge condition)

Plug the adapter DC into 5.5V socket in the bottom of the pen.One red indication light flashing continuously means under charging;

Two red indication lights keep on when charging over. Then charging is completed.

Remark: You can hear “Di” When press key.

IV Specification

a. Power adapter parameter

Input voltage: 100-240V AC 50-60Hz

Output voltage: 5.5VDC 500mA

Charge current: 300mA

Max power: 2.5 W

Negative ion output voltage: 5.5-7.5V

Negative ion current: 70-95uA

b. Size

Main body: 180mm* 30mm*30mm(L*W*H)

Adapter: 55mm*65mm*41mm(length*width*height)

c. Weight

Main body: 100g

Adapter: 80g

d. Battery specification

Rating voltage: 3.7V (inner lithium battery)

Rating capacity: 600mAh

e. Charging time: 3 hrs

f. Working duration: 1.5hrs

V Package Listing

a. 1 pc Main machine

b. 1 pc Power Adapter

c. 1 pc Chinese -English-Japanese Use manual

d. Without Package

VI Warming

- a.This machine is electronic product,please keep away from water and don't strike it;
- b.Don't use water or benzene products to rinse or clean it;
- c.Rotate left to remove the transparent head for inside clean;
- d. Don't use this product together with person who has infectious skin;
- e. Don't touch the machine with wet hands;
Stop use it when overheated,avoiding short-circuit or electric-shock;
- f. Don't try self-repaired,self-dismantled or self-modified at your peril;
- g.When use negative ion function ,people who has sensitive skin will feel slight electric shock(produced by roller massage head),is normal phenomenon;
- h.If this product impaired by falling off or striking, please stop use it to avoid electric shock or scratch skin;
- I.If there is malfunction,please turn over to professional technical person to repair or replace to adviod accident;
- g.You are not recommended to use the machine when charging;
- k.Unplug the power plug after use.

VII Attention

- a.People who suffer serious disease or infectious disease don't use this product
- b.Stop use when the power adapter is damaged
- c.Not for patients with heat disease or pacemaker
- d.Not for pregnant woman
- e.Not for menstrual period

