

# SM-271 Instruction

## PRINCIPLE:

This product uses Low-frequency circuit to raise the input voltage to 3000V~5000V. When the body touches the ozone bulb, high tension goes through inert gas inside the bulb, the human body, and earth to form a loop. When high tension goes through inert gas inside the bulb, the inert gas will ionized oxonium ion "O". Oxonium ion together with oxygen in the air regenerates ozone O<sub>3</sub>. That is  $O+O_2=O_3$ . So it can give the sterilization effect by using this product.

## EFFECT:

- 1、 Strong oxidation of ozone have the effect of disinfection
- 2、 In long-term use, it can help oxidation to resolve heavy metals, toxic, sulfide and other harmful substances accumulated in the body

## PROCEDURE

- 1、 Plug the power supply pin into the power supply jack
- 2、 Open the safeguard cover of the ozone bulb;
- 3、 Press the power supply switch (ON/OFF key), power pilot lamp(P) and low gear lamp(L) are on means that the mainframe is on working state;
- 4、 Gear selection key(VOL), press it the first time, change to Middle gear lamp(M) is on. Press it the second time, High gear lamp (H) is on. Press it the third time, it goes back to low gear lamp. This key can change the gears in turn.
- 5、 Press ON/OFF key again, this product will shut down.

**Note: Every time you press the key, you will hear a "B" sound.**

## ATTENTION:

1. This machine is electronic product and have fragile element, please keep away from water and don't strike it
2. Please use alcohol to clean the ozone bulb before and after use; Don't use benzene or other chemical solvent to clean it.
3. People with infectious skin should not share this product.
4. This product is strictly prohibited to repair, demount or alternate by oneself. If it causes problems in such way, the consequences are on one's own.



05. If the ozone bulb or power cord is broken, please do not use.
6. Please seize the power plug when you disconnect the power. Do not pull hard, or the power cord will be damaged and cause the accidents.
7. Do not use this product if your hand is wet. If the wire or plug is damaged or too hot, please stop using immediately, otherwise it may cause electric shock or short circuited.
8. Please do not use this product near water, such as bathroom, bathtub, washbasin and other high humidity places, otherwise it may lead to electric shock or fire.
9. It is normal phenomenon that you will hear some slight sound when the ozone bulb touches skin.
10. Should not let the ozone bulb stay more than 2 minutes in the same area.
11. When the product is working, the O<sub>3</sub> it made may cause a slight stimulation to respiratory tract. So we suggest you not to using it on the face for more than 20 minutes.
12. Eyes should be protected during use. This product should not touch eye balls.
13. Do not turn on the power until you use this product. Please shut down the power when you don't use it.
14. The bulb of this product shouldn't touch water.

Do not use this product in the following situation

- 1、 Body temperature is higher than 38° C
- 2、 People with acute diseases
- 3、 With heart disease or heart rate regulator.
- 4、 People has Extreme hypertension、 cerebral infarction and other terrible disease.
- 5、 Women must not use it in the abdomen during menstrual period or pregnancy;
- 6、 Do not allow children to use it by themselves.

### **ELECTRIC STANDARD**

1. Power input: AC110V/50HZ/60HZ or AC220V/50HZ/60HZ
2. Max power: 5W ± 10%
3. Dimensions: 185mm\*55mm\*65mm(Length\*Width\*Height)
4. Total Weight: 185g