

SR-901

I Before Operation

1. Please connect the USB charging line with the charging hole of the charging stand, then please make the charging stand connect with the power supply
2. The indication Lamp of charging stand will be red if connection is well, please make the main device Insert into the charging socket, the LCD of the main device will be off after full, 3 hours needed to be charged each time
3. Clean skin and put some essential on skin

II Operation

1. Long press "POWER", select the working mode and its intensity.
2. Long press "MODE" into mode selection, press "H/L" to select the intensity, when it is in lowest intensity, indicates turn off that mode.
3. Please move the device on skin from down to up, inner to outside, you will feel some warm on skin, it is more better to use some essential with this device, 2 times/week recommendation

III LED function

Mode1 Press LED SEL. ,Red light, wavelength: 620nm (+ - 5nm)
Make skin more bright, more elastic.

Mode2 Press LED SEL twice, Green light, wavelength:520nm(+5nm)
Enhance the absorption of beauty fluid and make the skin more clear and whiter.

Mode3 Press LED SEL three times,, Blue light, wavelength: 465nm (+ - 5nm)
Stimulate the production of collagen. Collagen is the essential protein used to repair damaged tissue and to replace old tissue. Best for fine lines, reduce pore size, swollen after surgery.

Mode4 Press LED SEL four times, Yellow light, wavelength: 590nm (+ - 5nm)

Detoxification, make the skin more healthy

Mode5 Press LED SEL five times, Pink light, wavelength: 700nm (+ - 5nm)
Skin whitening, suit for dark skin

Mode6 Press LED SEL, six times, Pink light flash.
Enhance the nutrition absorption.
Press LED SEL, seven times to turn off the light.

IV Parameters

Output voltage: DC5V 500mA

Charging time:3 hours

Consumption:MAX 3.4W

Standby time: 40 days

Capacity of battery: 750mAh

Working duration: 10 minutes

N.W:298g

Size: 165*45*46mm

Package size: 215*150*85mm