

Table of contents

Specification

Safety precautions

Important notice

Parts instruction

Quick user guide

Photon treatment

Warranty certificate

Specification

Photon-Ultrasonice waves skin care machine

Adaptor input power:AC100-240V

Adaptor output power:15V,400 m A

Dimension:17.5x4.7x4.4cm

Weight:150g

Ultrasonice frequency:3MHz+_5%

Safety Precaution

Note:Please read following information carefully before using your machine

This product is well-desinged precision instrument,in order to have your device functioning perfectly,pls read the below information thoroughly and follow al the instruments this will help you to ensure the warranty vailidity and to extend the life cycle of your product

- 1 Store your product always from the humidity and do not use it under water, otherwise the electronic components of your product be corroded or have short circuit
- 2 Never attempt to disassemble the product yourself, all the product components can only be maintained and repaired by professional technicians
- 3 Don't use any adaptor other than the one included in the package,using incompatible adaptor may result in unfunction and will avoid the warranty validity
- 4 Any people with the following condition must avoid using the products: patients with maligrant,women during her pregnancy or period ,patients with malignant,
- 5 Keep this product away from children ,children shall not play wit this product as the disassembled components may be damaged by the children accidently

Important notice

- 1 To ensure that this machine functions well,you can drop a tiny amount of water on the ultrasonic generator surface to see the water bobs,if the water bobs,it means that this machine is under normal function

- 2 This machine should be used under normal conditions,remember to turn off the power and unplug the cable when you stop caring skin,this will avoid over heating problem caused by the air vibration and to prevent any possible burning of shortening product lift cycle
- 3 Don't use the machine without putting any skin care product on you skin ,this will cause wave seriously reflection and result in discomfort or burning of your skin ,Gelatinous skin care is recommended
- 4 Avoid using this machine in one area for too long,move it around on your skin to avoid focusing effect caused by the wave energy,this effect may result in discomfort or burning of your skin
- 5 When you are using this machine,always adhering the ultrasonic wave generator completely to your skin,an incomplete contact will cause a big difference of impedance energy on the interface, this will also cause an over-heating problem and result in discomfort or burning of your skin of impedance and dispence energy on the interface,this will also cause an over-heating problem and result in a discomfort or burning of your skin
- 6 When you finish using the machine,wipe off the left on the ultrasonic wave generator,A soft cloth is recommended to prevent any possible scratch
- 7 Do not use water and any other chemical,cleanser or washing-up liquid to clean this product,this will cause it malfunction or shortened product life cycle

Please to allergy test before using for the first time


The simple allergy test:

Cleanse your wrist

Apply a portion of the special essence on the cleansed area, let it dry naturally(during the test,please don't touch water nor wash the essence away)if you experience allergy or any skin abnormalities,stop the test immediately and contact your doctor or nearest hospital straight away
If after 48 hours,no allergy or any other forms of skin abnormalities appear,you many enjoying this special unique essence and equipment regularly



Parts introduction



Operating Panel Instruction:

QUICK USER GUIDE

Method to care your face

1. Thoroughly remove make-up and impurities first. The ultrasonic waves skin-care treatment must be performed on the cleansed skin.
2. Connect the machine with the electric adapter. Turn power on and then the indicator light is on.



3. Apply appropriate portion essence on your face and avoid eye area. Please ensure the treatment area is well covered with the essence. (Note: We highly recommend you to use our MAYA designated essence for its best performance on conducting ultrasonic waves)
4. Apply appropriate portion essence on the ultrasonic waves generator.



Suggestion: You'd better use the designated essence to transfer ultrasonic waves.

Basic method to care your face



Step 1 it requires 8 minutes.
Using ultrasonic waves generator to massage skin equably.
 It requires 4 minutes. First, Using ultrasonic waves generator to care the right side of your face for 4 minutes (Please use a timer or your watch as a timer). Then working on the left side of your face for 4 minutes.
 Please operate the appliance in the order indicated in the diagram, working upward motion (about three times).



1. Along indicated route from jaw to the root of your ears.



2. Along indicated route from jaw to the front of your ears.



3. Along indicated route from the corner of mouth to the top of your ear




4. Along indicated route from the bottom of nose towards the area near your eyebrow.



6. Along indicated route from nose towards the area near temple.



5. Along indicated route caring both sides of your nose nosewing.



-6-

photon treatment

photon caring: Daily caring .

Select red light, nursing face , especially wrinkle area.
move detector every 5 seconds.

Recommended photon caring time is 10-20 minutes

photon caring: Dilute the spot, even skin color, whitening skin

Select blue light, nursing face , especially spot area.
move detector every 5 seconds.

Recommended photon caring time is 10-20 minutes



Lymphatic detoxing

photon caring: Acne (one time every night)

Lymphatic detoxing

1. Select ultrasonic function and intensity, perform Lymphatic detoxing according to next route.

Use blue light to care acne area 15-20 minutes.

Use green light to care 10 minutes.



Massage around acne.

- 0 -

[Guarantee Certificate]

Product model: _____

Purchase date: _____

Purchase place: _____

Sales Dept. _____

One year Repair Guarantee

In the unlikely event that your product develops a fault within the first one years of purchase:

(Please return your product along with your receipt and the original guarantee certificate and we will repair it free of charge provided you have operated it according to the Operating Instructions book)

What is not covered?

1. Breakdown or damages due to
 - (1) Not following the operating instructions.
 - (2) Incorrect or unauthorized repair
 - (3) Accidental damages during your own transportation, or you drop it on the floor or you deliberately damage or neglect the product after purchase.
 - (4) Fire, earthquake, weather such as flood, lightning, high wind.
2. Please keep your guarantee certificate safe, as no duplicated certificate will be issued.

-9-