

彩光 & 超音波嫩肤仪

Colour Photon & Ultrasonic beautiful skin instrument

中·英文使用说明书

Chinese & English Operation Instructions



Foreword:

The history of light-treated is very long, the ancient Greek, ancient Rome person have been known some simple diseases of human body, which can pass to take a sunbathe light to be eased or cure.

More than 100 year ago , Denmark doctor NIELS RYBERG GINSEN developed the first set optical instrument , which be used to cure skin tuberculosis, the curative effect was salience , so got the Nobel prize at that time .

The photon beautiful skin: the concept of the photon delicate skin originate the United States , translate from "Photorejuvenation", original meaning BE"light become rejuvenated in old age", from the typeface explain , be make use of strong pulse light instauration a skin of young and healthy. The photon beautiful skin is currently a kind of revolution hairdressing cure procedure in the world, it adopts patent of strong pulse light technique, can be not hurting the premise that the normal skin organizes to descend clearance face various dye and afferent sex spot a piece, stimulate gum to organize to increase to livinged at first in the meantime, instauration skin flexibility, make callous pore constringency, the crease reduces, making face skin the quality get whole promote, afresh shining the youth elegant appearance of health. It provides a safety and uninvolved method , can adapt the skin of different appearance

Light classification: :

1. The sunlight includes a lot of lights of different wave lengths, there are infrared light, visible light, (imp red, orange, yellow, green, black, orchid, purple ultraviolet light etc;
2. LCD figures seven colourful photon beautiful sk instruments chose to use seven kinds of colors in t natural rainbow light wave, is a skin reborn of medic treatment application technique;
3. Seven kinds of color lights is different to skin functi because of different wave-length:
 - (1) . Red light(R)-activate(wave-length 640 nm)
 - (2) . Blue light(B)-slow(wave-length 480 nm)
 - (3) . Green light(G)-parallelism(wave-length 518 nm)
 - (4) . Purple light(R+B)-drain(wave-length 640 nms)
 - (5) . The light(R+G) of Huang-protection(wave-leng 590 nm)
 - (6) . The month is only green (G+B)-right (wave-length 640 nm)
 - (7) . White light(R+B+G)-ducting(wave-length 510 nm)

Mostly technic guide line:

Output frequency:10W
 Output voltage:DC-24V
 Frequency:50-60HZ
 Packing:21cm X 18cm X 6cm.

Function elements:

- A. Pulse seven colourful kinetic energy light techniques through the clinical trial of several years and the thorough research of the medical expert, the following effect has already obtained a confirmation:
1. Can with function in cell membrane, regulate the fatty acids of cell, promote the creation of inside the cell, arouses the metabolism of cell, storing the energy that cell need thus;
 2. Increase cell membrane's to the oxygen and water to measure deeply, can make the cell absorb more oxygens;
 3. arouses natural defend a mechanism of body, increase blood plasma albumen etc. blood composition, so strengthen human body immunity dint;
 4. Help replication of collagen fiber;
 5. Promote the growth of new-born blood vessel;
 6. Regulate the nervous system of human body;
 7. Promote lymphoid system circulation;
 8. Speeding DNA inside the body and protein synthesizes, attaining a delicate skin and dispeling a wrinkly purpose thus.
- B. Ultrasonic: Flap to become bankrupt to living ultra sound of having the deep penetrability through each 1,000,000 time high frequency, ultrasonic function in human body behind absorb for the machine body, the voice can turn into thermal energy, making soft and

very small massage to the skin cell while warming up the function, promoting the blood circulation of the human body capillaries, accelerate metabolism to promote the nourishment material absorb, improve defense ability of the machine body, activate skin cell repress germ to breed.

Pulse seven colourful kinetic energy light character:

1. Pulse seven colourful kinetic energy photon release about the 400 nm-2000 nm wave-length;(can immediately visible light and little part of infrared light mix and become);
2. Pulse seven colourful kinetic energy lights has advantage more than the sun light(sunlight), because it doesn't contain ultraviolet light, will not tan and burn a skin, will not lead with the result that skin cancer, and don't bring injury to eyes;
3. pulse seven colourful kinetic energy lights also has advantage more than the laser(laser), because its wave-length ratio laser long, so pierce through the dint ratio laser geniality, but pulse seven colourful kinetic energies are unlike laser concentration only at 1:00, it can with geniality of the energy overlay more wide treatment scope;
4. pulse seven colourful kinetic energy light contain small quantity infrared light, will raise skin temperature slightly surname 1-3 degrees, very safe

The BIO biology character:

1. have no harm: Geniality all good reaction;
2. the reaction of slack: After ending operation treatment, the living creature effect produced by tiny electricity will not disappear immediately, but through a period of time after just disappear;

3. backlog effect: After many treatments, the living creature effect of living creature electric current will increase;
4. Living creature effect: The increment of the living creature tiny electricity, can cause some proliferation, come together or increment of the chemistry material reflect by going to change chemistry speed and mode in the body. For example after getting hurt of blaze sex a reflection, can accelerate wound to match more. The living creature tiny electricity anti-decrepitude, the protein synthesizes acceleration after the face treatment, the gum increases at first, promote rivet water ability and flexibility of skin;
5. The influence of the living creature film: Creation and maintenance certain stable chemistry set cent and selectivity deeply material function, living creature tiny electric power influence the osmosis ability of organism membrane, the n increasing the transparence, conductivity and energy of cell material conversion, promoting skin better absorb nourishment material.

Effect and usage scope:

- A. colourful light
1. red light(R):
 - a. living creature effect: Have cell to build up of incitement vitality, accelerate a blood circulation function;
 - b. physiotherapy function: Remove sallow, wrinkle, unmatue acne skin quality; Accelerate blood circulation function, promote activate of the skin organize;
 - c. apply scope: Any skin.
 2. green light(G):
 - a. living creature effect: It is a kind of color of nature, have to moderate, stabilize of the effect function lie between the orchid and yellow;
 - b. physiotherapy function: Regulate a callous gland body function, improve an oiliness skin, black-face role, blackhead, acne etc.;
 - c. apply scope: Oiliness skin.
 3. orchid light(B):
 - a. living creature effect: It is a kind of cold, pure and deep light, the usage will bring a fresh and cool felling, can be parallel with concentration energy, achievement ability with red only contrary;
 - b. physiotherapy function: Promote glair and bone gum synthesize at first, activate a skin, the constringency relaxs of skin;
 - c. apply scope: Any skin;
- B. Ultrasonic
1. beautify a skin: cancellation rimple and cancellation black spot, cancellation pimple;
 2. weight reduction result: clear away fat, partial weight reduction, whole body weight reduction;

The following circumstanc forbid to use:

A. Colourful light

1. pregnant woman;
2. using light sensitive medicine to the 560-1200 um :Tiny AN acids, four wreath vegetables;
3. just through insolation; (in two weeks)
4. diabetic;
5. have the sufferer of light allergy;
6. future the one who take allergic in the near;
7. the chronic dissolve a blood sufferer;

B. Ultrasonic

1. Body temperature more than 38 degrees persons;
2. Be easily suffer from polluting a sex paroxysm person;
3. Impatient disease sufferer;
4. Heart attack sufferer;
5. The extreme high blood pressure, brain obstruct etc. advanced case sufferer;
6. Afferent lump, whole body erythema lupus, cirrhosis;
7. Contact dermatitis, different dermatitis disease;
8. The issue of blood disease or external injury;
9. Have the women, pregnant woman of birth-control machine and menstruate inside womb women;

Instrumentation regulation

1. please careful reading this manual agowhen use the instrumentation;
2. hydrate to defend to bask over a long period of time, don't place instrument to be being subjected to heat and affecting by damp and cold of place;

3. The person of taboo can't use;
4. Don't put instrument to be falling off easily and touching easily of place;
5. If the instrument appears an abnormality, please don't tear open without authorization, should contact production factory house or supplier immediately;

Serve and maintain:

1. This instrument from purchase day to protect to fix for a year
2. Following circumstance are unavoidably to protect to fix service:
 - a. dismantle to pack to result in without authorization of breakdown;
 - b. rough it to result in of breakdown;
 - c. maintain to result in because of lacking of breakdown;
 - d. do not follow the guide of right instruction's result in breakdown;

Purchase location:

Purchase time: