

5 IN 1 BEAUTY INSTRUMENT

RU-315

Instruction Manual

Dear customer, thanks for buying our product!
For the best capability and safety, please read this instruction manual before any operation.
Please reserve the instruction manual for future reference.

Brief Introduction

RU-315 is a 5 in 1 multifunctional beauty instrument. It includes the functions of High frequency electrotherapy, Vacuum, Spray, Brushing and Galvanic induct and educe.

Main Functions:

1. High frequency electrotherapy. The instrument generates high frequency current and infrared ray to diminish inflammation and ease pain. The generated Ozone kill the bacteria, heal and astringe the acne. It has an effect on regulate sebum secretion, improve the texture of skin, accelerate the blood circulation and speed up metabolism, anti-wrinkle and anti-aging.
2. Vacuum. Refine the striper pore, acne, and pustule, clean the dirt and oil contamination inside the pore.
3. Spray. Convert the toning lotion, anti-inflammatory lotion or shrinkage water into steam and spray. It is more normative and hygeian.
4. Galvanic induct and educe. In accordance with opposites attract principle, it helps to induct the necessary essence which heal skin disease into deep skin and accelerate the absorption, meanwhile, it will educe the impurity.
5. Brushing. Use high speed rotated brush to clean the dead cells, make-up remnants, horny layer, acne etc to rejuvenate the skin. It is a safe, no side effect and efficient treatment. No scar left.

4. Welcome for any comment and suggestion.

Accessories:

Name	Qty
Mainframe	1 pc
Power line	1 pc
Galvanic probes	1 set
High frequency electrotherapy stick	1 pc
High frequency electrotherapy pipe	3 pcs
Vacuum glass pipe	3 pcs
Spray bottle	1 pc
Red & black pipe	1 set
Filter	1 pc
Brushing accessories	1 gun; 5 pcs brush heads

How to use?

1. High frequency electrotherapy

- 1) Connect the power line and the electropathy stick.
- 2) Plug the glass pipe into the stick.
- 3) Connect the power. Spin the knob to turn on and adjust the intensity.
- 4) Begin from the forehead, rotated and slide down. Treat about 5 minutes.
Dehydrated skin should treat with low intensity and short time; Oily skin with stronger intensity and longer time.

2. Vacuum & Spray

- 1) Select a glass pipe and connect to the rubber tube at one end and the filter at the other and then the filter to the vacuum outlet
- 2) Connect the power and press the “on off” button, press the hole on the pipe, it will generate attractive force.
- 3) If want to use the spray function, connect the sprayer to the tube, then connect the tube to the spray outlet. Hold it with thumb and long finger, and press the breather hole with forefinger. Air pressure will be engendered inside the bottle and spray out the toning lotion.

3. Galvanic induct and educe

- 1) Connect the galvanic accessories. Let the customer hold the hand-held galvanic wand, and the beautician hold the roller galvanic head. If operate by yourself, hold the wand with your left hand and the roller with right hand.
- 2) Press the “Power” button and spin the knob to adjust the intensity.
- 3) Select anode or cathode according to customers need. While the

instrument indicate anode, the wand is cathode and the roller anode. Embrocate liquid cosmetic essence on the skin, use a piece of cotton fleexe moistened by water or essence to wrap the roller and induct the nutrition. On the contrary, use cotton fleexe moistened by germicidal solution or toning lotion to wrap the roller and educe the impurity.

4) After operation, turn of the power, disinfect the accessories.

4. Brushing

- 1) Select a proper brushing head and connect to the brushing gun, then connect the gun to the brushing outlet.
- 2) Spin the knob to turn on and press the switch on the gun, the brush will begin to rotate.
- 3) Spin the knob to adjust the intensity and rotated frequency according to customers' skin condition.
- 4) Put the brushing head against the treating part and begin to work.
- 5) After finishing the treating, shut the machine off, clean and disinfect the brushing heads.

Attentions:

1. Keep the galvanic probes clean, wipe them with cotton gently. No water or uperization.
2. Never put the galvanic probes and high frequency electropathy stick against the eyeball.
3. During high frequency electrotherapy, adjust the current intensity from low to high. Do not let the electrotherapy stick to stay in the same position too long in case the skin get burn.

4. Make sure the power is off while connect or disconnect the electrode pipe. Keep the pipe clean and dry.
5. Keep the brushing gun and heads clean. While taking the brushing treatment, do not work to long in case to hurt the skin.
6. Person with angiocardopathy and skin disease could not take any of these treatments.
7. Pregnant women and patients with hemorrhagic disease please do not take any of these treatments.
8. Please pay attention not to connect the spray bottle with vacuum outlet, or the toning lotion will be vacuumed into the pump and damage the instrument.
9. Do not place the instrument in wet environment.

Maintenance Service:

1. The instrument is under warranty for one year except the accessories and the consumables. We would provide maintain service after the warranty period, but payable.
2. If you want to replenish or purchase the consumables and accessories, please contact the distributor for advice.
3. Breakdown cause by human as below is out of our warranty.
 - a. Breakdown cause by dismounting or transforming the instrument by customer himself/herself.
 - b. Breakdown cause by incorrect or unreasonable maintenance.
 - c. Breakdown cause by beating or dropping the instrument.
 - d. Breakdown cause by incorrect operation.
 - e. Damage cause by irresistible disaster like earthquake, fire disaster.