

INSTRUCTION MANUAL



INSTRUCTION MANUAL

1. Clean the part of skin then smear the special purpose jelly before using the appliance. The conductible jelly can be substituted by the cosmetic emulsion or the emulsion with vitamin (achieve better effect), the ways of massage refer to the latter chapter.

2. Press the power supply button A. The light NO.7 will be in the mode of continuous irradiation.

3. Press the output power button 3C. The light NO.6 will shine and the strength will be increased. Press the button 2B, then the strength will be slow down.

Ps: do not use sharp materials to damage the appliance or take apart the appliance.

USER`S MANUAL
ULTRASONIC BEAUTY&HEALTH INSTRUMENT

(1) Analysis of Body Massaging Location:

1. Upon applying the massager on thicker parts of subcutaneous tissues such as the waist hips, thigh, and calf, select the High Energy Output (2).
2. Upon applying the massager on parts like the face, where the muscles and fat appear thinner, select Low and Medium Energy Output (P.11 or P. 12 or 2).


(2) Cautions of Application:

1. before using, firstly clean the part to massaged.
2. Apply special GEL (or cream) on the parts to be slimmed or contracted or on the part with tiny wrinkles.
3. Regardless of whether harder or lighter massage on the skin, the same ultrasonic effect can be achieved (therefor, it is suggested to apply to a light massage pattern)
4. If it don't use, please put power off.
5. Please don't use at bathroom or wet place.
6. Finish working, use tissue or dry cloth to clean probe.
7. The ultrasonic massager has been indicated and operated by clinical experiments and based therapeutic apparatus. According to operating instruction for use.

(3) Two Available Massaging Methods:


1. Stroking Method:

Stroking is applied to rem over the fats in larger muscular parts, like the abdomen, hips, and thigh, stroke 2 to 3 seconds each time for a distance of approximately 10cm.



2. Rotating Method:

Rotating method is applied on the body parts with fewer muscles, like the face and arms. Its speed is similar to that of stroking.



(4) Product features

1. This product is equipped with the function selection of Continuous wave and Pulse wave.
2. The pulse is also known as a pulsation, which the emission method of instant termination after each emitting of ultrasonic waves.
3. Buzzer:
(a) At the start of operation, the buzzer function is activated by pressing any key. (b) Set the timer for 15 minutes and the buzzer function is activated after ending.

(5) Frequency of Daily Usage

1. This unit can be used at an interval of over 5 hours, up to 2-3 times daily. As an automatic setting unit, each application will last for 15 minutes after turning on the power.

2. Please note that too long of an emitting time without intervals or frequent emitting will cause muscle fatigue, which is an undesirable method and may emphasize tiny wrinkles.

3. Ultrasonic will enhance the internal skin to effectively contract the beauty care lotion of nutrition.

(6) In principle, people with the following conditions should be temporarily refrained from the emission of ultrasonic waves:

1. Suffering from body temperature over 38°C

2. Infectious disease patient.

3. Acute disease patient.

4. Cardiac illnesses patient.

5. Patient suffering from extreme hypertension and cerebral infarction.

6. People suffering from serious conditions in uterica, contagious dermatitis, and ectopic dermatitis.

7. Those who with angioma, herpes zoster, folliculitis acne, furunculus, carbuncle, general lupus erythematosus, general sclerosis, and Addison's Disease.

8. Bleeding illnesses or injuries.

9. Women implanted with intrauterine contraceptive devices.

10. Pregnant women.

11. Menstrual women, please don't use at the abdominal region.

(7) Specification

1. Frequency: 1 MHz (million times/second).

2. Each massage works up to 15 minutes and then device is turned off automatically.

3. Three level energy output (weak, medium, and strong): energy output.

4. Two output options available: continuous waves and pulse waves.

Marvlous Effect of Ultrasonic Massager

Basically, the following marvelous effects be accomplished by using the Ultrasonic massager:

- (1) Skin beauty
- (2) Body slimming
- (3) The user of marvelous witness

(1) Skin beauty by the ultrasonic massager

1. Removing wrinkles (Continuous emission waves)

The condition of beauty is the absence of wrinkles. Therefore, healthy skin means that your skin appears tender and elastic. When the elastic protein of your skin becomes aged, the skin loses its tender and elasticity and this is when wrinkles develop noticeably, these type of wrinkles are not tiny wrinkles, but remarkable ones appearing on the forehead and around the cheek, for some tiny wrinkles, owing to the reduced speed of the metabolism in your skin, the keratinized cells tend to remain. On the epidermis constantly, thus unable to carry out the replacement for forming of newly developed cells, resulting in coarse skin and tiny wrinkles. Under such circumstances, immediate use the ultrasonic massager will effectively improve wrinkles in bringing about a shiny face. The effect is remarkable.

2. Removing black spots (Continuous emission waves)

You need not worry about the sunlight so much because you are still young however, as you get older, gradually black spots become noticeable. Therefore, in using the ultrasonic massager for 15 minutes, the color of these spots will become lighter. For fine black spots, obvious effects can be seen just within a short period of time, displaying the power of ultrasonic waves in destroying melanin. Therefore, ultrasonic waves can make black spots become lighter with a bleaching effect, hence whitening your skin. Experiments have testified that in using the ultrasonic waves, the melanin would pass through blood vessels or lymphatic vessels to kidneys where it will be discharged from the body along with the urine.

3. Removing acne(Continuous emission waves)

Sebaceous glands are found in the deep part of hair pores from where sebum is secreted. These secreted fats together with the sweat secreted from the sweat glands tend to be mixed in forming the sebaceous membrane which protects your skin during the process of growing up from childhood to adulthood. There are times when hormones become unbalanced and the sebaceous glands controlled by hormones continuously secrete sebum. Furthermore, excessive secretion blocks hair pores and cause bacteria to get into your hair pores and from inflammations and this is the time in which acne develops. Therefore, it is necessary to wash away the oil from your face first. Then use the ultrasonic waves, which have the cleansing action, as well the ability to remove impurities from the hair pores, for the easy curing of acne.

(2) Slimming effect resulting from the ultrasonic massger

1. Removing fats

The major cause of fatness is over nutrition in which excessive nutrients are excreted out of the body while some are digested in and absorbed your body. These absorbed carbohydrates or proteins are then converted and stored inside the body in the form of fat. Therefore, fatness means excessive fats and people call this obesity. Of course, spots can burn fats and bring down weight, but most people fail to stick at it constantly and so the effect is less satisfactory. With the ultrasonic massger, you can remove fats. Instead of saying removing fat, it would more appropriate to say accepting the ultrasonic waves. The vibration of ultrasonic waves works to raise the temperature of some fat which will then be decomposed burned and finally turned into thermal energy for excretion out of the body.

2. Local slimming(continuous emission waves)

Ultrasonic waves can dissolve fats. For local fat parts, the application of ultrasonic waves can achieve the effect of local slimming. For parts which easily store fats such as the lower abdomen, hips, thighs, calves or double chin, apply, ultrasonic waves daily for 15 minutes once. IN TWO MONTHS, you will achieve a slim figure that envied by others.

3. Clinical experiments have indicated that ultrasonic waves could enhance the absorptive function of your skin.