



Bionic Golden Clip User Manual

Thank you for purchasing this instrument.

Please read the user manual carefully and use it correctly.

Be sure to read "notes" before using.

Please keep the instruction stored properly.

I. Product Overview

Bionic golden clip is a beauty instrument for eyes. Please read this manual carefully before using and use the beauty instrument in strict accordance with the instructions in this manual.

II. Product Introduction

Bionic Golden Clip

The principle of bionic golden clip: Activate the cell's metabolism through physical lifting technology, and realize the "kneading" to awaken skin active substances, stimulate collagen production, elastic fiber and hyaluronic acid to produce skin, make skin firmer and strengthen the body microcirculation, boost metabolism, eliminate facial edema, promote the discharge of deep toxins, mobilize cell activity, activate sleeping cells, relieve eye fatigue and make eye skin firmer, smoother and fuller.

Major Eye Problems



Eyestrain



Dark circles



Periocular edema



Crow's feet




Fat granule



Eye bags

Method of Use

Bionic Golden Clip



Frequency adjustment
【Press long to turn it off】

Bio current adjustment
【Press long to turn it off】

Instructions

1. Clean the operation area first before use (if there is makeup, first remove the makeup with makeup remover and then wash it with facial cleanser);
2. Connect the power cord. There is a round button on the top of the power supply to activate the golden clip button and the bio-current button. If you need to turn off the power, press the button for 2 seconds or longer;
3. The golden clip has three frequency adjustments, namely L (low speed), M (medium speed), and H (high speed). Similarly, the bio current has three gear adjustments, namely L (low gear), M (medium gear), and H (high gear). Both functions are displayed in the L position. If you need to adjust the gear size, press the round button to increase; when the micro current is working, the corresponding gear will always flash;
4. Apply the eye essence to the part to be operated during operation. During the operation, the essence is absorbed. It is recommended to brace the skin with right amount of force to keep the treatment surface tight and lifted and operate at a unified speed;
5. After using the instrument, wipe the remaining essence on the probe with a wet paper towel;
6. Cleaning and eye masks may not required after the treatment; (It is not recommended to wash the eyes with hot water. Warm or cold water is recommended).

Operating Method

Eyes: When operating, under the frequency of your choice, brace the skin with the device and keep the treated surface tight and lifted, from the inner corner of the eye to the outer corner of the eye, from the outer corner of the eye to the temple, and from the temple to the hairline. It is recommended to avoid the upper eyelid;(the upper eyelid is close to the eyeball therefore may result in discomfort).

Note: During operation, you will feel that the bionic golden clip is used to massage the eye skin in a "kneading" manner. At the same time, there will be a feeling of slight numbness. After the operation, the eyes may appear to be reddish, which is normal.

Note: The eye operation time is recommended to be around 10-20 minutes, and the maximum period is recommended not to go over 20 minutes. For normal skin it is recommended to receive the treatment once a day in the first five treatments. After five times, it can be operated once every 1-2 days. For those who has sensitive eyes, it is recommended to operate once every 2-3 days.

Equipment Maintenance

1. Ensure the power supply is stably adapted.
2. The instrument should not be placed in damp places or near water source, nor should it be directly exposed to sunlight.
3. The instrument should not be too close to a strong source of heat or radiation.
4. The instrument must be equipped with the original qualified parts provided by the manufacturer.
5. After each operation, the function handle and function head shall be cleaned and disinfected and properly kept.
6. When not using the instrument, please turn off the power, remove the plug, and place the instrument properly.

Note:

The instrument is a precision electronic product. It is forbidden to use the instrument in a humid environment, such as a bath room and other wet places. If it accidentally be soaked in water, please stop using it and return it to the factory for repair.

Warranty

1. One year warranty is guaranteed.
2. If the products fail in normal use within one year after purchase, we will provide free repair service.
3. The delivery cost of the product during the repair period shall be borne by the user.

Disclaimer

- ▲ When the fault occurs in the case of fire soaking in water and spontaneous combustion disaster.
- ▲ Failure occurs if the product is not used in accordance with the instruction manual.
- ▲ A fault that occurs when a consumer uses it improperly or carelessly.
- ▲ If the repair by the user is improper.
- ▲ Breakdown caused during transformation of this product by the user.
- ▲ Consumables are not covered by the warranty.