

Bionic Golden Clip & Nano Spray Gun

Cell Activation/Relieve Fatigue/Hydration



Nano Non-invasive
Mesotherapy
User Manual

- Please read the manual thoroughly before use.
- Please read "Dos and Don'ts" and "Contraindications" before use.
- Please store the user manual in place.

Overview

I. Product overview

Micro-particle non-invasive mesotherapy is a new non-invasive beauty instrument combined with bionic gold clip. Please read this manual thoroughly before using and strictly follow the instructions in this manual.

II. Product introduction


Non-invasive Micro Particle Injection Technology

The ultra-high-speed air flow passes through the particle generator and instantly decompose nutrients into hundreds of millions of fine particles, and then passes through the particle accelerator so that it accelerates nutrients to supersonic speed, instantly breaking through the skin barrier and sending nutrients to the skin.


Bionic Golden Clip

Working Principle: Activate the cell's metabolism through physical lifting technology, and realize the "kneading" to awaken skin active substances, stimulate collagen production, elastic fiber and hyaluronic acid to produce skin, make skin firmer and strengthen the body. Microcirculation increases metabolism, eliminates facial edema, promotes the discharge of deep toxins, mobilizes cell activity, activates sleeping cells, relieves eye fatigue and makes eye skin firmer, smoother and fuller.

Machine Display



Spray Gun



How To Use

How to use nano spray gun

1. Clean the intended treatment area before use (Remove makeup with makeup remover first and then facial cleanser);
2. Connect the power cord and press the round button to start;
3. Adjust the mode of the spray gun;

The screen shows the injection velocity mode; There are three modes for injection speed: A01, A02 and A03. The default mode is A02. It is recommended to switch to A01 before use, and adjust based on the customer's feedback. H00 is the manual mode, and the injection is done only after the output button is pressed for each use.

Click "OFF/ON" and blue light indicates "start"; No light means "pause"; Gear is the liquid output control valve.

4. Fill the liquid storage box when using, and it can be operated with moisturizing essence, whitening solution, hyaluronic acid and other products. Example: First, the liquid injector is used to draw 5ML of the original liquid to the syringe, and then the normal saline is used to dilute the original liquid (because the original liquid is relatively thick). After the dilution is done, the liquid storage box is injected through the injection port, and the device is ready to use after adjustment;
5. It is recommended to apply a hydrating mask or medical hydrating mask after use;

Note: It is recommended to clean the head of the gun with alcohol after each use; Inject purified water to clean the residual nutrients so as to avoid product blockage. (set the liquid output to the maximum)

Techniques

Cheeks: It is recommended to spray one spot at a time. Spray three times then move on to the next part; It is advised to stand the gun vertically at the cheeks to boost better product absorption; At the same time stroke skin upwards using hands to lift and tighten skin;

Note: It is recommended not to aim directly at the eyes, nostrils, ears and mouth (It is recommended to install positioning cap around the eyes and forehead);

After Treatment

1. It is recommended to apply a moisturizing mask every day within three days after use;
2. Always use sun protection when you go out.

Bionic Golden Clip

Bionic Golden Clip



Side A

Bio Current Adjustment
【Press long to turn off】



Side B

Frequency Adjustment
【Press long to turn off】

How To Use

How to use bionic clip

1. Clean the operation area first before use (if there is makeup, first remove the makeup with cleansing water and then wash it with facial cleanser);
2. Connect the power cord. There is a round button on the top of the power supply to activate the gold clip button and the bio-current.button. If you need to turn off the power, press the button for 2 seconds or longer;
3. The gold clip has three frequency adjustments L (low speed), M (medium speed), and H (high speed). Similarly, the bioelectric current has three gear adjustments L (low gear), M (mid gear), and H (high gear). Both functions are displayed in the L position. If you need to adjust the gear size, press the round button to increase; when the micro current is working, the corresponding gear will always flash;
4. Apply the eye essence to the part to be operated during operation. During the operation, the essence is absorbed. It is recommended to apply a little force to the skin and operate at a uniform speed.
5. After using the instrument, wipe the remaining media on the probe with a damp paper towel;
6. After the operation is completed, the eyes need not be cleaned, and the eye mask may not be applied after the completion; (It is not recommended to wash the eyes with hot water, warm or cold water is available)

OPERATING METHOD

Eyes: When operating, the frequency of the bionic gold clip is pulled from the inner corner of the eye to the outer corner of the eye, lifted from the outer corner of the eye to the temple, and the temple is pulled to the hairline. It is recommended to avoid the upper eyelid position;(the upper eyelid is close to the eyeball Poor comfort)

Note: During operation, you will feel that the bionic gold clip is used to massage the eye skin in a "pinch" manner. At the same time, there will be a feeling of micro-hemp. After the operation, the eyes may appear reddish. This is a normal phenomenon.

Note: The eye operation time is recommended to be around 10-20 minutes, and the longest is not recommended for more than 20 minutes. Normal skin is recommended to be used once every five times. After five times, it can be operated once every 1-2 days.

The eyes are sensitive. It is recommended to operate once every 2-3 days.

Major Eye Problems



Eyestrain



Dark circles



Periocular edema



Crow's feet



Fat granule



Eye bags

Am I a candidate ?



1. Yes, if you want to look younger



2. Yes, if you want to look better



3. Yes, if you want a quality life



4. Yes, if you want to more attractive



5. Yes, if you want his attention



6. Yes, if you want to look good even without

Equipment maintenance

Maintenance

1. Ensure the power supply is stable adaptation
2. The instrument should not be placed in damp places or water source accessories, and should not be directly exposed to sunlight
3. The instrument should not be too close to a strong source of heat or radiation
4. The instrument must use the original qualified parts provided by the manufacturer
5. After each operation, the function handle and function head shall be cleaned and disinfected and properly kept
6. In the case of not using the instrument, please turn off the power, remove the power, and place the instrument

Note

The instrument is a precision electronic product. It is forbidden to use the instrument in a humid environment, such as a bath room and other wet places. If you accidentally enter the water, please stop using it and return it to the factory for repair.

Warranty

1. The product is guaranteed for 3 months, and the company will be responsible for the maintenance of the products within the warranty period.
2. The starting date of the warranty period shall be subject to the date of the product invoice.
3. Please show the warranty card when repairing.
4. One of the following situations is not covered by the warranty.
 - Failure to follow the instructions.
 - Fault is caused by falling or collision.
 - Failure due to the customer's own modification and repair.
5. If the product has passed the warranty period or does not fall under free maintenance product, our maintenance service center will still be happy to serve you.