

USER MANUAL



NM-SN1S

Preface

Dear users:

We're pleased to present to you our NM-SN1S facial cleanser, which is a multifunctional cleaning and infusing beauty machine. By using advanced ultrasonic technology and fine hairbrush, we aim at a safe and, at the same time, effective result of cleaning pores, cleansing skin, and infusing nutrition. NM-SN1S facial cleanser is a beauty machine for professional use, which needs to be handled by professionals with proper training. Any improper use will result in adverse outcomes. Therefore we advise any personnel to read this guide thoroughly and follow the instructions strictly before operation.

We believe our quality product will bring in good returns on your end. Also, you can count on our worry-free after-sales service.

Thank you!

Contents

Part I

1. Brief Introduction	3
2. Advantages	3
3. Working Principles	4
4. Main Effects	4
5. Applicable Range	5
6. Inapplicable Range	5
7. Notes After Operation	5

Part II

1. Detailed Operation	6
2. How To Use	7
3. Technical Parameters	8
4. Packing List	8
5. Safety Precautions	9
6. Dos and Dont's	9
7. Troubleshooting	11
8. FAQs	11
9. Operational Diagrams	12

Part I

1. Brief Introduction

Facial cleanser is currently the most popular cleaning and infusing beauty machine that can effectively cleanse skin, remove blackheads and infuse nutrition to skin. The design of this beauty machine conforms to ergonomics that is comfortable, effort-saving and have clear effects. Facial cleanser is believed to replace the traditional washing with hands and its cleaning effect is 8 times cleaner than washing with hands. Moreover, it is convenient and easy to operate. It solves skin issues for beauty-lovers. This equipment is save and requires no injection, no medication, no operation, thus having no side-effects. It is used externally during the whole process and has an immediate effect. Based on the principle of ultrasonic vibration, facial cleanser uses hairbrush to deeply clean pores, remove clogs and make the metabolism of skin pores normal, thus reducing acne, blackheads, coarse pores and such skin problems.

2. Advantages

1. Facial cleanser is small, easy to carry and operate.
2. It is multifunctional to double cleanse skin and infuse nutrition to the skin.
3. It is easy to learn and operate with immediate effect, no side-effects.
4. The hairbrush is made of qualified materials and soft, thus having no hurt to skin and being eco-friendly.
5. Three modes can be switched freely according to your skin or feeling.
6. With top technology, solving 4 big facial problems, including cleaning& moisturizing, brightening, removing blackheads and reducing fine lines.
7. It is safer and more convenient with wireless charger.
8. Ultrasonic can infuse nutrition to the skin, at the same time, it will export the grease from pores, making moisture circulation in the pores, thus having a double effects to the skin.
9. It has special memorizing function with faster and more convenient operation.

3. Working Principles

3.1 Brush Cleansing

Based on the principle of high-frequency vibration of ultrasonic, vibrating the fine hair on the hairbrush, so as to achieve efficient and nonirritating cleaning and massage effect. It can help the face wash form high-density and even foam to thoroughly clean the face and achieve the beauty care.

1. Dry Skin Brush

Method: Clean the face skin with face wash.
It is suitable for dry skin and dry-mixed skin.

2. Oily Skin Brush

Method: Clean the face skin with face wash.
It is suitable for oily skin and oily-mixed skin.

3.2 Ultrasonic Lead In & Out

Based on the high-frequency vibration principle of ultrasonic, the high-frequency vibration reaches deeply into the 4~6 mm level of subcutaneous tissue. With its high energy and permeability, it vibrates tissues in every layers of your skin, exporting oil and dirt from pores. Also, it changes cells volumes so as to improve blood circulation and stimulates proliferation of cells, making it lustrous and elastic. Moreover, thanks to its high frequency vibration induced warm effect, it helps better absorption of nutrients into the skin at a higher speed.

4. Main Effects

1. Deep cleanse pores and booth skin absorption.
2. Reduce acne and blackheads.
3. Minimize enlarged pores and improve coarse skin.
4. Adjust skin sebum secretion and skin texture.
5. Deep channel in nutrients and help skin absorb nutrients better.
6. Improve uneven skin tone and dull and lusterless skin.
7. Stimulate cell regeneration and increase skin elasticity.

5. Applicable Range

1. Those with vigorous skin oil secretion, acne, blackheads and dull skin.
2. Those with oily skin and unbalanced water and oil production.
3. Those who don't clean for a long time and often wear makeup.
4. Those with poor skin absorption, metabolism and lack of nutrition.
5. Those with coarse pores, rough skin.
6. Those who love beauty and exquisite life.
7. Those with dull, sallow, uneven skin tone.
8. Those who work outside and bask in the sun.
9. Office workers who are in an airtight office and face their computers and mobile phones all day long.
10. Housewives who cook a lot.

6. Inapplicable Range

1. Those with severe heart disease.
2. Those with severe skin allergies, pustules and ulcers.
3. Those with skin diseases and infectious diseases, such as eczema and psoriasis.
4. Those with skin lesions, inflammation and other local signs of disease.
5. Those who have just had microneedling, intradermal injection, plastic surgery or implants.
6. It is prohibited for pregnant women and infants.

7. Notes After Operation




1. Keep skin hydrated and always use sun protection.
2. Apply at least 3 sheet masks a week.
3. Avoid using exfoliating products within a week.

Part II

1. Detailed Operation



2. How To Use

1. First charging should be 10 hours or above.
2. Insert the brush into the concave slot.
3. Wet skin and the brush with water;
4. Lather up the cleanser and apply on the skin or brush; Press for 2 seconds to turn it on. Default mode is cleansing mode. Click on/off button to switch energy level. Select the level that works for you and start cleansing from forehead (40sec), to nose and chin (40sec), and cheeks (each side for 20 sec). The device turns off automatically at 2 min. (Press the on/off button to turn off the device until it stops vibration);
5. 3 Levels of Cleaning
 -  Level 1 (blue) is applicable to sensitive skin
 -  Level 2 (purple) can be applied to cleanse pores and acne
 -  Level 3 (red) brings evident results of removing makeup residues
6. Press for 2 seconds to switch between cleansing and ultrasonic nutrition lead in and out and there's 3 levels of energy intensity;
7. The device is able to memorize the last setting;



Dry skin Brush

How To Use: Use with cleanser to clean facial skin. Each session lasts for 1-2 min and use 3-5 times per week. It applies to dry and combination skin.



Oily Skin Brush

How To Use: Use with cleanser to clean facial skin. Each session lasts for 1-2 min and use 5-7 times per week. It applies to oily and combination skin.

8. Red light stays on during charging and turns blue when it's fully charged;
9. Charging is needed when red light flashes;
10. Clean the device after treatment for next use.
11. Replace each brush every 3 months for better skin cleansing results.



Disinfection Light

Blue LED is used for sterilization of the brushes.



Clean the brush first and then place it on the stand to sterilize, best for 20-30 min.

3. Technical Parameters

Model:	NM-SN1S
Rated Voltage/Power (Main Device):	3.7V=4W
Waterproof Level:	IPX7
Charging Method:	Wireless
Rated Voltage/Power (Base):	5V=2.5W
Charging Time:	Around 8 Hours

4. Packing List

1 x Main Device	2 x Brush
1 x Power Cable	1 x Charging Base

5. Safety Precautions

1. Those with severe hypertension and heart disease.
2. Use carefully for pregnant women and patients with serious diseases. Please consult your doctor for advice.
3. Those with severe skin allergies, pustules and ulcers.
4. Those with skin diseases and infectious diseases, such as eczema and psoriasis.
5. Those with skin lesions, inflammation and other local signs of disease.
6. Do not operate around metal, plastic or silicone implants.

6. Dos and Dont's

1. Do not expose the instrument directly to the sunlight.
2. Pregnant women and sick patients should use this device with caution unless approved by a physician.
3. Use the device or train the operator of the device in strict accordance with the instructions in the manual.
4. Please do not disassemble the device by yourself. If the power cord and other non-replaceable parts are damaged, please contact customer service directly.
5. The device is waterproof and can be washed with water, but do not immerse the device in water for a long time.
6. When using this device , the operating part must be kept moist and avoid operation on dry skin .
7. Always use sunscreen and avoid direct sunlight after operation.
8. Do not use in vulnerable areas such as ears, eyeballs and throat.
9. During operation, the energy can be adjusted from low to high to avoid discomfort.
10. The operation should be close to the skin with even force to avoid skin discomfort.
11. Please wait after the repair phase of microneedling, intradermal injection, or plastic surgery to use the device.
12. Do not use the brush for more than 2 minutes when cleaning the face.

13. After each use of the brush, please put the brush down on the base for disinfection (the blue light on the base has the effect of disinfection and sterilization).
14. It is suggested that the operation time should be about 10 minutes when using this device for ultrasonic nutrition infusion. Do not place it on your eyeball during operation.
15. This device can be operated frequently.
16. Turn off the device after the operation, so as not to affect its service life.
17. After the operation of the device, wipe the device and handle with normal saline with tissue paper to effectively extend the service life (alcohol is not allowed).
18. Use in intervals for sensitive skin, for example, once every 2-3 days. A soft brush head is recommended.
19. If it has not been used for a long time (the battery of the product is out of power due to standby power consumption), when charging it again, the charging indicator does not show. Please wait for 6-10 minutes, and then it can be restored to normal charging. If the handle is heating up, please cut off the power immediately. Send it to after sales center for testing.
20. This product can only used with the brush head of this brand, or it will not bear all the consequences caused by using other non-official or certified brush head products.
21. This product is only used to clean the face. Do not use it for other purposes.
22. This product is for personal cleaning only. It is not recommended to use it for multiple users at the same time.
23. This product is only for indoor use. It should not be exposed to direct sunlight, and should be placed in a dry or ventilated place at room temperature.
24. If the battery has not been used for a long time, it should be charged once every three months to maintain the charging function of the battery.
25. When low power light flashes, it needs to be charged immediately.
26. It is recommended to remove the brush head every other week, rinse with warm water to eliminate the residue, put the protective cover upside down on the electrified charging base for blue light disinfection.

7. Troubleshooting

1. The device cannot be started?

- A. Make sure the battery is not drained, please charge it in time as needed.
B. If the above method cannot be handled, please contact the instrument dealer for assistance.

2. The device doesn't work with the cleansing function on?

- A. Make sure the device is in ultrasonic infusion function.
B. If the above method cannot be handled, please contact the instrument dealer for assistance.

3. Weakened vibration output?

- A. Too much pressure would affect it, please place the brush head gently on the face.
B. The gap between the brush head and the main body is too large or too small. Remove the brush head from the handle and then put it on the output shaft. Keep a certain distance between the brush head and the handle.

4. The device can only work for a short time?

- A. Make sure the charging base is connected properly to the main body. Keep the base clean and allow sufficient charging time for the instrument.
B. The battery may have reached the end of its service life.

8. FAQs

1. Q: Can I use this device every day?

- A: Yes. The operating head is made of high-quality materials, which is soft and does not hurt the skin. You can first soak the brush head before using. Brush gently on the skin, do not press too hard.

2. Q: How often do I change the brush?

- A: The brush wears out after using for some time under normal circumstances, and the ability to clean the oil and dirt of the skin will also be reduced. It is recommended to change the brush after 2~3 months for the best cleaning results.

3. Q: Can ultrasonic lead in function remove acne marks?

- A: Ultrasonic infusion is to lead skin brightening nutrients into the skin through ultrasonic high-frequency vibration, accelerate absorption, strengthen inhibition of melanin production and lighten spots, effectively improve dull skin, make the skin moisture and bright, and it has a special effect on pigmentation caused by acne significantly. No side effects, safe, reliable and easy to use. General outdoor activities are okay after treatment.

9. Operational Diagrams



STEP1

Remove makeup with makeup remover.



STEP2

Wet face. Lather up and apply foam evenly all over the face.



STEP3

Gently move the brush from nose the ear.



STEP4

Start from right cheek, then to the left cheek. Move on to forehead and nose.



STEP5

Rinse face clean with warm water.



STEP6

Apply sheet mask. Wait for 15 min.



STEP7

Wash face clean.



STEP8

Apply toner, emulsion, cream and sunscreen.